

Ovesco Research Update 38

November 2021 | Congress report UEG Week Virtual

OTSC® System

OTSC is superior and can be considered as the first-line treatment for larger peptic ulcers

J. Lau, Hong Kong, presented a poster on his multicentre RCT to compare OTSC to standard endoscopic treatment in patients with non-variceal upper gastrointestinal bleeding. Primary outcome was 30-day bleeding free probability. They finally recruited 190 patients and they were randomized into an OTSC arm (n = 93) and a standard arm (n = 97). The baseline characteristics of the patients were comparable with a mean age of 62.4 y vs 64.1 and most male patients (~80 %). The mean Glasgow Blatchford score was 10.9 and 11.0. 29.0% of patients in the OTSC arm and 35.0% in the standard arm were on anticoagulation or platelet aggregation inhibition. In both groups most bleedings occurred from a peptic ulcer (92.5% and 89.7%)

Analysis was by intention-to-treat. The 30-day bleeding free probability in OTSC group and standard treatment group was 96.8% ((90/93) vs 85.4% (83/97) respectively (HR, 0.21; 95% CI 0.06 to 0.73; P=0.006). Failure to control bleeding at first endoscopy (1 vs. 6; OR, 0.16, 95% CI 0.02 to 1.4) and 30-day recurrent bleeding (2 vs. 8; OR, 0.23, 95% CI 0.05 to 1.12) were lower in the OTSC group. In ulcers ≥ 10 mm in size, OTSC significantly reduced the rate of further bleeding (1 of 48 vs. 8 of 46; OR, 0.10, 95% CI 0.01 to 0.84). The 30-day re-intervention (2 vs. 8; p=0.138) and death did not differ significantly between the two groups (2 vs. 4; p=0.438). The cost to avert further bleeds was cheaper with use of OTSC. The ICER was – 647 USD.

Lau concluded that OTSC treatment can be considered as the first-line treatment for larger peptic ulcers.

OVER-THE-SCOPE CLIPS (OTSC) VERSUS STANDARD ENDOSCOPIC TREATMENT IN PATIENTS WITH ACUTE NON-VARICEAL UPPER GASTROINTESTINAL BLEEDING

Lau J, Tan C.H., Sun X, Song H, Li L, Li R, Li P, Feng J, Wang B, Leung W.K., Hartley I, Moss A.C., Suen B.Y., Yu Y, Chan F.K.L

STING II: OTSC first-line therapy is superior to standard therapy in acute nonvariceal**OGI bleeding with high risk of rebleeding**

B. Meier, Ludwigsburg, Germany presented the final data on the multicentre RCT STING II on the use of the OTSC System in severe, non-variceal upper GI bleeding. In the treatment of recurrent bleeding, OTSC has become established due to the significantly lower rebleeding rates compared to standard therapy (STING study). STING II was subsequently initiated to investigate whether OTSC treatment is superior to standard treatment also in acute, non-variceal upper GI bleeding with a high risk of rebleeding ("complete Rockall score" ≥ 7 points). Standard treatment consisted of at least 2 conventional clips or thermal procedure plus adrenaline injection. OTSC treatment consisted of primary OTSC application (injection allowed). According to the calculation of the number of cases (power 80%), n = 100 patients with endoscopically confirmed acute (non-variceal) upper GI bleeding and a "complete Rockall score" ≥ 7 points were included. Exclusion criteria included prior endoscopic

treatment within the last 4 weeks to rule out treatment of recurrent bleeding. Primary endpoint was successful haemostasis without recurrent bleeding within 7 days.

After randomisation, n = 52 patients could be evaluated in the standard arm and n = 48 patients in the OTSC arm. Both groups show good correspondence to the baseline. 42.3% of patients in the standard arm and 39.6% in the OTSC arm were on anticoagulation or platelet aggregation inhibition. The median age of patients was 79 years (51–96) in the standard arm and 78 years (42–92) in the OTSC arm. The median Rockall score was 8 points in both groups, and the predominant localisation was in the duodenum (46.1% standard vs. 60.4% OTSC). Peptic ulcer bleeding was predominant in both arms, with more than 80% in each. Approximately 60% of each bleed was active (Forrest Ia/Ib). Analysis of the data showed that the OTSC system was significantly superior to standard therapy at 91.7% (44/48) with 73.1% (38/52) in successful hemostasis (p = 0.019).

There was no persistent bleeding with the OTSC compared with 6 (11.5%) persistent bleeding in the standard arm (p = 0.027). All persistent bleeding in the standard arm could be successfully stopped with OTSC. In the OTSC arm, 8.3% (4/48) of rebleeds occurred within 7 days compared with 15.4% (8/52) in the standard arm (p = 0.362). Treatment time was comparable at 27 min (OTSC) and 28 min, respectively. There were no significant differences in other endpoints (blood transfusion, ICU stay, length of hospitalisation, mortality and need for second-line treatment).

Meier concluded that OTSC is superior to standard therapy for acute nonvariceal upper GI bleeding with a high risk of rebleeding. No additional time was required for OTSC therapy. Furthermore, the study once again highlighted the efficacy of OTSC salvage therapy after unsuccessful hemostasis.

Endoscopic Treatment of Non-variceal Upper GI-bleeding With High Risk of Recurrency - OTSC (Over-the-scope-clip) Versus Standard Therapy (STING2)

Meier B, Wannhoff A, Denzer U, Stathopoulos P, Schumacher B, Albers D, Hoffmeister A, Feisthammel J, Walter B, Meining A, Wedi E, Zachäus M, Pickartz T, Küllmer A, Schmidt A, Caca K

Retrospective study compares efficacy and complications of OTSC type-a vs. type-t clips in duodenal ulcer bleeding

M. Hollenbach, Leipzig, Germany presented a poster on a retrospective analysis comparing OTSC type-a versus type-t clips in the treatment of duodenal ulcer bleeding. For this purpose, data from 2009 – 2020 from 6 endoscopy centres was analysed and all patients treated with the OTSC System for duodenal ulcer bleeding during an emergency endoscopy were included in the analysis.

Finally, the data of 173 patients could be evaluated. Of these, 93 patients were treated with type-a and 80 patients with type-t. The analysis showed that the baseline characteristics age (71.2 y vs. 71.6 y, p = 0.255), gender (men: 69.9% vs. 67.5%, p = 0.735), use of anticoagulants (32.9% vs. 43.0%, p = 0.176) and Rockall score (7.2 vs. 7.4, p = 0.917) were comparable between the groups. However, there were some significant differences in the bleeding characteristics. Thus, the type-a group had significantly fewer ulcers with active bleeding (Forrest Ia/b) than the type-t group (51.1% vs. 62.5%, p=0.020). In addition, type-t was used more often as first-line treatment (95% vs. 77.8%, p=0.004). Analysis of the data showed that initial haemostasis (type-a: 93.5%, type-t: 90%, p=0.421) and bleeding-associated lethality were not significantly different (type-a: 3.2%, type-t: 7.8%, p=0.125). However, the OTSC-t group showed a significantly increase rate of rebleeding (37.2% vs. 6.5%, p<0.001).

Hollenbach concluded that therefore the type-a should be considered the standard of care for endoscopic treatment of duodenal ulcer bleeding, when an OTSC is used, because of the lower rate of rebleeding.

Treatment of duodenal ulcer bleeding – Efficacy of traumatic and atraumatic over-the-scope-clips (OTSC) – a multicenter retrospective analysis

M. Hollenbach, A. Schmidt, A. Decker, O. Möschler, C. Jung, N.-C. Mechie, T. Barhoom, A. Hegelein, R. Knoop, T. Blasberg, E. Wedi

Presentation of recently updated non-variceal upper GI bleeding (ESGE update 2021)

I. Gralnek, Haifa, Israel gave a talk on the recently updated ESGE Guideline – Endoscopic diagnosis and management of nonvariceal upper gastrointestinal hemorrhage. Besides several other changes regarding e.g., the pre-management with PPI, management of antiplatelet agents and the ideal timing of endoscopy Gralnek also spoke about the updated recommendations for the endoscopic therapy. In addition to the previous standard therapy with Endoclips now also the OTSC is recommended for the mechanical treatment of ulcer bleeding. Furthermore, the significance of the OTSC for the rescue therapy was highlighted. During the following Q&A session a variety of questions regarding OTSC were discussed.

Upper GI Bleeding: Whats' New in 2021? "Non-Variceal Bleeding"

Gralnek I

A. Meining, Würzburg, Germany, presented a case report of a 74-year-old patient with a history of gastrectomy with Roux-en-Y reconstruction for gastric carcinoma, right hemicolectomy for a tumor in the ascending colon that occurred several years later, and HIPEC therapy. Before being referred to the university hospital, the patient presented to a smaller center with acute pleural empyema on the left side with E. faecium and was treated with drainage. With drainage, the empyema disappeared. However, small air bubbles remained near the drainage, and there was persistent fluid formation with E. faecium. Gastroscopy with fluoroscopy failed to find leakage at the esophago-jejunal anastomosis. A colonoscopy was then performed, and a fairly large opening was found in addition to several diverticula. X-ray fluoroscopy revealed that a pleural-colonic fistula was present.

This unusual case was then closed with an optimized OTSC fistula closure technique developed by Meining. In this technique, he uses a needle knife to make an incision of the mucosal layer around the fistula opening. The size of the incision is approximately the same as the size of the OTSC cap, and the clip is placed in the incision. After placement of the clip, closure of the leak was demonstrated by injection of a contrast agent.

No leakage occurred later in the course, and the drainage could be extracted. The clip detached after 6 months and the last follow-up after 1 year showed no further complaints after the procedure.

Fistula from the colon to the pleura after surgery

Meining A

November 2021 | Congress report
Viszeralmedizin (Visceral medicine)

Joint Annual Meeting of the German Society for Gastroenterology, Digestive and Metabolic Diseases with Endoscopy section (DGVS) together with the German Society of General Surgery and Visceral Surgery (DGAV)

OTSC® System

Sting II: OTSC first-line treatment is superior to standard treatment for acute non-varicose upper GI bleeding with high risk of rebleeding.

B. Meier, Ludwigsburg, presented the final data of the multicentre randomised controlled STING II study on the use of the OTSC system in severe, non-varicose upper GI bleeding. In the treatment of recurrent bleeding, OTSC has become established due to the significantly lower rebleeding rates compared to standard treatment (STING study). STING II was consequently initiated to investigate whether OTSC treatment is superior to standard treatment also in acute, non-varicose upper GI bleeding with a high risk of rebleeding ("complete Rockall score" ≥ 7 points). Standard treatment consisted of at least 2 conventional clips or thermal procedure plus adrenaline injection. OTSC treatment consisted of primary OTSC application (injection allowed). The study included n = 100 patients with endoscopically proven acute (non-varicose) OGI haemorrhage and a "complete Rockall Score" ≥ 7 points according to case number calculation (power 80 %). Exclusion criteria included prior endoscopic treatment within the last 4 weeks to rule out treatment of recurrent bleeding. The primary endpoint was successful haemostasis without recurrent bleeding within 7 days.

After randomisation, n = 52 patients were evaluated in the standard arm and n = 48 patients in the OTSC arm. Both groups show good correspondence to the baseline. 42.3% of patients in the standard arm and 39.6% in the OTSC arm were on anticoagulation or platelet aggregation inhibition. The median age of patients was 79 years (51–96) in the standard arm and 78 years (42–92) in the OTSC arm. The median Rockall score was 8 points in both groups, and the predominant localisation was in the duodenum (46.1% standard vs. 60.4% OTSC). Predominantly peptic ulcer bleeding was present with over 80% in both arms. In each case approx. 60% of the bleedings were active (Forrest Ia/Ib). Analysis of the data shows that the OTSC System is significantly superior (p = 0.019) to standard treatment at 91.7% (44/48) to 73.1% (38/52) in successful haemostasis. There was no persistent bleeding with the OTSC compared with 6 (11.5%) persistent bleeding in the standard arm (p =

0.027). All persistent bleeding in the standard arm could be successfully stopped with OTSC. In the OTSC, 8.3% (4/48) of rebleeds occurred within 7 days compared with 15.4% (8/52) in the standard arm ($p = 0.362$). The treatment time was comparable at 27 min (OTSC) and 28 min. There were no significant differences with regard to the other endpoints (blood transfusion, stay in intensive care, length of hospitalisation, mortality and need for second-line treatment). Meier concluded that OTSC is superior to standard treatment for acute non-varicose upper GI bleeding with a high risk of rebleeding and that primary haemostasis by OTSC is highly effective (no persistent bleeding in the OTSC arm). No additional time is required for OTSC treatment. Furthermore, the subgroups that benefit from OTSC treatment still need to be better defined.

A 319 Die primäre OTSC-Therapie ist der endoskopischen Standardtherapie bei Hochrisikopatienten mit akuter oberer (nicht variköser) Gastrointestinalblutung überlegen (STING-2 Studie)

Meier B1, Wannhoff A1 1, Denzer U2, Schumacher B3, Albers D3, Feisthammel J4, Walter B5, Meining A6, Wedi E7, Zachäus M8, Pickartz T9, Küllmer A10, Schmidt A18, Caca K1

1Klinikum Ludwigsburg, 2Universitätsklinik Marburg, 3Elisabeth-Krankenhaus Essen 4Universitätsklinik Leipzig, 5Universitätsklinik Ulm, 6Universitätsklinik Würzburg, 7Sana Klinikum Offenbach, 8Helios Park-Klinikum Leipzig, 9Universitätsklinik Greifswald, 10Universitätsklinik Freiburg

Retrospective comparative study on OTSC type-a vs. type-t clips in duodenal ulcer bleeding

M. Hollenbach, Leipzig, reported on a retrospective analysis comparing OTSC type-a versus type-t clips in the treatment of duodenal ulcer bleeding. For this purpose, data from 2009 – 2020 from 6 endoscopy centres was analysed and all patients treated with the OTSC System for duodenal ulcer bleeding during an emergency endoscopy were included in the analysis.

Data from 173 patients has been evaluated so far. Data from a further 80 – 100 patients is to be included in the analysis after acquiring a centre at short notice. Of the patients evaluated so far, 93 were treated with type-a and 80 with type-t. The analysis showed that the baseline characteristics age (71.2 y vs. 71.6 y, $p = 0.255$), gender (men: 69.9% vs. 67.5%, $p = 0.735$), use of anticoagulants (32.9% vs. 43.0%, $p = 0.176$) and Rockall score (7.2 vs. 7.4, $p = 0.917$) were comparable between the groups. However, there were significant differences in the bleeding characteristics. The type-a group had significantly fewer ulcers with active bleeding (Forrest Ia/b) than the type-t group (51.1% vs. 62.5%, $p=0.020$). In addition, type-t was used more often as first-line treatment (95% vs. 77.8%, $p=0.004$). Analysis of the data to date showed that initial haemostasis (type-a: 93.5%, type-t: 90%, $p=0.421$) was comparable and bleeding-associated lethality was not significantly different (type-a: 3.2%, type-t: 7.8%, $p=0.125$). However, the OTSC-t group showed a significantly increased rate of rebleeding (37.2% vs. 6.5%, $p<0.001$).

Hollenbach concluded that the OTSC is highly effective in the haemostasis of bleeding duodenal ulcers in both first- and second-line treatment. Type-a should be considered the standard of care for endoscopic treatment of duodenal ulcer bleeding, if an OTSC is used, because of the significantly lower rate of rebleeding. However, it is important to wait for the other included cases before making a final assessment and to note the differences in the bleeding characteristics.

Analyse des traumatischen vs. atraumatischen over-the-scope-clip (OTSC) bei der Behandlung von duodenalen Ulkusblutungen

Hollenbach M1, Schmidt A2, Decker A2, Möschler O3, Jung C4, Mechie N.-C 5, Barhoorn T5, Hegelein A5, Knoop R6, Blasberg T5, Wedi E5

1Universitätsklinikum Leipzig, 2Universitätsklinikum Freiburg, 3Marienhospital Osnabrück, 4AUSL Romagna, 5Sana Klinikum Offenbach, 6Universitätsklinikum Göttingen

stentfix OTSC® System

First evaluations of the University Medical Centre Freiburg and the University Medical Centre Würzburg

show that the migration risk can be significantly reduced with stentfix OTSC

Lecture University Medical Centre Würzburg

K. Groneberg, Würzburg, presented initial data on the effectiveness of the stentfix OTSC system. The stentfix OTSC system has been examined to prevent stent migration in fully covered self-expanding metal stents (FC-SEMS). These are used in the endoscopic management of stenosis or leakage in the upper gastrointestinal tract and their main limitation is migration, which occurs in one third to one half of all patients, depending on the indication.

The data of 25 patients was retrospectively considered for the evaluation. In these patients, an FC-SEMS was fixed using stentfix OTSC within a period of 24 months following the introduction of the system in April 2019. The general feasibility and effectiveness of this new method were evaluated.

In 13 patients, stent implantation was due to stenosis (malignant = 4, benign = 9). In the other 12 cases, the stent was implanted to treat a fistula or leakage. The average duration of the stent implantation including clip application was 6.5 minutes and in regards to the clinical problem, a satisfactory result was achieved in all patients. In the follow-up examination, stent dislocation was detected in only 2 patients (8%). In 5 other patients (20%), the stentfix OTSC had detached from the mucosa without migration of the FC-SEMS. Fixation-associated side effects were not observed.

Groneberg concluded that based on this initial data, fixation of FC-SEMS in the oesophagus using stentfix OTSC is safe and easy to perform. In the short term, stent migration can be prevented in the vast majority of cases.

A 293 Effektivität der Fixierung von SEMS im Ösophagus mittels Stentfix OTSC

Groneberg K, Reimer S, Meining A, Universitätsklinikum Würzburg

Lecture University Medical Centre Freiburg

A first evaluation of the effectiveness of the stentfix OTSC was also presented by the University Medical Centre Freiburg. M. Schirmer investigated whether the use of the stentfix OTSC system significantly reduces the migration rate of oesophageal SEMS and performed a retrospective comparative monocentric analysis. The primary endpoint was stent migration. Univariable and multivariable logistic regression models were used to evaluate predictive factors for stent migration.

The intervention group consisted of 26 patients who underwent oesophageal SEMS insertion using the stentfix OTSC system between 02/2019 to 12/2020. As a control group, 51 cases of oesophageal SEMS inserts without fixation were analysed between 04/2017 and 01/2019. The main indications in the stentfix OTSC group were malignant stenosis (57.7%) and perforation/insufficiency (32%) which were treated with 50% partial and 50% full coverage stents. The application of the stentfix OTSC was technically successful in all cases and was accomplished without complications. In the control group without fixation, 71% had partially covered stents. The indication was malignant in 74.5% and benign in 25.5%. There was no difference in the other characteristics.

In both groups, data on the presence of stent migration could be collected in 92% of cases. This occurred significantly less frequently in the stentfix OTSC group (7.7%) than in the control group (33.3%, relative risk reduction (RRR): 76.9%, $p<0.001$). Univariable analyses showed that overstretching of the cardia tended to be associated with a higher risk of stent migration. The use of the stentfix OTSC led to a significant reduction in the risk of migration in this subgroup (8.3 vs. 50%, RRR: 83.4% $p = 0.024$). Stent removal after stentfix OTSC was performed in 11 patients (42%) without complications. The presence of malignant disease or the type of stent coating had no influence on stent migration.

Schirmer concluded that the use of the stentfix OTSC system can significantly reduce the risk of stent migration. This is especially true in patients with overstretching of the cardia.

A 294 Effektive Reduktion der Migrationsrate von Ösophagus-Stents durch ein neues Clipsystem (stentfix OTSC®)

Schiemer M, Schmidt A, Bettinger D, Mueller J, Schultheiß M,

Thimme R, Kuellmer A; Universitätsklinikum Freiburg

OTSG Xcavator™

A new and effective instrument for pancreatic necrosectomy, blood coagel and foreign body removal

M. Brand presented the first standardised evaluation of procedures with the new OTSG Xcavator. Between November 2020 and April 2021, 34 procedures with the OTSG Xcavator were recorded from 7 centres. Of the 34 procedures, n=24 were pancreatic necrosectomies (19 transgastric, 3 transduodenal).

In 67% (n=16) of the necrosectomies, access to the necrotic cavity was via a lumen apposing metal stent (LAMS). In one case the OTSG Xcavator could not be inserted into the cavity (transduodenal approach), in two other cases there was a dislocation of the LAMS. The technical and clinical success of necrosectomy was 97%, with a mean of 7 (2-19) pieces of necrosis removed. The mean procedure time was 47 min (15-120 min). Suction/irrigation was used to remove the necrotic material, and in some cases removal of the endoscope was necessary to clean the instrument.

Furthermore, n=4 cleanings of insufficiency cavities prior to Endovac treatment, n=4 removals of foreign bodies in oesophagus or stomach and n=2 removals of blood clots in case of ulcer bleeding in the duodenum were performed. Here, the technical and clinical success rate was 100%. Clinically relevant complications were not reported in any of the 34 cases.

Brand concluded that the OTSG Xcavator is a simple and safe tool for pancreatic necrosectomy. It is also useful and effective for foreign body removal and blood coagulation management in Forrest IIb situations.

OTSG Xcavator - Erste multizentrische Daten zum klinischen Einsatz

Brand M1, Bachmann J2, Hügler U3, Rahman I4, Schlag C2, Wedi E5, Braun G6, Möschler O7, Meining A1

1Universitätsklinikum Würzburg, 2Klinikum rechts der Isar der TUM, 3Kliniken der Stadt Köln - Holweide, 4University Hospital Southampton, 5Sana Klinikum Offenbach, 6Universitätsklinikum Augsburg, 7Marienhospital Osnabrück

FTRD® System

Report on STER and FTRD: procedure, areas of application and study data

In his lecture, J. Hochberger discussed STER (Submucosal Tunneling Endoscopic Resection) as a so-called exposure technique and presented an overview of the available study data. Two patient cases were displayed, a leiomyoma which could not be removed by atraumatic cranial salvage due to its tumour size. A second case showed a patient with a tumour originating from the propria muscle in the fundus. In both cases, the STER technique was successfully applied and incision closure was performed using the OTSC system. J. Hochberger concluded that the STER technique is still a limited technique due to the scarce and retrospective data, but could be a real alternative to surgery in the oesophagus.

In the second part of the lecture, endoscopic full-thickness resection with the FTRD System was presented as a so-called non-exposure technique, since the wall closure takes place before the resection. A review of data on the use of the FTRD System in the colon was shown, such as the WALL RESECT study with n=181 patients and the FTRD register with n=1178 patients. Both studies showed an R0 resection rate of about 80%. The presented Dutch FTRD register (n=367 patients) also achieved a similar R0 resection rate. When using the technique for early colon carcinomas, J. Hochberger concluded that eFTR is technically feasible and safe and the FTRD System allows an accurate histological classification of the tumour.

When the FTRD System is used for full-thickness resection of a tumour at the opening of the appendix, various past studies have shown an R0 resection rate between 64 and 85%. The appendicitis rate was between 11 and 17%. Another combination technique is hybrid FTRD, in which the lifting part is ablated using piecemeal EMR and the non-lifting part is resected using the FTRD System. An American abstract published in 2020 (Yuen et al.) with n=62 patients was able to achieve an R0 resection rate of 96% in this way. J. Hochberger concluded that full-thickness resection with the

FTRD System is an integral part of the colon and rectum repertoire, supported by clinical data. It also allows resections in difficult locations and in scar tissue. For use in early carcinomas, prospective studies would be desirable due to the limited amount of data available so far. The application possibilities in the stomach and small intestine already display very promising data.

STER und FTRD (Tumorsektion durch Tunnel oder Vollwand)

Hochberger J, Klinikum Berlin Friedrichshain

Successful treatment of complex early colorectal carcinomas with endoscopic full-thickness resection

A. Hoffmann presented a small case series (n=4) of patients with complex early colorectal carcinomas in whom the procedure of endoscopic transmural full-thickness resection was used.

The first case was a 75-year-old patient with an adenoma of approx. 3 cm in size on the right flexure, which could only be partially removed using the piecemeal technique (histology: severe dysplasia with focal pT1sm L0V0G1Rx early carcinoma). Despite scarring and difficult localisation, the recurrence was then successfully removed (histology: pT1sm1 L0V0G1R0).

In another case, a 69-year-old patient with a bioptic early carcinoma and a second suspicious polyp in the immediate vicinity was presented, in whom two endoscopic transmural full-thickness resections were successfully performed in parallel (histology: pTis L0V0G1/G2R0).

The remaining two cases were reports of successful endoscopic full-thickness resection for long-standing ulcerative colitis. During a routine outpatient colonoscopy, a 55-year-old female patient was found to have a conspicuous area that was difficult to define but clearly visible under virtual chromoendoscopy. The histology of the full-thickness resection showed colitis associated dysplasia with an early carcinoma pTis L0V0G2R0. In another case, a 56-year-old patient was found to have a biopsy-confirmed mucosal carcinoma in the rectum. In this case, full-thickness resection revealed a high-grade dysplasia that could be completely removed.

The authors conclude that endoscopic full-thickness resection with the FTRD can bypass surgery and that the resected tissue has good histological significance, especially in the case of early carcinomas.

A 338 Endoskopische Vollwandresektion als erfolgreiche Therapieoption bei komplexen kolorektalen Frühkarzinomen

Dorlöchter C, Hirsch A, Schmidt V, Conrad A, Hoffman A; Lehrkrankenhaus Klinikum Aschaffenburg Alzenau,

eFTR is becoming increasingly important in the resection of early colorectal carcinomas with "low risk" histology

P. Stathopoulos presented a unicentric, retrospective analysis of the effectiveness, safety and oncological follow-up of patients with early colorectal carcinomas treated with eFTR. From November 2016 to April 2021, 89 underwent eFTR, of which 31 patients (12 women, median age 76) had early CRC. All lesions (22 colon (71.0%), 9 rectum (29%)) were reached. The median histological lesion size was 14 mm (3-25 mm), median procedure time was 47.5 min (25-70 min). The technical success rate was 83.9% and R0 resection was achieved in 67.7% (n = 21). In 9/31 patients, eFTR was performed after questionably incomplete polypectomy of early adenocarcinoma, but no carcinoma remnants were seen in full-thickness resections. 4/31 patients had a "low-risk" situation, of which one patient underwent lateral oncological resection due to R1 status. In the remaining 18/31 patients, histology revealed a "high-risk" situation. Of these, 8 patients underwent oncological resection and 10 patients underwent clinical endoscopic follow-up. During follow-up (median 15.5 months, 0-47), one patient died from another cause, and no tumour recurrence was documented in the remaining patients. There were 3 minor rebleeds (9.4%); one ileus (3.1%), which was treated conservatively, and one perforation (3.1%), which was closed endoscopically with OTSC. Stathopoulos concluded that eFTR is a safe endoscopic procedure that allows accurate risk stratification and can spare patients with early CRC and "low-risk"

histology from oncological surgery. However, prospective studies with appropriate oncological follow-up are needed to evaluate the long-term effectiveness.

A 339 Endoskopische Vollwandresektion („endoscopic full-thickness resection“, eFTR) für das frühe kolorektale Karzinom (CRC) – eine retrospektive Analyse von 31 konsekutiven Fällen

Stathopoulos P, Zumblick M, Hoque S, Gress TM, Denkert C, Denzer UW; Universitätsklinik Marburg

Endoscopic clips including OTSC are highly effective standard treatment for defect closure

In his study, K. Kouladourous, Mannheim, presented the clip closure of bowel wall lesions in endoscopic mucosal resections (EMR). For this purpose, all EMR cases in the lower gastrointestinal tract, for the period 2007 – 2020 were retrospectively analysed and all patients with intrainterventional bowel wall injuries requiring treatment were identified. Periinterventional factors, including type and number of clips applied, and outcomes of treatment were assessed.

157 intrainterventional bowel wall injuries were identified from 3782 endoscopic mucosal resections. Endoscopic clip closure was possible in 148 cases (94%). 11 of these initially endoscopically treated patients required surgical resection, making the clinical success rate of endoscopic clip closure 92.6%. The size of the resected tumour and the defect, as well as the number and type of clips applied (through-the-scope clips or OTSC) had no influence on the success of endoscopic clip closure.

Kouladourous concluded that endoscopic clip closure of bowel wall lesions during endoscopic resection is a reliable treatment option and eliminates the need for surgery in more than 90% of cases. That's why this is, and should always be, the initial treatment of choice.

A 297 Resektionsbedingte Darmwandverletzungen: ist ein Loch gleich ein Loch? Eine retrospektive Analyse von 3782 endoskopischen Mukosaresektionen

Kouladourous K, Belle S, Kähler G; Universitätsklinikum Mannheim

Evaluation of GKV routine data underlines the position of endoscopic treatment for ulcer bleeding

D. Horenkamp-Sonntag from Techniker Krankenkasse, Hamburg, presented an evaluation of acute ulcer bleeding on the basis of GKV routine data from the period 2016-2020. These are available for 10.8 million TK insureds. The aim of the study was to evaluate the role of anticoagulant drugs (anticoagulation [AK] and acetylsalicylic acid [ASA]) and endoscopic haemostasis measures with clip insertion in acute ulcer bleeding in routine clinical practice. The data was tailored by using specific disease information (ICD K25-28, Z92) and endoscopic services (OPS codes) to focus on insured patients with ulcer disease of the stomach and/or duodenum.

The study was based on hospital cases with acute peptic ulcer without bleeding (n = 19,167) or with acute bleeding (n = 17,635), while 41,241 outpatient cases with ulcers in the upper GI tract were not included. The main findings of the interim analysis were:

1. Patients with ulcer bleeding were older than those without bleeding
2. Treatment with ASA showed only minor differences between the groups, while AK or the combination of ASA and AK was significantly more frequent in the bleeding population.
3. Clips were used in about 1/3 of the bleedings, with OTSC clips (in 7-12%) mainly in duodenal ulcers
4. Clips were used significantly (p<0.001) more often in patients with AK ± ASA than in those without coagulation, regardless of ulcer location. Based on hospital patient billing data, anticoagulation is a risk factor for acute ulcer bleeding in every 3rd patient. ASA use may be relevant for ulcer pathogenesis, but less so for bleeding. Clips are an important tool in haemostasis, especially under anticoagulants. The OTSC clip is currently preferred for duodenal ulcers.

A 317 Akute Ulkusblutung: Rolle der Gerinnungshemmung und Einsatz von Clips in der endoskopischen Blutstillung

Horenkamp-Sonntag D1, Koop H2, Skupnik C1

1Techniker Krankenkasse, Hamburg; 2Helios Klinikum Berlin-Buch

October 2021 | OTSC®: Systematic review confirms 85.2 % clinical success rate and 1 % adverse events in the management of iatrogenic GI perforations. The size of the defect ranged from 2 to 50 mm.

A systematic review of studies evaluating OTSC application for repair of iatrogenic perforations identified 12 studies with overall 191 patients. Pooled proportion analysis yielded a technical success rate of 89.1 %, and a clinical success rate of 85.2 %. The size of the defect was available in 7 studies and ranged from 2 to 50 mm. Complications related to OTSC application occurred in two patients (1 %).

C. Zhong et al., Affiliated Hospital of Southwest Medical University, Luzhou, China performed a systematic review aiming to evaluate the clinical safety and efficacy of the OTSC System for the management of iatrogenic GI defects. The use of endoscopic procedures such as endoscopic mucosal resection (EMR), endoscopic submucosal dissection (ESD), endoscopic ultrasound (EUS), endoscopic retrograde cholangiopancreatography (ERCP), and endoscopic full-thickness resection (EFTR) continues to increase worldwide and iatrogenic gastrointestinal (GI) tract perforation is a rare but severe complication of these procedures. The reported incidence in therapeutic endoscopy is between 0.6 % and 5.5 %. The application of the OTSC System is a promising minimally invasive approach for the management of iatrogenic GI defects avoiding surgical intervention.

To perform the systematic review, PubMed, Embase and Cochrane library were searched for eligible articles published between 01/2006 and 12/2018. Two independent reviewers selected the literature according to inclusion and exclusion criteria. Only studies reported in English and series with more than 3 cases were included. Studies describing chronic defects due to percutaneous endoscopic gastrostomy (PEG) extraction were excluded. The statistical analysis was conducted using Comprehensive Meta-Analysis software version 3.0. The overall success rate was estimated based on the pooled proportion with 95 % confidence intervals (CI). Statistical heterogeneity among studies was assessed with the I2 statistic, where an I2 value > 50 % was considered significantly heterogeneous.

Overall, 12 studies comprising 191 patients with iatrogenic GI defects were included. The major causes for iatrogenic GI defects were ESD (41.4 %) and EMR (16.2 %), while lesser frequent causes included EFTR (12 %), diagnostic endoscopy (7.9 %), EUS (6.8 %), endoscopic polypectomy (3.7 %), PEG placement or insertion complications (3.7 %), ERCP (2.1 %), and other causes (6.3 %). Localization of defects was in the esophagus (7.3 %), stomach (24.1 %), duodenum (38.2 %), jejunum (1.0 %), cecum (1.6 %), colon (17.2 %), and rectum (10.5 %). The size of the defect was available in 7 studies and ranged from 2 to 50 mm. Pooled proportion analysis yielded a technical success rate of 89.1 % (n = 182 patients; 95 % CI 81.6 % - 93.8 %, I2 = 41.06 %), and a clinical success rate of 85.2 % (n = 170 patients; 95 % CI 71.9 % - 92.8 %, I2 = 58.92 %). Complications related to OTSC application occurred in two patients (1 %). None of the included studies reported severe or fatal adverse events.

The authors concluded that endoscopic closure of iatrogenic GI defects with the OTSC System is a safe and effective approach. Further randomized controlled trials are necessary to compare the OTSC System to other treatment modalities.

Endoscopic management of iatrogenic gastrointestinal defects with the Over-The-Scope Clip (OTSC) system: an updated systematic review

Zhong C, Tan S, Ren Y, Luo X, Xu J, Fu X, Peng Y, Tang X.

October 2021 | OTSC® vs surgery in scope-induced duodenal perforations: Retrospective comparison study shows superiority of the OTSC regarding morbidity and mortality

The retrospective analysis included overall 20 patients with type I duodenal perforation. 8 patients were treated with OTSC and 12 by surgery. Baseline parameters and

perforation size were similar in both groups. 2 patients in the surgery group died vs none in the OTSC group ($p=0.48$). Median hospital stay was 2 days in the OTSC group vs. 22 days in the surgery group ($p=0.003$).

Dahale A. S. et al., Dr. D Y Patil Medical College and Hospital, Pimpri, Pune, India, conducted a study comparing OTSC and surgery in the management of scope-induced duodenal perforations. Retrospectively, data of 20 patients was analyzed. All included perforations had been detected and treated within 24 hours after the endoscopic procedure. Objectives of the study were spectrum, etiology, baseline parameters, perforation size, outcome, comorbidities, and duration of hospital stay. Of the 20 patients, 8 were treated with OTSC placement and 12 underwent surgery. Patient age in both groups was comparable. The majority of patients was female in both groups. Baseline parameters, etiology and comorbidities were similar in both groups. The median size of perforation was 1.5 cm in both, the OTSC and the surgery group.

All patients were treated with standard of care according to institutional protocols. Patients in the OTSC group started oral intake after 48 h of OTSC placement, while in the surgery group median time was 7 days. Median hospital stay was 2 days in the OTSC group vs 22 days in the surgery group ($p = 0.003$). Intravenous antibiotics were administered for 2 days in the OTSC group followed by oral antibiotics for 2 weeks, while in the surgery group, median duration of intravenous antibiotics was 19.5 days (range 7 – 45 days). There were no complications related to the OTSC. Seven patients from the surgery group had one or more post-surgical complications. Two patients had leak from the wound site and were managed conservatively. Two patients had delirium and spontaneous recovery within 2 days, one of whom also had left lower foot gangrene as thromboembolic phenomenon. One patient had local wound site infection, which was treated with wound care and antibiotics. Two patients in the surgical group died from sepsis and multiorgan failure, while there was no mortality in the OTSC group ($p = 0.48$).

The authors concluded, that OTSC is a feasible and better option in type I duodenal perforations and associated with a significantly shorter hospital stay, shorter application of intravenous antibiotics, less post-closure complications, earlier start of oral intake and lower mortality rate.

Management of scope-induced type I duodenal perforations: Over-the-scope clip versus surgery

Dahale AS, Srivastava S, Saluja SS, Sachdeva S, Dalal A, Varakanahalli S.

September 2021 | Hybrid EMR-FTRD®: In large colorectal lesions with non-lifting areas this hybrid-technique can close the therapeutic gap and avoid surgery

A retrospective cohort study of 17 patients undergoing hybrid EMR-FTRD for large colorectal lesions showed 100 % technical success and 76.4 % histologically confirmed R0-resection. No perforation or major bleeding occurred. Follow-up endoscopy showed recurrent adenomas in 2 patients.

M. Bauermeister et al., Sana Klinikum Lichtenberg, Berlin, Germany performed a retrospective evaluation of 17 patients undergoing hybrid EMR-FTRD for large colorectal lesions with non-lifting signs. Endpoints of the study were technical success, histological confirmation of R0-resection and adverse events.

All seventeen lesions were too large for en-bloc resection with the FTRD and all had positive lateral lifting signs. Therefore, a hybrid approach combining piecemeal EMR with FTRD was performed. After submucosal injection of saline and toluidine blue, the lateral parts of the adenoma were resected with piecemeal EMR. The non-lifting area of the adenoma was not touched. Afterwards, the FTRD was used for en-bloc resection of the remaining lesion. The specimen was retrieved and examined histologically. Endoscopy was finished by a final check of macroscopic success and exclusion of bleeding and perforation. A follow-up endoscopy was scheduled after 12 weeks for inspection including taking biopsies at the resection site.

Mean age of patients was 76 years (range 62-88 years). The lesions were located in the rectum ($n=6$), sigmoid colon ($n=5$), transverse colon ($n=2$), caecum ($n=3$), and at the appendiceal

orifice ($n=1$). The average size of the lesions was 29 mm (range 20-50 mm). Mean procedural duration was 69 min (range 35-160 min).

Technical success was achieved in all patients. 94.1 % of lesions (16/17) could be resected macroscopically complete with confirmed full-thickness resection. Histological work-up of the specimens showed R0-resection in 13 patients (76.4 %), unclear margins in 2 patients (11.8 %) and positive margins in 2 patients (11.8 %).

No immediate severe adverse events such as perforation or major bleeding occurred. Follow-up endoscopic examination was available in 12 patients of which 10 were without recurrent adenoma. One patient showed a stenosis in the sigmoid colon after resection and one recurrent adenoma was detected in 2 out of the 12 patients in follow-up endoscopy. The authors concluded that hybrid EMR-FTRD in the colorectum seems to be a safe and effective technique for large non-lifting lesions with positive lifting signs in the margins. This hybrid technique might close a therapeutic gap and thus avoid surgery. Further prospective evaluation of efficacy, safety and long-term outcome of this hybrid technique is necessary.

Hybrid resection of large colorectal adenomas combining EMR and FTRD

Bauermeister M, Mende M, Hornoff S, Faiss S

September 2021 | HemoPill® acute: multicentric retrospective study shows promising results in clinical routine

The use of the HemoPill acute demonstrated its ability to reliably detect or exclude bleeding in the upper gastrointestinal tract.

Between July 2019 and March 2020, 12 hospitals from Germany and Switzerland used the HemoPill acute in 61 patients and analysed it retrospectively. Indications for application were the clinical suspicion of upper gastrointestinal bleeding (UGIB), small bowel bleeding or recurrent bleeding after endoscopic hemostasis. Primary endpoints were technical success and detection/exclusion of bleeding. Secondary endpoints included adverse events and impact on clinical outcome.

The capsule was used in:

- 45 (73%) patients with suspected UGIB.
- 12 (20%) patients with suspected small bowel bleeding after negative esophagogastroduodenoscopy and tarry stool present.
- 4 (7%) patients to exclude rebleeding.

The technical success was 98%. 35/60 (58%) cases had a positive capsule finding (bleeding was detected), of which 20/35 (57%) had subsequent endoscopy showing bleeding. None of the 25 patients with a negative HemoPill acute finding bled again or showed signs of bleeding on further endoscopic examination. In addition, a negative result had an impact on further clinical management in 18/25 (72%) cases. In these 18 patients, emergency endoscopy was completely avoided in 10 cases and subsequent endoscopy in the remaining 8 cases. No serious adverse events occurred.

The authors evaluated the HemoPill acute as a promising tool for real-time detection of UGIB. Thus, on the one hand, it helps in risk stratification of patients, and on the other hand, it also detects the absence of fresh blood in hemodynamically stable patients who may be candidates for outpatient rather than inpatient treatment.

Telemetric capsule-based upper gastrointestinal tract – blood detection – first multicentric experience

Brunk T, Schmidt A, Hochberger J, Wedi E, Meier B, Braun G, Naser F, Schneider M, Kandler J, Bauerfeind P, Repp M, Weingart V, Brand M, Caca K, Wannhoff A, Messmann H, Karpyniec S, Kubisch I, Albert J, Neuhaus H, Schmitz L, Allescher HD, Meining A, Kuellmer A

September 2021 | FTRD®: A multicenter retrospective analysis evaluated the risk of appendicitis after EFTR of lesions involving the appendiceal orifice

A multicenter retrospective analysis of 50 patients undergoing EFTR for colonic lesions involving the appendiceal orifice showed an appendicitis rate of 14 %

(7/50 patients). Four patients could be managed conservatively, in 3 cases surgical appendectomy was necessary.

S. Schmidbaur et al., University Clinic Ulm, Ulm, Germany, conducted a multicenter retrospective analysis of patients treated with endoscopic full-thickness resection with the FTRD System for colonic lesions involving the appendiceal orifice. The study aimed to evaluate the occurrence of appendicitis.

The FTRD System is suitable for resection in difficult anatomical locations with high risk of perforation, such as para-diverticular or para-appendicular lesions. In addition, the prompt closure provided by the FTRD device has been described as minimizing the risk of peritoneal irritation during resection by shortening the contact time between bowel lumen and peritoneal cavity. However, as EFTR near the appendiceal orifice is associated with a subtotal appendectomy, cases of acute appendicitis after resection have been described. The objective of this study was to further evaluate the risk associated with EFTR at this challenging location.

Overall, 50 patients with colonic lesions near the appendiceal orifice and without prior appendectomy could be included in the study. EFTR was technically successful in 48 of these patients (96 %). In the two remaining cases, the clip could not be deployed due to high angulation of the colonoscopy; EMR could be successfully performed in these cases. Follow-up was terminated after a mean of 4 months after resection.

Seven patients (14 %) developed acute appendicitis during follow-up. Four appendicitis cases occurred during the first 10 days of post-interventional monitoring, the other three cases occurred after a latency period of < 1 month. In four patients, conservative management with intravenous hydration, antibiotics (3-5 days) and analgesics was sufficient. Three patients underwent surgical appendectomy. Three out of the four patients (75 %) with appendicitis during the first 10 days after the procedure were treated conservatively, and the other patient underwent surgery, whereas two out of the three patients (67 %) with later-onset appendicitis needed surgical therapy.

In one case of sessile serrated adenoma, post-interventional perforation of the cecum was observed and treated surgically. No other post-procedural complications (bleeding, severe pain, unexpected hospitalization) were reported.

The authors concluded that EFTR of lesions near or affecting the appendiceal orifice was associated with an acceptable complication rate. The risk of developing acute appendicitis was 14 %; however, 57 % of these cases could be treated conservatively owing to prompt detection. Further studies to determine risk factors for development of post-procedural appendicitis are mandatory.

Risk of appendicitis after endoscopic full-thickness resection of lesions involving the appendiceal orifice: a retrospective analysis

Schmidbaur S, Wannhoff A, Walter B, Meier B, Schäfer C, Meining A, Caca K

August 2021 | FTRD®: hybrid technique combining EMR and EFTR permits resection of large polyps

65 consecutive patients underwent EFTR of colorectal lesions, either via stand-alone EFTR (lesions <2 cm in size but not amenable to standard EMR; $n = 38$ lesions) or hybrid technique combining EMR and EFTR (lesions > 2cm; $n = 31$ lesions). Clinical success (91 %), technical success (83 %) and R0 resection rate (81 %) did not differ between the groups. Three adverse events occurred including two patients developing an acute appendicitis. S. Mahadev et al., Weill Cornell Medicine, New York, NY, USA performed a single-center retrospective study of consecutive patients with colorectal lesions unresectable by conventional EMR alone, that were treated via stand-alone EFTR with the FTRD System or using a hybrid technique combining the use of EMR and EFTR.

For lesions less than 2 cm in size but not amenable to EMR alone, EFTR was performed as a stand-alone procedure. For lesions larger than 2 cm, EMR was performed initially as part of a hybrid procedure with consecutive EFTR. The intention of the EMR procedure was to remove the lateral spreading component of the lesion in order to reduce tissue volume.

Consecutive EFTR would then remove fixed parts of the lesion and parts suspected of high-grade dysplasia or cancer. 38 stand-alone EFTR treatments and 31 hybrid treatments were performed in 65 patients (61 % male, average age 70 years, range 50-90 years). The most common indications were non-lifting polyps (43 %) and suspected high-grade dysplasia or carcinoma (38 %). The lesion size was significantly larger for hybrid EMR + EFTR (mean 39 mm, range 15-70 mm) than for stand-alone EFTR (mean 17 mm, range 7-25 mm, $p < 0.01$). Prior resection had been attempted in 20 cases (29 %), 13 in the stand-alone group (34 %) and 7 in the hybrid group (23 %).

Clinical success (91 %), technical success (83 %), and R0 resection rates (81 %) did not differ between stand-alone and hybrid groups. One reason for technical failure was inability to advance the FTRD through the colon to the target lesion, which occurred in 6 % (4 of 69 procedures). In 2 patients in whom the initial attempt to reach the target lesion was unsuccessful, a second attempt 1 to 3 months later was successful. Other reasons for technical failure were inability to pull tissue into the cap with grasping forceps (2 cases) tissue slippage out of the snare before resection (3 cases), and snare breakage (2 cases). In these cases the neoplastic tissue could be removed subsequently using hot snare and hot biopsy forceps. These events were considered technical failures for the purposes of the study, even if complete resection of neoplasia was achieved for clinical success.

96 % of patients were discharged home on the day of intervention. Three adverse events occurred, including 2 patients who developed acute appendicitis, and one patient with a transmural perforation due to improper deployment of the FTRD which could be successfully closed endoscopically. The authors concluded that a hybrid approach combining the use of EMR and EFTR maintains safety and efficacy while permitting the resection of significantly larger lesions than FTRD alone.

Outcomes of hybrid technique using endoscopic mucosal resection and endoscopic full-thickness resection for polyps not amenable to standard techniques (with video)

Mahadev S, Vareedayah AA, Yuen S, Yuen W, Koller KA, Haber GB.

July 2021 | Conference Report Digestive Disease Week

stentfix OTSC® System

Preventing migration of FCSEMS: stentfix OTSC shows advantages in comparison to no fixation and endoscopic suturing

D. Lew, et. al., Cedars-Sinai Medical Center, Los Angeles, California, US, presented findings from a retrospective cohort study conducted between January 2013 to October 2020 including 199 patients with FCSEMS placement in the GI tract. Multiple patients who required repeat procedures for stent adjustment and/or fixation, were counted as a separate event. The study included 438 procedures which were performed respectively in 60% (264) without fixation, 34% (150) with suturing, and 5% (23) with the recently available stentfix OTSC. Indications for stent placement included 44% benign stricture, 35% fistula/perforation, and 20% malignant stricture. Most stents (53%) were placed in the esophagus. In the stent migration rates at all weeks assessed, there was not a significant difference between stentfix OTSC and suturing ($p > 0.05$). Better results regarding the migration rate was achieved when comparing both stentfix OTSC and suturing to no fixation up to 8 weeks ($p = 0.02$). Beyond 8 weeks, there were no significant differences in all groups ($p > 0.05$). Median time to stent migration for no fixation was 3 weeks compared to 5 weeks with stentfix OTSC and suturing ($p = 0.005$).

Clinical success rate for the cohort occurred in 61% of the patients ($n = 121$). The total median procedure time for a fixation with stentfix OTSC was, with a difference of 24 min, significantly shorter than suturing (44 vs 68 minutes, $p = 0.002$). Rates of adverse events, including chest/abdominal pain and nausea/vomiting, were not significantly different but trended towards being lowest in the stentfix OTSC group at 9%, compared to 21% with no fixation, and 18% with suturing ($p > 0.05$). No perforations occurred. The authors found that stentfix OTSC and endoscopic

suturing seem to be equally effective in preventing stent migration when compared to no fixation. However, stentfix OTSC showed advantages by a decreased overall procedure time, less adverse events, and especially lower costs. Therefore, it may be preferred over endoscopic suturing. Anyway, further studies with larger sample size are needed.

Comparison of no stent fixation, full-thickness endoscopic suturing, and over-the-scope-clip (OTSC) in preventing migration of fully covered self expanding metal stents (FCSEMS)

Lew D1, Patel S1, Liu Q1, Gaddam S1, Gupta K1, Jamil LH2, Lo SK1, Park KH1

1Los Angeles, 2Royal Oak

How to prevent migration of esophageal stents: efficacy and safety of endoscopic fixation with suturing, stentfix OTSC and hemostatic clips

M. Coronel, et. al., The University of Texas MD Anderson Cancer Center, Houston, Texas, US, presented a single-center retrospective study reviewing patients who underwent esophageal FCSEMS with endoscopic stent fixation (ESF), performed from 9/1/2017 to 9/1/2020. A patient cohort was identified ($n = 21$) and matched (1:1 ratio) with a retrospective, consecutive patient cohort who underwent esophageal FCSEMS placement without ESF ($n = 21$). The primary outcome was to identify early (<30 days) or late (>30 days) stent migration rates in these two groups.

2 In the ESF group, endoscopic suturing was used in 67% (14) of the patients, stentfix OTSC in 23% (5), and hemostatic clips in 10% (2). The FCSEMS with ESF group showed no early stent migration and delayed migration in 3 patients (14%). These 3 patients had successful stent removal and in 1 patient another stent was placed with endoscopic suturing. The FCSEMS without ESF group showed early migration in 1 patient and delayed migration in 6 patients (33%) ($p = 0.14$). In 2 patients the stents could not be removed as they were noted to be impacted in the stomach.

While not statistically significant, migration rates were lower in the ESF group without causing additional complications when compared to the FCSEMS alone group. The authors concluded that endoscopic stent fixation appears to be safe and effective.

Efficacy and safety of endoscopic fixation to prevent migration of esophageal stents. A tertiary care center experience.

Coronel M1, Ge PS1, Kumar S1, Lum P1, Weston BR1, Ross WA1, Raju GS1, Lee J1, Coronel E1

1Houston

FTRD® System

Exploring the role of full-thickness resection with FTRD System in the management and staging of patients with suspected T1 colorectal carcinoma

A. Vareedayah, et. al, NYU Langone Health, New York, New York, USA, report on a single center retrospective study with prospectively collected data from patients who underwent FTR for colorectal lesions, focusing on those with histological evidence of carcinoma. Of a total of 64 patients, a subgroup of 12 patients (50 % male) with a mean age of 71, showed suspected malignancy in the biopsies taken. 9 underwent FTR alone and 3 as part of a hybrid approach for large, laterally spreading lesions with a combination of EMR followed by FTR (mean lesion size was 23 mm). 3 patients had known carcinoma prior to resection. Primary outcomes of this study were: technical success and clinical success (defined as macroscopically complete resection). Technical success was achieved in 9/12 cases; 2 failures were due to malfunction of the snare and could be resolved by resection of the lesion above the inserted clip with a separate snare. One case was inadvertent snare excision prior to clip deployment with post resection perforation closed by subsequent clip release. Clinical success was achieved in 100% of cases; R0 resection in 11 of 12 patients (1 patient showed invasion into the muscularis propria with a positive vertical margin). Of 4 patients who underwent further surgery, surgical pathology revealed one T1 stage with lymphatic invasion, two T2 and one T3 stages. On the second postoperative day after resection of a T1 adenocarcinoma at the appendiceal orifice, one patient presented to the

emergency department with nausea, vomiting and abdominal pain. CT showed peri-appendiceal inflammatory change that was resolved with IV antibiotics. No other adverse events were reported in this subgroup of patients. The authors conclude that these data suggest that FTR is technically feasible, effective, and safe for the removal of early-stage, low-risk colorectal cancer. No residual tumor was found in any of the patients who underwent surgery. Further studies and longer follow-up are needed to clarify the recurrence rates and long-term survival of these patients.

Final pathologic staging of patients with suspected T1 colorectal carcinoma with full thickness resection (FTR): A single center North American experience

Vareedayah AA1, Yuen PYS1, Ooka K1, Morales SJ1, Mahadev S2, Haber GB1

1New York, 2New York

An alternative to surgical resection: EFTR with FTRD System at an academic tertiary care cancer center setting

P. Ge, et. al, The University of Texas MD Anderson Cancer Center, Houston, US, Texas, presented on their initial clinical outcomes of EFTR at a major US-based academic cancer center. Demographics, procedural and technical characteristics, lesion characteristics, resection outcomes, and histopathologic diagnosis were recorded. Adverse events were recorded including perforation, 3 bleeding, pain/discomfort, intraluminal trauma during instrument insertion, and inadvertent injury to adjacent organs.

The EFTR was performed using the FTRD System for 3 gastric, 2 colonic, and 5 rectal lesions in a total of 10 patients (30 % female) with a mean age of 63.6 years. Indications included gastric neuroendocrine tumor (NET) (30 %), rectal NET (10 %), recurrent adenoma at EMR site (10 %), and prior polypectomy scars where the original pathology included adenocarcinoma or NET with positive margins (50 %). The gastroduodenal FTRD was utilized in all gastric cases and colonic FTRD was utilized in all colorectal cases. A 20 mm balloon was used to facilitate device passage through the upper esophageal sphincter in all 3 gastric cases and through the anal canal in 1 colorectal case. All 7 colorectal cases were notable for prior EMR scars and severe fibrosis. En bloc resection, R0 resection, and curative resection were achieved respectively in 100 %, 100 %, and 90 % cases. The lone non-curative resection was a rectal NET with negative resection margins but with lymphovascular invasion and perineural invasion on histopathology. Final pathology indicated NET in 4 patients, tubular adenoma in 1 patient, and clean polypectomy scar with no residual tumor in 5 patients. The mean specimen size was 26.5 mm (SD, 3.4 mm). Mean EFTR time was 9.9 min (SD, 6.6 min) and mean total procedure time was 51.4 min (SD, 16.0 min). Three minor adverse events were noted, including perineal pain, abdominal pain, and tenesmus, all of which were self-limited with spontaneous resolution within 1 day. No major adverse events occurred, and no patients required hospitalization. On short term follow-up (mean 1.7 months) no delayed adverse events could be noted.

The authors concluded that EFTR is an effective therapeutic option despite severe fibrosis and occupies a unique role in the endoscopic management of small fibrotic or subepithelial lesions which are otherwise unsuitable for conventional endoscopic resection techniques. Ongoing studies with long-term follow-up will seek to further validate these findings.

Exploring the role of endoscopic full thickness resection at an academic tertiary care cancer center setting

Ge PS1, Coronel M1, Tillman MM1, Badgwell B1, Bednarski BK1, Chang GJ1, Katz M1, Rodriguez-Bigas MA1, You YN1, Halperin DM1, Casanova DN1, Weston BR1, Bhutani MS1, Lee J1, Ross WA1, Coronel E1

1Houston

EFTR of upper gastrointestinal lesions with colonic FTRD - A retrospective observational of a 13 cases series

J. Nilsson, et. al, University of Alberta, Edmonton, Alberta, Canada presented outcomes from a retrospective observational case series consisting of 13 cases of duodenal (4) and gastric (9) lesion resections with the colonic FTRD System. Indications for EFTR were sub-epithelial tumor

(n=8), polyp (n=2) and scar-resection (n=3). The colonic FTRD could pass the upper esophageal sphincter or pylorus without dilatation and could be advanced to the lesion in 13/13 cases (100 %). One sub-epithelial lesion was too big for the cap and one scar could not be sucked into the cap. R0-resection rate for deployed clips was 10/11 (91 %). Technical success was achieved in 11/13 (85 %) of procedures. There were two superficial esophageal tears from FTRD insertion that required no therapy. No bleeding occurred in the postoperative period. This study further confirms acceptable efficacy and safety of EFTR in the upper GI use.

Endoscopic Full-Thickness Resection of upper gastrointestinal lesions using a colonic FTRD- A retrospective

observational case series of 13 FTRD cases

Nilsson JE2, Koch AD1, de Graaf W1

1Rotterdam, 2Edmonton

FTR in the UGIT: efficacy and safety of diagnostic FTRD and gastroduodenal FTRD

S. Yuen, et. al, NYU Langone Health, New York, New York, US, presented their retrospective single center, single endoscopist analysis of prospectively collected data on lesions resected with the FTRD Systems in the UGIT. The study included 6 patients (3 male) with a mean age of 67 years. Indications for FTR were gastric subepithelial tumors (4) or duodenal polyps (2). Three patients had endoscopic ultrasound prior to FTR. Mean lesion size was 15.6 mm in 5 patients and 40 mm in 1. Four patients had no prior attempts at removal and two were referred for FTR for polyp recurrence after prior (partly 4 multiple) polypectomy. For the 4 cm polyp recurrence a hybrid resection comprising EMR of the periphery and FTR of the central non-lifting component was used. Primary outcomes included technical success, clinical success, and R0 resection rate. Resections were performed using either the diagnostic FTRD or gastroduodenal FTRD, both with cap diameter of 19.5 mm and depth of 23mm. 5/6 patients underwent pre-dilation (Savary 20 mm dilator and/or 20 mm balloon) prior to passage of the FTRD through the upper esophageal sphincter. The FTRD comprising cap, snare and sheath was inserted and advanced to the resection site. Additional dilation of the pyloric channel was performed for duodenal access. A grasping forceps was used to gently pull the lesion into the cap, with adjuvant use of suction as required. The FTRD clip was then deployed, immediately followed by snare closure and electrosurgical excision of the entrapped tissue with a pure cut current. FTRD was advanced without difficulty after appropriate dilation or advancement over a balloon, and deployment was successfully performed. Clinical success, defined as macroscopically complete lesion resection, was achieved in 5/6 cases. In one antral leiomyoma, the deep margin was positive. Adverse events, including bleeding and perforation, did not occur. The authors concluded that FTR in the upper gastrointestinal tract was technically successful in 100 % of the cases. Pre-dilation and balloon assist are critical elements in successful deployment but were necessary as the major concern with this device was the large outer diameter with possible perforation or inability to pass the sphincters. UGIT FTR provides a simple, quick, and effective alternative for deep resection of subepithelial lesions and duodenal lesions notorious for complications of bleeding and perforation.

Efficacy and Safety of Full Thickness Resection (FTR) in the Upper Gastrointestinal Tract (UGIT): A Single Center North American Experience

Yuen PYS1, Vareedayah AA1, Morales SJ1, Ooka K1, Haber GB1

1New York

OTSC® System

Presentation of comparative data regarding OTSC vs standard therapy for the prevention of rebleeding in peptic ulcers

S. M. Chan, et. al, Chinese University of Hong Kong, Hong Kong, performed a multicenter randomized controlled trial from July 2017 till Oct 2020. Aim of the study was to compare the efficacy of the OTSC to standard endoscopic therapy in primary treatment of patients with peptic ulcer bleeding in

peptic ulcers ≥ 1.5 cm.

Inclusion criteria was a Forest la-IbB bleeding. The primary outcome was a clinical rebleeding within 30 days. 100 patients were enrolled and a crossover between the methods was allowed in case of failure.

There were two cases of successful crossover from failed standard therapy to OTSC and four unsuccessful crossovers after failed OTSC placement to the standard treatment. The overall rebleeding within 30 days was calculated as 14 % (7/50) in the OTSC arm and 16 % in the standard arm. The all-cause mortality was 4% (2/50) with OTSC vs 8% (4/50) in the standard arm (p=0.68). None of the patients required surgical intervention.

The authors state that, for large ulcers ≥ 1.5 cm, OTSC as primary hemostasis did not confer to an improvement in clinical outcomes in this study.

Several aspects of the study design remained unclear in the abstract presented such as the high percentage of patients (approx. 80 %) without active bleeding in the study population, also in reference to the primary endpoint being the re-bleeding rate at 30 days.

The Use of Over-the-scope-clip (OTSC) Versus Standard Therapy for the Prevention of Rebleeding in High Risk

Peptic Ulcers: A Randomized Controlled Trial

Chan SM1, Pittayanon R2, Wang H3, Chen J4, Kuo Y3, Teoh AY1, Yip H1, Tang RS1, Ng S1, Wong SH1, Mak JW1, Chan H1, Lau L1, Lui RN1, Wong M1, Rerknimitr R2, Chiu PW1, Ng

1Hong Kong, 2Bangkok, 3Taipei, 4Taipei

5

A Purely Endoscopic Management Approach for Type V Mirizzi Syndromes

S. Al Ghamdi, et. al, The Johns Hopkins Hospital, Baltimore, MD, present a case of a 94-year-old female who was transferred to their hospital with Type V Mirizzi Syndrome with many large stones in the cystic duct and gallbladder. In the first session stones were removed endoscopically.

Colonoscopy was performed one week later to evaluate the fistula. Several small gallstones were seen in the colon, along with a small fistula with intermittent extrusion of microvilli at the hepatic flexure. Contrast was injected into the defect, confirming the correct location of the CCF. The fistula tract was deepithelialized using argon plasma coagulation and a hemoclip was placed to mark the area. An OTSC clip (14/6 t) was successful deployed over the fistula and completely closed it, as later confirmed by a contrast injection. After the procedure, the patient's abdominal pain resolved completely. She returned five months later to undergo her final ERCP. Cholangioscopy revealed a residual stone at the CD takeoff that was successfully treated with EHL. An occlusion cholangiogram confirmed complete stone clearance.

This case demonstrates the use of ERCP and EHL for successful management of a large, impacted CD stone resulting in biliary obstruction and acute cholangitis. Endoscopic management of a concomitant CCF using an OTSC clip was also successful. This suggests the feasibility of a purely endoscopic management approach for Type V Mirizzi Syndrome.

A Purely Endoscopic Management Approach for Type V Mirizzi Syndromes

Ghamdi SA, Bejjani M, Ghandour B, Khashab MA

Baltimore

Endoscopically Directed Single Port Intra-gastric Fundoplication, Sleeve, and Myotomy: A Preclinical Study

H. Hernandez-Lara, et. al, Mayo Clinic Rochester, Rochester, MN, presented a video demonstrating the safety and feasibility of various peroral endoscope directed stapling procedures using a novel percutaneous endoscopic trocar and a laparoscopic stapler in 2 domestic pigs. OTSC clips were used to close the trocar insertion and the intra-gastric port site with no adverse events. The authors concluded that single port, endoscopically directed intra-gastric procedures such as fundoplication, sleeve, and myotomy using laparoscopic staplers are now feasible.

Endoscopically Directed Single Port Intra-gastric

Fundoplication, Sleeve, and Myotomy: A Preclinical Study

Hernandez-Lara H, Abu Dayyeh BK, Garcia Garcia de Paredes A, Rajan E, Storm AC.

Rochester

July 2021 | HemoPill® acute shows promising results in a non-fasting patient with a mid-gastrointestinal bleeding: a case report

Ovesco's first wireless sensor capsule for detection of upper and middle gastrointestinal (GI) bleeding was successfully used in a patient with a bleeding caused by Meckel's diverticulum and no patient preparation.

D. Wiedbrauck, S. Hollerbach and F. Wiedbrauck, AKH Celle, Germany reported on their experience with HemoPill acute to clarify an acute bleeding situation in a patient who had undergone EGD, colonoscopy and CT-angiography, which did not reveal the source of bleeding.

The cause of anaemia (hemoglobin 9.9 g/dL) and episodes of severe haematochezia were investigated in a 52-year-old male patient but could not be determined with the above procedures. Afterwards, the patient was stable and had eaten.

A few hours after the procedure, his hemoglobin level dropped severely to 6.4 g/dL and due to the food intake, a Video Capsule Endoscopy (VCE) could no longer be performed. Thus, they decided to use the HemoPill acute which detected blood after 6 h 58 min. Over time, it could be assumed that the bleeding would have to be in the ileum/terminal ileum. Afterwards, a VCE exposed a bleeding Meckel's diverticulum, and the patient underwent a segmental ileal partition.

The authors concluded that the HemoPill acute might serve as add-on pre-diagnostic tool to detect and localize bleedings in the upper and middle GI tract in non-fasting patients.

Detection of mid-gastrointestinal bleeding caused by Meckel's diverticulum using a novel telemetric sensor capsule in a non-fasting patient

Wiedbrauck D, Hollerbach S, Wiedbrauck F

July 2021 | BougieCap®: significant symptomatic improvement of esophageal strictures achieved by a single bougienage in patients with eosinophilic esophagitis

A cohort study of 50 patients with symptomatic eosinophilic esophagitis treated with a single dilation session using the BougieCap showed 100 % technical success, a drop of median symptom severity from 32 to 0 points, and no severe adverse events.

A. M. Schoepfer et al., Centre Hospitalier Universitaire Vaudois (CHUV) and University of Lausanne, Lausanne, Switzerland, presented a study evaluating technical feasibility, clinical efficacy, and safety of the BougieCap for esophageal stricture dilation in patients with eosinophilic esophagitis (EoE). The BougieCap is a bougienage instrument that allows optical and tactile feedback during endoscopic stricture dilation in the upper gastrointestinal tract. Fifty EoE-patients (30 % female, median age 41 years) with presence of esophageal strictures (esophageal diameter ≤ 14 mm) and stricture-related symptoms were included. Median disease duration in patients was 4 years, 50 % were being treated with swallowed topical corticosteroids, 10 % with proton pump inhibitors, 14 % with combined swallowed topical corticosteroids plus proton pump inhibitors, 14 % with elimination diet, and 12 % were without anti-eosinophil therapy. Symptoms were assessed before and 2 weeks after BougieCap treatment using the validated Eosinophilic Esophagitis Activity Index Patient Reported Outcomes (EeAI PRO) questionnaire (score ranges from 0 – 100 points). Endoscopic bougienage was technically successful in 100 %. All procedures were performed without introduction of a guidewire and without fluoroscopic guidance. Median esophageal diameter increased from 12 mm (IQR 12-13) to 16 mm (IQR 16-16, $p < 0.001$). Median symptom severity dropped from 32 points (IQR 27-41) to 0 points (IQR 0-10, $p < 0.001$) at 2 weeks post dilation. In one patient the BougieCap was temporarily lost after stricture dilation in the

hypopharynx but could be retrieved. No severe adverse events were reported. The authors concluded that in adults with EoE, endoscopic treatment of esophageal strictures using the BougieCap is technically feasible, safe and offers significant symptomatic improvement in the short term.

Technical feasibility, clinical effectiveness, and safety of esophageal stricture dilation using a novel endoscopic attachment cap in adults with eosinophilic esophagitis
Schoepfer AM, Henchoz S, Biedemann L, Schreiner P, Greuter T, Reinhard A, Senn J, Franke A, Burri E, Juillerat P, Simon HU, Straumann A, Sfroneeva E, Godat S

June 2021 | OTSG Xcavator™

A new dimension of endoscopic grasping: the OTSG Xcavator™

The removal of hard to grasp tissues and foreign bodies from the gastrointestinal tract often represents a challenge in endoscopy. In particular, pancreatic necrosectomy can be frustrating and tedious when using snares, graspers, etc. . Then as an endoscopist, you need a tool for reliable and effective removal of tissue. The OTSG Xcavator™ now represents a new dimension in endoscopic grasping. The extralarge grasping capacity of the OTSG Xcavator™ allows easy and effective removal of necrotic (pancreatic) tissue, blood clots, food boluses, and foreign bodies. Due to the external control, the working channel remains free for additional suction/irrigation capacity and allows the combination with further instruments (e.g. snares, graspers, etc.). Professor Meinert first presented the instrument at Endoscopy on Air 2020 and Dr. Brand published the first case study in the journal Endoscopy. Just recently, the OTSG Xcavator™ was also demonstrated at Endoskopie Live in Berlin by Professor Faiss, as well as by Dr. Möschler at the congress "Viszeralmedizin NRW".

June 2021 | Preliminary results from RCT study show that OTSC first-line treatment is superior to standard treatment for acute non-varicose upper GI bleeding with high risk of rebleeding

Summary of the presentation on the STING II study by Prof. K. Caca, Ludwigsburg at the DGE-BV 2021.

Subject of the Symposium on April 8th was telemetric haemorrhage detection with the HemoPill® acute followed by two presentations on hemostasis with the OTSC® System.

Prof. Dr. K Caca, Ludwigsburg: Primary haemostasis with the OTSC System is highly effective

Prof. K. Caca, Ludwigsburg, presented preliminary data on the RCT he coordinated (STING II) on the use of the OTSC System in severe, non-varicose upper GI bleeding. Despite good treatment options, the mortality of upper GI bleeding is still relatively high at 5–10%, especially in older and comorbid patients. Endoscopy is primarily successful in approx. 90% of cases with the current standard procedures, but the application of endoclips can be difficult, especially with a hard ulcer surface. In the treatment of recurrent bleeding, OTSC has become established due to the significantly lower rebleeding rates compared to standard treatment (STING study). STING II was subsequently initiated to investigate whether OTSC treatment is superior to standard treatment also in acute, non-varicose upper GI bleeding with a high risk of rebleeding ("complete Rockall score" ≥ 7 points). Standard treatment consisted of at least 2 conventional clips or thermal procedure plus adrenaline injection. OTSC treatment consisted of primary OTSC application (injection was allowed). According to the calculation of the number of cases (power 80%), $n = 100$ patients with endoscopically confirmed acute (non-varicose) upper GI bleeding and a "complete Rockall score" ≥ 7 points were included. Exclusion criteria included prior endoscopic treatment within the last 4 weeks to rule out treatment of recurrent bleeding.

The primary endpoint was successful haemostasis without recurrent bleeding within 7 days. Secondary endpoints were the need for second-line treatment, need for transfusion, hospitalisation duration, ICU stay and 30-day mortality.

After randomisation, $n = 52$ patients could be evaluated in the standard arm and $n = 48$ patients in the OTSC arm. Both groups show good correspondence to the baseline. 42.3% of patients in the standard arm and 39.6% in the OTSC arm were

on anticoagulation or platelet aggregation inhibition. The median age of patients was 79 years (51–96) in the standard arm and 78 years (42–92) in the OTSC arm. The median Rockall score was 8 points in both groups, and the predominant localisation was in the duodenum (46.1% standard vs. 60.4% OTSC). Mainly peptic ulcer bleeding was present with approx. 95% in both arms. Further details can be found in the video (presentation from 32:30 min).

Preliminary analysis of the data shows that the OTSC System is significantly superior ($p = 0.019$) to standard treatment at 91.7% (44/48) to 73.1% (38/52) in successful haemostasis (no rebleeding or persistent bleeding within 7 days). There was no persistent bleeding with the OTSC compared with 6 (11.5%) persistent bleeding in the standard arm ($p = 0.027$). All persistent bleeding in the standard arm could be successfully stopped with OTSC. In the OTSC, 8.3% (4/48) of rebleeds occurred within 7 days compared with 15.4% (8/52) in the standard arm ($p = 0.362$). The treatment time was comparable at 27 min (OTSC) and 28 min. Late rebleeding (day 8–30) occurred in 2 cases (4.2%) with OTSC versus none in the standard treatment arm ($p = 0.028$). This resulted in a total bleeding rate (incl. persistent bleeding) of 12.5% (6/48) with OTSC and 26.9% (14/52) with standard treatment ($p = 0.084$). There were no significant differences with regard to the other endpoints (blood transfusion, ICU stay, length of hospitalisation, mortality and need for second-line treatment). Prof. Caca summarised that OTSC is superior to standard treatment for acute non-varicose upper GI bleeding with a high risk of rebleeding and that primary haemostasis by OTSC is highly effective (no persistent bleeding in the OTSC arm). Furthermore, no additional time was required for OTSC treatment. The study is currently the largest RCT study regarding the first-line treatment of upper GI bleeding with the OTSC System. However, it should be noted that the data is currently still preliminary and not published.

STING II – randomised controlled data on first-line treatment of severe upper GI bleeding with the OTSC® System

(STING II - randomisiert kontrollierte Daten zur First-Line Therapie von schweren OGI-Blutungen mit dem OTSC® System)

K. Caca, Klinikum Ludwigsburg, Ludwigsburg

May 2021 | Symposium report DGE-BV 2021 – Part 2

OTSC® is recommended in current guidelines for the first-line treatment of large ulcers and the treatment of refractory ulcer bleeding

The 50th Conference of the German Society for Endoscopy and Imaging Procedures (DGEV) took place on April 8 – 10, 2021. Subject of the Symposium on April 8th was telemetric haemorrhage detection with the HemoPill® acute followed by two presentations on hemostasis with the OTSC® System.

Prof. Dr. Arthur Schmidt, Freiburg: OTSC System superior especially in large ulcers with fibrotic ulcer surface and in difficult anatomical location.

A. Schmidt, University medical centre Freiburg, gave an overview of the current data situation and guideline recommendations for OTSC treatment of ulcer bleeding in his lecture and presented study data on OTSC vs. TAE as well as OTSC vs. surgery. Schmidt reported that the OTSC System is mainly used for large ulcers with a fibrotic ulcer surface as well as for difficult anatomical locations (e.g. rear wall of the duodenal bulb) where other endoscopic procedures (e.g. injection in combination with endoclips) reach their limits. Compared to the treatment with endoclips, OTSC clips offer several advantages: Higher compression force (Kato et al., GIE 2012), better hold in the fibrotic tissue (Mönkemüller et al., Endoscopy 2014), and better visualisation of the bleeding source through the spacer cap (Mönkemüller et al., Endoscopy 2015). In large ulcers, the aim is not to close the entire ulcer with the OTSC, but to compress the vascular stump at the surface of the ulcer with the OTSC, which can achieve sustained haemostasis (Schmidt et al., Gastrointest Endosc Clin N Am 2020).

Many studies are available on OTSC treatment for bleeding peptic ulcers, some of them with a large number of cases. However, these are almost exclusively retrospective. A metaanalysis of this retrospective data (Weiland et al, Minim Invas Allied Technol 2019), shows a technical success rate of

93.0%, a clinical success rate of 87.5% and a rebleeding rate of 8.3%. The updated ESGE guideline recommendations (Gralnek et al., 2021) now also provide an optional recommendation for first-line treatment with the OTSC System for the treatment of large ulcers (> 2 cm, with visible vascular stump > 2 mm or located in high-risk vascularised regions).

In refractory bleeding, there is one prospective randomised trial to date (STING Study, Schmidt et al., Gas-troenterology 2018) comparing OTSC ($n = 33$) with previous standard treatment ($n = 33$). The primary end-point of the study was the occurrence of further bleeding (persistent bleeding or recurrent bleeding). This occurred in 57.6% of patients in the standard treatment group and only in 15.2% in the OTSC group ($p=0.001$). Thus, the OTSC resulted in a relative risk reduction of 73.6%.

As guideline recommendations for persistent ulcer bleeding, there is a DGVS guideline (Goetz et al., 2017) that recommends the use of OTSC clips or haemostasis sprays if bleeding persists after standard treatment (open recommendation, strong consensus). The ESGE guideline updated in 2021 also recommends the use of OTSC or haemostasis spray/powder in cases of persistent or refractory bleeding.

Schmidt also presented a retrospective multicentre study on OTSC vs. TAE in the treatment of refractory ulcer bleeding (Kuellmer et al., UEGW 2020, submitted for publication) and a retrospective monocentric study on OTSC vs. surgery (Kuellmer et al., ESGE Days 2021). In the comparative study OTSC vs. TAE, the patient groups were matched using propensity score matching so that 40 patients were available for each analysis. The analysis showed that OTSC is significantly superior to TAE in terms of both hospital mortality (5.0% vs. 22.5%) and length of ICU stay (5 days vs. 9 days). The clinical success (72.5% vs. 62.5%) and the rebleeding rate within 7 days (17.5% vs. 32.5%) were better with the OTSC, but the differences were not significant here. In summary, the study shows that TAE is associated with greater risks and a more severe disease course compared to OTSC treatment as a significantly more invasive procedure. The comparative study OTSC vs. surgery also showed significantly higher complication and mortality rates in the surgery group, as surgical care is a much more invasive procedure. These results are in line with the ESGE guideline (Gralnek et al., 2021) which recommends exhausting all endoscopic options before considering trans-arterial embolisation (TAE) in cases of persistent bleeding. Surgical treatment is only recommended if the TAE fails or if a TAE is not possible.

Superiority of the OTSC® System in salvage treatment of refractory ulcer bleeding: Technical aspects and current data A. Schmidt, University medical centre Freiburg, Freiburg

May 2021 | Symposium report DGE-BV 2021 – Part 1

HemoPill® & OTSC®: telemetric hemorrhage detection and latest data on hemostasis

The 50th Conference of the German Society for Endoscopy and Imaging Procedures (DGEV) took place on April 8 – 10, 2021. Subject of the Symposium on April 8th was telemetric haemorrhage detection with the HemoPill® acute followed by two presentations on hemostasis with the OTSC® System. HemoPill® acute is a suitable tool for risk stratification in upper GI bleeding in questions of emergency esophagogastrroduodenoscopy (EGD) and inpatient admission Dr. Thomas Brunk, Berlin: The HemoPill® acute reliably detects and excludes intraluminal blood/hematin T. Brunk, Vivantes Hospital Friedrichshain, Berlin, presented a retrospective observational study on HemoPill application in cases of suspected upper GI bleeding and a case series with HemoPill application and following double-balloon enteroscopy in cases of suspected small intestine bleeding. The HemoPill acute is a capsule with novel sensor technique for telemetric diagnosis of blood/hematin. An essential and unique advantage of the HemoPill is the possibility of application in patients without empty stomach. Measurements with the HemoPill allow a reliable discrimination between blood/hematin and other fluids.

Between 06/2019 and 03/2020, a retrospective observational

study with the HemoPill was performed. 61 cases of HemoPill application in 12 centers in Germany and Switzerland were evaluated. Indications for HemoPill application were suspected upper GI bleeding, small intestine bleeding, or recurrent GI bleeding. Primary endpoints of the study were correct detection and correct exclusion of haemorrhage, respectively, and technical success of the measurement. Secondary endpoints were complications and impacts on the further clinical course. 60 measurements could be included in the evaluation, one measurement was technically not utilisable, the technical success rate was 98 %. 35/60 measurements were positive, in 20 of these patients a source of bleeding was detected endoscopically, in 9 patients no source of bleeding was found by gastroscopy, in 6 patients no endoscopy was performed. The mean latency between HemoPill examination and endoscopic examination was 15h15min in the 9 patients without endoscopic findings of a bleeding source, which could also be a reason for the lack of detectability. In 25/60 cases the HemoPill measurement gave a negative test result, in 20 of these patients no bleeding source was found in the gastroscopy, in the 5 other patients of this group no endoscopy was performed, in none of these patients clinical signs of bleeding occurred until discharge. In summary, both, blood detection and haemorrhage exclusion corresponded to the results of endoscopy and the further clinical symptoms. Apart from one case with temporary mild swallowing symptoms no complications occurred. The results of HemoPill measurements had an impact on the further clinical course. A positive result led to a significant reduction of the latency to EGD (3h8min vs. 19h52min). A negative result led to a change in indication for EGD from urgent to elective in 10/25 cases, and to avoidance of control EGD in one case. The authors concluded that the HemoPill reliably detects and excludes upper GI bleeding. The application is comfortable and low in risk, the interpretation of the measurement is simple. Therefore, the HemoPill seems to be a suitable tool for risk stratification in upper GI bleeding. This interesting approach is currently subject of a prospective trial at the University Clinic Wuerzburg. Afterwards, T. Brunk presented a case series of 13 patients with melena and negative EGD who underwent HemoPill application. In 7/13 cases, blood was detected by the HemoPill and a double-balloon-enteroscopy was performed subsequently within 24h. In all 7 cases angiodysplasias were found and obliterated. In 3 of these 7 cases the lesion was still actively bleeding. In conclusion also active bleeding in the small intestine can be reliably detected by HemoPill application. For this, no preparation of the small intestine is necessary. The interpretation of the measurement results is fast and simple. A clear indication for HemoPill application in patients with melena and negative EGD should be the subject of further research.

HemoPill® acute: telemetric diagnosis of GI haemorrhage with novel sensor technique (HemoPill® acute: telemetrische Diagnostik von GI-Blutungen mit neuartiger Sensortechnik)

T.Brunk, Vivantes Klinikum Friedrichshain, Berlin

May 2021 | Congress report

50. Anniversary Congress of the DGE-BV (German Society for Endoscopy and Imaging Techniques)

Together with the professional societies DGD & DEGEA

Prof. Dr. Hagenmüller, Hamburg, in his lecture – Disasters and achievements in endoscopy – a look back and forward – which he gave at the opening of the congress, mentioned the clinical research of Prof. Caca & Prof. Schmidt regarding haemostasis with OTSC® and endoscopic full-thickness resection with the FTRD® System as one of the great accomplishments of DGE-BV members in the further development of endoscopy. There was also an evening symposium on HemoPill® & OTSC®: telemetric bleeding detection and latest data on haemostasis. In this context, the first data on the STING II study was presented, which investigated the OTSC first-line treatment of non-varicose upper GI bleeding with a high risk of rebleeding. Furthermore, several presentations dealt with the FTRD® System and OTSC®

System.

FTRD® System

Current update on the topic of full-thickness resection in the lower and upper GI tract

As part of the series of lectures on new methods and instruments in endoscopic treatment, K. Caca, Ludwigsburg, gave an update on full-thickness resection, which according to him has now become standard in clinics. EFTR with FTRD as a non-exposure technique has already found wide application to date. K. Caca supports this statement with a brief overview of the larger studies, such as the 2018 Wall Resect study with 181 patients and the two large national registers (German FTRD register with 1178 patients, Dutch register with 367 patients). All three studies show comparable results: technical success rate between 85–90%, R0 resection rate around 80%, serious complication rate between 2–2.5% and mortality of 0. The data also shows that the procedure can be used throughout the colorectum. In terms of indications, the distribution in the German FTRD register was about 2/3 difficult adenomas and 1/3 T1 carcinomas and submucosal tumours. However, no differences were found in terms of efficacy and safety with regard to size, location (colon vs. rectum) or pre-treatment. K. Caca shows that the data is identical to the Wall Resect study, with the special feature that in the registry centres from all over Germany have contributed data even with low case numbers and can use the procedure safely and effectively. In comparison, the ratio of indications is different in the Dutch registry with considerably more T1 carcinomas, but also here approximately the same results are achieved.

K. Caca then presented the hybrid FTRD. A combined technique of EMR at the edge and FTRD in the centre, e.g. when a lesion is too large or cannot be completely ablated with EMR, e.g. because there is a central non-lifting sign. In this regard, data from 10 patients has already been published by his team and updated again in 2019 (32 patients, n=4 left colon, n=26 right colon, examination time 40–140 minutes, no serious complications, 5 recurrences at follow-up). Another publication from the USA, with 62 advanced colorectal adenomas, also shows internationally comparable data. Here, 33 adenomas were resected with FTRD alone and 39 with hybrid technique with a technical success of 89% vs. 96%, an R0 resection rate of 97% vs. 96% and two serious complications (1 appendicitis and 1 perforation).

Another topic in K. Caca's update was T1 carcinomas. At the beginning, a video example was shown in which a curative resection was performed on a patient with a suspected carcinoma finding and thus a left hemicolectomy could be spared. He then presented the data of Küllmer et al. who studied 156 patients in two groups: Group 1 with post-resections of malignant polyps and Group 2 with primary non-lifting lesions with suspected carcinoma. In group 1, 76% had only a scar left in the re-resection and 22% had a carcinoma; 84% from group 1 could ultimately be classified as low-risk and 16% as high-risk lesions. In group 2, the results were reversed with 16.3% low-risk and 83.7% high-risk lesions. In the end, however, 16.3% of the lesions could be resected curatively according to the guidelines, which is why K. Caca concludes that it is worth trying if you are able to resect endoscopically.

The last aspect of the update was the full-thickness resection in the upper GI for which K. Caca also showed a video example – an EFTR with FTRD at the lower duodenal knee of a recurrence after surgical polypectomy. Here, too, the patient was spared surgery and the findings were successfully resected full-thickness with FTRD in . He briefly mentioned three clinical studies: 20 patients from Ludwigsburg, 10 patients from Rome and a publication from New York with 56 patients. Overall, the data in the duodenum shows more frequent bleeding than in the colon, which is why K. Caca also recommends from his own experience to always see that any source of bleeding in the duodenum is prophylactically coagulated. According to K. Caca, EFTR with FTRD is much less favourable in the stomach because the relevant lesions are usually larger and small lesions often do not have to be ablated. In the RESET study, he examined precisely these smaller findings, in which, however, a GIST or NET was found in only about 1/3 of the cases. However, an accurate diagnosis could be made in 100% of cases, which is why it is a good method when small findings are to be removed, e.g.

due to a patient request. For larger findings in the stomach, however, he chooses the GERD-X method, for example, for which data has also already been published.

In summary, K. Caca concludes that almost 100% of all benign GI lesions can be resected endoscopically today. However, the desired complete resection does not always have to be en bloc, which is why there are options such as hybrid FTRD. However, the limits should always be observed (e.g. submucosal infiltration). The FTRD-EFTR fills an important gap whilst the hybrid FTRD fills a small but difficult therapeutic gap. Sometimes, however, the full-thickness resection is only a large biopsy, which is then not therapeutic but diagnostic – but then supports the stratification and the definition of the further therapeutic procedure. In the end, however, K. Caca concludes that it always depends on the indication and localisation: in the colon and rectum, FTRD is to be regarded as the standard, in the duodenum it can be used particularly well for non-lifting and subepithelial findings, and in the stomach it is more suitable for smaller lesions.

Full-thickness resection update – (Update Vollwandresektion)

K. Caca, Ludwigsburg

Retrospective analysis of pain after full-thickness resection in the upper and lower GI tract

H. Heinrich, Zurich, reported on a retrospective analysis of pain in patients after full-thickness resection in the upper and lower GI tract. For this purpose, data from 107 patients (49 female, 58 male) who, between 2017 and 2021, received a full-thickness resection with the FTRD System was collected and analysed. Pain was assessed using the following factors: visual analogue scale in the endoscopy report, need for medication and hospitalisation.

Among the 107 patients evaluated (mean age 61.3 years), 34 (32%) had full-thickness resections in the upper GI tract and 73 (68%) in the lower GI tract. The R0 resection rate was n=81 (75%), with diagnoses varying and including GIST, NET, adenomas and others.

The evaluation of the data regarding pain after full-thickness resection showed that 80 patients (75%) had no pain at all. 27 patients (25%) reported pain, which presented as immediate pain in 19% (20/107) and delayed pain (> 48 h) in 6% (7/107). Ultimately, 8 of the patients were hospitalised with immediate pain due to age and comorbidities. Among the patients with delayed pain, two appendicitis and three perforations occurred, which had to be treated surgically. Two of the perforations were due to a passenger disorder in the form of constipation, the third perforation was the consequence of incorrect post-interventional feeding. Overall, hospitalisation was necessary in 12% of patients complaining of pain.

The authors concluded that delayed pain after full-thickness resection with the FTRD System is an indication for the occurrence of appendicitis or delayed perforation and also identified risk factors such as COPD and obesity. It was noted that consistent stool regulation by the patient and sufficient instruction on post-intervention nutrition are particularly important.

Pain after FTRD in the upper and lower GI tract – a single center experience

H. Heinrich¹, P. Bauerfeind¹

¹Triemli Hospital, Gastroenterology, Zurich, Switzerland

Report of perforated appendicitis as a late complication after EFTR

M. Raithel, Erlangen, presented the case of a 71-year-old female patient who had to undergo an emergency laparoscopy due to a late complication. The patient underwent EFTR with adenomectomy of the ostium appendicis 3 months previously. As the EFTR of the circular flat adenoma measuring approx. 22 mm was incomplete, the residual adenoma tissue was removed with a conventional polypectomy snare. The histopathological result showed a completely resected sessile-serrated adenoma without evidence of malignancy. The patient subsequently received prophylactic antibiotic treatment and could be discharged from hospital after 4 days following an uncomplicated course. However, 3 months after the procedure, the patient presented again with clinical symptoms of acute appendicitis. Emergency laparoscopic appendectomy revealed perforation of the appendix with local peritonitis. The defect localisation was laterally under the appendix closure due to the previously

conditioned FTRD closure. The postoperative course was unremarkable and the patient was able to leave the clinic 4 days after the emergency surgery.

M. Raithehl concluded that the case shows that there is an increased risk of complications with EFTR at the appendix region and therefore a distinction should be made between people with and without appendectomy. In contrast to the usual recording of complications (< 4 weeks after endoscopy) in EFTR, later complications would also have to be observed in the case of appendectomy.

Perforated appendicitis as a late complication after endoscopic full-thickness resection (EFTR) with OTSC closure (Perforierte Appendizitis als Spätkomplikation nach einer endoskopischen Vollwandresektion (eFTR) mit OTSC Verschluss)

Roßmeißl1, Y. Fessehaye-Seium2, H.-D. Allescher3, F. Heyder4, A. Brütting2, M. Raithehl1

1Malteser Waldkrankenhaus Erlangen, Medical clinic II, Erlangen, Germany

2Malteser Waldkrankenhaus Erlangen, Surgery, Erlangen, Germany

3Garmisch-Partenkirchen medical centre, Gastroenterology, Garmisch-Partenkirchen, Germany

4Gastroenterology specialist practice Höchststadt, Höchststadt/Aisch, Germany

Practical tips and tricks on the subject of FTRD

As part of the DEGEA Mini Symposia, I. Galla gave a brief insight into tips and tricks for FTRD from a nursing perspective. To start with, he briefly presented the field of application of the FTRD with regard to localisations and indications, then the instrument set in its composition as a procedure set (incl. the auxiliary instruments) and the individual components of the system. An essential part of the presentation was important information on the preparation of the procedure, e.g. clarification, resource planning, task distribution, endoscope selection with correct diameter, correct HF parameters, aids, important aspects regarding the assembly of the system, e.g. do not twist the clip, push the cap forwards enough, apply the cover correctly. Finally, I. Galla gave tips for after the procedure, such as correct disassembly or issuing the MRI pass, and referred to supporting application documents that can be requested via the Ovesco field service.

Mini Symposia (DEGEA) – Tips and tricks for the FTRD
Ingo Galla

OTSC® System

Retrospective comparative study conducted on OTSC type-a vs. type-t clips in duodenal ulcer bleeding

M. Hollenbach, Leipzig, reported on a retrospective analysis comparing OTSC type-a versus type-t clips in the treatment of duodenal ulcer bleeding. For this purpose, data from 2009 – 2020 from 6 endoscopy centres was analysed and all patients treated with the OTSC System for duodenal ulcer bleeding during an emergency endoscopy were included in the analysis.

Finally, the data of 173 patients could be evaluated. Of these, 93 patients were treated with type-a and 80 patients with type-t. The analysis showed that the baseline characteristics age (71.2 y vs. 71.6 y, $p = 0.255$), gender (men: 69.9% vs. 67.5%, $p = 0.735$), use of anticoagulants (32.9% vs. 43.0%, $p = 0.176$) and Rockall score (7.2 vs. 7.4, $p = 0.917$) were comparable between the groups. However, there were some significant differences in the bleeding characteristics, which then also influenced the discussion. Thus, the type-a group had significantly fewer ulcers with active bleeding (Forrest Ia/b) than the type-t group (51.1% vs. 62.5%, $p = 0.020$). In addition, type-t was used more often as first-line treatment (95% vs. 77.8%, $p = 0.004$). Analysis of the data showed that initial haemostasis (type-a: 93.5%, type-t: 90%, $p = 0.421$) and bleeding-associated lethality were not significantly different (type-a: 3.2%, type-t: 7.8%, $p = 0.125$). However, the OTSC-t group showed a higher rate of rebleeding (37.2% vs. 6.5%, $p < 0.001$).

The authors concluded that type-a should be considered the standard of care for endoscopic treatment of duodenal ulcer bleeding, if an OTSC is used, because of the lower rate of rebleeding. However, the significant differences within the bleeding characteristics were again pointed out during the discussion. No subgroup analysis was performed in the

evaluation to determine whether bleedings that initially bled actively (Forrest Ia/Ib) also bled more afterwards. Thus, no statement can be made as to whether bleeding that initially bled and was treated more frequently with the type-t clip influenced the higher number of rebleeds. Therefore, further research is needed in this regard.

Analysis of traumatic vs. atraumatic over-the-scope-clip (OTSC) in the treatment of duodenal ulcer bleeding (Analyse des traumatischen vs. atraumatischen over-the-scope-clip (OTSC) bei der Behandlung von duodenalen Ulkusblutungen)

M. Hollenbach1, A. Schmidt2, A. Decker2, O. Möschler3, C. Jung4, N.-C. Mechie5, T. Barhoom5, A. Hegelein5, R. Knoop6, T. Blasberg5, E. Wedi5

1University Hospital Leipzig, Department of Gastroenterology, Medical Clinic II, Leipzig, Germany

2University Medical Centre Freiburg, Medical Clinic II, Freiburg, Germany

3Marienhospital Osnabrück, Clinic for Internal Medicine and Gastroenterology, Osnabrück, Germany

4AUSL Romagna, Forlì-Cesena, Gastroenterology and Gastrointestinal Endoscopy, Forlì-Cesena, Italy

5Sana Klinikum Offenbach, Gastroenterology Clinic, Offenbach, Germany

6University Medical Centre Göttingen, Clinic for Gastroenterology, Göttingen, Germany

April 2021 | HemoPill® acute shows promising results in different use cases: two case studies

Ovesco's first wireless sensor capsule for immediate detection of upper and middle gastrointestinal (GI) bleeding was successfully used in patients with COVID-19 and for the assessment of GI-bleeding localisation.

PD Dr. A. Wannhoff and Prof. Dr. K. Caca, RKH Ludwigsburg clinic, Germany reported their experience with HemoPill acute in patients with suspected gastrointestinal bleeding. During the current SARS-CoV2 pandemic the HemoPill is used to reduce potential risk for endoscopy staff due to viral transmission.

HemoPill examination was performed in 2 patients. Patient #1 suffered from COVID-19, cardiac comorbidities as well as severe obesity. During the inpatient stay, the patient reported tarry stools and the haemoglobin levels dropped.

Patient #2 had pronounced anaemia, and gastrointestinal blood loss was suspected. Routine screening for SARS-CoV2 was performed by PCR test. Due to the delay in receiving the test result, a HemoPill investigation was carried out.

Use in patients with COVID-19

Wannhoff A, Caca K

Case Study, 26 January 2021

Case study #2:

PD Dr Edris Wedi, Sana Klinikum Offenbach, Germany, reports on his experience with HemoPill acute for the detection of acute bleeding or to assess the localisation of the source of bleeding in unclear cases to determine the required endoscopic procedure.

An 86-year-old patient with cardiac comorbidity was referred to the hospital for small intestine diagnostics for suspected moderate GI bleeding on ASA therapy. Externally, a gastro- and colonoscopy had already been performed without finding a source of bleeding. With a capsule endoscopy, the suspicion of lower small intestine bleeding was established. When the patient was admitted to their clinic, they had melaena and hemorrhagic anemia requiring transfusion with an Hb of 7.2 g/dl.

In the underlying case they wanted to evaluate whether the HemoPill acute can be used to assess the localisation of the source of bleeding, and whether the capsule can thus be helpful in selecting the endoscopic procedure.

Assessment of localisation of gastrointestinal bleeding

Wedi E

Case Study, 26 January 2021

April 2021 | FTRD®: international study shows favorable safety profile and efficacy of EFTR in the upper gastrointestinal tract

FTRD application in 56 patients with lesions in the upper GI tract showed a technical success rate of 93 %, a R0-

resection rate of 68 % and only mild to moderate adverse events in 21%.

K. Hajifathalian et al., Weill Cornell Medicine, Division of Gastroenterology and Hepatology, New York, NY, USA performed an international multicenter retrospective study aiming to evaluate efficacy and safety of the FTRD System for resection of lesions in the upper gastrointestinal tract.

56 patients who had undergone EFTR with the FTRD System for lesions of the upper GI tract in 13 centers were included. Patients had a mean age of 61 ± 14 years, ten patients were on antithrombotics before the procedure, one patient was on systemic steroids and one patient on immunosuppressive therapy. Four patients (7 %) had a diagnosis of end-stage renal disease, two patients (4%) a cirrhosis, ten patients (18 %) had diabetes, 7 patients (13 %) had congestive heart failure and coronary artery disease and nine patients (16 %) had a malignancy. Pre-procedural indications for EFTR were mesenchymal neoplasms including gastrointestinal stromal tumor (23.4 %), diagnostic full-thickness biopsy (18 %), initial resection of adenocarcinoma (5.9 %), initial resection of neuroendocrine tumor (5.9 %), resection of residual or recurrent adenoma (4.7 %), resection of residual or recurrent neuroendocrine tumor (4.7 %), resection of residual or recurrent mesenchymal neoplasms including GIST (3.5 %), and initial resection of adenoma (2.4 %). The majority of EFTRs was performed in the gastric antrum (38 %), followed by the gastric body (27 %), gastric cardia and fundus (19 %), duodenum (14 %) and esophagus (2 %). The average size of the lesions was 14 mm (range 3 – 33mm). Application of the FTRD was technically successful in 93 % of patients (52/56) leading to complete and partial resection in 43 (77%) and 9 (16 %) patients, respectively. The overall R0-resection rate was 68 % (38/56). A total of 12 (21 %) mild or moderate adverse events occurred. Follow-up endoscopy was reported for 31 patients (55%), on average 88 days after the initial procedure (IQR 68 – 138 days). 30 patients (97 %) did not have any residual or recurrent lesion on endoscopic examination and biopsy in follow-up, while in one patient (3 %) a residual adenoma was found.

The authors concluded that this international experience suggests that application of the FTRD in the upper GI tract is related to a high technical success rate, an acceptable R0-resection rate, and a low risk of adverse events or early recurrence.

Full-thickness resection device (FTRD) for treatment of upper gastrointestinal tract lesions: the first international experience

Hajifathalian K, Ichkhanian Y, Dawod Q, Meining A, Schmidt A, Glaser N, Vosoughi K, Diehl DL, Grimm IS, James T, Templeton AW, Samarasekera JB, El Hage Chehade N, Lee JG, Chang KJ, Mizrahi M, Barawi M, Irani S, Friedland S, Korp P, Aadam AA, Al-Haddad M, Kowalski TE, Smallfield G, Ginsberg GG, Fukami N, Lajin M, Kumta NA, Tang S, Naga Y, Amateau SK, Kasmin F, Goetz M, Seewald S, Kumbhari V, Ngamruengphong S, Mahdev S, Mukewar S, Sampath K, Carr-Locke DL, Khashab MA, Sharaiha RZ.

March 2021 | Analysis of large nationwide FTRD® registry shows 80 % R0-resection rate of difficult adenomas, early carcinomas and subepithelial tumors

Evaluation of the “German colonic FTRD registry” comprising 1178 colorectal FTRD procedures from 65 centers of different care levels demonstrated technical success in 88.2 % and R0-resection in 80.0 % of all cases. To date, this is by far the largest study of colorectal EFTR using FTRD.

B Meier et al., Ludwigsburg Hospital, Ludwigsburg, Germany, analyzed the German colonic FTRD registry to evaluate efficacy and safety of EFTR with the FTRD System in clinical routine.

Overall, 1178 colorectal FTRD procedures from 65 centers were included in the analysis. Median patient age was 69 years; lesions were located throughout colon (76.5 %) and rectum (23.5 %). 54.1 % of lesions were pretreated endoscopically. Main indication for EFTR was “difficult adenoma” in 67.1 %, defined as adenoma not amenable to conventional endoscopic resection, such as non-lifting lesions or lesions located at difficult anatomic locations. Other indications for EFTR were early carcinomas (18.4 %),

subepithelial tumors (6.8 %), no diagnostic EFTR (1.3 %). Mean lesion diameter was 15 mm, median procedure time was 35 minutes.

In the total cohort, 88.2 % of EFTR procedures were technically successful (complete macroscopic resection); full-thickness resection (visibility of all layers of the colonic wall including serosa within the resection specimen) was histologically confirmed in 89.9 % of all cases. R0-resection was achieved in 80.0 % of all cases. There was no significant difference for R0-resection between colon and rectum (78.9 % vs 83.6 %, $P = 0.11$) or between lesions < 20 mm and ≥ 20 mm (77.6 % vs 81.0 %, $P = 0.2$). Procedure-related adverse events occurred in 142 patients (12.1 %), thereof bleeding in 84 patients (which could be managed endoscopically in all but one case), perforation in 29 patients (11 managed endoscopically, 18 surgically), and inflammation/ infection in 15 patients (80 % successfully managed with antibiotics). In 3.1 % of all patients adverse events were classified as major events and 2.0 % of patients required surgical treatment due to adverse events. Endoscopic follow-up was available in 58.0 % and showed residual/recurrent lesions in 13.5 %, which could be managed endoscopically in most cases (77.2 %). In total, 96 patients of the total cohort (8.1 %) needed surgical therapy. In 62 patients (5.3 %), oncologic surgical resection was necessary, in 23 patients (2.0 %) surgery was due to adverse events and in 10 patients (0.8 %) surgery was required due to recurrent/residual adenoma, which could not be treated endoscopically.

The authors concluded that the study confirms favorable efficacy and safety of EFTR with the FTRD System in a real-world setting. In more than 90 % of patients with "difficult-to-resect" colorectal lesions not amenable to conventional endoscopic resection surgery could be avoided. The registry data also shows that EFTR with the FTRD has gained broad acceptance in the endoscopic community and has become part of the clinical routine.

Efficacy and Safety of Endoscopic Full-Thickness Resection in the Colorectum: Results From the German Colonic FTRD Registry

Meier B, Stritzke B, Kuellmer A, Zervoulakos P, Huebner GH, Repp M, Walter B, Meining A, Gutberlet K, Wiedbrauck T, Glitsch A, Lorenz A, Caca K, Schmidt A.

March 2021 | Ovesco HemoPill® acute was presented at the International Endoscopy Symposium

HemoPill® acute for gastrointestinal diagnosis and its future clinical routine – an Interview.

Dr. Sebastian Schostek, head of research and development of diagnostic systems at Ovesco Endoscopy AG talks with Prof. Alexander Meining from university hospital Würzburg about the clinical application of the HemoPill® acute. This device is a non-imaging swallowable capsule for the real-time detection of gastrointestinal bleeding.

Prof. Meining sees an important application of the HemoPill® acute for non-critically ill and stable patients, in order to be able to determine the right time for an endoscopy. With this tool, they get a quick answer whether they should do an emergency endoscopy or postpone it to a later date. Additionally, the patients do not need a preparation before application.

March 2021 | OTSC®: Systematic review of 3025 patients demonstrates efficacy in GI hemorrhage, perforation, anastomotic leak and stent fixation at an overall technical success rate of 94.4 %

A systematic review of studies with at least five patients examining OTSC application for any indication in the GI tract found 85 eligible articles with overall 3025 patients. Average clinical success rates were: 86.0 % for hemorrhage, 85.3 % for perforation, 55.8 % for chronic fistula, 72.6 % for anastomotic leaks, 92.8 % for defect closure following endoscopic resection, and 80.0 % for stent fixation. Procedural adverse events were reported in 2.1 %.

N. Bartell et al., University of Rochester Medical Center, Rochester NY, USA performed a systematic review aiming to

determine the clinical success and adverse event rates of OTSC application comprising all indications. A thorough search of the literature was conducted in the Pubmed database for eligible articles. Case reports and case series with less than 5 patients were excluded. Articles were included from January 1, 2007 to January 15, 2020. The following terms were used to perform the literature search: "over-the-scope-clip", "OTSC", "endoscopic fistula closure", "over-the-scope clip bleeding", "stent fixation", and "endoscopic perforation closure".

Technical success was defined as successful deployment of the OTSC clip to the targeted lesion or defect. Clinical success was defined as complete and durable resolution of the respective defect, hemorrhage and/or stent fixation.

A total of 85 articles with overall 3025 patients met the inclusion criteria. Overall technical success rate was 94.4 %, overall clinical success rate was 78.4 %. Per-indication clinical success rates were: 1120/1303 (86.0 %) for GI hemorrhage; 399/468 (85.3 %) for perforation; 347/622 (55.8 %) for fistulae; 284/391 (72.6 %) for anastomotic leaks; 205/221 (92.8 %) for defect closure following endoscopic resection (e.g. following EMR or ESD); and 16/20 (80.0 %) for stent fixation.

Adverse events related to OTSC deployment were only reported in 64 of 85 studies ($n = 1942$ patients), with an overall adverse event rate of 2.1 % (40/1942). Surgical intervention despite OTSC placement was required in 4.7 % of patients ($n = 143/3025$).

The authors concluded that this systematic review confirms the safety and efficacy of the OTSC System in the management of GI hemorrhage, perforations, anastomotic leaks, defects created by endoscopic resections and for stent fixation.

Clinical efficacy of the over-the-scope clip device: A systematic review

Bartell N, Bittner K, Kaul V, Kothari TH, Kothari S

February 2021 | EFTR with the FTRD® System offers alternative to radical surgery for smaller complex colorectal lesion

A prospective multicenter study analyzing the data of the Dutch colorectal eFTR registry comprising 367 EFTR procedures shows a technical success rate of 83.9 % and a R0 resection rate of 82.4 %. Overall adverse events occurred in 9.3 %.

L.W. Zwager et al., University of Amsterdam, Amsterdam, The Netherlands, presented a prospective observational multicenter study aiming to investigate the technical and clinical success, and safety of EFTR with the FTRD System in complex colorectal lesions in current clinical practice. For the study, the Dutch colorectal eFTR registry was analyzed. The registry comprises nationwide gathered prospective data from all consecutive EFTR procedures in 20 hospitals (5 academic and 15 non-academic).

Between July 2015 and October 2018, a total of 367 EFTR procedures were carried out in 362 patients (63 % men, mean age 69 years). Indications for EFTR were: difficult polyps (non-lifting sign and/or difficult location; $n = 133$), primary resection of suspected T1 colorectal cancer ($n = 71$), relapse T1 colorectal cancer ($n = 150$), and subepithelial tumor ($n = 13$). Technical success was achieved in 83.9 % of all procedures. In 21 procedures (5.7 %), EFTR could not be performed because the lesion could not be reached or could not be suctioned into the cap. In the remaining 346 procedures, full-thickness resection was achieved in 83.2 % (288/346) and R0-resection was achieved in 82.4 % (285/346). The median diameter of resected specimens was 23 mm.

Overall adverse events occurred in 34 procedures (9.3 %), 10 patients (2.7 %) required emergency surgery for 2 immediate and 5 delayed perforations (equivalent to an anastomotic leak) and 3 cases of appendicitis (after resection of the basis of the appendix). The authors concluded that EFTR with the FTRD System is an effective and safe en-bloc resection technique for complex colorectal lesions. The technique has the potential to obviate the need for radical surgery of lesions considered incurable with conventional endoscopic resection techniques.

Endoscopic full-thickness resection (eFTR) of colorectal lesions: results from the Dutch colorectal eFTR registry

Zwager LW, Bastiaansen BAJ, Bronzwaer MES, van der Spek

BW, Heine GDN, Haasnoot KJC, van der Sluis H, Perk LE, Boonstra JJ, Rietdijk ST, Wolters HJ, Weusten BLAM, Gilissen LPL, Ten Hove WR, Nagengast WB, Bekkering FC, Schwartz MP, Terhaar Sive Droste JS, Vlug MS, Houben MHMG, Rando Munoz FJ, Seerden TCJ, Beaumont H, de Ridder R, Dekker E, Fockens P; Dutch eFTR Group.

February 2021 | OTSC® placement in the colon in routine endoscopy: administrative data from large health insurance company shows high safety and efficacy.

Colonic OTSC placement was mostly performed in patients undergoing polypectomy (62.4 %), but also for closure of iatrogenic perforations (15.8 %), colonic bleeding (10.1 %) and not further specified other underlying diseases (11.7 %). Surgical interventions within 10 days after OTSC placement were performed in 1.6 % of patients after polypectomy, 13.8 % of patients after iatrogenic perforation, and 7.8 % of patients with colonic bleeding. Only one case of surgical intervention during further follow-up could potentially be traced back to an endoscopic procedure complication.

D. Horenkamp-Sonntag et al., Techniker Health Insurance (Techniker Krankenkasse), Hamburg, Germany, performed a study investigating OTSC use in the colon in routine endoscopy using their own administrative data as a large statutory German health insurer. As application of OTSC during colonoscopy generates a specific code in the German diagnosis-related group (DRG) reimbursement system, respective patient data can be identified and analysed. The health insurance database comprises approximately 10 million insured patients, it contains all procedures performed in hospitals as well as in outpatient care. The results of the study thereby can be considered as representing German routine clinical practice in terms of indications, effectiveness, and complications of endoscopy of the cases treated with OTSC in the colon.

Overall 505 patients undergoing colonoscopy with simultaneous OTSC application, but not undergoing upper GI endoscopy on the same day, were identified. According to indications for OTSC clipping, the cohort was divided into 4 groups: patients with iatrogenic perforations ($n = 80$; Group A), patients undergoing polypectomy ($n = 315$; Group B), patients with colonic bleeding ($n = 51$; Group C), and patients with various underlying diseases ($n = 59$; Group D).

13.8 % of patients of group A ($n = 11$) underwent an operative procedure during the short follow-up period of 10 days after clipping, in 9 of these cases surgery was carried out within 24 hours after clipping, in 8 patients laparoscopic suturing of the defect was the only operative procedure, in 3 patients surgery included colonic resection. In Group B early colonic surgery within 10 days after clipping was necessary in 5 patients (1.6 %). Surgical interventions included resection in 3 cases (all of which had colorectal cancer), laparoscopic suturing in one case, and diagnostic laparoscopy without any further procedure in one patient. In Group C, 7.8 % of patients ($n = 4$) had to undergo surgical resections (persistent bleeding $n = 1$, colorectal cancer $n = 2$), while in 6 patients early repeat colonoscopy was performed for recurrent bleeding. During further follow-up (days 11 – 30), 17 patients underwent resection for colonic neoplasms ($n = 12$) or persistent or recurrent bleeding ($n = 4$). Only in one patient (Group A), a plausible cause for surgery could not be identified and might represent delayed insufficiency of the endoscopic closure site.

The authors concluded that OTSC placement in the colon is safe and effective in clinical routine.

Use of over-the-scope clips in the colon in clinical practice: results from a German administrative database

Horenkamp-Sonntag D, Liebentraut J, Engel S, Skupnik C, Albers D, Schumacher B, Koop H.

January 2021 | Market approval for LiftUp® in Saudi Arabia

Ovesco is proud to announce that LiftUp® has received market approval in Saudi Arabia. It is an injection agent for endoscopic resection techniques such as ESD, EMR and

POEM. Due to its thermoreversible properties it becomes a gel at body temperature and forms a stable and long-lasting cushion in the submucosa.

Preclinical evaluation shows that LiftUp is an effective alternative to common injection solutions¹ and could fasten up endoscopic resection because of its long-lasting and pronounced lifting effect². Together with our procedural concepts EMR+ and ESD+ it delivers an optimized endoscopic resection technique. Further information on LiftUp can be found here: [Flyer LiftUp](#)

1 Endoscopic submucosal dissection with a novel high viscosity injection solution (LiftUp) in an ex vivo model: A prospective randomized study

Wedi E, Koehler P, Hochberger J, Maiss J, Milenovic S, Gromski M, Ho C, Conrad G, Baulain U, Ellenrieder V, Jung C. *Endosc. Int. Open.* 2019 07(05). doi: 10.1055/a-0874-1844

2 Preclinical evaluation of a novel thermally sensitive copolymer (LiftUp) for endoscopic resection.

Wedi E, Ho C, Conrad G, Weiland T, Freidinger S, Wehrmann M, Hochberger J

January 2021 | Coag Dissector: preclinical evaluation of Cold-ESD shows that electrocautery-free dissection with the Coag Dissector is feasible and safe

Cold-Endoscopic Submucosal Dissection (C-ESD) was performed in live porcine models. 12 C-ESDs were performed: 4 in the esophagus, 4 in the colon and 4 in the rectum. 5 / 8 (62 %) colorectal dissections were performed completely free of electrocautery. All esophageal cases needed at least an initial incision with HybridKnife to gain submucosal access, the subsequent dissection however could be completed with cold technique in 2 / 4 cases in the esophagus. Two perforations occurred in the colon, both during circumferential incision with the biopsy forceps, there were no perforations with the dissector.

A. Parra-Blanco and M. Fraile-López, Nottingham University Hospitals NHS Trust, Nottingham, UK, conducted a preclinical proof-of-concept study aiming to evaluate the feasibility and safety of cold-endoscopic submucosal dissection (C-ESD) in an in vivo porcine model. In C-ESD, the steps of endoscopic submucosal dissection are in principle performed without electrocautery. Circumferential or partial incision (depending on the technique applied and the location and the organ) is performed with a biopsy forceps. Submucosal dissection is carried out using the Coag Dissector (Ovesco, Tuebingen), an endoscopic Maryland dissector which allows blunt dissection of the submucosa when opening and closing the clamp. Besides, large vessels can be pre-coagulated with the Coag Dissector.

Overall, 12 dissections were performed: 4 in the esophagus, 4 in the colon and 4 in the rectum. Feasibility of C-ESD was tested with different ESD techniques: Conventional, pocket-creation method and tunneling. Usual traction methods were applied. Median dissection times, including circumferential incision and dissection for each location were 62, 64.5 and 88 minutes, respectively.

5 / 8 (62 %) colorectal dissections were performed completely free of electrocautery. All esophageal cases needed at least an initial incision with Hybrid-Knife to gain submucosal access, the subsequent dissection however could be completed with cold technique in 2 / 4 cases in the esophagus.

Coagulation with the Coag Dissector was applied prophylactically in vessels larger than 1 mm and in cases of minor bleeding. No major bleeding was encountered. Two perforations occurred in the colon, both during circumferential incision with the biopsy forceps. The authors concluded that full C-ESD is feasible in the colorectum, whereas a small hot incision is needed in the esophagus. However, in 50 % of the colonic cases, there were perforations caused by the biopsy forceps making the circumferential incision. Therefore, precautions in feasibility studies of colonic C-ESD in humans should be taken and some modifications in the design of the biopsy forceps would be desirable.

Is it time for Cold-Endoscopic Submucosal Dissection? A feasibility study in an esophageal and colorectal live porcine model

Parra-Blanco A, Fraile-López M.

January 2021 | OTSC® in first-line treatment of NVUGIB: RCT shows significant benefits of OTSC in rates of rebleeding, severe complications and blood transfusions

53 patients with NVUGIB were randomized to OTSC hemostasis or standard endoscopic treatment. Rebleeding occurred in 4 % in OTSC group vs 28.6 % in standard group (p = 0.017), severe complications did not occur in OTSC group, but in 14.3 % of patients with standard endoscopic treatment, 0.04 units of red blood cell transfusions were needed in patients treated with OTSC post-randomization vs. 0.68 units in standard therapy group.

D.M. Jensen et al., CURE Digestive Diseases Research Core Center, Los Angeles, CA, USA, presented a randomized controlled trial (RCT) on initial endoscopic treatment of severe nonvariceal upper gastrointestinal bleeding (NVUGIB) with OTSC compared to standard endoscopic hemostasis.

53 patients with bleeding ulcers or Dieulafoy's lesions and major or lesser stigmata of hemorrhage were randomized to OTSC (n = 25) or endoscopic standard therapy (n = 28). Patients and their healthcare providers were blinded to treatments and made all postrandomization management decisions. Major stigmata of recent hemorrhage were defined as active spurting or pulsatile arterial bleeding, non-bleeding visible vessel, or adherent clot. Lesser stigmata of recent hemorrhage were defined as oozing bleeding without clot or visible vessel or flat spot with arterial blood flow underneath, detected by Doppler probe. Standard endoscopic hemostasis was performed with hemoclips or multipolar electrocoagulation. After randomization all patients received high dose intravenous infusions of PPI (proton pump inhibitors) for 72 hours, followed by twice daily oral PPI for 27 days in patients with peptic ulcer bleedings.

Baseline patient characteristics and endoscopic findings were similar between the OTSC and the standard therapy group. However, patient outcomes showed significant differences in OTSC vs standard group in rates of rebleeding (4 % vs. 28.6 %; p = 0.017; relative risk 0.10, 95 % CI 0.01 – 0.91; number needed to treat 4); severe complications (0 % vs. 14.3 %); and post-randomization units of red cell transfusions (0.04 vs. 0.68). All rebleeds occurred in patients with major stigmata of hemorrhage and none in patients with lesser stigmata of hemorrhage.

The authors concluded that OTSC significantly reduces rates of rebleeding, severe complications, and post-randomization red cell concentrate transfusions. Patients with major stigmata of hemorrhage benefit significantly from hemostasis with OTSC. Based upon these RCT results and those of the RCT by Schmidt A et al. (*Gastroenterology* 2018), current guidelines for standard endoscopic hemostasis of severe NVUGI hemorrhage should be reevaluated and updated.

Randomized Controlled Trial of Over-the-Scope-Clip as Initial Treatment of Severe Non-Variceal Upper Gastrointestinal Bleeding

Jensen DM, Kovacs TOG, Ghassemi KA, Kaneshiro M, Gornbein J.

Ovesco Research Update 37

December 2020 | Conference of UEGW 2020

Conference Report of the United European Gastroenterology Week (UEGW) 2020:

- **HemoPill® acute detects small bowel bleeding rapidly and reliable**
- **OTSC® is superior over TAE as salvage therapy in refractory peptic ulcer bleeding in terms of ICU stay and in-hospital mortality**
- **EFTR with the FTRD® is a safe and effective method for treating malignant or difficult colonic lesions**

The 28th United European Gastroenterology Week (UEGW) was taking place digitally on October 11-13, 2020. Ovesco technology and procedures were presented in talks, posters and a live-broadcasted endoscopy.

HemoPill

First application experience with the telemetric capsule HemoPill acute shows fast and reliable haemorrhage proof

T. Brunk et al., Vivantes Hospital in the Friedrichshain, Berlin, Germany, presented a study on the application of the HemoPill acute in case of suspected small bowel bleeding. 13 patients (5 female, 8 male, age 28 – 84 years, Glasgow-Blatchford-Score 6 – 12 (median 10, SD 2)) with acute GI bleeding and without findings in esophagogastroduodenoscopy (EGD) were included in the study.

The HemoPill acute is a sensor-based telemetric capsule for detection of haemorrhage in the upper gastrointestinal tract. It determines photometrically an intensity ratio of violet and red light and calculates the so-called HemoPill indicator (HI). The measurement is independent from food components. It is visualized in real time via radio on a mobile HemoPill Receiver. HI of 1.0 – 1.5 indicates fresh blood/ haematin, other gastrointestinal content inclusively bilirubin leads to HI of 0.5 – 0.9.

Application of the HemoPill acute was performed via swallowing in 9 cases and via endoscopic placement into the duodenum in 4 cases. The application was technically successful in all cases, no complications occurred. In 7/13 cases the HI was ≥ 1.0 (median 1.4; SD 0.7) and indicated bleeding of the small bowel. In these patients double-balloon-enteroscopy was performed within 24 h (median 22 h). Angiodysplasia was detected as probable bleeding source and treated with Argon-Plasma-coagulation. In all other 6 cases with HI < 1 standard elective evaluation of the colon and small intestine was performed. Signs of active bleeding were not found. In 1 case video-capsule-enteroscopy provided evidence of a not-bleeding jejunal ulcer.

The authors concluded that based on this preliminary data the HemoPill acute rapidly and reliably detected active small bowel bleeding in patients with suspected acute GI bleeding and negative EGD. Further comparative investigation is needed to better define the value of this promising new system.

HemoPill acute: preliminary results using a sensor based telemetric capsule system in patients with EGD negative acute mid GI-bleeding

Brunk T, Tauchmann C, Berger AW, Hochberger J, Berlin, Germany.

OTSC System

OTSC superior to TAE in refractory peptic ulcer bleeding – study shows significantly lower in-hospital mortality and shorter ICU stay

A. Kuellmer et al., Medical Center University of Freiburg, Freiburg, Germany, presented a retrospective multicenter study comparing OTSC vs TAE (transarterial angiographic embolization) as salvage therapy for refractory peptic ulcer bleeding.

Primary endpoint of the study was clinical success defined as the combined endpoint of successful hemostasis and no re-bleeding within 7 days. Secondary endpoints were adverse events, length of hospital stay, days on intensive care unit (ICU), number of blood transfusions and mortality. Statistical analysis was performed for the total cohort and a matched cohort after adjustment of differences in baseline characteristics with propensity score matching (PSM).

Overall, 128 patients with peptic ulcer bleeding refractory to standard endoscopic therapy were included in the study. 66 patients were treated with OTSC, 62 patients with TAE.

Between the two comparison groups there were no significant differences regarding age, Charlson comorbidity index, Rockall score, Helicobacter pylori status, ongoing anticoagulation, NSAID intake, primary hemostasis rate in first line therapy and number of endoscopic treatment attempts before salvage therapy. Also, in both groups, the proportion of patients with ulcer size > 20mm was similar (27.3 % vs. 33.9 %, p = 0.48). Most ulcers were in the duodenal bulb (65 % in OTSC group; 85.5 % in the TAE group; p = 0.014). The OTSC group included significantly less Forrest Ia bleedings (19.7 % vs. 38.7 %, p = 0.02) and significantly more Forrest Ib bleedings (63.6 % vs. 43.5 %; p = 0.03). PSM was performed to control for these biases and resulting in treatment groups of n = 40 each, with no significant

differences in ulcer localization and bleeding characteristics. Clinical success was higher in the OTSC group but did not reach statistical significance (72.5 % vs. 62.5 %; $p = 0.474$) while TAE patients stayed significantly longer in ICU (4.9 vs. 9.2 days, $p = 0.009$) and in-hospital mortality was significantly higher in the TAE group (5.0 % vs. 22.5 %, $p = 0.048$). The 7-day re-bleeding rate was higher in the TAE group (17.5 % vs. 32.5 %; $p = 0.196$). Also, severe adverse events occurred more often in the TAE group (3.0 % vs. 7.5 %, $p = 0.308$). The authors concluded that OTSC treatment for refractory peptic ulcer bleeding shows at least similar efficacy compared to TAE, but significantly lower mortality rates and significantly shorter ICU stay.

OTSC vs TAE as salvage therapy for refractory peptic ulcer bleeding

Kueller A¹, Mangold T¹, Bettinger D¹, Maruschke L¹, Wannhoff A², Caca K², Edris W³, Jung C⁴, Kleemann T⁵, Schulz T⁵, Thimme R¹, Schmidt A¹. ¹Freiburg, Ludwigsburg, Offenbach³, Goettingen⁴, Cottbus⁵, Germany.

The OTSC is established as sole second line treatment for recurrent peptic ulcer bleeding

K. Caca, Hospital of Ludwigsburg and University Heidelberg, Germany, gave a lecture on new tools for the management of upper GI bleeding.

First-line endoscopic hemostatic therapy achieves very high success rates of about 90 %, independent from which hemostasis technique is used: injection techniques, thermal coagulation or conventional clips. However, the mortality rate of acute upper GI bleeding is 5.8 % (Lau et al. Lancet 2012). Predictors for adverse outcome are the patient's age and comorbidities, and especially re-bleeding (Chiu et al. Clin Gastroenterol Hepatol 2009). In case re-bleeding occurs, re-endoscopy is still superior to surgery concerning complications and associated with a success rate of about 75 % (Lau et al., NEJM 1999). But if re-endoscopy is not successful and surgical salvage therapy is necessary, the mortality rate is quite high, it rises to 29 % (Jairath et al. B J Surg 2012).

Over the past few years, new endoscopic hemostasis techniques have been developed to increase success rates of secondary endoscopic hemostasis. These are the so-called "topic substances" (Hemospray, EndoClot and PuraStat) and the OTSC.

There are no randomized studies evaluating topic substances; observational studies show that they can stop haemorrhage even in diffuse bleedings in nearly 95 %, but quite high rebleeding rates of about 30 – 40 % occur. So, these techniques are established as rescue therapy or "add-on" therapy especially for diffuse bleedings.

The OTSC however is established as the sole second line treatment for recurrent peptic ulcer bleeding. It overcomes typical problems of through-the-scope clips, which can hardly grasp a centric vessel in a large ulcer with fibrotic base and cannot be well applied in ulcers at the posterior duodenal wall or the duodenal knee because of the tangential position and the narrow distance. The OTSC allows a better visualization of the ulcer due to suction into the cap of the device and/or the possibility of grasping by an OTSC Anchor. Besides, OTSC has the advantage of higher compression force (as high as a surgical seal) and thereby better hold in fibrotic tissue.

The STING trial (Schmidt A et al. Gastroenterology 2018) is a randomized controlled trial in patients with recurrent peptic ulcer bleeding ($n=66$, 33 OTSC, 33 standard endoscopic hemostasis). It showed a relative risk-reduction of 73.6 % in the OTSC group compared to the group with standard endoscopic therapy. 14 patients crossed over to the OTSC group after failed standard endoscopy (10 patients) or rebleeding after primary successful standard therapy (4 patients), so no significant differences could be found in rates of angiographic embolization, surgery, transfusion requirements and mortality.

Regarding the role of OTSC in first line therapy, there are no prospective trials so far; a retrospective study (FLET Rock study, Wedi E et al., Surg Endosc 2018) evaluating OTSC "first line" vs. a matched control group, showed that observed mortality, re-bleeding, and mortality after re-bleeding were significantly lower with OTSC as first line therapy.

A study prospectively evaluating OTSC as first line therapy in 100 patients with acute NVUGIB and Rockall Score ≥ 7

(STING II) is ongoing.

New tools for treatment of upper GI bleeding. Therapy update: Non-variceal upper GI bleeding

Caca K, Ludwigsburg, Heidelberg, Germany.

FTRD System

EFTR is a safe and effective method for treating malignant or difficult colonic lesions

S. Sferazza et al., Santa Chiara Hospital, Trento, Italy, presented a prospective cohort study evaluating endoscopic full-thickness resection with the FTRD System for colorectal lesions unsuitable or difficult to remove both with standard and advanced endoscopic resection techniques.

20 subsequent patients of two tertiary referral centers were included in the study. 70 % of patients were male, median age of the patients was 71.5 years. Indications for EFTR were malignant features (40 %), recurrence or non-lifting sign (50 %) and intra-diverticular or intra-appendicular location (10 %). 16 % of procedures were performed in an outpatient setting. No differences were found in patient or procedural characteristics between inpatient or outpatient settings. 25 % of lesions were located in the rectum, 15 % in the sigmoid, 15 % in the descending colon, 20 % in the ascending colon, 20 % in the cecum and 5 % in a surgical anastomosis. Technical success was achieved in 95 % of the lesions. Median procedural time was 15 minutes (IQR 15 – 20). Mean lesion size was 19 mm (9 – 40 mm). No immediate peri-procedural complications occurred. One patient developed acute appendicitis requiring surgery. Median follow-up was 5 months (IQR 0-17 months). During follow-up, 2 recurrences were found, both were < 10 mm in size and removed endoscopically.

The authors concluded that endoscopic full-thickness resection is a safe and effective method for treating malignant or difficult colonic lesions. The procedure is not time-consuming and can be performed by expert endoscopists with minimal prior experience with the device.

Endoscopic full-thickness resection for the management of difficult colorectal lesions: a prospective cohort study

Sferazza S¹, Vieceli F¹, Maida MF², Michielan A¹, Conte E³, Di Mitri R³, de Pretis G¹, ¹Trento, ²Caltanissetta, ³Palermo, Italy.

Appendicitis as a delayed complication after EFTR

M. Raitheil et al., Maltese Wood Hospital, Erlangen, Germany, presented a case report of perforated appendicitis as a delayed complication 3 months after successful EFTR.

A 71-year-old woman with family history positive for colorectal carcinoma and > 20 mm hyperplastic polyp at the appendix opted for EFTR for resection after careful clarification and discussion of endoscopic or surgical resection. The endoscopic resection was performed with the FTRD after chromo-endoscopic visualization, magnification with near focus view and marking of the lesion. The lesion was mobile, could be well grasped and the clip safely placed around it. However, because of incomplete EFTR, further polyp tissue had to be removed with a standard polypectomy snare to achieve complete resection. At the end of the procedure the lesion was macroscopically completely resected and a clear coagulation zone present. The histological preparation showed a R0 resection of a sessile serrated adenoma at and into the appendiceal orifice, which was completely resected. The patient recovered quickly without systemic signs of inflammation and was discharged from hospital 4 days later. At day 92 after EFTR the patient presented with right sided abdominal discomfort, fever, nausea and vomiting. Also, laboratory and sonographic examinations were consistent with acute appendicitis. Emergency laparoscopic appendectomy had to be performed, showing a perforated appendix, local hemorrhagic peritonitis with complete adhesion and scarring of the appendix and the cecum. The postop course was uneventful, the patient left the hospital at day 4 and is doing well.

Perforated appendicitis as a delayed complication 3 months after successful endoscopic full-thickness resection (EFTR) with the over-the-scope clip device (OTSC): Case report and assessment of recent studies involving the appendix

Raitheil M¹, Roßmeißl A¹, Fezehaye J¹, Allescher HD², Braun A¹, Bruehlting A¹, Albrecht H³, ¹Erlangen, ²Garmisch-

Partenkirchen, ³Nuernberg, Germany.

Live EFTR of a cecal non-lifting lesion

E. Dekker, Academic Medical Center, Amsterdam, the Netherlands, moderated the UEG Live Demo II Session with I.M. Gralnek, Afula, Israel, and R. Bisschops, Leuven, Belgium.

The case presented was a 75-year-old female patient with a small non-lifting lesion in the cecum detected during screening colonoscopy. The lesion was suspected for early cancer. The colonoscopy images presented showed a small sessile polyp with non-lifting sign 1 cm in size. For this type of lesion endoscopic full-thickness resection is a very suitable therapeutic option.

The procedure was performed by B. A. Bastiaansen, P. Fockens and team in Amsterdam (NL) and broadcasted live. The lesion presented itself as very vulnerable and irregular lesion with partly granular structure and in other parts complete loss of mucosal structure. The endoscopist estimated a deep submucosal invasion and stated that ESD in this case would be impossible, because the risk of adverse events would be very high, and it would not be radical enough while EFTR would yield a safe full-thickness specimen with good assessability by the pathologist.

Advancement of the endoscope with mounted FTRD to the lesion was successful, the grasper was used to grasp as much tissue as possible, gentle suction was applied and the lesion was pulled slowly into the cap. The endoscopist emphasized the importance of having a good visualization of the white ring and watching the white ring moving forward when setting the clip application mechanism in motion to ensure that the clip has been deployed and no perforation is caused by the subsequent resection. After attainment of a good imaging, the clip was released, the lesion resected, and the endoscope was pulled out.

The specimen was about 3.5 cm with a lesion a bit larger than expected, the lateral margin was nice and clear at the ventral side and quite tight but sufficient at the dorsal side.

Evaluating the final endoscopic result, all four edges of the clip were clearly visible, tissue assessment of the margins was difficult but did not hint any conspicuities. Histology of the specimen showed a resection in R0.

UEG Live Demo Live Endoscopy II

Dekker E, Amsterdam, the Netherlands

December 2020 | OTSC® superior to TAE in refractory peptic ulcer bleeding – study shows significantly lower in-hospital mortality and shorter ICU stay

Compared to patients treated with transcatheter arterial embolization (TAE) for refractory peptic ulcer bleedings, OTSC showed higher efficacy (72.5% vs. 62.5% ($p = 0.474$)) and significantly reduced length of stay in ICU (4.9 vs. 9.2 days $p = 0.009$) as well as a significantly reduced in-house mortality (5.0% vs. 22.5% ($p = 0.048$)).

A. Küllmer et al., Medical Center University of Freiburg, Freiburg, Germany, presented a retrospective multicenter study comparing OTSC vs TAE (transarterial angiographic embolization) as salvage/bailout therapy for refractory peptic ulcer bleeding at this year's virtual UEG Week. Results of the study were also presented at this year's DGVS/DGAV conference (conference of the German Society of Gastroenterology, and Digestive and Metabolic Diseases with Endoscopy Section (DGVS) and the German Society of General Surgery and Visceral Surgery (DGAV)).

Primary endpoint of the study was clinical success defined as the combined endpoint of successful hemostasis and no re-bleeding within 7 days. Secondary endpoints were adverse events, length of hospital stay, days on intensive care unit (ICU), number of blood transfusions and mortality. Statistical analysis was performed for the total cohort and a matched cohort after adjustment of differences in baseline characteristics with propensity score matching (PSM). Overall, 128 patients with peptic ulcer bleeding refractory to standard endoscopic therapy were included in the study. 66 patients were treated with OTSC, 62 patients with TAE. Between the two comparison groups there were no significant differences regarding age, Charlson comorbidity index, Rockall score, Helicobacter pylori status, ongoing anticoagulation, NSAID intake, primary hemostasis rate in

first line therapy and number of endoscopic treatment attempts before salvage therapy. Also, in both groups, the proportion of patients with ulcer size > 20mm was similar (27.3 % vs. 33.9 %, $p = 0.48$). Most ulcers were in the duodenal bulb (65 % in OTSC group; 85.5 % in the TAE group; $p = 0.014$). The OTSC group included significantly less Forrest Ia bleedings (19.7 % vs. 38.7 %, $p = 0.02$) and significantly more Forrest Ib bleedings (63.6 % vs. 43.5 %, $p = 0.03$). PSM was performed to control for these biases and resulting in treatment groups of $n = 40$ each, with no significant differences in ulcer localization and bleeding characteristics. Clinical success was higher in the OTSC group but did not reach statistical significance (72.5 % vs. 62.5 %, $p = 0.474$) while TAE patients stayed significantly longer in ICU (4.9 vs. 9.2 days, $p = 0.009$) and in-hospital mortality was significantly higher in the TAE group (5.0 % vs. 22.5 %, $p = 0.048$). The 7-day rebleeding rate was higher in the TAE group (17.5 % vs. 32.5 %, $p = 0.196$). Also, severe adverse events occurred more often in the TAE group (3.0 % vs. 7.5 %, $p = 0.308$). The authors concluded that OTSC treatment for refractory peptic ulcer bleeding shows at least similar efficacy compared to TAE, but significantly lower mortality rates and significantly shorter ICU stay.

OTSC vs TAE as salvage therapy for refractory peptic ulcer bleeding

Kueller M, A, Mangold T, Bettinger D, Maruschke L, Wannhoff A, Caca K, Edris W, Jung C, Kleemann T, Schulz T, Thimme R, Schmidt A.

United European Gastroenterology Week Virtual 2020, Congress lecture, 2020 October 11

November 2020 | Conference of DGVS and DGAV 2020

Conference of the German Society of Gastroenterology, and Digestive and Metabolic Diseases with Endoscopy Section (DGVS) and the German Society of General Surgery and Visceral Surgery (DGAV):

- HemoPill® – Swallowable sensor capsule proves small bowel bleeding fast and reliable.
- FTRD® System – Endoscopic full-thickness resection with the FTRD closes an important gap in the endoscopic treatment of mucosal and submucosal lesions.
- OTSC® System – The OTSC is a well-established method of treatment for acute perforations.

The 75th annual conference of the German Society of Gastroenterology, and Digestive and Metabolic Diseases with Endoscopy Section (DGVS) took place together with the 14th autumn conference of the German Society of General Surgery and Visceral Surgery (DGAV) on September 16 – 19, 2020, originally planned to be held in Leipzig, Germany, due to the SARS-CoV-2 pandemic held in digital form, .

HemoPill

First application experience with the telemetric capsule HemoPill acute in 13 patients with acute GI bleeding and negative EGD

T. Brunk et al., Vivantes Hospital in the Friedrichshain, Berlin, presented a study on the application of the HemoPill acute in case of suspected small bowel bleeding. 13 patients (5 female, 8 male, age 28 – 84 years, Glasgow-Blatchford-Score 6 – 12 (median 10, SD 2)) with acute GI bleeding and without findings in esophagogastroduodenoscopy (EGD) were included in the study.

The HemoPill acute is a sensor-based telemetric capsule for detection of acute bleedings in the upper gastrointestinal tract. It photometrically determines an intensity ratio of violet and red light and calculates the so-called HemoPill indicator (HI). The measurement is independent from food components. It is visualized in real time via radio on a mobile HemoPill receiver. HI of 1.0 – 1.5 indicates fresh blood/haematin, other gastrointestinal content inclusively bilirubin leads to HI 0.5 – 0.9.

Application of the HemoPill acute was performed via swallowing in 9 cases and via endoscopic placement into the duodenum in 4 cases. The application was technically successful in all cases, no complications occurred. In 7/13

cases the HI was ≥ 1.0 (median 1.4; SD 0.7) and indicated bleeding of the small bowel. In these patients double-balloon-enteroscopy was performed within 24 h (median 22 h). Angiodysplasia was detected in all 7 cases, 3 angiodysplastic lesions showed active oozing. All angiodysplastic lesions were treated with Argon-Plasma-Coagulation. Post-interventional bleeding did not occur in any case. In one of the 6 cases with HI < 1 a subsequent video-capsule-enteroscopy provided evidence of a non-bleeding jejunal ulcer.

The authors concluded that based on this preliminary data the HemoPill acute detects small bowel bleeding fast and reliable. The diagnostic-therapeutic sequence of HemoPill acute and double-balloon-enteroscopy could constitute an interesting, accelerated treatment pathway for small bowel bleeding. Further comparative investigation is needed to closer define the role of and indications for the HemoPill acute.

HemoPill acute – First application experience with a sensor-based telemetric capsule for fast detection of bleeding of the small bowel after negative EGD
(HemoPill acute – Erste Anwendungserfahrungen mit einer sensorbasierten, telemetrischen Kapsel zum raschen Nachweis von Dünndarmblutungen nach negativer ÖGD)

Brunk T¹, Tauchmann C¹, Berger AW¹, Hochberger J¹, ¹Berlin.

FTRD System

Evaluation of the German FTRD Register (1178 cases) shows a R0-resection rate of 80.0 % and major complications in 3.1 %

B. Meier, Hospital of Ludwigsburg, presented an evaluation of the German FTRD register. It comprised data on 1178 FTRD applications in 65 centers. Indications for endoscopic full-thickness resection with the FTRD System were: "difficult adenomas" (i.e. non-lifting sign, difficult anatomic localisation) in 67.1 % of cases, early carcinomas in 18.4 %, subepithelial tumors in 6.8 % and resection for diagnostic purposes in 1.3 %. Average lesion size was 15 x 15 mm. 54.1 % of the lesions were pre-treated endoscopically.

The rate for technical success (resection macroscopically complete) was 88.2 %, R0-resection could be confirmed in 80.0 % of all cases. No significant differences were found comparing the R0-resection-rates of smaller vs larger lesions neither in colonic vs rectal lesions. Complications occurred in 12.1 % of cases, in 3.1 % complications were defined as major. Consecutive surgical intervention was performed in 2.0 % of interventions. During follow-up, residual/recurrent lesions were found in 13.5 %, these could be treated endoscopically in 77.2 %.

The authors concluded that this study represents the currently largest study on endoscopic full-thickness resection with the FTRD. It confirms the good effectivity and safety of the FTRD in the real world setting for the resection of difficult lesions in the lower gastrointestinal tract.

Effectivity and safety of endoscopic full-thickness resection in the lower gastrointestinal tract: Results of the German FTRD Register.

(Effektivität und Sicherheit der endoskopischen Vollwandresektion im unteren Gastrointestinaltrakt: Ergebnisse des deutschen FTRD Registers.)

Meier B¹, Stritzke B², Kueller M³, Zervoulakos P⁴, Huebner GH⁵, Repp M⁶, Walter B⁷, Meining A⁸, Gutberlet K⁹, Wiedbrauck T¹⁰, Glitsch A¹¹, Lorenz A¹², Caca K¹, Schmidt A³.
¹Ludwigsburg, ²Tuebingen, ³Freiburg, ⁴Koeln, ⁵Halle a. d. Saale, ⁶Altenburg, ⁷Ulm, ⁸Wuerzburg, ⁹Delmenhorst, ¹⁰Duisburg, ¹¹Greifswald, ¹²Berlin.

EFTR is a safe and fast resection method for small well-differentiated rectal NET

M. Brand, University Hospital Wuerzburg, presented a study which compares trans-anal endoscopic microsurgery (TEM) and endoscopic full-thickness resection (EFTR) for treatment of small well-differentiated neuroendocrine rectal tumors. Neuroendocrine tumors (NET) of the rectum are rare and mostly discovered during colonoscopy. For well-differentiated tumors without risk factors (G1/2, Ki67 < 10 %, < 10 mm, pT1, LO) local resection is the treatment of choice. It can be performed by transanal endoscopic microsurgery (TEM) or endoscopic procedures (EMR, ESD, EFTR). While complete resection of submucosal lesions by EMR is often not possible, EFTR completely resects also tumors with deep submucosal

infiltration.

For the study, two patient collectives with transanal resection of rectal NET were analysed retrospectively. These were overall 28 patients. Patient collective 1 included 13 patients which had undergone TEM under general anaesthesia, patient collective 2 included 15 patients which had undergone endoscopic full-thickness resection with the FTRD in sedation. In the EFTR collective R0 resection was achieved in all interventions, in the TEM collective one resection was incomplete (R1). In the TEM group the size of the resectates and the diameter of the lesion were tending to be larger (resectate area 2.9 ± 2.2 vs 2.4 ± 1.5 cm²; lesion diameter 6.7 ± 4.2 vs 4.6 ± 2.2 mm). The procedural time was significantly shorter in the EFTR group (19.2 ± 4.2 min vs 48.9 ± 29.1 min, $p < 0.01$). The technical success rate was 100 % in both groups, no peri-procedural complications occurred. In all follow-up examinations of EFTR patients the OTSC Clip had spontaneously detached.

The authors concluded that EFTR represents a safe and fast resection method for small well-differentiated rectal NET and shows various advantages over TEM.

Economic analysis of the endoscopic full-thickness resection (Ökonomische Analyse der endoskopischen Vollwandresektion).

Rathmayer M, Schmidt A, Schepp W, Heinlein W, Albert JG, Gundling F. Munich, Freiburg, Stuttgart.

Live-broadcast of a lecture by K. Caca on basis and limitations of endoscopic resection techniques

K.Caca, Hospital of Ludwigsburg, lectured on basis and limitations of endoscopic resection techniques. At the beginning, he presented the techniques of snare polypectomy, endoscopic mucosal resection (EMR) and endoscopic submucosal dissection (ESD) with the aid of different studies and meta-analyses. He pointed out that the success story of endoscopic polyp treatment began with these resection techniques, which relevantly impacts public health in prophylactic but also therapeutic ways. However, the three above-mentioned conventional polypectomy techniques have weaknesses.

Based on the CARE study by Pohl H et al. (Gastroenterology 2013), a prospective study in two centers with more than 1700 patients, K. Caca demonstrated, in which cases the conventional resection techniques reach their limits. The study showed a high rate of incomplete resections in case of large polyps (diameter 15 – 20 mm, resected incompletely in 23.3 %), in case of sessile serrated adenoma (SSA/P, 31.0 % incomplete resection), in case of resection in piecemeal technique (20.4 % incomplete resection) and in case of difficult-to-resect lesions (17.8 % incomplete resection). A current meta-analysis including 32 studies with overall more than 2700 patients by Djinnbachian R et al. (Gastroenterology Online ahead of print 2020) showed similar findings.

This important therapeutic gap in endoscopic resection of mucosal/submucosal lesions is closed by the technique of endoscopic full-thickness resection (EFTR). EFTR is applied when the concerned lesion is not limited to the mucosa (negative lifting-sign caused by scarring or malign infiltration), when it has an endoscopically-difficult-to-reach localisation (appendiceal orifice, diverticulum, curve/fold), or is a subepithelial lesion (mostly tumors of the upper GI tract). EFTR can be performed with different techniques, these can be categorized in two groups according to the procedural approach: either the resection and thereby de facto a perforation of the gastrointestinal wall is performed at the beginning and closed subsequently (for example in STER submucosal tunnelling endoscopic resection) or closure of the wall is performed primarily by wall duplication before resection (so-called non-exposure techniques i.e. FTRD System or GERDX-suturing System).

By now, there is a multitude of data evaluating the FTRD System, retrospective and prospective, monocentric and multicentric, as well as register data from Germany and the Netherlands. The technical success rate is in the overall view between 75 % and 97 %, the R0-resection rate between 75 % and 90 %.

Data of the German colonic FTRD Register is accessible online now (Meier B et al., AJG 2020, 1178 cases), and evaluates cases from medical centers with different care levels. However, results are comparable to earlier studies which shows that the technique is well learnable and performable anywhere. It is an interesting fact that the study

did not find differences in R0-resection rates that were dependant on the lesion size (lesion size < 20 mm R=-resection rate 78 % vs lesion size ≥ 20 mm R0-resection rate 81 %).

In a video K. Caca showed a full-thickness resection with the FTRD in a lesion with suspected carcinoma in a 84 year-old patient. Previously performed surface biopsies had not yielded clear results. Full-thickness resection with the FTRD was technically unproblematic, macroscopically diagnosed as R0, histological examination showed a tubular adenoma with HGien with focal transition into a moderately differentiated adenocarcinoma. R0 resection could be confirmed histologically, the tumor was classified as pT1 L0 V0 R0, depth of infiltration was < 1 mm. Hence, EFTR had spared the patient a surgical intervention.

A study by Kuellmer A et al. (GIE 2019) analysed data from 159 FTRD applications in lesions with non-lifting sign. The cohort was divided into two groups: in group 1 the lesions were pre-treated endoscopically, yet could not be confirmed histologically as R0-resection (n=64). In group 2 the existence of a carcinoma was primarily suspected due to tumor characteristics (n=92). In group 1, the full-thickness specimen was identified histologically in 76.6 % as scar, in 21.9 % as carcinoma (thereof 84.1 % as low-risk and 15.9 % as high-risk). In group 2 histological diagnosis was a carcinoma in 100 % of cases, thereof 16.3 % low-risk and 83.7 % high-risk carcinomas. In all patients with low-risk carcinoma the FTRD treatment spared a surgical intervention.

After that, K. Caca lectured on lesions, which cannot be completely resected by EMR due to negative-lifting or localisation, but are also too large for en-bloc-EFTR. In these lesions, a hybrid technique is applied, in which the lateral margins are removed as far as possible with piecemeal-EMR and the non-lifting region is resected with the FTRD. A first study on this procedure was published by Meier et al. (Surg Endoscopy 2017), it evaluates 10 cases with hybrid EMR-EFTR. Currently, a working group at the Hospital of Ludwigsburg performs a study in 32 patients who underwent hybrid EMR-EFTR. Four treated lesions were located in the left colon, 28 in the right colon. Procedural time was 40 – 140 minutes. No SAEs occurred. During follow-up 3 relapse lesions were found, all could be treated with FTRD again.

Yuen et al presented data on hybrid EMR-EFTR at the DDW 2020 online. The study included 62 advanced colorectal adenomas, 33 of these were resected with FTRD only and 29 with hybrid technique. Median size of the lesions was 19 mm (7-40 mm) in group 1 and 36 mm (15 – 60 mm) in group 2. The technical success rate was 89 % in group 1 and 96 % in group 2. R0 resection was confirmed in 97 % of interventions of group 1 and 96 % of interventions of group 2. Overall two major complications occurred, one appendicitis and one perforation.

K. Caca summarized that almost 100 % of all benign GI lesions nowadays can be resected endoscopically, which was not the case 20 – 30 years ago. He emphasized the immense importance of always aiming for a complete resection of the lesion, which is performed with forceps in lesions < 3 mm, and with cold snare in lesions of 4 – 10 mm size. When en-bloc resection is not possible, piecemeal EMR can be performed. It is important to respect the limits of the endoscopic technique: In case of signs for obvious deep submucosal infiltration, endoscopic resection is not indicated, but surgical resection must be performed. The FTRD has closed an important therapeutic gap in endoscopic resection of mucosal/submucosal lesions. Hybrid EMR-EFTR closes a small but annoying therapeutic gap.

K.Caca ended his lecture with a slogan known from surgery, which he extended in his own way: When in doubt, take it out! But take it out endoscopically!

Basis and limitations of endoscopic resection techniques (Basis und Grenzen endoskopischer Resektionstechniken).

Caca K, Ludwigsburg

OTSC System

Live-broadcast of a lecture by J. G. Albert on endoscopic perforations

J. G. Albert, Robert-Bosch-Hospital Stuttgart, lectured on endoscopic perforations. First of all he gave an overview of the pathophysiology of transmural gastrointestinal defects: perforations can develop from ischaemic causes (intestinal

obstruction, necrosis, vascular occlusion), infections can play a role (appendicitis, diverticulitis), an inflammation (fulminant colitis ulcerosa) or neoplasia can be the cause, or the defect can result from mechanical injury (trauma, iatrogenic). It is important to keep this in mind during endoscopy, because in case of corresponding findings (i.e. colitis ulcerosa) a defect can already be present or easier to cause.

Generally we distinguish leakage (post-operative anastomotic insufficiency), perforations (acute discontinuity) and fistulas (conjunction between i.e. bowel and cavity caused by inflammation/neoplasia or secondary after perforation/leakage).

Perforations during endoscopy are rare, they occur during EGD in about 0.03 % of cases, during (ileo-) colonoscopy in 0.5 %. In colonoscopy screenings, the occurrence of perforations is also dependent on the indication. In a study by Pox C et al., (Gastroenterol 2012) with data on colonoscopy screenings in Germany a perforation occurred in 0.022 % of cases. If a polypectomy was performed during colonoscopy, the ratio was 0.046 % while in colonoscopy without polypectomy the ratio was 0.012 %. Besides, patient age and sex play a role for the degree of perforation risk. The older a patient is, the higher is the risk, furthermore the perforation risk is significantly higher in female patients than in male patients.

Other studies, however, report significant higher perforation rates (Heldwein W, Endoscopy 2005; Heresbach, Endoscopy 2010), in these too, the risk is dependent on the performed procedure, for snare polypectomy a perforation rate of 1.1 % is reported, for EMR up to 7 %. For ESD, the perforation rate is even substantially higher, retrospective meta-analyses on resections of gastric neoplasias (Facciorusso A 2014, Lian J 2012) show a 4-fold higher risk of ESD in comparison to EMR. The decision whether the therapy of a perforation is carried out surgically or endoscopically should always be an interdisciplinary decision on/with the patient. It is a fact that a tight suture cannot only be reached surgically, but just as good with endoscopic clips (OTSC, conventional clips). The therapy mustn't focus only on closure of the defect/the leakage, but must also take into account that an infection can have occurred, which also needs treatment (antibiotics, effective drainage). Primary fast diagnostic endoscopy by an experienced endoscopist is reasonable, but subsequently an interdisciplinary decision-making should take place with participation of an endoscopist, a surgeon and a radiologist. Especially in cases of delayed diagnosis of a perforation this interdisciplinary decision-making is very important, in which the appropriate approach for effective drainage and defect closure is discussed. However, emphasis should be laid on a fast course of action, because the mortality rate significantly increases in case of treatment delay > 24h.

The most often applied treatment procedures for endoscopic perforation closure are cSEMS and OTSC. A multicenter retrospective study (Farnik H, Plos One 2015) shows that the stent is more frequently used in larger perforations (median leakage size treated with cSEMS 12.6 mm vs with OTSC 7.1 mm), and also in case of already existing infection the cSEMS was more often applied than the OTSC. The OTSC is favoured in acute perforations as well as in endoscopically iatrogenic perforations, while the cSEMS is rather applied in case of older leakages and in surgically iatrogenic caused defects.

Endoscopic perforation: How do we prevent them? How do we treat them? (Endoskopische Perforation: Wie vermeiden? Wie behandeln?)

Albert JG, Stuttgart

November 2020 | EMR+: pilot study evaluating data from the first in human clinical cases shows 100 % R0-resection rate and no severe adverse events

The novel EMR+ technique uses an external additional working channel (AWC®) allowing for tissue retraction and triangulation. This pilot study describes the early safety and efficacy data from the first 6 in human clinical cases. En-bloc resection and R0-status was achieved in all patients, no severe adverse events occurred.

A. Sportes et al., Department of Gastroenterology, Institut Arnault Tzanck, Saint Laurent du Var, France, presented a pilot study evaluating the first use of the novel EMR+

technique in humans. EMR+ is a grasp and snare technique utilizing an over the scope additional working channel (AWC, Ovesco Endoscopy), which allows the introduction of an additional grasping tool for traction or counter-traction. Positioning of the AWC can be customized to the clinical scenario, depending on the position of the lesion and optimal angle of exit for the second tool. Thereby EMR+ allows for en-bloc resection of mucosal lesions ≥ 2 cm suspected of harboring high-grade dysplasia or early adenocarcinoma.

EMR+ was performed in 6 patients in 2 endoscopy centers (University Hospital of Goettingen, Germany and Institut Arnault Tzanck, Saint Laurent du Var, France) for lesions of the upper (n = 2) and lower (n = 4) gastrointestinal tract. The median age of the patients was 76 years, 3 patients were male, 3 female. The average estimated lesion size was 30.8 mm. Mean procedure time was 25.5 min. All resections were technically successful (en-bloc resection confirmed by endoscopist) and clinically successful (histologically confirmed R0-status in all cases). No severe adverse events such as major bleeding or perforation occurred in any of the procedures. In 2 cases, there was intra-procedural minor bleeding which could be managed via a hemoclip application. There were no post-interventional adverse events. During a follow-up of 6 months no further endoscopic or surgical treatment was subsequently needed.

The authors concluded that the EMR+ technique showed a good safety and efficacy profile in this pilot study. The technique may offer an alternative to available techniques to remove complex, large mucosal-based lesions in the GI tract. **Novel modified endoscopic mucosal resection of large GI lesions (> 20 mm) using an external working channel (AWC) may improve R0 resection rate: initial clinical experience**

Sportes A, Jung CFM, Gromski MA, Koehler P, Seif Amir Hosseini A, Kauffmann P, Ellenrieder V, Wedi E. BMC Gastroenterology 2020; 20:195. <https://doi.org/10.1186/s12876-020-01344-6>

Oktober 2020 | OTSC® closure of the entry site in POEM and G-POEM procedures is related to less clip-related adverse events when compared to closure with conventional clips

In 46 POEM and 26 G-POEM procedures, OTSC closure of the entry site was related to 13 % and 0 % clip-related adverse events, whereas closure with conventional clips led to 21.7 % and 7.7 % clip-related adverse events. Only a single OTSC clip was used in all procedures compared to 5 and 6 conventional clips in POEM and G-POEM procedures, respectively.

Ó.V. Hernández Mondragón et al., Endoscopy Department, National Medical Center Century XXI, Mexico City, Mexico, published a study comparing OTSC and conventional clips for entry site closure in per-oral endoscopy myotomy (POEM) and gastric peroral endoscopy myotomy (G-POEM) procedures.

46 POEM and 26 G-POEM procedures were included in the study. In the POEM group, 23 patients received entry site closure (ESC) with OTSC and in the remaining 23 patients ESC was performed with conventional clips (CC). In the G-POEM group, ESC closure was performed in 13 patients with OTSC and in 13 with CCs.

In POEM procedures, no differences were found between the OTSC and CC arm regarding demographic data and pre-interventional clinical characteristics. Besides, there were no differences between the OTSC and the CC arm regarding tunnel length, myotomy length, incision size and procedure time. However, in the CC arm compared with the OTSC arm, there were longer clip placement times (6.5 ± 1.1 vs 3.2 ± 0.7 minutes; p = 0.01), higher number of clips used (5 vs 1; p = 0.01) and more clip related adverse events (21.7 % vs 13 %; p = 0.01). Clinical POEM success was 100 % for both arms.

In G-POEM procedures no demographic or pre-interventional clinical differences were found between the OTSC and the CC arm. Also procedural characteristics were similar. However, longer clip placement times (6.9 ± 1.4 vs 3.1 ± 1.1; p = 0.02), a higher number of clips used (6 vs 1; p = 0.01) and more clip related adverse events (7.7 % vs 0 %; p = 0.02) were reported in the CC arm compared with OTSC group. G-POEM had a 100 % clinical efficacy for both arms.

The authors concluded that OTSC represents a safe and

effective alternative for entry site closure in POEM and G-POEM procedures. Further studies are needed to recommend OTSC as the first option for submucosal tunnel closure in these cases.

Conventional clips vs over-the-scope-clips for the closure of the entry site in POEM and G-POEM procedures

Hernández Mondragón ÓV, Gutiérrez-Aguilar R, García Contreras LF, Palos-Cuellar R, Blanco Velasco G, Monroy Teniza ZA

Rev Esp Enferm Dig. 2020; 112 (5): 338-342.

Oktober 2020 | ESD+: preclinical evaluation of the novel technique for endoscopic en-bloc resection shows superiority over conventional ESD regarding procedural time and safety particularly in flat gastric lesions in retrograde positions

The novel ESD+ technique uses the external additional working channel (AWC[®]), which allows for use of additional endoscopic instruments and thereby facilitates grasping and mobilization of the lesion. A preclinical study evaluated the technique in the EASIE-R1 Simulator (explanted pig stomach, 64 lesions). Especially in retrograde lesions ESD+ was significantly faster with a median of 22.5 vs 34.0 min in 3 cm retrograde lesions ($p = 0.002$) and 34.5 vs 41.0 min ($p = 0.011$) in 4 cm retrograde lesions. Muscularis damage occurred in one case with ESD+ technique and in 6 cases with conventional ESD (3.13 % vs 18.75 %, $p = 0.045$).

R. F. Knoop et al., Department of Gastroenterology and Gastrointestinal Oncology, University Medical Center, Georg-August-University, Goettingen, Germany, conducted a prospective study aiming to systematically evaluate the novel ESD+ technique and to compare it to conventional ESD. The ESD+ technique is characterized by the use of the external additional working channel (AWC), which allows bimanual working with a standard endoscope. Thereby for example the combination of grasper and ESD knife or ESD coagulation dissector is made possible. In contrast to a dual channel endoscope with a close and fixed distance between the channels, the use of the AWC enables triangulation and thereby intraluminal tissue traction, flexibility and a better overview.

The study was conducted in the EASIE-R simulator, a pre-clinical ex-vivo animal model using porcine stomachs. Standardized lesions of 3 cm or 4 cm were set in antegrade or retrograde positions. Overall 64 procedures were performed (32 with ESD+ and 32 with conventional ESD, thereof 16 with antegrade and 16 with retrograde position, respectively). The ESD+ technique as well as the conventional ESD technique were both reliable and showed en-bloc resection rates of 100 %. Across all groups, ESD+ was associated with significantly shorter procedural time compared to conventional ESD (24.5 vs 32.5 min, $p = 0.025$). Especially in retrograde lesions ESD+ was significantly faster with a median of 22.5 vs 34.0 min in 3 cm retrograde lesions ($p = 0.002$) and 34.5 vs 41.0 min ($p = 0.011$) in 4 cm retrograde lesions. No perforations occurred, neither with ESD+ nor with conventional ESD. Muscularis damage occurred in one case with ESD+ technique and in 6 cases with conventional ESD (3.13 % vs 18.75 %, $p = 0.045$).

The authors concluded that this preclinical study suggests that ESD+ allows faster and safer resection of flat gastric lesions compared to conventional ESD. This may be particularly relevant for difficult lesions in retrograde positions. **Endoscopic submucosal dissection with an additional working channel (ESD+): a novel technique to improve procedure time and safety of ESD**

Knoop RF, Wedi E, Petzold G, Bremer SCB, Amanzada A, Ellenrieder V, Neese A, Kunsch S.
Surg Endosc. 2020 Jul 16. doi: 10.1007/s00464-020-07808-w.

Oktober 2020 | 2020 Latest ESGE Position Statement on diagnosis and management of iatrogenic endoscopic perforations recommends OTSC[®] for defects in size > 10 mm

and < 20 – 25 mm

OTSC application is recommended for esophageal iatrogenic perforation > 10 mm, for gastric iatrogenic perforations > 10 mm, for early recognized (< 12 h) non-periampullary duodenal iatrogenic perforations, for immediately diagnosed jejunal and ileal iatrogenic perforations and for colorectal iatrogenic perforations > 10 mm

The revised position paper from the European Society of Gastrointestinal Endoscopy (ESGE) on inadvertent iatrogenic perforations occurring during gastrointestinal endoscopy was published in August 2020 by G. A. Paspatis et al., Venizelion General Hospital, Heraklion, Crete-Greece. It updates the 2014 ESGE recommendations on iatrogenic perforations.

The development process included meetings and online discussions of the project committee, which was divided into subgroups each in charge of a series of key questions. Searches were performed on Medline and the Cochrane Central Register of Controlled Trials. Each subgroup developed draft proposals that were presented to the entire group for general discussion. A draft prepared by the project committee chair was sent to all group members for review and besides reviewed by two external reviewers and then sent to the ESGE national societies and individual members. After agreement on a final version, the manuscript was submitted for publication.

The ESGE recommendations include the implementation of a written policy regarding the management of iatrogenic perforation at each endoscopy center, a comprehensive report in case of endoscopically identified perforation, as well as careful evaluation and CT scan in case of specific post-endoscopy symptoms. The therapy recommended for iatrogenic perforation is immediate endoscopic closure if possible. Through-the-scope clips may only be used for small defects (< 10 mm), the OTSC is recommended for larger perforations (> 10 mm but < 20 – 25 mm). The authors state that the wide application of OTSC clips in everyday practice has enhanced the efficacy of endoscopic closure, even in defects larger than 30 mm. Further treatment recommendations include switch to carbon dioxide endoscopic insufflation, diversion of digestive luminal content and decompression of tension pneumoperitoneum or pneumothorax. Further management should be based on the estimated success of the endoscopic closure and on the general condition of the patient. Hospitalization and surgical consultation are recommended in case of not possible or failed endoscopic closure and in patients whose clinical condition is deteriorating.

The position paper gives detailed recommendations for iatrogenic esophageal perforation, gastric perforation, periampullary and biliopancreatic ductal perforation, duodenal and small bowel perforation and colorectal perforation. ESGE recommends first-step endoscopic treatment for esophageal perforation with conventional clips in case of < 10 mm size, OTSC for perforations > 10 mm, stents can be used for larger defects > 20 mm. For gastric perforations, ESGE recommends conventional clips for perforations ≤ 10 mm, and OTSC clips or omental patching or the combined technique using endoloop and conventional clips for perforations > 10 mm. ESGE suggests nonsurgical management in the majority of ERCP-related periampullary or biliopancreatic ductal perforations. In case of severe leakage or sepsis, surgery is required. For duodenal iatrogenic perforation, ESGE recommends endoscopic treatment with OTSC clips or a combination of conventional clips and endoloop or stents if the perforation is recognized within 12 hours after the procedure. In case of failed endoscopic treatment or late diagnosis (> 12 h) of the perforation and clinical deterioration, management should be surgical. Also jejunal and ileal iatrogenic perforation should be treated endoscopically when recognized immediately and surgically in the case of failed endoscopic closure or postprocedural recognition with severe clinical status. For colorectal perforations, ESGE recommends the use of conventional clips for lesions < 10 mm and OTSC clips for defects > 10 mm.

In summary, OTSC clips are recommended in the ESGE position statement on iatrogenic endoscopic perforations for defects > 10 mm with esophageal, gastric, non-periampullary duodenal (when recognized early), jejunal

and ileal (when immediately diagnosed) and colorectal location.

Diagnosis and management of iatrogenic endoscopic perforations: European Society of Gastrointestinal Endoscopy (ESGE) Position Statement – Update 2020
Paspatis GA, Arvanitakis M, Dumonceau J, Barthet M, Saunders B, Turino SY, Dhillon A, Fragaki M, Gonzalez J, Repici A, van Wanrooij RLJ, van Hooft JE.
Endoscopy 2020; 52: 792–810. <https://doi.org/10.1055/a-1222-3191>

Ovesco Research Update 36

September 2020 | EFTR with the FTRD[®] should be considered as an alternative to surgery in gastrointestinal lesions that are not amenable to conventional endoscopic resection

A meta-analysis comprising 18 studies with overall 730 patients with difficult gastrointestinal lesions showed a R0 resection rate of 82 %, a full-thickness resection rate of 83 %, and an enbloc resection rate of 95 %. The pooled estimate for perforation was < 0.1 % and for bleeding 2 %.

O. I. Brewer Gutierrez et al., Division of Gastroenterology and Hepatology, Johns Hopkins Hospital, Baltimore, Maryland, USA, published a meta-analysis aiming to summarize the current scientific evidence on efficacy, safety and clinical outcomes of patients with gastrointestinal neoplasia treated with the FTRD. The primary outcome of interest was the rate of histologically confirmed R0 resection, secondary outcomes were enbloc resection rate, rate of full thickness resection, adverse events, and post-EFTR surgery. The authors performed a comprehensive literature search and used the random-effects model to calculate pooled estimates and generate forest plots.

The meta-analysis included 18 studies with overall 730 patients and 733 lesions for analysis. Indications for EFTR were difficult residual colorectal adenoma, adenoma at a diverticulum or appendiceal orifice, and early cancer ($n = 634$), colorectal subepithelial lesions ($n = 42$), upper gastrointestinal lesions ($n = 51$), and other colonic lesions ($n = 6$). Median size of the lesions was 13.5 mm. 22 EFTR attempts failed. The pooled overall R0 resection rate was 82 % (95 % CI: 75 – 89). The pooled overall FTR rate was 83 % (95 % CI: 77 – 89). The pooled overall enbloc resection rate was 95 % (95 % CI: 92 – 96). The pooled estimates for perforation and bleeding were < 0.1 % and 2 %, respectively. Following EFTR, a total of 110 patients underwent surgery for any reason (pooled rate 7 %, 95 % CI: 2 – 14). The pooled rates for post-EFTR surgery due to invasive cancer, for non-curative endoscopic resection and for adverse events were 4 %, < 0.1 %, and < 0.1 %, respectively. No mortality related to EFTR was observed.

The authors concluded that EFTR with the FTRD is safe and effective for difficult mucosal and submucosal gastrointestinal lesions with high rates of complete resection and an acceptable rate of adverse events. These findings emphasize the importance of optimizing and standardizing the EFTR technique, of attaining widespread implementation of this procedure and thereby sparing the need of surgery in patients with difficult mucosal and submucosal gastrointestinal lesions.

Endoscopic full-thickness resection using a clip non-exposed method for gastrointestinal tract lesions: a meta-analysis

Brewer Gutierrez OI, Akshintala VS, Ichkhanian Y, Brewer GG, Hanada Y, Truskey MP, Agarwal A, Hajjiyeva G, Kumbhari V, Kalloo AN, Khashab MA, Ngamruengphong S
Endosc Int Open. 2020 Mar;8(3):E313-E325. doi: 10.1055/a-1073-7593. Epub 2020 Feb 21.

September 2020 | OTSC[®] clipping for recurrent peptic ulcer bleeding is cost-effective as compared to a repeat of standard treatment

A cost analysis based on mathematic modelling determined the average cost of repeating standard treatment in patients with persistent/recurrent peptic ulcer bleeding (\$ 6,578) and the average cost of second-line OTSC treatment after failed standard therapy,

which was lower (at average \$ 6,298). The average cost of standard treatment was based on the cost of one hemoclip and a gold probe. However, many patients treated with standard endoscopy receive more than one hemoclip, which would make the cost of repeat standard treatment even higher.

J. X. Yu et al., Division of Gastroenterology and Hepatology, University of Michigan, Ann Arbor, Michigan, USA, evaluated by mathematic modelling the cost-effectiveness of OTSC clips for the treatment of peptic ulcer bleeding as first-line and second-line therapy. A decision tree was used to model the cost, effectiveness, and rates of persistent/recurrent bleeding of OTSC versus standard treatment (ST) for the management of peptic ulcer bleeding. Three possible treatment strategies were modeled for a patient with peptic ulcer bleeding: OTSC clipping after first-line standard therapy, a repeat of standard therapy after first-line standard therapy, and standard therapy after first-line OTSC clipping. It was assumed that if the second-line therapy remained unsuccessful, the patient would undergo interventional radiology or surgery.

The average cost of standard treatment was based on the cost of one hemoclip and a gold probe. The cost of an esophagogastroduodenoscopy and for hospitalization were estimated from Medicare and Medicaid data, the probability of persistent/recurrent bleeding with standard therapy and OTSC therapy was obtained through review of the literature. Results showed that the first treatment strategy of first-line ST followed by OTSC was the most cost-effective strategy, costing \$ 6,298 per patient and resulting in 0.0686 QALYs. The second strategy with ST followed by ST cost \$ 6,576 and resulted in 0.0659 QALYs. The third strategy with OTSC as first-line treatment and second-line ST cost \$ 6,490 and resulted in 0.0687 QALYs. As the probability of persistent/recurrent bleeding following ST increases, OTSC treatment becomes preferred in more scenarios. Specifically, for medium- or high-risk ulcer with Rockall ≥ 4 , OTSC as first line therapy followed by ST is preferred. The authors concluded that in patients with persistent/recurrent peptic ulcer bleeding after standard endoscopic treatment a repeat of the standard approach is not cost-effective. Instead, OTSC clipping should be used. For first-line therapy, OTSC clipping should be considered for patients with intermediate or high-risk ulcers as determined by the Rockall score.

Clipping Over the Scope for Recurrent Peptic Ulcer Bleeding is Cost-Effective as Compared to Standard Therapy: An Initial Assessment

Yu JX, Russell WA, Asokkumar R, Kaltenbach T, Soetniko R

Gastrointest Endoscopy Clin N Am 30 (2020) 91-97. <https://doi.org/10.1016/j.giec.2019.09.004>

July 2020 | stentfix OTSC® System: case series shows the new stentfix device to prevent effectively esophageal stent migration

stentfix OTSC was used in five complex cases who already had experienced stent migration and were considered unfit for surgery. Placement of stentfix OTSC resulted in immediate technical success in all patients. Stent migration or fixation-related adverse events did not occur in any patient. M Conio et al., Gastroenterology Department, General Hospital, Sanremo, Italy, published a case series describing the use of stentfix OTSC in five consecutive patients. All patients had been treated previously with a fully covered self-expandable metal stent (FCSEMS) and experienced FCSEMS migration. All patients were considered unfit for surgery and treated with the stentfix OTSC for fixation of the stent to prevent further migration. Four patients were male, one patient female. Median age of the patients was 58 years (range 53 – 87 years). Indications for FCSEMS placement were esophageal perforation (n=1), trachea-esophageal fistula (n=1), broncho-esophageal fistula (n=1), advanced esophagogastric adenocarcinoma (n=1) and colorectal anastomotic stricture (n=1). In three patients one stentfix OTSC was placed, in the remaining two patients two stentfix OTSC were placed. Placement of stentfix OTSC resulted in immediate technical success in all patients. Stent migration or fixation-related adverse events did not occur in any

patient. Stent and stentfix OTSC dwell time was 60 days (n=1), 90 days (n=1), 105 days (n=1), 150 days (n=1), and unknown (n=1). The authors concluded that the stentfix OTSC System represents an effective, easy to perform, and safe endoscopic technique for endoluminal fixation of fully covered self-expandable metal stents. Further experience is, however, necessary.

A newly designed over-the-scope-clip device to prevent fully covered metal stents migration: a pilot study

Conio M, Savarese MF, Baron TH, De Ceglie A. Techniques and Innovations in Gastrointestinal Endoscopy 2020 doi: <https://doi.org/10.1016/j.tige.2020.05.004>

July 2020 | OTSC® management of non-acute full-thickness gastrointestinal defects represents an effective and safe alternative to a potentially morbid surgical intervention

92 patients with 117 non-acute full-thickness gastrointestinal defects (65 fistulas, 52 leaks) underwent OTSC therapy. With a median follow-up time of 5.5 months, overall defect closure success rate was 66.1 %. Only 14.9 % of patients underwent operative management. No complications related to endoscopic intervention occurred and no patients required urgent surgical intervention.

D.J Morrell et al., Department of Surgery, Division of Minimally Invasive and Bariatric Surgery, Penn State Health Milton S. Hershey Medical Center, Hershey, PA, USA, conducted a study evaluating endoscopic management of non-acute full-thickness gastrointestinal defects (FTGID) using the OTSC (Ovesco Endoscopy AG, Tuebingen, Germany).

Non-acute FTGID, when compared to acute perforations, are typically more difficult to manage due to a delayed diagnosis and/or chronicity. The open surgical management of these defects is often challenging and related to substantial morbidity, requiring washout and drainage of infection, takedown or recreation of surgical anastomoses, deviation stoma formation and/or parenteral nutrition.

For the study, all patients undergoing OTSC management of FTGID between 2013 and 2019 were identified. Acute perforations and defects requiring endoscopic suturing were excluded. The endoscopic methodology to approach FTGID included treating underlying infection, optimizing nutrition, removing foreign bodies, de-epithelizing/ablating mucosa, relieving downstream obstructions and collecting GI effluent. Success was strictly defined as complete FTGID closure, patients with multiple FTGID were only counted as successfully managed if all FTGID were closed at the end of follow-up.

Overall, 92 patients (mean age 54.6 \pm 15.4 years, 64.1 % female) with 117 FTGID (65 fistulas, 52 leaks) met the inclusion criteria. The majority of defects was located in the upper GI tract (64.1 %). The estimated defect diameter was 4.6 mm for leaks and 6.8 mm for fistulas. Overall prior intervention rate was 40.2 % with a mean of 0.56 prior interventions. The OTSC was used for treatment in all cases. Median follow-up time was 5.5 months (interquartile range 2.2 – 14.8). Long-term complete FTGID closure was achieved in 66.1 % (55.0 % fistulas vs 79.6 % leaks, p = 0.007). Additional closure attempts were required in 22.2 % of FTGID. Patients with failed FTGID closure in comparison to patients with successful closure were more likely to have had a history of radiation treatment at the defect site (12.9 % vs 1.8 %, p = 0.0329). Besides, a history of previous smoking was also more common in patients with unsuccessful closure (67.7 % vs 32.1 %, p = 0.0014).

No complications related to endoscopic intervention occurred. There were 4 mortalities from causes unrelated to the FTGID. Only a small subset of patients (14.9 %) ultimately underwent surgical management, no patients required urgent operative intervention.

The authors concluded that OTSC management of delayed or chronic full-thickness gastrointestinal defects is safe and effective and constitutes an attractive alternative to potentially morbid surgical management.

Over-the-scope clip management of non-acute, full-thickness gastrointestinal defects

Morrell DJ, Winder JS, Johri A, Docimo S, Juza RM, Witte SR, Alli VV, Pauli EM.

Surg Endosc (2019). <https://doi.org/10.1007/s00464-019-07030-3>

Ovesco Research Update 35

June 2020 | First clinical cases with LiftUp® reveal promising results, proving it as valuable, safe and fast

The novel submucosal injection agent LiftUp was successfully used in a first clinical case series with 10 patients. Resection of colorectal polyps (median lesion size of 20.5 mm) was macroscopically complete in 100 %. R0 resection was achieved in 78 %. No adverse events were reported.

M. Schaefer et al. reported on the first experience with the novel agent LiftUp for submucosal injection. 10 consecutive patients (median age 74.5 years; 6 male, 4 female) with colorectal polyps ≥ 15 mm were included. Lesions were localized in the rectum (n=5), caecum (n=2), sigmoid (n=1), transverse colon (n=1), and ascending colon (n=1) and the median lesion size was 20.5 mm (range: 15-30 mm). To assess endoscopic resection after LiftUp injection, EMR (8 patients) or ESD (2 patients) were performed. Median volume of LiftUp was 5.25 ml (range 2-12 ml) and median procedure time was 2.25 minutes for EMR and 5.5 minutes for ESD. In all 10 cases macroscopic complete resection was achieved and the R0 resection rate was 78 % (7/9). No adverse events occurred.

The authors concluded that endoscopic resection with LiftUp is feasible, safe, and fast. However, larger studies are needed for better evaluation of efficacy and safety.

First human in-vivo experience with novel agent for submucosal injection (LiftUp)

Moritz Schaefer, Benjamin Meier, Andreas Wannhoff, Karel Caca.

ePoster submitted at DDW 2020; No. 3351549.

June 2020 | BougieCap: case study of post-esophagectomy stricture in the high cervical esophagus successfully dilated to 12 mm with the BougieCap

The patient had experienced esophageal perforation due to bougienage with a through-the-scope balloon, which led to emergency thoracotomy with esophagectomy and gastric pull-up. A 5 mm stricture developing at the anastomotic site 15 months later was safely and effectively dilated with the BougieCap

H. Ramrakhiani and G. Triadafilopoulos, Stanford University School of Medicine, Stanford, CA, USA, presented a case study on stricture dilation with the BougieCap.

A 62-year old man had previously presented with chest pain, nausea, vomiting, and odynophagia. Diagnostic endoscopy with circumferential biopsies showed evidence of acute ischemic injury of the esophageal wall. Four weeks later, the ischemic area had turned into a stricture, extending 10 cm long from mid to distal esophagus and as narrow as 4 mm. A balloon dilation was undertaken with a through-the-scope balloon, which was complicated by an esophageal perforation. Emergency thoracotomy was necessary, associated with distal esophagectomy, gastric pull-up and placement of a feeding jejunostomy. 15 months later the patient presented with progressive difficulty in swallowing associated with 16 lb weight loss. Barium swallow examination revealed a tight stricture at the anastomosis with a narrowed lumen diameter of 5 mm. The stricture was located just distal to the upper esophageal sphincter. This time, the decision for bougienage with the BougieCap was made. Over the course of a few weeks, the patient underwent three upper endoscopies with BougieCap dilation of the stricture to 12 mm size with good results. At a follow-up examination 6 months later, the patient reported being able to eat a normal diet and had gained 10 lb.

The authors discussed that currently available bougienage techniques for esophageal strictures (balloons and mechanical dilators) have limitations when used for the treatment of high cervical esophageal strictures or eosinophilic esophagitis-related tight esophageal strictures. The transparent BougieCap allows for visualization of mucosal disruption in real time during stricture dilation. Thus

the authors conclude that the device may help to avoid overstretching and reduce the risk of complications. As the BougieCap is a relatively new device, prospective studies with larger numbers of patients and with long-term follow-up are still outstanding.

Negotiating Dire Straits with a BougieCap

Ramrakhiani H, Triadafilopoulos G.

Digestive Diseases and Sciences; <https://doi.org/10.1007/s10620-020-06236-z>

June 2020 | Conference Report of DDW 2020 Conference report of the Digestive Disease Week (DDW) 2020

- **FTRD® System: EFTR enables new therapeutic options for the treatment of lesions and is a feasible and safe modality, which leads to a high R0 resection rate**
- **First clinical cases with LiftUp® reveal promising results, proving it as valuable, safe and fast**
- **OTSC® shows superiority & cost-effectiveness over standard therapy for first-line UGIB and shows high success rates in the management of non-acute, full-thickness gastrointestinal defects**

The conference of the Digestive Disease Week (DDW) 2020 scheduled for May 2-5, 2020 in Chicago, had to be cancelled due to the current situation. However, many scientific publications on Ovesco products were submitted this year again. The most important abstracts on Ovesco products are summarized in the report attached below.

FTRD System

The FTRD System is a safe and useful method for resection of colonic lesions including the appendix. Y. Ichkhanian et al. presented the first multicenter international study assessing the outcome of using FTRD for endoscopic full-thickness resection of appendiceal lesions. To study the rate of appendicitis after lesion resection, 56 consecutive patients between 11/2016 and 11/2019 underwent resection of colonic lesions involving the appendix. Most lesions were of Is 25 (45 %) and Ila 20 (36 %) Paris class and had a mean circumferential appendiceal orifice involvement of 66±29 %. The technical success rate was 95 % (53 patients) and in 47 cases successful EFTR was achieved. R0 resection was achieved in 92 % (49/53), antibiotics were administered in 81 % (43/53) for a median of 5 days. Follow-up colonoscopy showed no residual lesions in 19/53 patients. In 9 (19 %) out of 47 successful EFTR cases appendicitis was diagnosed. Three patients were treated with antibiotics alone, 6 underwent appendectomy.

The authors concluded that the FTRD is safe and feasible for the resection of appendiceal lesions, with appendicitis occurring in about 20 % of cases.

Endoscopic full-thickness resection of polyps involving the appendiceal orifice: first multicenter international study

Yervant Ichkhanian, Mohammed Barawi, Bachir Ghandour, Talal Seoud, Shyam Thakkar, Truptesh H. Kothari, Krystle Bittner, Asad Ullah, Vivek Kaul, Shivangi Kothari, Edris Wedi, Patrick Aepli, Thomas E. Kowalski, Brianna Shinn, Reem Z. Shariha, Qais M. Dawod, Srihari Mahadev, Kaveh Hajifathalian, Gianluca Andrisani, Heinz Albrecht, Andrew Q. Giap, Shou J. Tang, Yehia M. Naga, Ayla S. Turan, Shai Friedland, Benjamin Tharian, Shayan S. Irani, Andrew S. Ross, Laith H. Jamil, Daniel Lew, Andrew Nett, Omid Sanaei, Nasim Parsa, Jad Farha, Manol Jovani, Vivek Kumbhari, Anthony N. Kallou, Mouen A. Khashab

EFTR with the FTRD in the colorectum is safe and more efficient when compared to conventional endoscopic resection

B. Meier et al. report from the German colonic FTRD registry in which the use of the FTRD for resection of colorectal lesions is analysed. Data from 1178 FTRD procedures provided by various German endoscopy centers were collected between 09/2015 and 10/2019. Whereas the R0 resection rate (overall 80 %) was higher in the rectum (83.6 %) than in the colon, the full-thickness resection rate (overall 89.9 %) in the rectum (83.3 %) was lower when

compared to the colon. In 92 % of cases (1086/1178) follow-up data were available, endoscopic follow-up in 58 % (683/1178) with a mean follow-up time of 22 weeks (median 14 weeks, range 1-202 weeks). Dislocation of the FTRD Clip appeared in 69 % and in 31 % it was still in situ. Procedure-related adverse events were reported in 12.1 %, of which 2 % required surgical treatment.

The authors summarized that FTRD shows safe, effective, and fast results in endoscopic full-thickness resection of colorectal lesions.

Efficacy and safety of endoscopic full-thickness resection with the FTRD in the colorectum: evaluation of a large PMCF analysis

Benjamin Meier, Arthur Schmidt, Armin Kuellmer, Julius Mueller, Andreas Wannhoff, Klaus Weiss, Georg Huebner, Michael Repp, Alexander Meining, Benjamin Walter, Karel Caca.

EFTR is a more reliable in the therapy of colorectal neoplasia with high risk of incomplete removal when compared to endoscopic submucosal dissection (ESD). S. Suchanek et al. compared endoscopic full-thickness resection (EFTR) with the endoscopic submucosal dissection (ESD) in colorectal neoplasia (size ≤3 cm) therapy and primarily assessed R0 resection. The Bicentric prospective randomized study included 35 patients (71 % men, mean age 66 years) with T1 cancers (46 %), non-lifting adenomas (29 %) and local residual neoplasia (23 %). The 19 lesions (54 %) treated by EFTR were mostly localized in the right colon, while the 16 lesions treated by ESD were mostly localized in rectum. Technical success in the EFTR group was achieved in 90 % and in the ESD group in 94 % (p-value 0.999). R0 resection was successful in 17 patients (90 %) treated with EFTR and in 13 (81 %) patients treated with ESD. In 5/11 patients (46 %; EFTR) and in 6/10 patients (60 %; ESD) curative resection of malignant lesions was performed. Perforation occurred in 4 cases (11 %; all in ESD) of which 2 (6 %) needed surgical therapy and 2 (6 %) were treated endoscopically.

The authors concluded that both EFTR and ESD are technically effective methods for the treatment of colorectal neoplasia with high risk of incomplete removal, but EFTR seems to be the more safe method comparing to ESD. **Endoscopic full thickness resection versus endoscopic submucosal dissection in colorectal neoplasia therapy – bicentric prospective randomized study**

Stepan Suchanek, Premysl Falt, Ondrej Ngo, Nagyija Brogyuk, Renata Chloupkova, Ondrej Urban, Ondrej Majek, Miroslav Zavoral

EFTR+ EMR as hybrid technique enables new possibilities in resection of large colorectal lesions

W. Yuen et al. reported on the comparison of using the FTRD alone and using a hybrid EFTR + EMR technique in consecutive patients with lesions unresectable by conventional EMR alone. Many large colorectal lesions can be targeted by combining full-thickness resection (FTR) and endoscopic mucosal resection (EMR) by treating the laterally spreading components with EMR and invasive, non-lifting portions with FTR. Out of the 62 patients (total) 33 were treated with FTR alone and 29 were treated with hybrid EMR + FTR. The mean lesion size for the hybrid group was 36 mm (15-60 mm) and 19 mm (7-40 mm) for the FTR group. Technical success for FTR was achieved in 55 out of 62 patients (89 %), of which 53 (96 %) had R0 resection margins. No difference in R0 resection rate appeared when comparing hybrid FTR + EMR (23/24, 96 %) to FTR (30/31, 97 %) alone. In 2 cases adverse events occurred: one patient developed acute appendicitis and one patient suffered an inadvertent perforation.

The authors concluded that FTRD is a safe and effective method to resect large and complex colorectal lesions with high technical and clinical success rates. A hybrid EMR + FTR technique can expand the pool of resectable lesions.

Hybrid technique using endoscopic mucosal resection and endoscopic full-thickness resection for large colorectal lesions: initial north American experience

William Yuen, Sri Hari Mahadev, Poi Yu Sofia Yuen, Kristen Koller, Ashley A. Vareedayah, Gregory B. Haber

FTRD is an effective and safe method for lesions in the GI tract not amenable to traditional endoscopic

resection

A. Agnihotri et al. published a systematic review and meta-analysis to assess the efficacy and safety of the endoscopic full-thickness resection device (FTRD). Only studies with more than 9 patients were searched in bibliographic databases (PubMed, Ovid, Scopus, EMBASE, Web of Science, Cochrane and Google scholar) and abstracts published in major gastroenterology conferences within the last two years were evaluated. Overall, 15 studies with 1752 patients (mean age ranged from 59.7 to 70 years; 60 % male) with the following indications were included: recurrence or incomplete resection after prior polypectomy (50.7 %), primary non-lifting adenoma (17.8 %), known intramucosal cancer (9.9 %), submucosal lesions (9.2 %) and lesions involving appendiceal orifice (6.6 %). Pooled technical success was achieved in 88.3 % of cases (95 % CI: 0.86-0.89, I2=0 %) and most adverse events (overall pooled rate 12.3 %) were minor.

The authors confirmed previous study results by classifying the FTRD as effective and overall safe modality for endoscopic full-thickness resection.

Safety and efficacy of endoscopic full-thickness resection device (FTRD) in the gastrointestinal tract: a systematic review and meta-analysis

Abhishek Agnihotri, Mohammad Bilal, Shailendra Singh, Laura C. Horton, Corey S. Miller, Brian Chan, Douglas K. Pleskow, Mandeep Sawhney, Tyler M. Berzin, Jonah Cohen

FTRD in therapy for lesions in the upper GI tract is associated with a high clinical success rate and low risk of recurrence

K. Hajifathalian et al. presented an international multicenter retrospective study in which 56 patients from 13 centers (average age 61±14 years) with mesenchymal neoplasms including gastrointestinal stromal tumour (n=23, 41 %), adenomas (n=7, 13 %), and hamartomas (n=6, 11 %) were included. The average lesion size was 14 mm (range 3 to 33 mm), of which 47 (84 %) were localized in the stomach, and 8 (14 %) in the duodenum. The technical success rate with the FTRD was 93 %. In 43 patients (77 %) complete resection and in 9 patients (16 %) partial resection was achieved. R0 resection was successfully performed in 38 (68 %) cases and 12 (21 %) cases of adverse events were reported. Follow-up endoscopy of 31 cases (55 %) showed no residual or recurrent lesion in 30 patients (97 %), in one patient (3 %) residual adenoma was found. Based on the results, the authors concluded that FTRD in the upper GI tract shows high technical success rates, as well as an acceptable risk of complication and histologically complete resection.

Full-thickness resection device (FTRD) for treatment of upper gastrointestinal tract lesions: the first international experience

Kaveh Hajifathalian, Yervant Ichkhanian, Qais M. Dawod, Alexander Meining, Arthur Schmidt, Nicolas Glaser, Kia Vosoughi, David L. Diehl, Ian S. Grimm, Theodore W. James, Adam W. Templeton, Jason B. Samarasekera, Nabil E. Chehade, John Lee, Kenneth J. Chang, Meir Mizrahi, Mohammed Barawi, Shayan S. Irani, Shai Friedland, Mohammad A. Al-Haddad, Thomas E. Kowalski, George Smallfield, Gregory G. Ginsberg, Norio Fukami, Shou J. Tang, Franklin Kasmin, Martin Goetz, Stefan Seewald, Vivek Kumbhari, Saowane Ngamruengphong, Saurabh Mukewar, Kartik Sampath, Srihari Mahadev, David L. Carr-Locke, Mouen A. Khashab, Reem Z. Shariha

RESECT+

First case-series to use LiftUp for human in-vivo study

M. Schaefer et al. reported on the first experience with the novel agent LiftUp for submucosal injection. 10 consecutive patients (median age 74.5 years; 6 male, 4 female) with colorectal polyps ≥15 mm were included. Lesions were localized in the rectum (n=5), caecum (n=2), sigmoid (n=1), transverse colon (n=1), and ascending colon (n=1) and the median lesion size was 20.5 mm (range: 15-30 mm). To assess endoscopic resection after LiftUp injection, EMR (8 patients) or ESD (2 patients) were performed. Median volume of LiftUp was 5.25 ml (range 2-12 ml) and median procedure time was 2.25 minutes for EMR and 5.5 minutes for ESD. In all 10 cases

macroscopic complete resection was achieved and the R0 resection rate was 78 % (7/9). No adverse events occurred. The authors concluded that endoscopic resection with LiftUp is feasible, safe, and fast. However, larger studies are needed for better evaluation of efficacy and safety.

First human in-vivo experience with novel agent for submucosal injection (LiftUp)

Moritz Schaefer, Benjamin Meier, Andreas Wannhoff, Karel Caca

OTSC System

OTSC as first-line therapy for severe non-variceal UGI hemorrhage is recommended as the treatment of choice when compared to standard endoscopic hemostasis

D. Jensen et al. presented two randomized controlled trials (RCT) on the first-line treatment of UGIB where Doppler probe assisted (DEP) or OTSC hemostasis is compared with STANDARD endoscopic hemostasis. Individual data from the two successive RCT's were used to analyze the rates of recurrent bleeding and severe complications for patients with bleeding peptic ulcers (PUB's) or Dieulafoy's lesions (DL). 69 patients were treated with DEP, 24 with OTSC and 101 with STANDARD endoscopic treatment. All cases were randomized concerning the presence of stigmata of hemorrhage – major SRH (spurting arterial bleeding, non-bleeding visible vessel, or adherent clot) or lesser SRH (oozing bleeding or flat spot). The rates of 30-day rebleeding for OTSC were 4.17 % and more than double for DEP in 10.14 %. STANDARD endoscopic hemostasis had a rebleeding rate of 25.7 %, which undermines the significant superiority of DEP/OTSC.

The authors concluded that regarding the higher rates of rebleeding and severe complications for the STANDARD treatment group, the DEP or OTSC are the new initial treatments of choice for severe NVUGIB from PUB's or DL's.

What treatment to use for severe non-variceal UGI hemorrhage: STANDARD endoscopic hemostasis, Doppler guided hemostasis, or over-the-scope-clip (OTSC)?

Dennis M. Jensen, Thomas O. Kovacs, Rome Jutabha, Kevin A. Ghassemi, Marc Kaneshiro, Mary Ellen Jensen, Jeffrey Gornbein

OTSC as first-line therapy for NVUGIB is associated with lower costs and higher effectiveness in QALY when compared to standard therapy

J. Yu et al. evaluated costs and effectiveness in quality adjusted life years (QALY) for first-line therapy of severe NVUGIB with the OTSC System, the doppler probe assisted hemostasis (DEP) and standard therapy based on two randomized controlled trials (RCT). The three therapies were compared for lesions with major SRH (spurting bleeding, visible vessel, adherent clot) and lesser SRH (oozing bleeding, flat spot). 98 % of lesions included in the study were peptic ulcer disease (PUB) or Dieulafoy's (DL), while 2 % were Mallory Weiss tears. Cost research included procedure and hospitalization costs, as well as equipment costs. QALY was evaluated based on literature review. Analysis revealed that for high risk-stigmata OTSC costs less and shows better effectiveness in QALY results. Main reason for this is, that due to low rebleeding rates, fewer reinterventions are necessary. For lesser risk stigmata, DEP and OTSC showed same effectiveness in QALY, while DEP therapy (\$6147) costs \$18 less than OTSC (\$6165). Standard therapy had higher costs and lower QALY rates when compared to OTSC/DEP for both stigmata.

The authors concluded that the OTSC System is the therapy of choice for severe NVUGIB, especially for high-risk stigmata. For lesser-risk stigmata, DEP is preferred as it has a small advantage in cost-effectiveness.

Clipping over the-scope is cost effective for first line therapy of severe non-variceal UGI bleeding lesions with major stigmata

Jessica X. Yu, W. A. Russell, Dennis M. Jensen, Roy M. Soetikno

OTSC application is a valuable tool as part of combination of endoscopic and cystoscopic treatment for closure of colovesicular fistula

S. Siegal et al. reported on a male patient with diverticulitis presented with fecaluria, pneumaturia and urine per rectum and who was unable to undergo surgery. The workup revealed a colovesicular fistula which was treated with a combined approach of endoscopic and cystoscopic closure. For an easier fistula tract identification, initial cystoscopy was performed. Wire access through the bladder was placed to guide the therapeutic endoscope to the fistula in an endoscopic rendezvous maneuver. After successful deployment of the over-the-scope clip, the bladder mucosa was ablated by cystoscopic laser. Fluoroscopy showed that complete closure of the fistula was achieved.

The authors summarized that combined endoscopic and cystoscopic maneuvers are safe and reliable and are useful for patients to avoid surgical procedures.

Combined endoscopic and cystoscopic closure of colovesicular fistula

Steve R. Siegal, John Knoedler, Jeffrey S. Scow, Eric Pauli

OTSC is an effective alternative with high success rates to surgical intervention for the management of non-acute, full-thickness gastrointestinal defects

D. Morrell et al. published a systematic review and meta-analysis to investigate the use of OTSC for non-acute full-thickness gastrointestinal defects (FTGID). Databank research including MEDLINE and the Cochrane Library was performed for 43 studies in which the OTSC system was the primary treatment modality. These data were analyzed using a general variance-based approach and successful management was defined as radiologic

6 or clinical evidence of complete defect closure after conclusion of the follow-up period. 691 patients, of which most had gastrocutaneous fistulae (74), enterocutaneous fistulae (45), gastrogastroic fistulae (39), and gastric sleeve leaks (29) were evaluated. The pooled success rate for OTSC treated full-thickness GI defects was 67 % (95 % CI: 59-75 %). Analysis of subgroup studies of fistulas (n=410) and leaks (n=210) showed 52 % and 77 % pooled success. As conclusion, the authors see the OTSC as promising endoscopic treatment modality for non-acute GI defects. High rates of successful management as well as high effectiveness is expected when using the OTSC for defects usually treated with potentially morbid surgical interventions.

Over-the-scope clip management of non-acute, full-thickness gastrointestinal defects: a systematic review and meta-analysis

David Morrell, Christopher Hollenbeak, Eric Pauli

OTSC is a safe endoscopic method or hemostasis in acute variceal bleeding especially as a rescue treatment*

G. Sirinet et al. presented a study to evaluate the use of the OTSC System as first-line treatment or rescue therapy for endoscopic hemostasis in patients with severe variceal bleeding. Out of 21 consecutive patients (14 male) with a mean age of 38, 14 patients (rescue treatment group) had previously been treated unsuccessfully with conventional endoscopic methods and 7 patients (first-line therapy group) have not had endoscopic treatment before. In the rescue treatment group, 5 patients had anastomotic varices, 5 had fibrotic esophageal varices, and 4 had fundal varices. In first-line therapy group, 2 patients had anastomotic varices, 2 had fibrotic esophageal varices and 3 had fundal varices. Immediate hemostasis was achieved in all patients. Three patients with fundal varices suffered rebleeding, 2 in primary hemostasis group, and 1 in rescue therapy group, which could be treated successfully with synthetic glue injection. No adverse events related to OTSC application occurred. The authors concluded that the OTSC seems to be a safe and useful method for hemostasis in acute variceal bleeding, especially as rescue treatment.

*Note: Usage of the OTSC System Set is contraindicated for the treatment of bleedings of esophageal varices.

Is over-the-scope-clip an effective rescue therapy for serious variceal bleeding?

Goktug Sirin, Ali Erkan Duman, Hasan Yilmaz, Altay Celebi, Sadettin Huelague

May 2020 | OTSC®: 100 % clinical success rate without any complications in patients with

perforated peptic ulcer

106 patients with perforated peptic ulcer were treated with OTSC (n = 26) or conservative therapy (n = 80). In the OTSC group, technical and clinical success was achieved in 100 % of patients without any complications, the median operation time was 10 min. In the control group, the clinical success rate was 57.5 %, the mortality rate was 13.8 % and subsequent surgeries were required in 30 % of patients.

J. Wei et al., Endoscopic Center, The First Affiliated Hospital of Fujian Medical University, Fuzhou, China, conducted a retrospective study on 106 patients with perforated peptic ulcer.

26 of those patients were treated with OTSC and 80 were treated with pharmacotherapies as the initial choice. No significant differences in age (p = 0.073), gender composition (p = 0.182), Boey score (p = 0.847) or MPI (Mannheim Peritonitis Index) score (p = 0.113) were noted between the two treatment groups. The mean lesion size of perforation was 5.0 ± 1.0 mm in the OTSC group and lesions were located to duodenal (n = 18), gastric (n = 7), and anastomotic (n = 1) regions. In the control group, the sites of perforation were not specified in up to 70 % of patients who were examined radiologically. The sites of perforation were located in the duodenal bulb (n = 17) and stomach (n = 7) in the patients who underwent laparotomy. The technical success rate in the OTSC group was 100 %. None of the patients experienced any complications associated with OTSC placement. The mean procedure time was 10.0 ± 2.5 min. The clinical success rate was 100 % (26 / 26) in the OTSC group and 57.5 % (46 / 80) in the control group. Subsequent surgery for diagnosis and/or treatment of recurrent ulcer after discharge from hospital including abdominal laparotomy and peritoneal lavage were not required for any of the patients in the OTSC group but were required for 30.0 % (24 / 80) of patients in the control group (p < 0.001). All patients in the OTSC group were discharged from the hospital. However, the mortality rate was 13.8 % (11 / 80) in the control group. The causes of death were uncontrolled sepsis followed by advanced multiple organ failure (8 / 11), heart failure (2 / 11), and gastrointestinal bleeding (1 / 11). The time to oral feeding was significantly shorter in the OTSC group (3.5 days, IQR 2.0 – 5.25) compared with the control group (7.0 days, IQR 5.0 – 9.0; p < 0.001). However, no significant difference was noted in the length of hospital stay (p = 0.439) or antibiotic use (p = 0.237).

The authors discussed that the advantage of OTSC therapy for perforated peptic ulcer therapy is that the procedure enables clear location of the perforation, accurate evaluation of lesion sizes, and reliable assessment of the patients' response to therapy.

The authors concluded that the OTSC-based endoscopic treatment of perforated peptic ulcers is associated with a very high clinical success rate and low complication rate.

Over-the-scope-clip applications for perforated peptic ulcer

Wei J, Xie XP, Lian TT, Yang ZY, Pan YF, Lin ZL, Zheng GW, Zhuang Z.

Surgical Endoscopy 2019; 33:4122-4127. <https://doi.org/10.1007/s00464-019-06717-x>

May 2020 | Ovesco FTRD® is recommended in the “ASGE guideline for endoscopic full-thickness resection and submucosal tunnel endoscopic resection” for the treatment of subepithelial tumors and epithelial neoplasia with associated fibrosis

The document reviews current techniques and devices used for EFTR, clinical applications and outcomes.

The Technology Committee of the American Society for Gastrointestinal Endoscopy (ASGE) prepared an ASGE guideline for endoscopic full-thickness resection and submucosal tunnel endoscopic resection, aiming to review current techniques and devices used for EFTR and to review clinical applications and outcomes. The document also contains a number of links to videos illustrating different techniques of exposed and non-exposed EFTR.

The authors state that with the development of reliable endoscopic closure techniques and tools, endoscopic full-thickness resection (EFTR) is emerging as a therapeutic alternative to surgical resection for the treatment of subepithelial tumors that involve the muscularis propria and for epithelial neoplasia with significant fibrosis.

EFTR techniques are classified in two groups: "exposed" and "non-exposed" EFTR. In exposed EFTR, the full-thickness resection is undertaken first, with subsequent closure of the defect, the temporary exposure of the peritoneal cavity to the intestinal lumen being the basis for the term "exposed". Further classification of the exposed techniques is made into tunneled and non-tunneled techniques. In non-exposed EFTR, a secure serosa-to-serosa apposition is achieved before full-thickness resection of the isolated lesion.

Ovesco's FTRD as the first FDA-cleared device for non-exposed EFTR in the colon is described with detailed technical data, a video illustration and a summary of some studies evaluating the device. A few studies directly comparing different resection and closure techniques are also presented.

The authors summarize that exposed and non-exposed EFTR techniques are emerging as less invasive alternatives to surgical resection of subepithelial tumors that involve the muscularis propria or mucosal neoplasia with associated fibrosis. However, the degree to which endoscopic resection can be used in place of oncologic surgical resection remains to be determined by additional data on variable lesion types, sizes, locations, and stages evaluated in future studies. Continued refinement of techniques and dedicated tools may expand the use of EFTR and make it accessible to more endoscopists.

ASGE guideline for endoscopic full-thickness resection and submucosal tunnel endoscopic resection

Aslanian HR, Sethi A, Bhutani MS, Goodman AJ, Krishnan K, Lichtenstein DR, Melson J, Navaneethan U, Pannala R, Parsi MA, Schulman AR, Sullivan SA, Thosani N, Trikudanathan G, Trindade AJ, Watson RR, Maple JT for the ASGE Technology Committee
VideoGIE. 2019 Jun 29;4(8):343-350.
doi:10.1016/j.vgie.2019.03.010.

April 2020 | OTSC® Proctology: rectovaginal fistula repair with the OTSC Proctology System is feasible, safe and effective

16 patients with rectovaginal fistulas, 13 of those having undergone previous surgical interventions, were treated with the OTSC Proctology. The healing rate at a median follow-up of 10.2 months was 43.7 %. Complications were pain in 4 patients and spontaneous clip detachment in 2 patients. No patients reported postoperative incontinence.

Y Tong et al., Michallon University Hospital, Grenoble, France, published a non-randomized prospective monocentric study evaluating feasibility, safety, and effectiveness of the OTSC Proctology for rectovaginal fistula repair. The OTSC Proctology is a minimally invasive technique avoiding sphincter damage, and has the advantage of maintenance of compression force on the internal fistula opening.

16 consecutive patients with a median age of 40.1 years were enrolled in the study. The causes of rectovaginal fistula were adverse events following proctectomy or pelvic radiotherapy (n = 6), obstetric trauma (n = 5), inflammatory bowel disease (n = 3), and unknown (n = 2). 13 patients had undergone previous surgical interventions before the clipping procedure, thereof 7 patients up to three surgical interventions and 6 patients more than 3 surgical interventions. The fistula tract was drained by loose setons in 6 patients and 11 patients had a temporary diverting stoma at time of the clipping procedure.

The placement of the OTSC Proctology was technically successful in all cases. Healing was defined as an absence of symptoms of leakage reported by the patient, and confirmed by clinical examination, radiological contrast control through the lower limb of the colostomy if present, and examination under general or epidural anesthesia. Successful primary healing of the fistula was observed in 7 patients (43.7 %) at 8 months and this result persisted at a

median follow-up of 10.2 months (range 8 – 36 months). Short-term postoperative complications were rare: one patient reported pain due to a perianal abscess that required an immediate removal of the clip and drainage. Two cases of spontaneous clip detachment were observed at 1-month follow-up. Chronic pain and ulceration around the clip led to surgical removal of the clip in three further patients. Most recurrences appeared shortly after the procedure, with a median delay of 45 days (range 16 – 217). No patients reported postoperative incontinence.

The authors concluded that the OTSC Proctology is feasible, safe and effective for rectovaginal fistula repair. Further trials are needed, incorporating comparisons to well-established techniques, a longer follow-up period, and a larger cohort.

Short-term outcomes of the over-the-scope clip proctology system for rectovaginal fistula repair: a prospective study

Tong Y, Trilling B, Sage P, Girard E, Faucheron J.
Techniques in Coloproctology 2019; 23:245-249.
https://doi.org/10.1007/s10151-019-01948-5

March 2020 | EMR+ enables a grasp-and-snare technique and consequently facilitates en-bloc resection of larger lesions compared to conventional EMR

In lesions > 2 cm, EMR+ outdoes its advantages: comparative trial demonstrates en-bloc resection rate of 86.36 % in 3 cm lesions with EMR+ vs. 18.18 % with conventional EMR

For larger sessile or laterally spreading gastrointestinal polyps ≥ 2 cm, en-bloc resection with the conventional EMR technique is technically very difficult. EMR with a piecemeal technique can be performed but this is associated with higher recurrence rates. The classical EMR technique can be improved with a new technique using an external additional working channel (AWC) that is mounted on a standard endoscope similar to the setup used with the FTRD® System. The technique is called EMR+.

Knoop RF et al., Department of Gastroenterology and Gastrointestinal Oncology, University Medical Center, Goettingen, Germany, evaluated the EMR+ technique prospectively in flat lesions comparing it with a conventional EMR. The study was conducted in an ex vivo animal model with porcine stomachs placed into the EASIE-R simulator. Prior to intervention, standardized lesions were set by coagulation dots, measuring 1 cm (n = 12 per group), 2 cm (n = 22 per group), 3 cm (n = 22 per group), or 4 cm (n = 20 per group).

Overall, 152 procedures were performed. In lesions of 1 cm, both EMR and EMR+ were very reliable with en-bloc resection rates (RO) of 100 %. In 2 cm lesions, EMR+ en-bloc resection rate was significantly higher (p = 0.02), conventional EMR already dropped to 54.55 % en-bloc resection rate, while EMR+ still achieved an en-bloc resection rate of 95.44 %. Conventional EMR did not provide sufficient resection rates for lesions of 3 cm or even 4 cm (18.18 % and 0 %). EMR+ still achieved very satisfying results in 3 cm lesions (86.36 %, p < 0.01) but also relevantly decreased at 4 cm (60.00 %, p < 0.01). From 3 cm on, EMR+ was significantly faster than conventional EMR. In 3 cm lesions, procedural time was 5 min with EMR+ vs. 12.5 min with conventional EMR, p < 0.01. In 4 cm lesions, procedural time was 5.5 min with EMR+ vs. 15 min with conventional EMR, p < 0.01.

In lesions of all sizes, the resection area was significantly larger in the EMR+ groups. In 1 cm lesions, the median resection area was 4.44 cm² for EMR+ vs. 3.14 cm² for conventional EMR, p = 0.012. In 2 cm lesions, it was 5.94 cm² for EMR+ vs. 3.30 cm² for conventional EMR, p < 0.01. At 3 cm it was 9.62 cm² for EMR+ vs. 1.50 cm² for conventional EMR, p < 0.01. In 4 cm lesions, EMR+ median resection area was 13.37 cm² vs. 4.02 cm² for conventional EMR, p = 0.03. A perforation rate of 15 % (3/20) was observed in the 4 cm-group treated with EMR+.

The authors summarized that EMR+ enables a grasp-and-snare technique and consequently facilitates en-bloc resection of larger lesions compared to conventional EMR. In lesions > 2 cm, EMR+ outdoes its advantages, especially concerning en-bloc resection rate. At 3 cm, EMR+ reaches

its best discriminatory power. At 4 cm, also EMR+ comes to its inherent limits and the risk of perforations rises. Then, ESD or surgery should be considered. The authors concluded that EMR+ could help to close a therapeutic gap in interventional endoscopy with manageable technical complexity, time and costs.

Endoscopic mucosal resection with an additional working channel (EMR+) in a porcine ex vivo model: a novel technique to improve en-bloc resection rate of snare polypectomy

Knoop RF, Wedi E, Petzold G, Bremer SCB, Amanzada A, Ellenrieder V, Neesee A, Kunsch S.
Endoscopy International Open 2020; 08: E99-E104.
https://doi.org/10.1055/a-0996-8050

Ovesco Research Update 34

March 2020 | OTSC® beyond standard use: new pictorial publication which explains the use of OTSC to salvage in challenging clinical situations

Because of its tissue apposition capabilities and easy en face application, OTSC is used as a single therapy or as part of a multimodal combination approach to treat complex defects of the luminal gastrointestinal tract

Moenkemüller K et al., Department of Gastroenterology, Otto-von-Guericke University, Magdeburg, Germany, published an atlas of examples of challenging clinical situations treated with the OTSC, aiming to describe its use beyond traditional indications. Special emphasis is laid on key technical aspects and potential limitations of the device in the context of complex endoluminal gastrointestinal disorders.

Diagnostic pictorial material of complex clinical conditions, pictorials of the endoscopic view before, during and after therapy, graphics showing the steps of OTSC application, and a supplementary video provide an overview on the topic as well as detailed information on the different clinical pictures and therapeutic steps.

Topics treated in the publication include OTSC application for bleeding ulcers, located in difficult positions, for prophylaxis and management of perforations associated with endoscopic resection of submucosal tumors, for closure of complex gastrointestinal wall defects such as gastro-pleural and entero-enteric fistulas, for esophageal stent fixation, for acquisition of a full-thickness specimen for diagnostic purpose and many more.

The authors emphasize, that although using the OTSC clip is generally straightforward, it is essential that the therapeutic endoscopist has a clear understanding of the advantages and potential limitations of this device.

The Use of the Over the Scope Clips Beyond Its Standard Use: A Pictorial Description

Moenkemüller K, Martínez-Alcalá A, Schmidt AR, Kratt T.
Gastrointest Endoscopy Clin N Am 30 (2020) 41-74.
https://doi.org/10.1016/j.giec.2019.09.003

March 2020 | OTSC® is associated with better procedural risk profile when compared to standard surgical correction in patients with non-resolving enterocutaneous fistula

In 7/10 patients with enterocutaneous fistula, OTSC application resulted in complete fistula closure. While NSQIP predicted rates of mortality, any complication, and median length of stay were 1.1 %, 34.5 %, and 9.5 days, respectively; patients successfully closed with OTSC experienced 0 mortalities, 0 complications, and had a median length of stay of 4 days.

J. Roy et al., Department of Surgery, University of South Alabama, Mobile, AL, USA, described their technical approach for OTSC closure of enterocutaneous fistulas (ECF) and the outcomes of ECF patients at their tertiary care hospital treated with OTSC between 07/2015 and 10/2017.

Enterocutaneous fistulas represent one of the most protracted and difficult problems in colorectal surgery with substantial morbidity and mortality. Closure rates without operative intervention in the era of advanced wound care and parenteral nutrition are markedly variable in reports, with most studies demonstrating closure rates in 30 – 35 % of cases. Conventional surgical intervention consists of extensive lysis of adhesions, bowel resection, and complex abdominal wall reconstruction and represents a potentially perilous endeavor. Mortality rates for ECF patients in different studies vary considerably but are as high as 10 – 20 %. Besides, ECF treatment is related to substantial costs for health care systems.

For the actual study, 10 patients (5 female, 5 male, average age 54.5 years, range 32 – 76) met the inclusion criteria. 7 patients had acute ECF (≤ 30 days old), 3 suffered from chronic ECF (> 30 days old). None of the patients had undergone prior abdominal or pelvic radiation. Median NSQIP-predicted risks for surgical correction of the fistula were 1.1 % for mortality, 34.5 % for periprocedural complications and 9.5 days for the length of stay.

Technical success was achieved in 9/10 OTSC applications. The overall clinical success rate for fistula closure was 70 %. Closure of acute fistulas resulted in an 86 % success rate (6/7), while closure of chronic fistulas resulted in a 33 % (1/3) success rate. There were no mortalities or operative complications associated with OTSC application, the median length of stay was 4 days.

The authors concluded that the OTSC is an effective adjunctive measure to improving rates of successful ECF closure. OTSC treatment of non-resolving ECFs can substantially decrease morbidity, hospitalization time, need for major surgical intervention, and possibly mortality.

Endoscopic technique for closure of enterocutaneous fistulas

Roy J, Sims K, Rider P, Grimm L, Hunter J, Richards W. *Surg Endosc.* 2019 Oct;33(10):3464-3468. <https://doi.org/10.1007/s00464-018-06646-1>

February 2020 | Conference Report of the American College of Gastroenterology 2019 Annual Scientific Meeting-ACG 2019

- **OTSC®:** Besides a RCT showing OTSC's superiority in NVUGIB treatment, many case reports were presented documenting successful OTSC use in critical clinical cases
- **FTRD®:** Successful R0 resection of difficult adenomas, subepithelial lesions and early carcinomas located in the upper and lower GI tract by EFTR with the FTRD System
- **BougieCap:** The BougieCap allows for direct visual control and optical feedback during bougienage procedures

The American College of Gastroenterology 2019 Annual Scientific Meeting (ACG) 2019 was held on October 25-30, 2019, in San Antonio, TX USA. Ovesco products were presented during various plenary and poster sessions throughout the conference. A summary of those presentations can be found below:

OTSC System

Besides a RCT showing OTSC's superiority in NVUGIB treatment, many case reports were presented documenting successful OTSC use in critical clinical cases

OTSC treatment of NVUGIB is associated with 85 % lower rebleeding rate when compared to standard hemostasis methods

Jensen DM et al., CURE Digestive Diseases Research Center, Los Angeles, CA, United States presented the results from their recent Randomized Controlled Trial (RCT) that was conducted at 2 academic medical centers which evaluated the outcomes of patients treated with OTSC for initial endoscopic treatment of severe non-variceal upper gastrointestinal bleeding (NVUGIB) to patients treated with standard hemostasis RX methods (e.g. hemoclips and/or multipolar probe with epinephrine pre-injection). 49 patients who met clinical and esophagogastroduodenoscopy (EGD) criteria for PUB's or UGI Dieulafoy's lesions (with major SRH or spots with arterial blood flow) were randomized in a 1:1 allocation to either standard RX (as described in above RCT) or OTSC. All patients received high dose PPI's after randomization and were followed prospectively for 30 days. The proportion of patients with 30-day rebleeding was significantly higher with standard RX – 28 % (7/25) – compared to OTSC RX – 4.2 % (1/24). The OTSC rebleeding rate was 85 % lower than standard group (relative risk 0.149 with 95 % confidence intervals of 0.006, 0.863). The number needed to treat (NNT) was 4.2. In patients with severe UGI bleeding from PUB's or Dieulafoy's lesions, primary endoscopic treatment with OTSC resulted in significantly lower rates of rebleeding and severe complications and fewer RBC transfusions than standard endoscopic hemostasis. Results appear to relate OTSC's superior ability to obliterate arterial blood flow underneath SRH and thereby reduce lesion rebleeding.

Randomized Controlled Trial (RCT) of Over-the-Scope Clip (OTSC) as Initial Endoscopic Treatment of Severe Non-Variceal Upper Gastrointestinal Bleeding (NVUGIB)

Jensen DM^{1,2,3}, Kovacs TOG^{1,2}, Ghassemi KA^{1,2}, Kaneshiro M^{1,2}, and Gombein, J⁴

¹CURE Digestive Diseases Research Center, Los Angeles, CA, United States. ²David Geffen School of Medicine at UCLA and Ronald Reagan UCLA Medical Center, Los Angeles, CA, United States. ³West Los Angeles Veterans Administration Medical Center, Los Angeles, CA, United States. ⁴Department of Medicine – GIM - DOMSTATS, UCLA, Los Angeles, CA, United States

The OTSC as valuable tool for salvage hemostasis in duodenal ulcer

Meir Mizrahi, MD et al, University of South Alabama College of Medicine, Mobile, AL, presented a case where traditional hemostasis methods failed, and an OTSC was used as salvage/rescue hemostasis therapy. The patient was a 45-year-old Asian male with a history of cirrhosis secondary to alcohol presented with jaundice and melena with a hemoglobin of 5.9 gm/dl. EGD showed severe duodenitis. A side viewing duodenoscope showed a large, cratered duodenal bulb ulcer with a visible vessel. Epinephrine was injected and 5 hemoclips were placed. Despite this, there was still active oozing of blood. All the hemoclips were removed with snare. OTSC System was loaded onto a standard gastroscope and advanced to the lesion. Target tissue was sucked into the cap and OTSC clip was released by turning the hand wheel. The clip closed itself and firmly anchored the tissue resulting in complete hemostasis. Proton pump inhibitors were continued and follow up endoscopy at 1 month showed healed ulcer site with clip in place. While OTSC have been well established in leaks, perforations and fistulas, their role in UGI hemorrhage is evolving. Several randomized trials comparing OTSC with conventional endoscopic hemostasis therapies are ongoing and results are awaiting, however our case reiterates that OTSC platform should be available and be part of endoscopic hemostasis armamentarium.

Dominant Claw Control: Role of OTSC Clip as Salvage Hemostasis Therapy

Gilad Shapira, MD¹, Manoj Kumar, MD², Yazan Fahmawi, MD², Meir Mizrahi, MD³

¹University of South Alabama College of Medicine, Mobile,

AL; ²University of South Alabama, Mobile, AL; ³University of South Alabama College of Medicine, Spanish Fort, AL

Clip retention time of OTSC is median 8.9 months, clip location and indication for placement are key factors in time to clip loss

Xiao Jing Wang, MD et al, Mayo Clinic, Rochester, MN presented their findings from a retrospective review of OTSC usage from 1/1/2011 to 1/1/2019. The primary endpoint was clip retention time (time when the clip was placed to time the clip was last seen or was physically removed). Secondary endpoints were rate of complications and factors that influence clip retention time and clinical success. A total of 358 OTSC were placed in 299 patients (55.8 % male, n = 184) with mean age 64 + 15.3 years. Clip placement was clinically successful in 81 % of cases (n = 265). Median time to clip loss as assessed by Kaplan-Meier curve was 267 days (95 % CI 196 - 406). Adverse events occurred in 34 clip placements (10.7 %) and included bleeding (n = 8, 24.2 %), perforation (n = 5, 15.2 %), ulceration (n = 3, 9.1 %), migration (n = 15, 45.5 %), and infection (n = 1, 3.0 %). Univariate analysis revealed that age, gender, type of OTSC clip used, and performing endoscopy did not impact clip retention rate. A size 14 clip had marginally higher rates of retention compared to others. Clips placed in the rectum or esophagus were lost earlier compared to those in the small/large intestine and stomach. Clips placed for fistulas and those associated with adverse event were associated with earlier clip loss. In evaluation of clinical success, clips placed for fistulas had a relative risk of clinical failure of 2.08 (95 % CI 1.29 - 3.36) compared to all other indications. OTSC are retained for a median of 8.9 months with clip location and indication for placement being key factors in time to clip loss.

Clip It and Forget It? A Natural History of 358 Over-the-Scope Clips

Xiao Jing Wang, MD¹, Daniel Gonzalez Izundegui, MD¹, Sneha Singh, MBBS¹, Gerardo Calderon, MD², Oscar Garcia Valencia, MD³, Andrew C. Storm, MD⁴

¹Mayo Clinic, Rochester, MN; ²Indiana University Hospital, Rochester, MN; ³University of Miami, Rochester, MN; ⁴Mayo Clinic College of Medicine and Science, Rochester, MN

Post-operative duodenal-enteric fistula in patient with multiple gunshot wounds successfully closed with OTSC

David H. Kruchko, et al, Advocate Lutheran General Hospital, Chicago, IL presented a case where an OTSC was used successfully for immediate closure of a duodenal-enteric fistula caused by a gunshot wound, which saved the patient from undergoing additional surgery. A 41-year-old man presented with multiple gunshot wounds to the abdomen. An emergent exploratory laparotomy was performed, with partial small bowel resection for gunshot wounds to the jejunum and ileum. His postoperative course was complicated by an intra-abdominal abscess. A sinogram revealed a duodenal-enteric fistula with a fluid collection, requiring Jackson Pratt drain placement and antibiotic therapy. Given the lack of improvement, the patient was taken for EGD. A 10 mm fistulous tract was identified at the junction of the 2nd and 3rd portion of the duodenum. A 12/6t OTSC was deployed over a pediatric colonoscope. The scope was passed to the level of the fistula, with fluoroscopy for localization. Using suction, the fistula was pulled into the cap of the OTSC device and the clip was successfully deployed, closing the fistula. Clip placement was confirmed via fluoroscopy. The next day, upper gastrointestinal stomach X-Ray was completed and did not show extravasation to suggest a leak at the duodenum. His diet was advanced and after a prolonged hospital course, he was discharged.

Endoscopic Closure of Post-Operative Duodenal-Enteric Fistula Following Repair of Abdominal Gunshot Wound

David H. Kruchko, DO¹, Natasha Shah, MD², Ryan T. Hoff, DO², Asif Lakha, MD²

¹Advocate Lutheran General Hospital, Chicago, IL; ²Advocate Lutheran General, Hospital, Park Ridge, IL

Boerhaave's Syndrome successfully repaired by placement of 2 OTSC clips and esophageal stent

Eric Then, MD et al, St. Barnabas Hospital, Bronx, NY,

presented a case of Boerhaave's Syndrome successfully repaired by using 2 OTSC clips and placement of an esophageal stent. Patient is a 65-year-old male with a medical history of chronic obstructive pulmonary disease who presented the emergency department with chest pain, shortness of breath and productive cough after an episode of retching 2 days prior. Notable workup in the emergency department included a chest x-ray which showed left lower lobe pneumonia, pneumothorax, and a left pleural effusion. A chest tube was then placed and a computed tomography angiography of the chest was also done. This showed pneumomediastinum with air around the mid to distal esophagus. A barium esophagram was subsequently done which showed contrast leak at the distal esophagus. An upper endoscopy was then performed and showed a large linear perforation near the gastroesophageal junction. This was closed with 2 OTSC clips in addition to placement of an overlapping partially covered stent. In order to ensure complete closure an esophagram was done with radioisotope. This showed complete closure of the perforation, with no evidence of leakage of the radiotracer. After the procedure the patient was able to tolerate oral feedings with no clinical sequelae. During his hospital stay the patient completed treatment with antibiotics, his chest tube was successfully removed and he was subsequently discharged home. Surgical repair is currently the treatment of choice in BS. Present day, endoscopic measures are challenging this consensus. More recently cases of BS have been successfully treated with the use of endoscopic suturing with esophageal stent placement. Our case was successfully treated by using 2 OTSC and placement of an esophageal stent.

Successful Endoscopic Repair of Boerhaave's Syndrome: Changing the Narrative

Eric O. Then, MD¹, Michell Lopez, MD¹, Andrew Ofosu, MD², Vijay Gayam, MD³, Vijay S. Are, MD⁴, Vinaya Gaduputi, MD, FACC⁵

¹St. Barnabas Hospital, Bronx, NY; ²Brooklyn Hospital Center, Brooklyn, NY; ³Interfaith Medical Center, Brooklyn, NY; ⁴Stormont Vail Health System, Topeka, KS; ⁵SBH Health System, Bronx, NY

Staple-Line Leak Following Sleeve Gastrectomy successfully closed with OTSC

Augustine Tawadros, MD et al, Robert Wood Johnson Medical School, Rutgers University, New Brunswick, NJ, presented a case of endoscopic management with an OTSC for a staple-line leak following gastrectomy. A 38-year-old woman, with limited past medical history, presented 6 weeks after SG with a staple-line leak. Her immediate post-op period had been complicated by a peri-splenic abscess, which was managed with IV antibiotics (ABX) and percutaneous drainage (PD). Upper GI series (UGIS) at this time did not reveal a SLL. After clinical and radiographic resolution, the drain was removed. Several days later, the patient presented with nausea, vomiting, and non-bloody watery diarrhea. She was found to be febrile and tachycardic; labs, including a test for *Clostridium difficile*, were negative. A repeat CT scan confirmed an upper abdominal fluid collection. Subsequent UGIS revealed a SLL in the proximal stomach with a suture line fistula. After several weeks of NPO status, ABX and PD failed to lead to leak closure, the patient underwent endoscopy with successful placement of a gastric OTSC clip. Repeat UGIS confirmed fistula closure. After several weeks of additional conservative management, the PD was removed. Several months later, the patient remains stable.

Endoscopic Management With an Over-the-Scope Gastric Clip of a Staple-Line Leak Following Sleeve Gastrectomy

Augustine Tawadros, MD¹, Kevin Skole, MD², Wai Y. Chau, MD³

¹Robert Wood Johnson Medical School, Rutgers University, New Brunswick, NJ; ²University Medical Center of Princeton at Plainsboro, Plainsboro Township, NJ; ³University Medical Center of Princeton at Plainsboro, Plainsboro, NJ

Successful closure of Boerhaave with use of OTSC alone

Muaataz Azawi, MD et al, Elmhurst Medical Center, Astoria, NY, reported a case of a 57-year-old male who presented with trauma in setting of alcohol intoxication. In the ED

patient was found to have a left hemopneumothorax found in CT scan. The patient was admitted for lactatemia and chest tube management. Afterwards patient developed sepsis due to empyema which prompted him to undergo for VATs procedure. Patient remained septic and CT with contrast revealed an esophageal leak and Boerhaave's syndrome was suspected. GI was consulted and patient underwent EGD. A 8 mm perforation was found in the left latero-posterior portion of the esophagus consistent with Boerhaave syndrome. To repair the defect, the tissue edges were approximated using suction and one 11 mm OTSC was successfully placed (Bear claw). Thereafter patient underwent g-tube and j-tube placement. Next, the patient was discharged home and was seen in the clinic, where the g-tube and j-tube was removed.

The current management of BS includes conservative, endoscopic, and surgical treatments. The survival rates for each treatment are 75 %, 100 %, and 81 %, respectively. The estimated mortality rate of patients with BS is 20 % – 40 %. Endoscopy can be utilized to confirm the diagnosis of BS, particularly in equivocal cases. Endoscopic treatment with stenting can be effective and less invasive. In a systematic review that, endoscopic stenting had a success rate of 81 % but endoscopic intervention was required in 17 % and surgical reintervention in 10 % patients. More recently, the use of OTSC has been reported in treatment of boerhaave's alone or in combination with local therapies or esophageal stent. These endoscopic interventions are typically performed within several days of the perforation. Here we report successful closure of boerhaave with use of OTSC alone after 14 days from perforation.

Boerhaave's Syndrome Was Treated With Over-the-Scope Clip OTSC Closure (Bear Claw Closure): A Rare Case Report

Muaataz Azawi, MD¹, Soohwan Chun, MD², Raghav Bansal, MD³

¹Elmhurst Hospital Center, Astoria, NY; ²Elmhurst Hospital Center, Elmhurst, NY; ³Icahn School of Medicine at Mount Sinai - Elmhurst, Elmhurst, NY

Gastric perforation by foreign body successfully managed by removal of foreign body with endoscopic forceps and perforation closure with OTSC

Andy Tien, MD et al, Kaiser Permanente Los Angeles Medical Center, Los Angeles, CA, presented a case of non-surgical endoscopic management of a sealed nail perforation in the gastric antrum.

A 48-year-old homeless male with a history of multiple abdominal surgeries due to gunshots and stab wounds complicated by abdominal hernias presented with 1 week of abdominal pain after drinking an open bottle of soda. Vital signs were stable on presentation. Labs were unremarkable. Computed tomography (CT) revealed a needle-like foreign body in the antrum of the stomach that appeared to project beyond the lumen without signs of pneumoperitoneum free perforation or fluid. An upper endoscopy was performed to further characterize the foreign body. Upon entering the stomach, a 2 mm round metallic circular object was visualized along the lesser curvature of the antrum, suggestive of a nail embedded into the gastric antrum. General surgery was consulted given concern for perforation, however, the patient was deemed a high-risk surgical candidate due to hostile abdominal scarring. The decision was made for endoscopic repair of the nail perforation. On repeat endoscopy, the nail was again visualized in the lesser curve of the antrum; rat-tooth forceps were used to remove a 2.8 cm rusty nail which was safely pulled through an overtube. The scope was reinserted with an OTSC attached, which was then successfully deployed at the site of perforation. Patient tolerated the procedure well with no signs of peritonitis and follow-up CT four days post-procedure showed the clip in place. There are limited cases in the literature of endoscopic intervention in a sealed nail perforation in the stomach in an adult. We found that an OTSC clip was effective in sealing a small perforation without evidence of peritonitis on follow-up examinations. Coordination with the surgical service is important in endoscopic repair of perforations. Overall, endoscopic removal of a sealed nail perforation is a viable option for high-risk patients.

Non-Surgical Endoscopic Management of a Sealed Nail Perforation in the Gastric Antrum

Andy Tien, MD, MS, Elizabeth Dong, MD, Marianne Fahmy, MD

Kaiser Permanente Los Angeles Medical Center, Los Angeles, CA

Successful endoscopic closure of pyelo-duodenal fistula using the OTSC System

Erica Duh, MD et al, University of California Irvine Medical Center, Orange, CA, presented a case of successful endoscopic closure of a pyelo-duodenal fistula using the OTSC System. A 62-year-old female presented with sharp, severe abdominal pain. She denied back pain, hematuria, fever, and dysuria. She had no significant past medical history and did not take any medications. CT abdomen and pelvis demonstrated a large right-sided staghorn calculus extending into the renal pelvis. The patient underwent percutaneous nephrolithotomy, which was complicated by perforation of the renal pelvis and abdominal compartment syndrome. Subsequently, a proximal ureteral stricture was refractory to stenting and balloon dilation. CT urography demonstrated a PDF. After placement of a metallic ureteral stent and nephrostomy tube, a nephrostogram was negative for extravasation, so the patient underwent EGD to evaluate for the resolution of the PDF. During the procedure, iohexol contrast and methylene blue were inserted via nephrostomy tube, and fluoroscopic imaging showed contrast in the right kidney and bladder without extravasation. However, on endoscopy, methylene blue extravasation led to identification of two mucosal defects in the anterior duodenal sweep. The larger defect was closed using the helix tissue retractor and OTSC System (12/6 t-type). The smaller defect was closed with 2 hemostasis clips. Follow-up renal scan and CT imaging showed no evidence of contrast extravasation. The OTSC is an innovative endoscopic technique most often used for hemostasis and for treatment of endoscopic or surgical complications (e.g., fistula, leakage, or perforation). Most reports document successful OTSC use for fistula closure including gastrocutaneous and esophageal fistulas. Our case is, to our knowledge, the second ever documented PDF closure using the OTSC System.

Successful Endoscopic Closure of a Pyeloduodenal Fistula Using the Over-the-Scope Clip System

Erica Duh, MD¹, Meredith Clary, MD², Jason Samarasena, MD, FACC¹, Ralph Clayman, MD³, Kenneth J. Chang, MD, FACC²

¹University of California Irvine Medical Center, Orange, CA;

²University of California Irvine, Orange, CA; ³University of California Irvine, Irvine, CA

Duodeno-pleural fistula secondary to transarterial chemoembolization successfully treated with OTSC

Peter Natov, BSc et al, Center and Tufts University School of Medicine, Boston, MA presented a case of endoscopic closure of a duodeno-pleural fistula secondary to transarterial chemoembolization. A 49-year-old woman with a history of chronic hepatitis B cirrhosis and HCC presented to the emergency room with right-sided chest pain. She was diagnosed with HCC 16 years ago in China and underwent partial hepatectomy and chemotherapy. Two years prior to current presentation, MRI demonstrated a 2.0 x 2.7 x 3.3 cm liver lesion consistent with recurrent HCC. TACE was performed. She subsequently reported abdominal pain. Endoscopy discovered a duodenal ulcer, and imaging noted a developing phlegmon and free air. Management with bowel rest and parenteral nutrition led to clinical improvement without the need for surgery. On this presentation, she described new onset pleuritic chest pain and dyspnea. Imaging revealed a large loculated right pleural effusion with gas foci and visible extraluminal enteric contrast surrounding the duodenal bulb, suggestive of an empyema secondary to a duodeno-pleural fistula. Pleural fluid analysis showed a pH < 6.8, high amylase and lipase levels, and cultures of multiple enteric organisms. Fistula closure was attempted endoscopically. A 4 mm opening was found in the duodenal bulb. The mucosa was denuded with argon beam plasma coagulation and an 11/6t OTSC successfully closed the defect. Extravasated contrast was not seen on imaging two days later.

Endoscopic Closure of a Duodeno-Pleural Fistula Secondary to Transarterial Chemoembolization

Peter S. Natov, BSc1, Sukeerti Kesar, MD2, Robert F.

Yacavone, MD3

¹Tufts Medical Center and Tufts University School of Medicine, Boston, MA; ²Tufts Medical Center, Boston, MA; ³Tufts University Medical Center, Boston, MA

Gastro-colo-cutaneous fistula with OTSC closed with hemoclips and OTSC

Hassan Zia, MD et al, University of Oklahoma Health Sciences Center, Oklahoma City, OK, presented a case of endoscopic closure of a gastro-colo-cutaneous fistula with an OTSC System. A 55-year-old male with Stage IV gastroesophageal junction neuroendocrine carcinoma with ongoing chemoradiotherapy had a PEG tube placed for dysphagia. There were no post procedural complications. Two months post procedure, patient presented with pancytopenia and a two-week history of persistent diarrhea, malodorous peristomal drainage and weight loss. Computerized tomography (CT) imaging with intravenous (IV) contrast confirmed the presence of a gastro-colo-cutaneous fistula with internal bolster located in the transverse colon. Due to immunocompromised status, endoscopic closure was pursued with a colonoscopy and esophagogastroduodenoscopy (EGD). A pediatric colonoscope with distal cap attachment was used, an internal bolster with colo-cutaneous fistula with no cologastric fistula was seen. The site was tattooed with 6 mL of spot (carbon black). The gastrostomy tube was removed, and three hemostatic clips were used for closure. No air leak was noted over the exterior abdominal wall. This was followed by EGD, the residual gastrocolic fistula was closed by a OTSC 10/3t. A CT scan with oral contrast demonstrated no residual leak/fistula. Total parenteral nutrition (TPN) was initiated. Patient continued to do well, and post-procedural hospital stay was uneventful.

Endoscopic Closure of GastroColoCutaneous Fistula: A Rare Complication of Percutaneous Endoscopic Gastrostomy Tube Placement

Hassan Zia, MD, Maham Hayat, MD, Ralph Guild, MD
University of Oklahoma Health Sciences Center, Oklahoma City, OK

Two patients, in which the OTSC System was successfully used as salvage therapy for refractory NVUGIB after failed conventional management

Shivantha Amamath, MD et al, Staten Island University Hospital, Northwell Health, Staten Island, NY, presented a case series of OTSC as salvage therapy after failed conventional management.

85-year-old female with Pulmonary Hypertension and Atrial Fibrillation was admitted for hypovolemic shock from hematemesis. Endoscopy revealed a large posterior duodenal bulb ulcer with adherent clot for which epinephrine injection was used. The ulcer's size and location prompted angiography and coil embolization of the gastroduodenal artery (GDA). Patient was in hypovolemic shock again due to a massive bleed. She was high risk for surgical intervention due to pulmonary hypertension, and repeat endoscopy was performed for hemostasis. An 11/6t-type OTSC clip was successfully deployed at a large vessel oozing blood, after removal of a large adherent clot by snare. No further intervention was needed afterward. 76-year-old male with extensive cardiac history was admitted with burns and inhalation injury from a house fire, requiring tracheostomy and gastrostomy placement. During hospitalization, patient had multiple episodes of melena and required blood transfusion. Endoscopy revealed two large anterior and posterior duodenal bulb ulcers, with a visible vessel, managed by electrocautery. Due to persistent bleeding, he underwent coil embolization of the GDA. Two days later, patient had recurrent bleeding. Repeat endoscopy with OTSC was used for successful hemostasis of a large visible vessel in the posterior duodenal bulb. Traditional endoscopic management carries a 20% rebleeding risk and successful secondary hemostasis drops from 90% to 75%. Anatomic variation in blood supply to the duodenal bulb, either from branches of celiac axis or superior mesenteric artery may be the reason why our patients failed angiographic coil embolization of the GDA. To our knowledge, we report the first two cases where OTSC was used for successful salvage of NVUGIB that failed initial endoscopic and angiographic interventions.

Over-the-Scope Clip as Salvage Therapy for Refractory Non-Variceal Upper Gastrointestinal Hemorrhage After**Failed Conventional Management: A Case Series**

Shivantha Amamath, MD, BSc (Hons) Biology, MRSB¹, Jobin Philipose, MD², Jeffrey Abergel, MD¹, Hafiz Khan, MD²

¹Staten Island University Hospital, Northwell Health, Staten Island, NY; ²Staten Island University Hospital, Staten Island, NY

FTRD System**Successful R0 resection of difficult adenomas, subepithelial lesions and early carcinomas located in the upper and lower GI tract by EFTR with the FTRD System****Esophagus dilation was performed to enable insertion of colonic FTRD in the stomach for resection of suspected GIST**

Meir Mizrahi, MD et al, University of South Alabama College of Medicine, Mobile, AL, presented a case of endoscopic full-thickness resection (EFTR) of stomach lesion using colonic FTRD using a novel technique for upper esophagus dilation to assist in the insertion of colonic FTRD through the mouth. A 49-year-old female was referred for endoscopic resection of submucosal lesion suspected for gastric stromal tumor (GIST) prior to bariatric gastric bypass surgery. EGD revealed the lesion in the gastric body. Endoscopic ultrasound (EUS) confirmed it as well demarcated hypoechoic lesion measuring 0.8 cm x 0.8 cm arising from the intersection of submucosal layer and muscularis layer without adjacent lymphadenopathy. FTRD under general anesthesia was planned. EFTR technique: To assist the insertion of colonic FTRD, upper esophagus was dilated with Savary dilator size 51 French for 1 minute, then with 55 French for 1 minute and then with 60 French for 5 minutes. Pediatric colonoscope installed with FTRD was gently advanced into the stomach lumen. Circumferential marking of the lesion was done using FTRD marking probe. Favorable positioning was achieved. FTRD grasper was used to grasp and mobilize the tissue into the cap. Clip was applied and snare was hooked to high frequency generator performing high frequency snare resection. Endoscope was removed along with resection specimen. Histopathology revealed fibrotic nodule with complete resection (R0). Endoscopy at 3-month revealed no recurrence.

Overcoming the Challenge of Full Thickness Resection of Gastric Lesions Using Colonic Full Thickness Resection Device (FTRD): Novel Technique

Yazan Fahmawi, MD¹, Manoj Kumar, MD¹, Meir Mizrahi, MD²

¹University of South Alabama, Mobile, AL; ²University of South Alabama College of Medicine, Spanish Fort, AL

Meta-analysis including 480 patients with FTRD treatment demonstrates pooled R0 resection rate of 81.12% and technical success rate of 90.39%

Meir Mizrahi, MD et al, University of South Alabama College of Medicine, Mobile, AL, presented a meta-analysis findings from a study conducted to confirm efficacy and safety of the FTRD System. A comprehensive literature review was performed. Studies with < 5 cases were excluded. Rates of histologic complete resection (R0), technical success, and complications were extracted. Pooled estimates and the 95% CI were calculated depending on heterogeneity. Heterogeneity was assessed using I² statistics. Nine studies (5 retrospective and 4 prospective) including 480 patients with 489 lesions were included. Indications were difficult adenomas (355), early carcinoma (82), and subepithelial lesions (SEL) (47). Locations of the lesions: stomach (5), duodenum (20), proximal colon (220), distal colon (98), and rectum (141). 57% of patients were male with median of age 66.5 years (range 20 - 92). Mean size of the lesions was 15.36 mm (range 1 - 40). Median procedure time was 46.56 minutes (range 3 - 190). Mean hospital stay was 3.23 days (range 0 - 11). Pooled overall R0 resection rate and technical success rate were 81.12% (95% CI: 77.63 - 84.6; 156.7%) and 90.39% (95% CI: 87.76 - 93.00; 110.2%), respectively. The pooled R0 resection rate was 81.59% (95% CI: 76.96 - 85.63) for difficult adenomas, 83.47% (95% CI: 73.19 - 91.00) for early carcinomas, and 81.59% (95% CI: 68.25 - 91.05) for SELs. For locations, R0 resection rate were 78.42% (95% CI: 72.17 - 85.50),

76.97% (95% CI: 62.56 - 87.88), and 80.59% (95% CI: 67.56 - 90.09) for proximal, distal colon, and rectum, respectively. The complications were as follows: minor bleeding 2.89 (95% CI: 1.60 - 4.76), major bleeding 0.88 (95% CI: 0.26 - 2.15), perforation 2.23 (95% CI: 1.12 - 3.95), postpolypectomy syndrome 1.26% (95% CI: 0.48 - 2.68), local trauma 0.9% (95% CI: 0.27 - 2.19), and other complications 1.91% (95% CI: 0.90 - 3.54). Of 50 patients with peri-appendicular lesion, the rate of appendicitis was 14.29% (95% CI: 6.38 - 26.23). The rate of complications that required surgery was 1.55 (95% CI: 0.66 - 3.06). No procedural death related was found. There was no correlation between mean lesion's size and complications rates.

Efficacy and Safety of Endoscopic Full-Thickness Resection Device (FTRD) for Gastrointestinal Lesions**Full-Thickness Resection: A Meta-Analysis Study**

Yazan Fahmawi, MD¹, Manoj Kumar, MD¹, Gilad Shapira, MD², Nicholas Ludvik, MD², Meir Mizrahi, MD³

¹University of South Alabama, Mobile, AL; ²University of South Alabama College of Medicine, Mobile, AL; ³University of South Alabama College of Medicine, Spanish Fort, AL

Single Center Experience with FTRD System shows difficulties in endoscope advance / suction of fibrotic lesions, but high R0 resection rate

Meir Mizrahi, MD et al, University of South Alabama College of Medicine, Mobile, AL, presented data from a retrospective case series of 17 patients who underwent EFTR with FTRD their institution. Complete resection (R0), technical success, and clinical outcomes were reviewed. The indications for EFTR were recurrent/residual adenoma (12 cases), primary non-lifting adenoma (2 cases), submucosal lesions resection (2 cases), and early carcinoma (1 case). Eleven patients were males and six were females. The median age was 60 years (range 32 - 84). Locations of the lesions were: (3) in the cecum, (8) in the ascending colon, (2) in the transverse colon, (2) in the hepatic flexure, (1) in the rectum, and (1) in the stomach. Technical success rate was 13/17 (76%). Three of the unsuccessful attempts were due to inability to advance the endoscope to the cecum in one case and inability to pull up fibrotic lesions in the other two cases. The fourth unsuccessful case was due to device malfunction as the lesion was snared but the clip was not deployed. The mean maximum diameter of lesions after resection was 18.3 mm (range 13 - 22). R0 resection rate was 14/14 (100%). Microscopic full-thickness resection was achieved in 43% (6/14). The median procedure time was 67 minutes (range 31 - 140).

Endoscopic Full-Thickness Resection of Gastrointestinal Lesions With the Over-The-Scope Device: A Single Center Experience

Yazan Fahmawi, MD¹, Manoj Kumar, MD¹, Baret Bercier, MD², Rufaat Mando, MD³, Meir Mizrahi, MD⁴

¹University of South Alabama, Mobile, AL; ²University of South Alabama College of Medicine, Mobile, AL; ³University of South Alabama College of Medicine, Mando, AL; ⁴University of South Alabama College of Medicine, Spanish Fort, AL

Successful R0 resection of residual sigmoid carcinoid tumor by FTRD

Pratik Patel, MD et al, East Carolina University - Vidant Medical Center, Greenville, North Carolina presented a technically and histologically successful full-thickness resection (FTR) of a sigmoid goblet cell carcinoid tumor. A 61-year-old male presents for screening colonoscopy and was found to have an 8 mm submucosal mass in the distal sigmoid. Pathology of biopsies showed hyperplastic colonic mucosa with marked reparative stromal and epithelial changes and mixed features of prolapse. He was referred for EUS of the submucosal lesion, which revealed hypoechoic submucosal lesion in distal sigmoid measuring 7.4 x 4.7 mm. FNA was attempted but unsuccessful as the mass was small and mobile. Therefore, decision was made to perform endoscopic mucosal resection (EMR). Pathology on EMR of sigmoid lesion revealed goblet cell carcinoma. Deep and lateral margins were positive. Chromogranin A level was mildly elevated. He was referred to surgical oncology who recommended 3-month surveillance colonoscopy given low Ki index. CT Scan of Abdomen and Pelvis as well as Octreotide scan were both negative for

distant metastases. At 3-month follow up colonoscopy, the decision was made to use Full-Thickness-Resection-Device (FTRD) of the sigmoid lesion given positive margins. Pathology of FTRD was consistent with neural sheath tumor and deep and lateral margins were negative. Repeat colonoscopy revealed well healed scar at site of FTRD and cold biopsies around the scar were benign. He will follow up in 3 years for repeat colonoscopy.

Full-Thickness-Resection-Device for the Removal of a Submucosal Sigmoid Carcinoid Tumor

Pratik Patel, MD¹, Ethan Phan, MD¹, Narasimha Swamy Gollol Raju, MD¹, Prashant Mudireddy, MD²

¹East Carolina University, Greenville, NC; ²Brody School of Medicine, Greenville, NC

Successful R0 resection of duodenal neuroendocrine tumor using the FTRD Device

Andrew Groff, BS et al, Penn State University Milton S. Hershey Medical Center, Hershey, PA, presented a successful R0 resection of a T2, grade 2 neuroendocrine tumor using the FTRD. A 68-year-old-male with history of Barrett's esophagus was found to have a low-grade neuroendocrine tumor of the duodenal bulb on surveillance upper endoscopic screening. Computed tomography (CT) imaging of the abdomen revealed a hypervascular duodenal bulb mass, without obvious extension to the wall, lymphadenopathy, or evidence of metastasis. Endoscopic ultrasound (EUS) with fine needle aspiration (FNA), revealed a 1.2 cm mass arising from the submucosal layer with pathology indicating a low-grade neuroendocrine tumor. Due to the small tumor size, EFTR by FTRD was planned. Upon advancing the endoscope with the FTRD System attached, a non-bleeding neuroendocrine tumor was visualized. The mass was carefully pulled into the cap and after banding, ligation and cutting; there was a full-thickness en-bloc resection with all layers of the duodenum apparent. Pathology of the specimen documented a grade 2 neuroendocrine tumor with clean margins, confirming R0 resection. The patient was discharged the same day without associated complications. The FTRD System offers a novel endoscopic treatment approach to tumors extending beyond the submucosa that require full-thickness resection and would otherwise require higher risk surgical intervention. The technology allows for safe tissue closure, valid histologic evaluation of the en-bloc specimen, and minimal thermal injury, thus reducing the morbidity and cost otherwise associated with surgery.

A Minimally Invasive Approach to Resection of a Duodenal Neuroendocrine Tumor Using the Full-Thickness-Resection-Device

Leonard Walsh, MD¹, Andrew Groff, BS², Matthew Moyer, MD, MS¹

¹Penn State University Milton S. Hershey Medical Center, Hershey, PA; ²Penn State University Milton S. Hershey Medical Center and College of Medicine, Hershey, PA

Successful full-thickness resection of a recurrent colonic adenoma using the FTRD

Pradeep Yarra, MD et al, University of Kentucky, Lexington, KY presented a case of complete adenoma resection using the FTRD System. A 66-year-old male with a past medical history of CAD s/p CABG, hypertension, hyperlipidemia, liver cirrhosis and diabetes mellitus presented for surveillance colonoscopy. Three months prior to this he had a piece meal EMR of a 22 mm flat Tubulo-villous Adenoma in the proximal ascending colon. After this, a colonoscopy was passed under direct vision. EFTR device was passed and when it reached the proximal ascending colon, the previously placed tattoo was seen. A postpolypectomy scar was found between the tattoos. This was a non-lifting lesion and was granular lateral spreading and 22 mm in size. Thermal marking was placed and the FTRD clip was now passed. FTRD clip was deployed and using a snare within cap technique a full-thickness resection of 22 mm was performed.

EFTR as a Means for Complete Resection for Adenomas
Pradeep Yarra, MD, Praneeth Kudaravalli, MD, Moamen Gabr, MD, MSc
University of Kentucky, Lexington, KY; University of Kentucky College of Medicine, Lexington, KY

Meta-analysis including 559 patients with colorectal lesions resected by EFTR shows R0-resection rate of 81 % and overall combined rate of any complication of 9.2 %

Chandra S. Dasari, MD et al, Kansas City VA Medical Center, Kansas City, KS presented findings from their FTRD meta-analysis. An electronic database search was conducted in PubMed/MEDLINE, Embase, Google Scholar and Cochrane databases to identify studies that used EFTR for colorectal lesions using the Full-Thickness-Resection-Device (FTRD). The primary outcome was the rate of technical success and R0 resection. Secondary outcomes included complications (bleeding, perforation and post polypectomy syndrome) and the total procedure time. Pooled rates were reported with 95 % Confidence Interval (CI) with heterogeneity (I²). We used a random effects model to calculate pooled rates and used the I-square statistic to quantify heterogeneity. A total of 7 studies were included in the final analysis - 559 patients with 563 lesions removed using FTRD; mean age 67 years, males 59.3 % and mean lesion size of 15.7 mm. Mean follow up of the studies was 4.0 months. The colorectal lesions resected by EFTR included: 32 % cancers, 0.01 % neuroendocrine tumors and 67 % adenomas. The pooled overall technical success was 91 % (95 % CI 88 % - 93 %, I² = 0 %) with a R0 resection rate of 81 % (95 % CI 74 % - 86 %, I² = 62 %). The pooled rates of immediate bleeding, perforation, and post polypectomy syndrome were 5 %, 3 % and 2 % respectively. The overall combined rate of any complication (bleeding, perforation and post polypectomy syndrome) was 9.2 % (6.1 % - 13.5 %). The mean total procedure time was 47.76 min (95 % CI 40 - 55, I² = 84 %). EFTR using FTRD System appears to be an effective and safe technique for the resection of non-lifting colorectal lesions, with a > 90 % technical success rate, R0 resection rate of > 80 % and an overall complication rate of < 10 %.

Efficacy and Safety of Endoscopic Full-Thickness Resection (EFTR) of Colorectal Neoplasms Using the Full-Thickness-Resection-Device (FTRD): A Systematic Review and Pooled Analysis

Chandra S. Dasari, MD¹, Venkat Notalapati, MD¹, Abhiram Duvvuri, MD², Viveksandeep Thoguluva Chandrasekar, MD³, Afeerah Malik³, Divyanshu Kohli, MD³, Prateek Sharma, MD, FACP⁴

¹Kansas City VA Medical Center, Kansas City, KS;

²University of Kansas Medical Center, Kansas City, KS;

³Kansas City VA Medical Center, Overland Park, KS;

⁴Kansas City VA Medical Center, Kansas City, MO

BougieCap

The BougieCap allows for direct visual control and optical feedback during bougienage procedures

Patient with prolonged history of esophageal stricture and complicated treatment experienced significant improvement of symptoms and no complications with BougieCap treatment

George Triadafilopoulos, MD, FACP et al, Stanford University School of Medicine, Mountain View, CA presented a case successful esophageal stricture dilation using the BougieCap. A 62-year-old male presented with dysphagia. He had a prior history of peptic stricture at the gastroesophageal junction, treated with multiple balloon dilations that was complicated by an esophageal perforation. Patient underwent an emergent thoracotomy, esophagectomy, and gastric pull-up with placement of a feeding jejunostomy. Four months after surgery, he noted recurrence of dysphagia, with food getting lodged in the area below his throat and upper chest. Barium swallow revealed a tight stricture at the anastomosis that narrowed the lumen down to 5 mm in size. Upper endoscopy confirmed a stricture located 18 cm from the incisors. Using an endoscopic BougieCap the stricture was dilated serially from 7 to 12 mm in size with good results. Expected mucosal disruption was identified at the end of the procedure. On follow up, patient reported significant improvement in his symptoms without any complications. Use of Savary-Gilliard Dilators for stricture dilation provides only haptic control and often requires the use of fluoroscopic guidance. Our case highlights the use of a novel device - the BougieCap - for treatment of esophageal strictures which allows for direct visual control and optical feedback during the procedure.

Ability to visualize the tension on the tissue while passing the stricture with the cap avoids overstretching and may help reduce the risk of complications, allowing for a more effective procedure. Moreover, by eliminating the need to switch instruments, the use of BougieCap can reduce procedure times.

The Bougie Cap: A Novel Device for Treatment of Esophageal Strictures

Hannah Ramrakhiani, George Triadafilopoulos, MD, FACP
Palo Alto High School, Palo Alto, CA; Stanford University School of Medicine, Mountain View, CO

For questions and further information:

Ovesco Endoscopy USA, Inc.

120 Quade Drive

Cary, NC 27513

U.S.A.

customerservice@ovesco-usa.com

January 2020 | EFTR is a fast, safe and effective option for rectal NETs < 20 mm

In 40 patients with rectal neuroendocrine tumors (NETs), EFTR resulted in 100 % macroscopically and histologically complete resection. Median procedure time was 18.5 minutes. No major adverse events occurred.

B. Meier et al., Department of Gastroenterology, Ludwigsburg Hospital, Ludwigsburg, Germany, conducted the first larger study evaluating endoscopic full-thickness resection (EFTR) of rectal neuroendocrine tumors (NETs) with the FTRD® device. Rectal NETs are subepithelial tumors which bear potential for malignancy in about a quarter of cases. In consequence, depending on tumor size and risk factors, endoscopic or surgical resection is necessary. Hitherto, the optimal endoscopic approach is not defined.

For the study, the German FTRD registry, an online database comprising data on FTRD employment in 31 German endoscopy centers, was searched for cases of rectal NETs. 40 cases of rectal NETs were identified and analyzed retrospectively. The median lesion size was 8 mm (range 3 – 25 mm); small lesions (< 10 mm) accounted for 67.5 % (27/40), medium size (10 – 20 mm) for 30 % (12/40), and large lesions (> 20 mm) for 2.5 % (1/40). 15 % (6/40) were recurrent NETs and had been pretreated with EMR (4/6) or multiple forceps biopsies (2/6).

Resection of the lesion using the FTRD device was successful in all cases. The median procedure time was 18.5 minutes (range 7 – 60 minutes). Resection was macroscopically and histologically complete in all cases. Full-thickness resection was achieved in 38/40 cases (95.0 %). Histological examination of the resected specimen revealed in 70 % (28/40) low grade NETs without lymphovascular infiltration (G1, L0, V0) and without other risk factors. In 30 % (12/40) histology showed granulation tissue or scarring.

No major adverse events occurred. Procedure-related minor adverse events were observed in 12.5 % (5/40). In four patients, minor periprocedural bleeding was seen and managed endoscopically by coagulation and/or injection. In one case, an FTRD snared rupture occurred during resection. The procedure was completed using a conventional resection snare. Endoscopic follow-up was available in 80 % (32/40) and performed over median 12 weeks (range 1 – 49 weeks). In no case, evidence of residual or recurrent tumor was found.

The authors concluded that EFTR with the FTRD is safe and effective for resection of smaller rectal NETs. It is associated with short procedural time, low complication rate and simplicity of the technique. EFTR with FTRD seems to be preferable to EMR/ESD for small and medium size rectal NETs, however, prospective comparative trials are needed to clearly define the role of EFTR in the treatment of rectal NETs.

Full-thickness resection of neuroendocrine tumors in the rectum

JMeier B, Albrecht H, Wiedbrauck T, Schmidt A, Caca K.
Endoscopy. 2020 Jan; 52(1):68-72.
<https://doi.org/10.1055/a-1008-9077>

January 2020 | Systematic review and meta-analysis of OTSC® therapy for digestive bleeding, covering 851 cases shows consistently high rates of definitive hemostasis, technical and clinical success

In 851 patients from 21 studies, definitive hemostasis was reached in 87.8 %, technical success in 97.8 %, and primary clinical success in 96.6 %.

V.T. Chandrasekar et al., Department of Gastroenterology and Hepatology, University of Kansas School of Medicine, Kansas City, Kansas, United States, conducted a comprehensive electronic database search for articles using OTSC for hemostasis aiming to determine technical and clinical success rates in achieving hemostasis as well as the rate of re-bleeding. All articles describing the use of OTSC for GI bleeding were reviewed. Case reports and smaller case series of fewer than five patients were excluded. The primary outcome was the rate of definitive hemostasis after primary technical and clinical success and without re-bleeding during follow-up.

A total of 21 studies comprising overall 851 patients met the inclusion criteria. 62.2 % of patients were males, median patient age was 69.7 years. 687 patients (80.7 %) were treated for upper GI bleeding and the remaining 164 patients (19.3 %) were treated for lower GI bleeding. OTSC was utilized as first-line treatment in 645 patients (75.8 %) and as second-line treatment in 206 patients (24.2 %). The definitive hemostasis rate was 87.8 % (95 % CI: 83.7 % – 92 %) during a median follow-up of 56 days. The rate of definitive hemostasis was 86.6 % (95 % CI: 81.9 % – 91.3 %) for upper GI bleeding and 89.5 % (95 % CI: 85.2 % – 93.8 %) for lower GI bleeding. The technical success rate was 97.8 % (95 % CI: 96.7 % – 98.9 %) and the primary clinical success rate was 96.6 % (95 % CI: 95.1 % – 98.2 %). Re-bleeding occurred in 10.3 % of patients (95 % CI: 6.5 % – 14.1 %). The failure rate of OTSC as first-line treatment was 9 % (95 % CI: 5.2 % – 12.8 %) and 26 % (95 % CI: 16.1 % – 36.0 %) when used as second line treatment. Only two adverse events requiring intervention were reported in 851 patients. In one patient with bleeding duodenal ulcer perforation occurred during OTSC placement which required surgery and another patient developed duodenal obstruction 1.8 months after OTSC placement; the obstruction could be resolved by three balloon dilations.

The authors concluded that this systematic review evaluating OTSC treatment of gastrointestinal bleeding showed high rates of definitive hemostasis, technical success, and primary clinical success, along with low re-bleeding rates. More randomized-controlled trials were desirable.

Efficacy and safety of Over-The-Scope-Clips for gastrointestinal bleeding: a systematic review and meta-analysis.

Chandrasekar VT, Desai M, Aziz M, Patel HK, Gorrepati VS, Jegadeesan R, Rai T, Sathyanurthy A, Murino A, Hassan C, Repici A, Sharma P..

Endoscopy 2019; 51: 941-949. <https://doi.org/10.1055/a-0994-4361>

January 2020 | The OTSC Proctology® is a safe, effective, minimally invasive, and sphincter-sparing instrument for the treatment of complex cryptoglandular fistulas

In 96 patients with complex anorectal fistulas, OTSC Proctology application lead to an overall healing rate of 65 % and a success rate of 79 % for first-line fistula therapy.

Prosst RL, Proctological Institute Stuttgart, Germany, and Joos AK, Center for Rectal and Colonic Diseases Mannheim, Germany, published follow-up data on 96 consecutive patients who underwent treatment of complex anorectal fistulas with the OTSC Proctology.

A total of 100 OTSC Proctology procedures were performed in 96 patients (median age 50 years, range 20-80 years, 64 male, 32 female) with 55 transsphincteric, 38 suprasphincteric, 2 extrasphincteric and 5 rectovaginal fistulas. All but 11 fistulas (8 Crohn's disease, 3 ulcerative

colitis) were of cryptoglandular origin. In 73 cases, the clip was used as first-line therapy, whereas in 27 cases, OTSC Proctology was applied in recurrent fistulas. Minimally invasive surgery with the OTSC Proctology was performed according to a standardized technique: the fistula tract was debrided using a special fistula brush, and the clip was applied on the internal fistula opening. The median operation time was 32 min (range 17 – 66 min). There were no major intraoperative technical problems. The transanal clip procedure was fast and easy to perform in experienced hands. In 93 % of procedures, the clip reliably closed the internal fistula opening already after the first application. Clip malplacement with inadequate fistula closure requiring clip removal and application of a new clip occurred in only 7 patients. All patients found the postoperative pain to be well controlled with standard pain medication. Patients were discharged home on the third and fourth postoperative days. 6-months follow-up data was available for 99 clip applications. On examination 6 months after surgery, 65 % of patients had no clinical signs or symptoms of their previous fistula and were considered to be healed. The healing rate for first-line fistula therapy was 79 %, whereas in recurrent fistulas, the success rate was 26 %. The healing rates of transsphincteric, suprasphincteric, extrasphincteric, and rectovaginal fistulas were 61, 74, 100, and 20 %, respectively. OTSC Proctology was successful in 45 % of fistulas associated with inflammatory bowel disease.

The authors discussed that with an overall healing rate of 65 % and a success rate of 79 % for first-line fistula treatment, the clip achieves results at least as good as established surgical techniques like the advancement flap procedure and the fistulectomy with primary sphincter reconstruction. Moreover, this is done without traumatizing the sphincter muscle and without the risk of fecal incontinence, even if surgery fails. When compared to another sphincter-preserving minimally-invasive technique, the fistula plug, which renders a maximal healing rate of 50-60 %, OTSC Proctology obviously achieves better results. The authors concluded that the OTSC Proctology provides convincing results as first-line treatment for complex cryptoglandular fistulas. It is a safe, effective, minimally invasive, and sphincter-sparing procedure with postoperative pain comparable to other types of fistula surgery.

Short-term outcomes of a novel endoscopic clipping device for closure of the internal opening in 100 anorectal fistulas.

Prosst RL, Joos AK.

Tech Coloproctol 2016; 20:753-758.

December 2019 | Conference Report: United European Gastroenterology Week (UEGW) 2019

- OTSC® shows superiority over conventional therapy in high-risk ulcer bleeding and variceal bleeding
- First clinical cases with stentfix OTSC® yield promising results, proving the tools as valuable, safe and effective
- EFTR is a feasible and safe treatment of T1 colorectal cancer, which delivers optimal histology for risk assessment and leads to a high R0 resection rate

The 27th United European Gastroenterology Week (UEGW) was held on October 19-23, 2019, in Barcelona, Spain. Several workshops, talks and posters presented original research with Ovesco technology and procedures.

OTSC System

OTSC as first-line therapy for high-risk GI ulcer bleeding is associated with shorter procedure time and less rebleeding when compared to case-match controls with conventional therapy

R. Oleas et al., Instituto Ecuatoriano de Enfermedades Digestivas (IECED), Guayaquil, Ecuador, presented a case-match control study assessing the OTSC as first-line therapy in comparison to a combined therapy with conventional hemoclips and adrenaline injection in the management of high-risk bleeding peptic ulcers. The following bleeding ulcers were considered as high-risk ulcers: those located in a major arterial territory, those having an endoscopically visible large-caliber artery > 2mm,

and those with a fibrotic ulcer base and high-risk endoscopic stigmata (Forrest classification types I and II). 95 consecutive patients (mean age 60.9 ± 19.1, 32.6 % female) presenting with high-risk ulcer GI-bleeding between 05/2014 and 09/2018 were included, 46 received an OTSC as primary therapy and 49 matched cases received the combined conventional therapy. Most lesions were gastric ulcers (71.6 %). 6 cases of rebleeding occurred: 2 in the OTSC group and 4 in the combined therapy group (p=0.444). The median procedure time was 11 min (10-15) for OTSC and 20 min (15-40) for combined therapy (p<0.001).

The authors concluded that the OTSC is safe and effective for first line single therapy of high-risk bleeding peptic ulcers. It should be considered the treatment of choice in patients with high-risk bleeding peptic ulcers.

Over-the-scope clip as first-line therapy in the management of high-risk bleeding peptic ulcers: a case-match control study

Robles-Medrande C, Alcivar-Vasquez J, Oleas R, Baquerizo-Burgos J, Olmos JI, Rubio-Cordova M, Puga-Tejada M, Pitanga-Lukashok H, Guayaquil, Ecuador.

OTSC is a safe and effective alternative to sclerotherapy in the management of actively bleeding fibrotic esophageal varices

R. Oleas et al., Instituto Ecuatoriano de Enfermedades Digestivas (IECED), Guayaquil, Ecuador, reported on the use of the OTSC in the management of actively bleeding fibrotic esophageal varices. Repetitive endoscopic band ligations are associated with fibrosis of the esophageal wall. Fibrosis of the esophageal wall, however, impairs band ligation in case of rebleeding. Out of 95 patients presenting with actively bleeding esophageal varices between 09/2016 and 01/2019, 5 patients could not be treated with band ligation due to fibrosis of the esophageal wall. Thus, the OTSC was deployed over the bleeding fibrotic varix as rescue therapy. The cirrhosis etiology was alcohol in one patient and NASH in 4 patients. 2 patients had received 2 previous band ligations, 2 patients had received 3 previous band ligations and 1 patient had received a single previous band ligation. 4 patients were staged as Child-Pugh B and 1 patient as Child-Pugh C. During follow-up endoscopy, 2 patients were submitted to further band ligation for variceal eradication. Neither re-bleeding nor mortality did occur. The authors concluded that the OTSC is a safe and effective alternative in the management of fibrotic esophageal varices. Randomized controlled trials are required to validate these data.

Over-the-scope clip as a rescue therapy for fibrotic bleeding esophageal varices: A single-center experience

Oleas R, Alcivar-Vasquez J, Alvarado-Escobar H, Puga-Tejada M, Robles-Medrande C, Guayaquil, Ecuador.

Endoscopic clipping is a therapeutical option for complications after esophago-gastric surgery

A. Meining, University Hospital Wuerzburg, Germany, held a talk on therapeutic options for complications after esophago-gastric surgery. When diagnosis of complications is made early and perforation/leakage is small, endoscopic clipping can be a valuable treatment option. The ESGE Guidelines recommend the following endoscopic treatment options for early diagnosed gastric perforations, as long as the patient shows no signs of sepsis: Endoclips for small defects, OTSC for large defects if available, when OTSC is not available, endoloop in combination with clips is recommended. If endoscopic treatment fails or the patient develops signs of sepsis/peritonitis, surgical management is recommended. A. Meining presented a clinical case of a patient with leak after oesophagectomy, closure of the leak was attempted with stents, hemoclips, glueing, and OTSC. Mucosal incision prior to placing another OTSC was performed. This led to a successful treatment as it reduced the tension on the tissue and therefore closed the leak.

Crash Course: Complications after esophago-gastric surgery: What to do? Endoscopic clipping is still an option

Meining A, Wuerzburg, Germany.

OTSC for prophylactic closure of the mucosal defect after ESD for SNADET-treatment is associated with

higher rates of complete closure, lower rates of delayed adverse events and shorter procedure time when compared to conventional clips and laparoscopic hand-suturing

O. Dohi et al., Kyoto Prefectural University of Medicine, Kyoto, Japan, presented a study comprising consecutive patients with superficial non-ampullary duodenal epithelial tumor (SNADET), who were treated by Endoscopic submucosal dissection (ESD) and prophylactic defect closure using conventional clips, laparoscopic hand-saw suturing, or OTSC. Inclusion criteria were SNADETs larger than 10 mm without submucosal invasion with clinically low risk for lymph node metastasis. ESD was performed using a needle-type knife (Flush Knife) or a scissor-type knife (Clutch Cutter). The primary end point was to evaluate the en-bloc resection rate and intraoperative complication rate (perforation and bleeding) for comparison of Flush Knife and Clutch Cutter. The secondary endpoint was to evaluate the rate of delayed complications for comparison of the three prophylactic closure techniques.

A total of 77 lesions in 75 patients (male/female 54/23; median age 67 years) were resected by ESD. 37 ESD procedures were performed using the Flush Knife and 40 procedures using the Clutch Cutter. There were 14, 13, and 50 cases using conventional clip, laparoscopic hand-sewn suturing and OTSC for prophylactic closure of the mucosal defect after ESD, respectively. R0-resection rate was 83.8 % using the Flush Knife and 97.5 % using the Clutch Cutter, the difference did not reach statistical significance. The rates of complete closure were 78.6 % with conventional clip, 92.3 % with hand-sewn laparoscopic suture and 98.0 % with OTSC ($p=0.13$). The procedure time with OTSC was significantly shorter than that with the conventional clip and the laparoscopic suturing ($p<0.01$). Adverse events such as delayed bleeding and delayed perforation occurred in 15.4 % and 7.7 % with conventional clip, in 2.3 % and 0 % with laparoscopic suturing, and in 0 % and 2.3 % with OTSC, respectively ($p=0.56$ and $p=1$). The authors concluded that ESD using Clutch Cutter and prophylactic defect closure with OTSC seems to be the best alternative for minimally invasive treatment of SNADET.

Clinical impact of endoscopic submucosal dissection using clutch cutter with over-the-scope clip closure for superficial non-ampullary duodenal epithelial tumor

Dohi O, Yoshida N, Ishida T, Takayama S, Nakano T, Majima A, Inoue K, Kamada K, Takagi T, Konishi H, Naito Y, Itoh Y, Kyoto, Japan.

100 % technical and clinical success rate with OTSC for EUS related iatrogenic duodenal perforation

H. El Bacha et al., Ibn Sina Hospital, University Mohamed V, Rabat, Morocco, presented a retrospective study of patients with EUS related iatrogenic perforation. The study included all consecutive patients with EUS related perforation observed in the gastroenterology and hepatology unit between 2011 and 08/2018, who were diagnosed immediately and received conservative endoscopic management. Endoscopic management consisted of immediate OTSC-clipping. Patients with primary surgical management and conservative non-interventional management, as well as perforation resulting from an endoscopic intervention were excluded. 13 perforations in 8504 EUS procedures had occurred (0.15 %). One was a large duodenal tear requiring immediate surgery, one other was misdiagnosed and early discharged, but readmitted 24 h later with peritonitis which lead to emergency surgery. These two patients were excluded from the study. A total of 11 patients were included (all women, mean age 75, range 68-88). Perforations were located in the superior flexure of the duodenum in 9/11 cases (81 %), in the descending part of the duodenum in 1/11 cases (9 %), and in the inferior duodenal flexure in 1/11 cases (9 %). Defect size ranged from 10-15 mm. OTSC clipping was technically and clinically successful in all cases. 3/11 patients (27 %) had a stay in the intensive care unit for less than 72 h, total hospital stay ranged from 3- 22 days.

The authors concluded that duodenal perforation is a potentially serious adverse event of diagnostic EUS, but conservative endoscopic treatment with OTSC represents a feasible, efficient and safe treatment that can prevent surgery in most instances.

Over the scope clips for EUS duodenal perforation

El Bacha H, Rabat, Morocco, Prat H, Paris, France.

stentfix OTSC

The stentfix OTSC System is a viable approach to avoid migration of a fully covered self-expandable metallic stent in the treatment of trachea-esophageal fistula

M. Conio et al., General Hospital, Sanremo (IM), Italy, presented a clinical case, where the stentfix OTSC System was used. In the treatment of tracheoesophageal fistulas, the placement of fully covered self-expandable metallic stents (FC SEMS) is currently one of the most commonly used strategies. The main limitation of FC SEMS is migration, occurring in at least one-third of patients and requiring further endoscopic interventions. For the anchorage of esophageal stents, external or internal fixation by clips is employed. The stentfix OTSC is a new designed version of the OTSC with a rounded clip preloaded on an applicator cap especially designed for stent fixation. It is intended for use in the gastrointestinal tract for the fixation of metallic stents.

A 53-year old man with locally advanced pharyngeal squamous cell carcinoma was treated with pharyngolaryngectomy and adjuvant radiotherapy and received endoscopic percutaneous gastrostomy (PEG) and tracheostomy. One year thereafter, he developed complex hypopharyngeal stricture with severe dysphagia. Multiple mechanical and pneumatic dilatation sessions led to the occurrence of a trachea-esophageal fistula. For fistula treatment, a 16 mm FC SEMS was placed. Two weeks later, dysphagia occurred due to stent migration. The migrated stent was removed and another stent was placed and anchored with the stentfix OTSC System to treat the fistula. For treating the hypopharyngeal stenosis a similar second stent was placed. No adverse events occurred. Three months later, the patient can swallow a semiliquid diet.

The authors concluded, that the stentfix OTSC System is a viable approach to avoid stent migration. Its efficacy needs to be fully explored.

A new endoscopic device to prevent fully covered metal stent migration

Conio M, Savarese MF, Crespi M, De Ceglie A, Sanremo, Italy.

OTSC Proctology

The OTSC Proctology is a valuable tool within the great variety of approaches for the treatment of cryptoglandular fistulas

J. Warusavitarne, St. Mark's Hospital Harrow, London, UK, held a talk about the treatment of cryptoglandular fistulas. There is a whole toolbox of endoscopic and surgical approaches and even stem cell therapy available for the treatment of cryptoglandular fistulas. The different techniques are used to either lay open the track, disconnect the track from the gut, excise or obliterate the track, fill the track, or correct the immunopathology. Important factors to consider when treating a fistula are the questions, whether the fistula is intersphincteric or transsphincteric, how many tracts the fistula does have, the size of the internal opening, the degree of muscle damage from repeated sepsis, if there is faeces or air coming from the opening and the amount of discharge. The approach has to be tailored to the patient after careful consideration of all influencing details. The OTSC Proctology was presented as valuable tool within the great variety of approaches for the treatment of cryptoglandular fistulas with high success rates reported in clinical case series.

Cryptoglandular fistula: What is new?

Warusavitarne J

FTRD® System

Dutch nationwide prospective cohort study shows 88.8 % R0 resection rate and 2.4 % major adverse events in EFTR of T1 colorectal cancer

L. Zwager, Amsterdam UMC, University of Amsterdam, Amsterdam, the Netherlands, presented a prospective multicenter cohort study evaluating technical and clinical success rates and safety of EFTR for T1 colorectal cancers (T1 CRC). Consecutive patients of 21 Dutch hospitals were included between 09/2015 and 04/2019. Inclusion comprised all scheduled T1 related procedures, both, primary treatment for lesions with optical diagnosis of T1

CRC and secondary treatment after previous (potentially) incomplete resection of T1 CRC. Technical success was defined as en-bloc resection with no macroscopic evidence of residual lesion judged by the endoscopist. Clinical success was defined as R0 resection with tumor-free lateral and deep resection margins and possibility of discrimination between high-risk vs. low-risk T1 CRC. A lesion was defined as high-risk if one of the following risk-factors was present: poor differentiation, lymphatic or vascular invasion, deep submucosal invasion ($\geq 1000 \mu\text{m}$) or incomplete resection (R1/Rx resection). Besides, adverse events were evaluated. 247 procedures were included. Indications for EFTR were primary resection for suspected T1 in 81 cases and re-resection after previous incomplete resection of T1 CRC in 166 cases. Technical success was achieved in 211/247 cases (85.4 %), histopathology could not be obtained in 15/247 cases (6.1 %) because the lesion could either not be reached or not be retracted into the cap. In the remaining 232 cases amenable to EFTR, the R0 resection rate was 88.8 % (206/232). Final histopathology confirmed residual adenocarcinoma in 33.2 % ($n=77/232$). Discrimination between high-risk and low-risk carcinoma was achieved in 97.4 % (75/77). Low-risk T1 CRC was identified in 22.1 % ($n=17/77$) and high-risk T1 CRC in 75.3 % ($n=58/77$). Additional surgery was performed in 41.4 % (24/58) of the high-risk cases, of which 21/24 (87.5 %) had no residual cancer or lymph node metastasis. Endoscopic surveillance strategy was initiated in 46.6 % (27/58). The overall adverse event rate was 8.5 % ($n=21/247$), with emergency surgery in 2.4 % (6/247) for 2 immediate and 4 delayed perforations. In conclusion, the authors rated EFTR a feasible and safe treatment for T1 CRC, both as primary treatment and secondary treatment after previous incomplete resection. EFTR delivers optimal histology for risk assessment and leads to a high R0 resection rate, avoiding surgery in most cases.

Endoscopic full-thickness resection is feasible for T1 colorectal cancers – a Dutch nationwide prospective cohort study

Zwager L, Amsterdam, Bastiaansen BAJ, Amsterdam, van der Spek BW, Alkmaar, Heine GDH, Alkmaar, Bronzwaer MES, Amsterdam, Haasnoot KJC, Alkmaar, van der Sluis H, Zwolle, Perk L, the Hague, Boonstra JJ, Leiden, Rietdijk ST, Amsterdam, Schwartz MP, Amersfoort, Wolters HJ, Groningen, Weusten BLAM, Nieuwegein, Gijzen LPL, Eindhoven, ten Hove WR, Leiden, Nagengast WB, Groningen, Bekkering FC, Capelle aan den IJssel, Terhaar sive Droste JS, Hertogenbosch, Fockens P, Amsterdam, Dekker E, Amsterdam, the Netherlands.

EFTR of adenomas involving the appendiceal orifice is associated with an intermediate risk for developing appendicitis and consecutive risk for appendectomy

B. Walter et al., University Hospital Ulm, Ulm, Germany, presented the results of a retrospective analysis of 38 patients from 2 centers with adenomas involving the appendiceal orifice treated with the colonic FTRD. Objective of the study was to evaluate post treatment complications in acute and long term follow up (appendicitis, mucocele). All patients had received prophylactic antibiotic treatment for an average of 3.7 days in a row started preinterventionally. No acute severe events were reported for the EFTR procedure. Within follow-up, symptoms of appendicitis occurred in 9 patients (23.7 %). In 3 cases, conservative treatment was sufficient, 5 patients were transferred to appendectomy. No development of a mucocele was reported during long term follow-up (median follow-up time 11 months, range 6-32 months).

The authors concluded that EFTR of adenomas involving the appendiceal orifice is associated with an intermediate risk for developing appendicitis and consecutive risk for appendectomy. Patients must be thoroughly informed about the risks.

Intermediate risk of appendicitis following full-thickness resection of adenomas arising from the appendiceal orifice – a retrospective analysis

Walter B, Ulm, Germany, Wannhoff A, Ludwigsburg, Germany, Schmidbaur S, Ulm, Germany, Meier B, Ludwigsburg, Germany, Meining A, Würzburg, Germany, Caca K, Ludwigsburg, Germany.

Spanish multicenter study proves EFTR as safe and

feasible for various colorectal lesions

Uchima H et al., University Hospital Doctor Josep Trueta de Girona, Girona, Spain, from the Endoscopic Resection Working Group of the Spanish Society of Digestive Endoscopy, presented a multicenter study evaluating efficacy and safety of the FTRD System in colorectal lesions. Consecutive patients of 10 Spanish endoscopy centers treated with the colonic FTRD were assessed. 71 EFTR were scheduled. Indications for EFTR were: recurrent lesions with non-lifting sign (46.47 %), not pretreated lesions with non-lifting sign (23.94 %), residual lesions with non-lifting sign (11.26 %), appendicular lesions (2.8 %), suspected T1 lesions (7 %), suspicious scar (4.2 %) and subepithelial lesion (4.2 %). In 3 patients passage of the sigmoid was not possible with the FTRD cap. In the other 68 patients, technical success rate was 85.2 %, the en-bloc resection rate was 83.8 %. Final histology revealed LGD adenoma (40 %), HGD adenoma (23 %), intramucosal adenocarcinoma (4.47 %), SSP (5.87 %), T1sm1 (2.9 %) and others (2.8 %). In one case, the OTSC clip was not deployed resulting in an intraprocedural perforation. There were 2 cases of delayed perforation and 1 case of delayed bleeding. 10 patients underwent surgery: 3 for perforation, 1 for intraappendicular lesion and 6 for advanced adenocarcinoma. During follow-up, 3 residual/recurrent adenomas were treated endoscopically, all of which showed benign histology.

The authors concluded that EFTR is a safe and feasible technique for colorectal lesions. Evaluation of the insertion with a long cap (e.g. "proVE" cap) and traction of the lesion prior to EFTR is highly recommended.

Safety and feasibility of endoscopic full-thickness resection in colorectum using over the scope clip. A multicenter Spanish experience

Uchima H, Girona, Barcelona and Badalona, Barquero D, Barcelona, Esteban Lopez-Jamar JM, Madrid, Espinos JC, Terrassa, Marin-Gabriel J, Madrid, Roson P, Malaga, Fernandez Cadenas F, Oviedo, Palacio Galan MA, Oviedo, Puig I, Manresa, Rodriguez Sanchez J, Ciudad Real, Fraile Lopez M, Oviedo, Ortega Alonso A, Malaga, Arribas-Anta J, Madrid, Mel L, Madrid, Sabado F, Madrid, Garcia Lledo, Madrid, Fernandez-Simon A, Barcelona, Mata A, Barcelona, Albeniz E, Pamplona, Spain.

Cases of intended EFTR, where part of the lesion was trapped inside the OTSC, are evaluated

H. Uchima et al., University Hospital Doctor Josep Trueta de Girona, Girona, Spain, reported on 10 EFTR cases, where the lesion was trapped in the OTSC clip and could not be resected properly. Locations of EFTR were appendix (n=1), stump (n=1), right colon (n=1), transverse colon (n=2), left colon (n=2), sigma (n=2), and rectosigmoid junction (n=1). Mean diameter of the lesions was 19 mm. Indications for EFTR were: recurrent lesion with non-lifting sign (6 cases), native lesion with non-lifting sign (n=1), residual lesion with non-lifting sign (n=2), and appendicular lesion (n=1). In 8 cases, partial resection of the lesion was performed, in 2 cases only biopsies were taken. Final histology revealed LGD in 4 cases, HGD in 2 cases, intramucosal adenocarcinoma in 1 case, SSP in 2 cases and advanced adenocarcinoma > sm2 in 2 cases. 3 patients underwent surgery (appendicular lesion and advanced adenocarcinomas), 3 residual lesions were treated endoscopically and in 4 cases the scar showed no residual tissue.

The authors concluded that in some cases of intended EFTR, where part of the lesion was trapped inside the OTSC, residual lesion may be easily treated endoscopically or must not be treated at all when histology shows scar without residual dysplasia.

EFTR with OTSC in colorectum: what happens when the lesion is trapped in the over-the-scope-clip and is not resected

Uchima H, Girona and Barcelona, Barquero D, Barcelona, Esteban Lopez-Jamar JM, Madrid, Espinos J, Terrassa, Marin-Gabriel J, Madrid, Roson P, Malaga, Fernandez Cadenas F, Oviedo, Palacio Galan MA, Oviedo, Puig I, Manresa, Rodriguez Sanchez J, Ciudad Real, Fraile Lopez M, Oviedo, Ortega Alonso A, Malaga, Arribas-Anta J, Madrid, Mel L, Madrid, Sabado F, Madrid, Garcia Lledo J, Madrid, Fernandez-Simon A, Barcelona, Mata A,

Barcelona, Eduardo A, Pamplona, Spain.

For questions and further information:

Ovesco Endoscopy AG
Scientific Information Service
Friedrich-Miescher-Straße 9
D-72076 Tübingen
science@ovesco.com

Ovesco Research Update 33**December 2019 | New randomized-controlled trial: OTSC® in first line has superior ability to reduce rebleeding in NVUGIB when compared to standard treatment****A randomized controlled trial presented at the American College of Gastroenterology (ACG) meeting shows the rebleeding rate after first-line OTSC treatment of NVUGIB to be 85 % lower than after standard treatment.**

At this year's ACG Annual Scientific Meeting, taking place on October 25-30 in San Antonio, TX, USA, Dr. D. Jensen UCLA School of Medicine, Los Angeles, CA, presented a RCT evaluating OTSC as first-line endoscopic treatment of severe non-variceal upper gastrointestinal bleeding (NVUGIB). The RCT was conducted at two academic medical centers. Patients with peptic ulcer bleeding or Dieulafoy's lesions, who met emergency endoscopy criteria were included in the study and randomized in a 1:1 allocation to either standard treatment (conventional clips and/or multipolar electrical probe with epinephrine pre-injection) or OTSC. All patients received high dose PPIs after randomization and were followed for 30 days.

Overall 49 patients could be included in the study, 24 patients were allocated in the OTSC arm and 25 patients in the standard treatment arm. The two groups were well matched regarding demographics, risk factors, lesion type and stigmata of hemorrhage. During the 30 day follow up, rebleeding occurred in one patient in the OTSC group (1/24, 4.2 %) and 7 patients in the standard treatment group (7/25, 28 %; p=0.024). The rebleeding rate was 85 % lower in the OTSC group than in the Standard treatment group (relative risk 0.149, 95 % confidence intervals of 0.006 and 0.863). The number needed to treat (NNT) was 4.2. Severe complications did not occur in the OTSC group (0/24, 0 %), but did occur in 4 patients in the standard treatment group (4/25, 16 %, p=0.041). Angiographic embolization was not necessary in the OTSC group, but was necessary in 2 patients (8 %) of the standard treatment group (p=0.157). The mean number of red blood cell transfusions (\pm standard deviation) after randomization was 0.4 ± 0.2 in the OTSC group and 0.76 ± 1.64 in the standard treatment group, a Kruskal-Wallis test performed for statistical comparison.

In summary, patients with severe NVUGIB, who received primary endoscopic treatment with the OTSC, experienced significantly lower rates of rebleeding, significantly less severe complications, and needed fewer red blood cell transfusions than patients treated with standard endoscopic hemostasis. The authors concluded that the OTSC has a superior ability to obliterate arterial blood flow underneath stigmata of hemorrhage and thereby reduce lesion rebleeding.

Randomized Controlled Trial (RCT) of Over-the-Scope Clip (OTSC) as Initial Endoscopic Treatment of Severe Non-Variceal Upper Gastrointestinal Bleeding (NVUGIB).

Jensen DM, Kovacs TOG, Ghassemi KA, Kaneshiro M, Gornbein J. Los Angeles, CA, USA.

ACG 2019, Annual Scientific meeting & Postgraduate Course, Oct 25-30, 2019, San Antonio, TX, USA.

December 2019 | Conference Report of DGVS and DGAV 2019

The 74th annual conference of the German Society of Gastroenterology, and Digestive and Metabolic Diseases with Endoscopy Section (DGVS) took place together with the 13th autumn conference of the German Society of General Surgery and Visceral Surgery (DGAV) on October 2-5, 2019 in Wiesbaden, Germany.

Ovesco products were presented in talks, posters, meetings, video forums, workshops and hands-on training

sessions.

FTRD® System

Endoscopic full-thickness resection with the FTRD System is highly cost-effective when compared to surgical and endoscopic alternatives.

Additional costs per R0 resection in relation to EFTR are +27,000 € for the surgical alternative and +4,800 € for EMR/ESD

A. Kuellmer, University Hospital Freiburg, presented a cost-effectiveness study of the endoscopic full-thickness resection (EFTR) of difficult-to-remove lesions with the FTRD System in comparison to available alternative therapies. The data for analysis was drawn from the patients of the WALL-RESECT study and a hypothetically designed control group. The control group was generated by a systematic survey of the treating physicians of all participating centers of the study (University Hospital Freiburg, Hospital Ludwigsburg, Evangelical Hospital Duesseldorf, Elisabeth-Hospital Essen, University Hospital Ulm, University Hospital Wuerzburg, Hospital Augsburg, Hospital Dortmund, Hospital Krefeld, and Hospital Boeblingen-Sindelfingen) and comprised various alternative therapies (ESD, EMR, surgical resection). Costs and revenues were determined from the DRG-Report 2019. The defined effectiveness parameter was the primary endpoint of the WALL-RESECT study, the R0 resection. For the control group, the R0 resection rates published in current literature were used. ICER (Incremental Cost-Effectiveness Ratio) and ACER (Average Cost-Effectiveness Ratio) served as cost effectiveness criteria. Data was presented from the payers' point of view and from the care providers' point of view.

In the control group, 106 patients underwent surgical resection (59 %), 45 patients underwent EMR (25 %) and 29 patients underwent ESD (16 %). Case costs for the different procedures were 2,758 € for EFTR, 8,895 € for surgical resection, 1,712 € for endoscopic alternatives (EMR/ESD), and 5,828 € for the whole group of all alternative therapies ("casemix-alternative"). From the care providers' point of view, the ACER for EFTR is 3,582 €, for surgical resection 8,924 €, for EMR/ESD 2,985 € and for the casemix-alternative 7,195 €. The ICER in relation to EFTR, meaning the additional costs per R0 resection, are +27,066 € for the surgical alternative, +4,788 € for the endoscopic alternative and +76,748 € for the casemix alternative.

The authors concluded that endoscopic full-thickness resection with the FTRD System is highly cost-effective when compared to surgical and also endoscopic alternative therapies.

Endoscopic full-thickness resection in the colorectum with the FTRD System – a cost-effectiveness analysis. (Die endoskopische Vollwandresektion im Kolorektum mit dem FTRD System – eine Kosteneffektivitätsanalyse.)

Kuellmer A, Behn J, Beyna T, Schumacher B, Meining A, Messmann H, Neuhaus H, Albers D, Birk M, Probst A, Faehndrich M, Frieling T, Goetz M, Caca K, Thimme R, Schmidt A. Freiburg, Ludwigsburg, Duesseldorf, Essen, Ulm, Wuerzburg, Augsburg, Dortmund, Krefeld, Boeblingen-Sindelfingen.

Costs for EFTR are significantly lower by 5,094 € when compared to surgical segment resection

M. Rathmayer, Inspiring-health GmbH, Munich, presented first results of a study, which compares the costs for open-surgical and laparoscopic procedures with EFTR. In addition, an analysis of cost-covering by the G-DRG System was performed. Performance numbers were queried at the German Federal Statistical Office (DeStatis). Healthcare costs will be represented by the respective DRG revenues of 2019. For calculation of the actual costs, cost data of hospitals shall be used, which participate in the cost calculation of the Institute for Remuneration System in Hospitals (Institut fuer Entgeltsystem im Krankenhaus - InEK). The proof of cost coverage shall be provided by the comparison of the actual costs of the hospitals with the standard deviation of the respective DRG as published by InEK.

2016, EFTR was performed in 166 hospitals 763 times (OPS 2016: 5-452.25; from 2017: 5-452.65); in 2017, the

FTRD System was used in 234 hospitals 1,096 times. In comparison, in 2016 and 2017, 5,922 and 8,184 open surgical and laparoscopic segment resections of the colon with anastomosis were performed (OPS 2019: 5-455.01 and 5-455.02). EFTR is represented in the G-DRG System by DRG G47B. The average costs are 3,058 € and are remuneration with 3,069 €. Open-Surgical and laparoscopic segment resection has an average cost of 8,138 € and a remuneration of 8,164 € (DRG G18C).

In summary, the number of deployments of the FTRD System increased by 44 % from 2016 to 2017 to 1,096 cases in 2017. The costs for the health care system are significantly lower by 5,094 € when compared to surgical segment resection. A cost-cost comparison of cases in hospitals, which participate in the InEK-calculation, will show, whether this advantage can also be proved on case level. In addition, the cost coverage in the G-DRG-System 2019 in these hospitals will be evaluated and compared with the cost data from 2018.

Economic analysis of the endoscopic full-thickness resection (Ökonomische Analyse der endoskopischen Vollwandresektion).

Rathmayer M, Schmidt A, Schepp W, Heinlein W, Albert JG, Gundling F. Munich, Freiburg, Stuttgart.

Pooled analysis from 5 hospitals shows 77.9 % R0 resection rate and 1.3 % rate of major complications for EFTR of various lesions in the whole GI-tract

B. Lewerenz, Hospital Bogenhausen, Munich, presented a study comprising data of 232 patients from 5 different hospitals in Southern Germany (focal care and maximum care hospitals), in whom an EFTR was performed between 11/2014 and 08/2019. Technical success rate, resection status, and procedural complications were retrospectively analysed.

EFTR was performed in 232 patients on the basis of the following indications: 69.4 % adenomas with negative lifting sign, thereof 22.8 % relapse adenomas, 28.0 % residual adenomas after uncomplete resection and 18.5 % native adenomas. In 21.1 % of cases, a proven colonic or gastric adenocarcinoma existed or was suspected to a high degree, in 0.9 % of cases there was a DALM (Dysplasia associated lesion or mass) and in 8.2 % a submucosal tumor or NET (neuro-endocrine tumor). In one case, biopsy was performed on suspicion of M. Hirschsprung. Localisation of EFTR was in 25.4 % in the ascending colon, in 25.0 % in the rectum, 14.7 % in the sigmoidal colon, in 11.6 % in the transverse colon, in 8.2 % in the descending colon, in 8.2 % in the cecum, in 3.9 % at the appendiceal orifice, in 1.7 % at a colonic surgical anastomosis, in 0.9 % in the stomach and in 0.4 % in the duodenum. The median size of the resected lesions was 16.22 mm (range 3 – 50 mm). EFTR could not be performed in 8 cases (3.4 %) due to impossibility of endoscope advancement with FTRD cap and in 29 cases (12.5 %) resection was technically not successful (n=1 marking of the lesion not possible, n=7 tissue cannot be fixed, n=8 suction in FTRD cap not sufficient because of severe scarring, n=1 clip malfunction, n=3 tissue avulsion by grasper, n=9 macroscopically diagnosed residual adenoma after EFTR). Thus, the technical success-rate was 84.8 %, the R0 resection rate was 77.9 %. Post-interventionally, in 3 cases (1.3 %) major complications occurred (1 perforation, 2 acute appendicitis).

The authors concluded that EFTR as minimally invasive method is highly successful in the resection of various lesions in the whole gastrointestinal tract with an overall very low procedural complication rate.

Technical success, resection status, and procedural complication rate of the endoscopic full thickness resection with the “full-thickness Resection Device” (FTRD) – results of a pooled analysis from various hospitals. (Technischer Erfolg, Resektionsstatus sowie prozedurale Komplikationsrate der endoskopischen Vollwandresektion mittels „Full-Thickness Resection Device“ (FTRD) – Ergebnisse einer gepoolten Analyse aus verschiedenen Kliniken).

Krutzenbichler I, Dollhopf M, Lewerenz B, Lilje T, Fuchs M, Herrmann S, Vitali F, Rath T, Kaiser C, Eigler A, Diepolder H, Schepp W, Gundling F. Munich, Kaufbeuren, Erlangen.

The FTRD System is a safe and effective instrument for resection of adenomas with high-grade dysplasia and

of small early-stage adenocarcinomas

P. Stathopoulos, University Hospital Marburg, presented monocentric data supplying information on the deployment of the FTRD in the colon and rectum. Between 11/2016 and 08/2019, 50 patients (32 % female, median age 66 ± 13 years) were treated with the FTRD System. 23 lesion (60.5 %) were located in the colon and 15 lesions (29.5 %) in the rectum. Median size of the lesions was 18.9 mm (range 3-50 mm). Indication for EFTR was in 21 cases a relapse/residual adenoma, in 10 cases an adenoma with negative lifting-sign, in 12 cases an incompletely resected early-stage carcinoma or adenoma with high-grade dysplasia, in 6 cases a rectal NET, in 1 case a submucosal tumor, and in 1 case the suspicion of M. Hirschsprung. All lesions could be reached with FTRD cap. In 2 lesions, a hybrid technique was performed (EMR +EFTR) due to the big size of the lesions (lesion size 40 and 50 mm, respectively). The median procedural time was 50 min (10 – 110 min). Resection was technically successful in 94 % (47/50, in 2 cases insufficient suction into FTRD cap due to severe scarring, in one case cutting of the loop at the side of clip). The R0 resection rate was 76 % (38/50). In 11 patients, histological examination proved the existence of an adenocarcinoma (R0 in 6/11 cases, 55 %). 2 of the adenocarcinomas were diagnosed as sm2 stage, 8 of the adenocarcinomas as sm3 stage and one as T2 stage, 4 patients underwent oncological surgical revision. In 7 cases, complications occurred: 4 bleedings which could be managed conservatively, 2 perforations (1 treated surgically, 1 closed endoscopically with OTSC) and 1 Subileus (treated with balloon dilation and decompression tube). In 34 patients, follow-up data was available, in 2 cases a residual/relapse adenoma was found, both could be curatively resected by de novo EFTR.

The authors concluded that EFTR represents a safe and effective method, which usefully complements the present endoscopic set of instruments. According to the authors, EFTR plays a particular role in the resection of adenomas with high-grade dysplasia and of small early-stage adenocarcinomas.

Endoscopic full-thickness resection (eFTR) of colorectal lesions – a retrospective analysis of 50 consecutive cases. (Endoskopische Vollwandresektion („endoscopic full-thickness resection“, eFTR) für kolorektale Läsionen – eine retrospektive Analyse von 50 konsekutiven Fällen.)

Stathopoulos P, Gallmeier E, Gress TM, Denkert C, Denzer U. Marburg.

Medium risk of appendicitis after EFTR of adenomas near the appendiceal orifice

S. Schmidbaur, University Hospital Ulm, reported on a study examining how high the risk for post-interventional appendicitis is after full-thickness resection at the appendiceal orifice. Data of patients of 3 centers (Hospital Ludwigsburg, Hospitals of the Neumarkt district, and University Hospital Ulm), in which an EFTR of adenomas near the appendiceal orifice had been performed between 2014 and 2018, was retrospectively evaluated. Patients which had undergone appendectomy prior to EFTR were excluded from the study. 50 patients (16 male, median age 65.8 years, range 46-83 years) fulfilled the inclusion criteria. All patients received peri-interventional antibiotics for median 3.7 days. Intra-interventionally, no major complications occurred. During in-hospital and further follow-up (4 months), 7 patients developed signs of appendicitis (14 %). In 4 patients, the signs arose within 10 days after the intervention, in the remaining 3 patients about 1 month after intervention. In 4 cases, conservative therapy was possible, in 3 patients, surgical appendectomy was necessary. In summary, 57 % of patients with appendicitis could be treated conservatively. Development of mucocele was not observed in any patient.

The authors concluded that a medium risk of appendicitis must be presumed for EFTR of adenomas near the appendiceal orifice. Thus, it is necessary to sufficiently inform the patients about the risk and the eventual necessity for surgical appendectomy. Possibly, the risk can be minimized by longer-time post-interventional antibiotics. Further studies on this would be desirable.

Moderately increased risk of appendicitis after eFTR of adenomas near the appendiceal orifice – a

retrospective analysis. (Moderat erhöhtes Appendizitisrisiko nach eFTR von Adenomen nahe des Appendixabgangs – eine retrospektive Analyse.)

Walter B, Schmidbaur S, Wannhoff A, Meier B, Meinig A, Caca K. Ulm, Ludwigsburg, Wuerzburg.

Application-oriented lectures at the conference underline the establishment of the instrument in everyday clinical practice

Various other lectures, held e.g. by B. Schuhmacher, Elisabeth Hospital Essen, K. Caca, Hospital Ludwigsburg, and S. Hollerbach, AKH Celle, explained application-oriented and practically the deployment of the FTRD in everyday clinical practice. By means of clinical examples, photos, videos and illustrations, the correct deployment and avoidable mistakes were illustrated and appropriate indications were explained as well as therapy alternatives. Besides, data from actual studies evaluating the deployment of the FTRD in the upper GI tract was presented.

Full-thickness resection in the colorectum. Why do we need endoscopic FTRD? (Vollwandresektion im Kolorektum. Warum brauchen wir die endoskopische FTRD?)

Hollerbach S, Celle.

Endoscopic Full Thickness Resection (EFTR).

Caca K, Ludwigsburg.

Endoscopic Full-Thickness Resection in the duodenum (Endoskopische Vollwandresektion im Duodenum).

Caca K, Ludwigsburg.

Submucosal tumors in the upper gastrointestinal tract.

Tunnel techniques and full-thickness resection – the solution for many clinical conditions? (Subepitheliale Tumoren im oberen Gastrointestinaltrakt. Tunneltechniken und Vollwandresektion – die Lösung für vieles?)

Schumacher B, Essen.

Tumor resection by tunneling and by full-thickness resection in the upper GI-tract (Tumorsektion durch Tunnel und per Vollwandresektion im oberen GI Trakt).

Caca K, Ludwigsburg.

RESECT+

Deployment of an additional working channel leads to the development of new techniques from EMR and ESD, which show significant advantages regarding resectable lesion size, procedural time and complication rate

ESD+ is in ex-vivo model related to significantly shorter procedural time (24.5 vs. 32.5 min) and causes significantly less muscular damage (1/32 vs. 6/32, p=0.04) than conventional ESD

S. Kunsch, University Medicine Goettingen, presented a study at the ex-vivo model, which compares the new resection technique ESD + with the conventional ESD technique. The study was performed in the ex-vivo model of a porcine stomach in the EASIE-R-Simulator. Lesions of 3 cm and 4 cm were manually placed and afterwards resected by ESD or ESD+, respectively. Overall 64 interventions (32 ESD and 32 ESD+; each (n=16) 3 cm and (n=16) 4 cm lesions, thereof again half (n=8) antegrade and the other half retrograde) were performed. En-bloc resection was reached in all resections. Procedural time was significantly shorter with ESD+ than with conventional ESD (median 24.5 min vs. 32.5 min). Particularly in 3 cm and 4 cm lesions resected in retroflexion, this advantage could be seen. Perforations did not occur, however, in the ESD+ group, there was significantly less muscular damage than in the ESD group (1/32 vs. 6/32; 3.2 % vs. 18.7 %; p=0.04). In conclusion, the authors rated ESD+ a promising new resection method, which shows significant advantages concerning procedural time and complication risk when compared to conventional ESD.

Endoscopic submucosal dissection with an additional working channel (ESD plus) as new method improving the traditional ESD in the es vivo model. (Endoskopische Submukosadisektion mit einem zusätzlichen Arbeitskanal (ESD plus) als neue Methode zur Verbesserung der traditionellen ESD im Ex-vivo-Model.)

Knoop RF, Wedi E, Ellenrieder V, Neesse A, Kunsch S, Goettingen.

Lesions of 3 cm or 4 cm size can be resected en-bloc by EMR+ (86.36 % and 60 % en-bloc resection rate, respectively), this is insufficient up to not possible by conventional EMR (18.8 % and 0 % en-bloc resection rate, respectively)

R.F. Knoop, University Medicine Goettingen, presented a study evaluating the EMR+ technique in the ex-vivo model. EMR+ in contrast to EMR is performed using an additional working channel (AWC), which enables the additional use of a grasping instrument. The study was performed prospectively in the ex-vivo model of a porcine stomach in the EASIE-R-Simulator. Lesions (overall n=152) of 1 cm, 2 cm, 3 cm, and 4 cm in size were manually placed by coagulation. In 1 cm-lesions, EMR and EMR+ achieved the same en-bloc resection rate of 100 %. In 2 cm-lesions, the en-bloc resection rate dropped to 54.55 % with EMR, while EMR+ still achieved 95.44 %. Conventional EMR did not achieve sufficient results in 3 cm and 4 cm lesions (18.8 %; 0 %). However, EMR+ still achieved an en-bloc resection rate of 86.36 % in 3 cm lesions and dropped to 60 % in 4 cm lesions.

In conclusion, the authors rated the EMR+ technique a promising new resection technique with the potential to improve the safe en-bloc resection rate of lesions ≥ 2 cm with reasonable technical and financial effort.

EMR+: Presentation of a new endoscopic resection method in the ex-vivo model. (EMR+: Vorstellung einer neuen endoskopischen Resektionsmethode im Ex-vivo Model.)

Knoop RF, Wedi E, Ellenrieder V, Neesse A, Kunsch S, Goettingen.

OTSC® System

The clinical application of the OTSC System is trained in workshops and hands-on training sessions

Hands-on training of hemostasis techniques in small groups attracted lively interest

On Wednesday, hands-on trainings sessions of endoscopic-interventional techniques took place. Training on two different topics was offered: first hemostasis and clips, and second stent placement. In the beginning, a short lecture gave an overview over the respective topic, then hands-on training of the respective procedures was performed in small groups. Within this part, also possible complications of the hemostasis procedures were presented and acceptable therapeutic solutions worked out. The mainly practical instruction in small groups on the basis of examples from everyday endoscopic practice found widespread appeal among the very interested endoscopists.

DEGEA offers instruction on leakage closure and variceal bleeding at the training model

The German Society for Endoscopy Professionals (Deutsche Gesellschaft für Endoskopiefachberufe – DEGEA) offered several workshops on Saturday, including one workshop with the topic: "What is when... difficult situations and prevention of complications". In the theoretical part, S. Loeffler, Ovesco Endoscopy AG, Tuebingen, spoke about perforation closure, foreign body removal and stent fixation. In the practical part, various techniques were trained under instruction by different tutors at the training model; besides others, the OTSC System and endoclips were employed for leakage closure and variceal bleeding. The program was well attended by experienced colleagues and beginners.

Endoscopic closure techniques are highly effective for the management of acute perforations

A. Schmidt, University Hospital Freiburg, reported on perforation management in endoscopy. First, he defined the term perforation and differentiated it from leakage/insufficiency and fistula development on the basis of clinical characteristics. Then he explained that a paradigm change has taken place within the last years in the management of free abdominal air due to the introduction of highly effective endoscopic closure techniques. The presence of extraluminal air does not automatically imply the need for surgical therapy any more. Based on the recommendations of the ESGE Position Statement 2014, he explained general measures in case of

acute perforation and criteria for decision making between endoscopic and surgical treatment. He emphasized, how important it is for therapeutic success, to keep the time between discovery of the perforation and start of therapy as short as possible. After that, he presented common endoscopic closure techniques. These included through-the-scope clips, OTSC clips, stents (SEPS, SEMS) and sponge/drainage for endoscopic vacuum therapy. Endoscopic suturing techniques were not part of the lecture. Depending on defect size, localization and local expertise, the different endoscopic closure techniques are employed for perforation management. Herein, inter-disciplinary peri-interventional patient care is of great importance. In cases of late diagnosis or not certain/ not possible closure, or clinical worsening, surgical therapy is still the treatment of choice.

Perforations in endoscopy – do I still need my surgeon? (Perforationen bei der Endoskopie – brauche ich meinen Chirurgen noch?)

Schmidt A, Freiburg.

For questions and further information:

Ovesco Endoscopy AG
Scientific Information Service
Friedrich-Miescher-Straße 9
D-72076 Tübingen
science@ovesco.com

November 2019 | EFTR is a useful tool for staging of early colorectal cancer and simultaneously can be sufficient treatment for low-risk lesions

In 156 patients with histologic evidence of adenocarcinoma, EFTR was technically successful in 92.3 % and R0 resection was achieved in 71.8 %. Exact risk stratification was possible in 99.3 %. Severe procedure-related adverse events were recorded in 3.9 % of patients.

A. Kuellmer et al., University Medical Center Freiburg, Germany, conducted a multicenter, retrospective study aiming to evaluate efficacy, safety and clinical value of EFTR for early colorectal cancer. Current international guidelines recommend endoscopic resection for T1 colorectal cancer with low-risk histology features and oncologic resection for those at high risk of lymphatic metastasis. Exact risk stratification is therefore crucial to avoid under-treatment as well as over-treatment.

Records of 1234 patients undergoing EFTR for various indications at 96 centers were screened for eligibility. A total of 156 patients with histologic evidence of adenocarcinoma were identified. This cohort included 64 cases undergoing EFTR after incomplete resection of a malignant polyp (group 1) and 92 non-lifting lesions (group 2).

Technical success was achieved in 144/156 patients (92.3 %). Mean procedural time was 42 minutes. R0 resection was achieved in 112/156 patients (71.8 %). A subgroup-analysis was performed assessing the R0 resection rate for subgroup 1 (previous incompletely resected malignant polyps, n=64) and subgroup 2 (non-lifting lesions, n=92). In group 1, a R0 resection rate of 87.5 % and in group 2 a R0 resection rate of 60.9 % had been reached (p<0.001).

Severe procedure-related adverse events (all perforations) were recorded in 6 cases (3.9 % of patients). One perforation was diagnosed immediately during the intervention and closed successfully with an over-the-scope clip. Five secondary perforations were observed after 1-8 days and required surgical intervention. Discrimination between high-risk versus low-risk tumor was successful in 155/156 cases (99.3 %). In group 1, 84.1 % were identified as low-risk lesions, whereas 16.3 % in group 2 had low-risk features. In total, 53 patients (34 %) underwent oncologic resection due to high-risk features. 98 patients (62 %) were followed endoscopically.

The authors concluded that in colorectal cancer, EFTR is technically feasible, effective and safe. It allows exact histological risk stratification and can avoid surgery for low-risk lesions. Prospective studies are required to further define indications for EFTR in malignant colorectal lesions and to evaluate long-term outcome.

Endoscopic full-thickness resection for early colorectal cancer.

Kuellmer A, Mueller J, Caca K, Aepli P, Albers D, Schumacher B, Glitsch A, Schaefer C, Wallstabe I, Hofmann C, Erhardt A, Meier B, Bettinger D, Thimme R, Schmidt A.

Gastrointest Endosc 2019; 89(6):1180-89.

November 2019 | Resection techniques EFTR, EMR and ESD chosen as one of the top 3 endoscopic advances of 2018

The Gastrointestinal Endoscopy Editorial Board chooses resection techniques namely EMR/submucosal dissection and full-thickness resection as one of the top 3 endoscopic advances of 2018. Votes from each individual of the 9-member board were tallied to identify a consensus list of 10 areas of endoscopic research they considered a top priority based on the criteria of significance, novelty, impact on national health and impact on global health. With 7 votes, gastrointestinal resection techniques including EFTR were listed as number 3 of this list.

To identify major new advances in gastrointestinal endoscopy in 2018, the American Society for Gastrointestinal Endoscopy's (ASGE) Gastrointestinal Endoscopy (GIE) Editorial Board reviewed original research articles pertaining to GI endoscopy published during 2018 in Gastrointestinal Endoscopy and 10 other leading medical and gastroenterology journals. These journals were selected on the basis of impact factor and included medical journals (Annals of Internal Medicine, Journal of the American Medical Association, Lancet, New England Journal of Medicine), gastroenterology journals (American Journal of Gastroenterology, Clinical Gastroenterology and Hepatology, Gastroenterology, Gut, Hepatology), and specific endoscopic journals (endoscopy, GIE). Reviewing the titles, each member of the 9-member editorial board independently identified 10 areas of endoscopic research they considered a top priority based on the criteria of significance, novelty, impact on national health, and impact on global health. The list from each member was compiled, and the votes were then tallied to identify the consensus "Top 10 topics in GI endoscopy in 2018". The list is arranged in order of priority based on the number of votes for each topic: adenoma detection (9 votes), bariatric endoscopy (8 votes), EMR/submucosal dissection/full-thickness resection (7 votes), artificial intelligence (7 votes), expandable metal stents for palliation of biliary obstruction (6 votes), pancreatic therapy with lumen-apposing metal stents (6 votes), endoscopic reprocessing (6 votes), Barrett's esophagus (6 votes), interventional EUS (5 votes), and GI bleeding (5 votes).

Acknowledging there may still be some bias and subjectivity, the board trusts that the list of the top 10 advances will guide educators by defining new endoscopic techniques that need to be propagated in clinical practice, and focus investigators on priority areas for research.

Gastrointestinal Endoscopy Editorial Board top 10 topics: advances in GI endoscopy in 2018.

Cohen J, Desilets DJ, Hwang JH, Baig KRKK, Leung FW, Maranki JI, Okolo PI, Swanstrom LI, Chak A.

Gastrointestinal Endoscopy 2019; 90(1):35-43; DOI: <https://doi.org/10.1016/j.gie.2019.03.020>

October 2019 | EMR+: preclinical evaluation of the novel technique for endoscopic en-bloc resection proves it to be effective, safe, and easy and fast to perform

The novel EMR+ technique uses an external additional working channel (AWC®) and a poloxamer agent for submucosal injection (LiftUp®). A preclinical study evaluated the technique in the EASIE-R1 Simulator (explanted pig stomach, 22 lesions) and subsequently in 13 in vivo resections in 3 domestic pigs. In vivo, resection was macroscopically complete in 92.3 %. Median procedure time was 5 min and median size of the resected specimens was 35 mm x 35 mm x 11 mm. Adverse events did not occur except from one minor

periprocedural bleeding managed by coagulation.

B Meier et al., Ludwigsburg Hospital, Ludwigsburg, Germany, presented a preclinical study evaluating a modified EMR technique (EMR+) using an additional working channel (AWC, Ovesco Endoscopy) and a new agent for submucosal injection (LiftUp, Ovesco Endoscopy). The EMR+ technique is a modified grasp and snare technique and allows for en-bloc resection of larger polyps, which could be related with lower recurrence rates when compared to piecemeal endoscopic mucosal resection (EMR) and with lower complication rates when compared to endoscopic submucosal dissection (ESD).

After the technique was established, 22 ex vivo resections were performed (EASIE-R1 Simulator, explanted pig stomach). Primary aim of this setting was to develop the concept of EMR+ (e.g. positioning of instruments, injection technique, resection technique). The procedure time was measured from the beginning of submucosal injection to the extraction of the resected specimen and lasted in median 7 min (range 5-11 min). Median size of the resected specimens was 30 mm x 26 mm x 11 mm.

Afterwards, the technique was transferred and evaluated in 13 EMR+ resections in in vivo porcine models (stomachs of 3 domestic pigs). In vivo, median procedure time was 5 min (range 3-12) and median size of resected specimens was 35 mm x 35 mm x 11 mm. Major adverse events did not occur. During one resection (7.7 %), minor periprocedural bleeding occurred which could be managed by coagulation. The authors concluded that EMR+ showed to be effective and safe in the porcine model and allowed for en-bloc resection of lesions up to 40 mm. Further evaluation of the novel technique in clinical comparative trials is needed.

Novel technique for endoscopic en-bloc resection (EMR+) – Evaluation in a porcine model.

Meier B, Wannhoff A, Klinger C, Caca K.
World J Gastroenterol 2019; 25(28): 3764-3774

September 2019 | OTSC® Proctology: Pilot RCT shows efficacy and safety in achieving permanent closure of cryptoglandular anal fistulas

When compared to fistulectomy and primary sphincter reconstruction, the OTSC Proctology showed comparable results (93 % primary healing rate) but a significantly reduced number of hospitalization days and medications after treatment.

Mascagni D et al., Department of Surgical Sciences, Sapienza University, Rome, Italy, published a randomized controlled trial comparing the OTSC Proctology with standard fistulectomy and primary sphincter reconstruction for the treatment of anorectal low trans-sphincter fistula. Between 02/2012 and 03/2013, patients affected by trans-sphincter anal fistula were consecutively enrolled in the trial. Patients were randomized to receive fistulectomy with primary sphincter reconstruction of OTSC Proctology. Demographic characteristics, comorbidities, previous anorectal treatments, and recurrent fistula data were acquired. Post-operative therapy data and Wexner scores (30 and 60 dd) were acquired during follow-up.

30 consecutive patients were included in the study. 15 patients underwent the OTSC Proctology procedure and 15 patients underwent standard fistulectomy. The primary healing rate was 93.3 % in the OTSC group, in one patient a relapse fistula developed after 30 days which had to be re-treated surgically. No relapse fistula occurred in the fistulectomy group. The mean postoperative stay was 1.3 days in the OTSC group and 3.6 days in the fistulectomy group. The mean medications required for complete healing was 3.2 (range 2-5) in the OTSC group and 8.9 (range 7-12) in the fistulectomy group. The Wexner Score at 30 days was 0 (range 0-0) in the OTSC group and 3 (range 0-6) in the fistulectomy group. After 60 days, the Wexner score was 0 in both groups. After 6 months, all patients underwent MRI examination, which showed complete healing in all patients. In the 3 years of follow-up, none of the patients showed recurrence of the fistula or incontinence.

The authors concluded that the OTSC Proctology is an effective and safe tool for achieving permanent closure of the internal opening of anal fistulas, with excellent results in terms of pain, postoperative incontinence, healing time, and

days of hospitalization.

OTSC® Proctology vs. fistulectomy and primary sphincter reconstruction as a treatment for low trans-sphincter anal fistula in a randomized controlled pilot trial.

Mascagni D, Pironi D, Grimaldi G, Romani AM, LA Torre G, Eberspacher C, Palma R, Sorrenti S, Pontone S.
Minerva Chir. 2018 Feb 1. doi:10.23736/S0026-4733.18.07617-4. [Epub ahead of print]

Ovesco Research Update 32**September 2019 | The BougieCap provides visual control for effective and safe dilation of esophageal strictures and prevents complications from over-dilation**

In 50 patients with complex benign esophageal stenoses, endoscopic bougienage with the BougieCap was successful in 96 %. Symptoms of dysphagia decreased significantly after bougienage (59.0 points at Day 0 vs. 28.6 points at Day 14; p<0.001). No severe adverse events occurred.

B Walter et al., University Clinic Ulm, Germany, presented a multicenter study (Ulm, Southampton and Essen) evaluating dilation of benign esophageal stenoses with the BougieCap. The BougieCap allows, in contrast to Savary bougies, direct visual control of the process without the need for x-ray.

50 patients (25 f, 25 m, median age 67.1 ± 16.8 years) with benign stenosis of the esophagus and clinically apparent symptoms of dysphagia were included. Genesis of the stenosis was peptic (n=23), radiation (n=13), anastomotic (n=6), caustic ingestion (n=4), Post-ESD (n=2), EoE (n=1) and unknown (n=1). Dilation was successful in 96 % of cases (48/50). In two cases with narrow and longer strictures (one peptic, one post-radiochemotherapy) located in the lower esophagus (length 40 mm, diameter 2 mm and 3 mm, respectively), bougienage failed due to high resistance, which caused buckling of the endoscope in the pharynx. Bougienage was aborted and balloon dilation was performed successfully. The mean number of sequential bougienages was 2.3 ± 0.7. A stabilizing guidewire was used in 10 cases, 8 of which were with a pediatric scope. Severe complications did not occur. In two cases, a BougieCap was lost in the stomach; no clinical discomfort of complications resulted. The BougieCaps evacuated spontaneously with the stool. Mean dysphagia-associated symptoms, using the DHI (Dysphagia Handicap Index), decreased from 59.0 (moderate dysphagia) to 28.6 points (mild dysphagia) over the short-term follow-up period of 14 days (p<0.001).

The authors concluded that endoscopic treatment of benign esophageal stenoses with the BougieCap allows direct visual control of the dilation process and of beginning mucosal lacerations. Thus, in contrast to the conventional blind method, over-dilation and re-traumatization are reduced and the dilation process can be performed with better adaptation to the stenosis. Usage of a guidewire is reasonable and necessary in special cases (i.e. very high-grade stenosis, usage of a pediatric gastroscope).

The BougieCap – a new method for endoscopic treatment of complex benign esophageal stenosis: results from a multicenter study.

Walter B, Schmidbaur S, Rahman I, Albers D, Schumacher B, Meining A.
Endoscopy 2019; DOI: <https://doi.org/10.1055/a-0959-1535>.

August 2019 | gastroduodenal FTRD®: Full-thickness resection of gastric SETs enables definitive histological diagnosis, better than conventional biopsy

Gastric EFTR in 29 patients with gastric SET led to en-bloc resection in 89.7 %, 76 % R0 resection and final histological diagnosis in all resected specimens. The only periprocedural complication of minor bleeding was managed endoscopically in all cases.

Dr B. Meier et al., Department of Gastroenterology, Ludwigsburg Hospital, Ludwigsburg, Germany, initiated

the RESET trial in March 2017, a prospective observational multicenter pilot trial evaluating efficacy and safety of clip-assisted endoscopic full-thickness resection (EFTR) for gastric subepithelial tumors (SET). Gastric SET up to 15 mm were included in the trial. Primary endpoint was complete en-bloc resection. Secondary endpoints were R0 resection, full-thickness resection, adverse events and recurrence rate at 3-months follow-up.

29 patients underwent gastric EFTR. Histology after conventional biopsy prior to EFTR was only able to define histological tumor type in 31.2 %. Complete en-bloc resection with the gFTRD System was achieved in 89.7 %. Histology of the specimen defined the histological tumor type in 100 %. 76 % of all SET was resected histologically complete (R0), requiring no further diagnostic intervention; full-thickness-resection specimens were obtained in 65.5 %. In the remaining cases resection reached deep into the wall. In 31 % minor bleeding was observed and managed endoscopically in the same session. Follow-up was available in 79.3 % and showed OTSC clip detachment in 78.3 % and the OTSC clip in position in 21.7 %. No signs of residual or recurrent tumors were observed after 3 months.

The authors concluded that EFTR of gastric SET with the gFTRD is safe and effective. It provides a definitive diagnosis (including sufficient risk stratification in case of GIST or NET) in contrast to conventional biopsy. R0 resection is achieved in most cases and might obviate the need for further surveillance endoscopies.

Endoscopic full-thickness resection of gastric subepithelial tumors with the gFTRD-system: a prospective pilot study (RESET trial).

Meier B, Schmidt A, Glaser N, Meining A, Walter B, Wannhoff A, Riecken B, Caca K.
Surgical Endoscopy <https://doi.org/10.1007/s00464-019-06839-2>

August 2019 | Conference Report of the Digestive Disease Week 2019

The 50th Digestive Disease Week (DDW) 2019 took place on May 18 – 21, 2019 in San Diego, CA, USA.

Ovesco products were presented in talks, posters, state-of-the-art-lectures and debates.

National and international multicenter studies confirm that the FTRD® allows for fast, safe, and effective resection of difficult colorectal and upper GI lesions

FTRD System

U.S. multicenter FTRD study shows R0 resection rate of 79 % in patients with difficult colonic lesions
Y. Ichkhanian et al., Baltimore, Maryland, USA, performed a retrospective study at 26 U.S. tertiary-care centers assessing feasibility, efficacy and safety of the FTRD for the resection of colonic lesions. Patients who underwent EFTR using the FTRD for lower GI lesions between 10/17 and 01/19 were included. Outcomes were technical success (defined as en-bloc resection) and R0 histologic margin. Chi square test was used to assess the association between lesion type, size, and location with the two outcomes.

A total of 95 patients (mean age 65.5 years, 38.9 % female) underwent resection of colonic lesions using the FTRD. Inclusion criteria were the following lesions of the lower gastrointestinal tract: difficult adenomas (defined as non-lifting, recurrent, residual or involving appendiceal orifice/ diverticular opening) (n=63), low-risk adenocarcinomas (n=21), and sub-epithelial tumors (n=11). Lesion location was in the proximal colon in 58 (61 %) patients, followed by distal colon in 17 (18 %), and rectum in 20 (21 %). 65 patients (70 %) received propofol sedation, 21 patients (23 %) received general anesthesia. Prophylactic antibiotics were used in 32 patients (33.8 %). Mean pre-resection lesion diameter and total procedure time were 15.5 ± 6.4 mm (range 3-30 mm) and 59.7 ± 31.8 min, respectively. 71 patients (76.3 %) were discharged post-procedurally. Technical success was achieved in 80 (84.2 %) patients. R0 resection was achieved in 75 patients (79 %). A total of 5 (5.3 %) adverse events occurred, 3 AEs were mild (3.6 %) (1 minor bleeding and 2

iatrogenic strictures), and 2 AEs (2.1 %) were severe (1 appendicitis, 1 perforation).

The authors concluded that EFTR is a feasible, safe and effective technique for EFTR of difficult colonic lesions. Surgical management can be avoided in the vast majority of cases.

Non-exposure full-thickness resection of colonic lesions in the U.S.: The FTRD experience.

Ichkhanian Y¹, Vosoughi K¹, Sharaiha RZ², Hajifathalian K², Tokar JL³, Templeton AW⁴, James TW⁵, Grimm IS⁵, Mizrahi M⁶, Samarasekera JB⁷, Chahade NEH⁷, Lee J⁷, Chang KJ⁷, Barawi M⁸, Irani SS⁹, Friedland S¹⁰, Korc P¹¹, Aziz AAdam A¹², Al-Haddad MA¹³, Kowalski TE¹⁴, Novikov AA¹⁴, Diehl DL¹⁵, Smallfield G¹⁶, Ginsberg GG¹⁷, Oza V¹⁸, Pannu D¹⁸, Fukami N¹⁹, Pohl H²⁰, Lajin M²¹, Kumta NA²², Tang SJ²³, Amateau SK²⁴, Ngamruengphong S¹, Kumbhari V¹, Brewer Gutierrez OI¹, Khashab MA¹.

¹Baltimore, ²New York, ³Philadelphia, ⁴Seattle, ⁵Chapel Hill, ⁶South Alabama, ⁷Orange, ⁸Detroit, ⁹Indianapolis, ¹⁰Stanford, ¹¹Newport Beach, ¹²Chicago, ¹³Indianapolis, ¹⁴Philadelphia, ¹⁵Danville, ¹⁶Richmond, ¹⁷Philadelphia, ¹⁸Florence, ¹⁹Scottsdale, ²⁰New York, ²¹Michigan, ²²New York, ²³Jackson, ²⁴Minneapolis

EFTR of colonic lesions with evidence of adenocarcinoma allows exact histological risk stratification and avoidance of surgery in low-risk lesions

A. Küllmer et al., University of Freiburg, Freiburg, Germany, presented a retrospective study at 96 German centers which evaluated EFTR with the FTRD for risk stratification of early colorectal cancer. 156 patients with evidence of adenocarcinoma were included. 64 patients underwent EFTR after incomplete resection of a malignant polyp (group 1) and 92 patients had non-lifting lesions. Technical success was achieved in 144/156 (92.3 %), mean procedural time was 42 minutes. R0 resection was achieved in 112/156 (71.8 %). Subgroup analysis showed a R0 resection rate of 87.5 % in Group 1 and 60.9 % in Group 2 ($p < 0.001$). Severe procedure-related adverse events were recorded in 3.9 %. Discrimination between high- vs. low-risk tumor was successful in 155/156 cases (99.3 %). 84.1 % of Group 1 were identified as low-risk lesions, whereas 16.3 % in group 2 had low-risk features. In total 53 patients (34 %) underwent oncologic resection due to high risk features whereas 98 patients (62 %) were followed endoscopically.

The authors concluded that EFTR in early colorectal cancer is technically feasible and safe. It allows exact histological risk stratification to avoid surgery for low-risk lesions.

Endoscopic full-thickness resection in colorectal cancer: experience.

Küllmer A¹, Mueller J¹, Caca K², Aepi P³, Dakkak D⁴, Schumacher B⁴, Glitsch A⁵, Schäfer C⁶, Wallstabe I⁷, Hofmann C⁸, Erhardt A⁹, Meier B², Bettinger D¹, Thimme R¹, Schmidt A¹

¹Freiburg, ²Ludwigsburg, ³Luzern, ⁴Essen, ⁵Greifswald, ⁶Neumarkt, ⁷Leipzig, ⁸Mainz, ⁹Wuppertal

EFTR with the FTRD is a fast, safe and effective option for rectal neuroendocrine tumors

B. Meier et al., Ludwigsburg hospital, Ludwigsburg, Germany, presented a study evaluating EFTR resection of rectal neuroendocrine tumors (NET). Rectal NET are subepithelial tumors with potential for malignancy. Prevalence is rare but increasing over the last decades. Between 09/2015 and 05/2017, data of 40 cases of rectal NET were collected. Lesions were located in the lower (13/40), middle (24/40) and upper (3/40) rectum and had a median size of 8 mm (SD 4.43, range 3-25 mm). In 15 % (6/40) rectal NET were reported as recurrent NET.

All lesions could be reached and resected with FTRD. The median procedure time was 18.5 minutes (range 7-60 min). All lesions could be resected macroscopically and histologically complete. Histology after EFTR showed low-grade NET (G1) without lymphovascular infiltration (L0, V0) and without other risk factors in 70 % (28/40). In 30 % (12/40) histology revealed granulation tissue or scarring. Procedure-related adverse events were observed in 12.5 % (5/40). In four cases (10 %) minor periprocedural

bleeding was observed and managed endoscopically (coagulation and/or injection). In one case (2.5 %) a technical problem was observed (rupture of the FTRD snare). Endoscopic follow-up was available in 80 % (32/40) and conducted after a median time of twelve weeks (range 1-49 weeks) after resection. In 72 % (23/32) the OTSC had spontaneously detached and in 28 % (9/32) the clip was still in situ. No evidence of a residual or recurrent lesion could be found.

The authors concluded that EFTR with the FTRD is a fast, safe and effective option for rectal NET < 20 mm. IT should be considered as first-line therapy for rectal NET < 20 mm without risk factors.

Full thickness resection of neuroendocrine tumors in the rectum.

Meier B¹, Albrecht H², Wiedbrauck TH³, Schmidt A⁴, Caca K¹

¹Ludwigsburg, ²Neumarkt, ³Duisburg, ⁴Freiburg

86 % R0 resection rate achieved in large colorectal lesions resected with FTRD alone or hybrid technique

A. Vareedayah et al., NYU Langone Health, New York, New York, USA, performed a retrospective study on consecutive patients treated with the FTRD device at 9 North American centers. A total of 64 patients was included. The primary indication for EFTR was non-lifting adenoma (37 patients, 58 %). Other indications included suspected high-grade dysplasia/cancer, lesions in the appendiceal orifice and subepithelial lesions. The mean size of the resected lesions was 24 mm. The FTRD could not be advanced to the lesion in 4 patients. 56/60 lesions (93 %) were successfully removed by FTRD alone (when < 2 cm) or a combination of FTRD and EMR (when > 2 cm). R0 (including clear vertical margin for hybrid procedures) was achieved in 48/56 patients (86 %). Complications occurred in 3 patients. Perforation in one patient was treated endoscopically. Rectal bleeding in one patient did not require intervention. One patient required surgery for appendicitis 72 h after the procedure.

The authors concluded that EFTR is a safe and effective method for resection of large colorectal lesions.

Initial north American experience with endoscopic full-thickness resection of colorectal lesions: a multicenter retrospective cohort study.

Vareedayah A¹, Yuen PYS¹, Skinner M¹, Koller K¹, Alkaade S², Diehl DL³, Al-Haddad MA⁴, Templeton AW⁵, Hwang JH⁶, Stavropoulos SN¹, Cohen J¹, Mendoza Ladd AH⁷, Grajales-Figueroa G⁸, Mahadev S⁹, Haber GB¹

¹New York, ²St. Louis, ³Danville, ⁴Indianapolis, ⁵Seattle, ⁶Redwood City, ⁷El Paso, ⁸Mexico City, ⁹New York

First international experience of FTRD application in the upper GI tract shows high technical and clinical success rates

K. Hajifathalian et al., Cornell, NYC, New York, United States, presented an international multicenter retrospective study including patients who had an endoscopic resection of an upper GI tract lesion using the FTRD System between 08/2017 and 11/2018. 54 patients from 11 centers with endoscopic resection of an upper GI tract lesion using FTRD were included. The most common lesions were mesenchymal neoplasms (n=22), followed by ectopic pancreas or scar tissue (n=12), adenomatous polyps (n=7), hamartomatous polyps (n=6), adenoma with high-grade dysplasia or adenocarcinoma (n=3), carcinoid tumor (n=3), and hyperplastic polyp (n=1). 1 lesion was located in the esophagus, 10 in the cardia/fundus, 15 in the stomach body, 20 in the antrum, and 8 in the duodenum. The average size of lesions was 14 mm (SD 8 mm). Deployment of FTRD was technically successful in 92 % of the patients (n=50) leading to complete and partial resection of the target lesion in 41 (76 %) and 9 (16 %) patients, respectively. Histological margin of resection was reported in 47 patients, and in these FTRD led to R0 resection in 74 % (n=35), R1 resection in 21 % (n=10), and Rx in 4 % (n=2). Thus a R0 resection could be achieved in 35 patients (65 %) out of the 54 patients the procedure was attempted. Out of 4 patients with adenoma with high-grade dysplasia or adenocarcinoma, 3 had an R0 resection. FTRD was complicated by intra-procedural minor bleeding in 5 patients (11 %) and major bleeding in 4 patients (9 %),

all of which were controlled endoscopically. 9 patients (17 %) were taking antithrombotic medication at the time of the procedure. Previous biopsy or attempted resection by hot snare, EMR or ESD had been performed in 32 patients (59 %).

The authors concluded that these results suggest a high technical and clinical success rate with low risk of early recurrence and an acceptable complication rate for FTRD in the upper GI tract.

Full-thickness resection device (FTRD) for treatment of upper gastrointestinal tract lesions: the first international experience.

Hajifathalian K¹, Dawod QM¹, Issa D¹, Meining A³, Schmidt A⁴, Vosoughi K², Ichkhanian Y², Ngamruengphong S², Kumbhari V², Samarasekera JB⁵, El Hage Chehade N⁵, Tang SJ⁶, Kasmin F¹, Templeton AW⁸, Fukami N⁹, Goetz M¹⁰, Sampath K¹, Glaser N⁴, Mahadev SH¹, Mukewar S¹, Call-Locke DL¹, Hwang JH¹¹, Sharaiha RZ¹, Khashab MA²

¹New York City, ²Baltimore, ³Ulm, ⁴Freiburg, ⁵Irvine,

⁶Jackson, ⁷New York City, ⁸Seattle, ⁹Scottsdale,

¹⁰Tuebingen, ¹¹Stanford

Hybrid EMR-EFTR is associated with 76 % negative vertical margins and low complication rate

P.Y.S. Yuen et al. presented a comparative study, comparing consecutive patients treated with hybrid EMR-EFTR due to large colorectal lesions (n=17) with patients who underwent EFTR alone due to colorectal lesions (n=14). Mean lesion size in the EFTR group was 16.8 mm, mean lesion size in the hybrid group was 33.6 mm. Mean procedure time was 67 minutes in the EFTR group and 100 minutes in the hybrid group. Histology confirmed R0 resection in 93 % (13/14 patients) in the EFTR group and negative vertical margins in 76 % (13/17 patients) of lesions in the hybrid group. Two of four in the hybrid group with positive vertical margins were due to technical difficulty (snare malfunction and lesion not entrapped by snare) and subsequently removed with EMR. Complications were relatively scarce. In one patient the snare was inadvertently closed, with tissue resection prior to clip deployment resulting in a wall defect. One patient, who did not take the antibiotic medication as prescribed, developed appendicitis, which required surgery 72 hours after the procedure. The authors concluded that Hybrid EMR-EFTR for colorectal lesions is a safe and effective method for resection of lesions that are otherwise too large for EFTR alone. There were no adverse events related to deployment of the clip into tissue with EMR defect. This approach is an alternative to ESD or surgery.

A novel hybrid technique using endoscopic mucosal resection (EMR) and endoscopic full-thickness resection (EFTR) for large colorectal neoplasms unresectable by EMR alone

Yuen PYS¹, Vareedayah AA¹, Skinner M¹, Hoerter NA¹, Koller K¹, Mahadev SH², Haber GB¹

^{1,2}New York City

OTSC® treatment of high-risk peptic ulcer bleeding is proven to be safe, effective and cost-effective when compared to standard treatment

OTSC System

Treatment of recurrent peptic ulcer bleeding with the OTSC is associated with on average \$2160 savings per patient when compared to standard TTS clipping

J. X. Yu et al., University of Michigan, Ann Arbor, Michigan, USA, presented a study evaluating the cost-effectiveness of OTSC clips as compared to through-the-scope clips. A decision tree was used to model the costs, effectiveness (quality-adjusted life years) and rates of persistent/recurrent bleeding were compared in OTSC versus standard therapy for the treatment of peptic ulcer bleeding. The costs were estimated based on 2016 CMS reimbursement rates. Cost-effectiveness of the modalities was determined by the incremental cost-effectiveness ratio. The initial procedure cost was estimated to be the cost of an EGD with hemostasis and the cost of either an OTSC or 2 TTS was estimated using actual costs from a large health care system in the US. The authors assumed that the patients who were successfully treated incurred

the cost of an admission for gastrointestinal bleeding with less than major comorbidity or complication. Patients who did not have clinical success with either the OTSC or standard therapy would incur the cost of an admission for gastrointestinal bleeding with major comorbidity or complication. The primary outcome of interest was the total cost. Sensitivity analyses were performed to ensure the robustness of the results.

The total cost to treat a patient with recurrent bleeding was \$8368.56 using the over the scope clip and \$10,528.55 using TTS. Thus, the use of OTSC clips, on average, resulted in \$2160.00 savings per patient. The findings can be regarded as robust as sensitivity analyses showed that OTSC remains cost effective if the rate of further bleeding after OTSC remains lower than 55 % or remains higher than 17 % with standard therapy using TTS clips.

The authors concluded that a strategy to treat recurrent peptic ulcer bleeding using the OTSC is associated with both a higher efficacy and a lower cost. Gastroenterologists should consider using the over the scope clip rather than standard therapy when the risk of rebleeding after standard therapy is higher than 55 %.

Over the scope clips for recurrent peptic ulcer bleeding is cost effective as compared to the through the scope clips

Yu J X¹, Kaltenbach T R², Keyes E², Soetniko R M²
¹Ann Arbor, ²San Francisco

State-of-the-art lecture on the OTSC for severe upper GI bleeding

J. Hochberger, Vivantes hospital Friedrichshain, Academic Teaching Hospital of Humboldt University Charité, Berlin, Germany, held a state-of-the-art lecture on the OTSC treatment of severe upper GI bleeding. The OTSC is a Nitinol macro clip which provides controlled mechanical tissue compression. It has the working principle of a surgical stapler, therefore, tissue perfusion is preserved after clip application. The application of the clip is easy and similar to a band ligation device. The lesion is targeted, the OTSC cap is brought in connection to the tissue, the target tissue is suctioned into the cap and the OTSC clip is released by turning the hand wheel. In acute bleeding, suction is often enough to pull the target tissue into the application cap. For fibrotic tissue or tangential application, a forceps or OTSC Anchor can be used for transferring the target tissue into the cap. For this maneuver, the OTSC Anchor is positioned and tissue is fixated with the anchor, the OTSC cap is aligned to the lesion by pulling the Anchor and advancing the endoscope. Thereby, the tip of the OTSC Anchor can be mobilized into the cap, the anchor spikes may remain external. In the next step, the clip is released. After clip application the OTSC Anchor is detached from the tissue. The OTSC bench data was collected on hemostasis (GIE 2012; 75: 152-9). OTSC showed a persistent pressure increase after application in comparison to the sloping pressure curve achieved with conventional clips. With OTSCs, a significantly lower number of clips was needed for effective hemostasis and a significantly shorter time to effective hemostasis was needed with OTSC vs 2 TTS clips. In summary, the OTSC is easy to apply, application is fast, one single OTSC is sufficient in most cases. The OTSC provides a strong and reliable mechanical closure with maintained tissue perfusion. Early clip loss is rare with OTSC. Its special clinical strength are chronic peptic ulcers with fibrotic base. Limitations for the OTSC are rare, application is limited when access to the bleeding vessel with the OTSC and clip housing is insufficient. This can be the case when there is a stenosis between endoscope and target area, then prior dilation is necessary. Another reason for insufficient access is a lateral position of the bleeding source, in this case, traction into the housing is necessary, use with a side viewing endoscope is not possible. Another limitation for OTSC application is diffuse tumor bleeding, in such a case, spray, injection and multiple conventional clips must be used.

Different types of over-the-scope clips have been developed. The OTSC t has teeth with small spikes, it provides compression plus anchoring. The OTSC a has round teeth, it provides mainly compression. The OTSC gc has prolonged teeth with spikes for gastric wall closure.

There are three different hood sizes (8.5-11 mm, 10.5-12 mm, 11.5-14 mm) and 2 different cap depths (3 mm and 6 mm) on the market. The OTSC t (traumatic) is used for stomach and chronic duodenal ulcers at the level of the bulb. The OTSC a (atraumatic) is used for the small intestine and colon parts with thin wall.

Recent studies with the OTSC are

- A Meta-Analysis of 20 studies (n = 510 patients) regarding OTSC hemostasis showing a high technical and clinical success rate of 93.0 % and 87.5 % (Weiland T et al., MinInVTh 2019)
- The multicenter FLETRock trial (n=118) of the OTSC as first line therapy revealing a total success rate of 92.5 %. The re-bleeding rate and re-bleeding associated mortality determined in comparison to the prognostic Rockall score was significantly reduced (Wedi E et al., Surg.End. 2018)
- The randomized controlled STING trial (n=66) proving that OTSC clipping is significantly superior to former standard therapy techniques in the treatment of severe recurrent UGIB (Schmidt A et al., Gastro 2018)
- A large multicenter cohort study (n=286) on OTSC first-line treatment revealing superior technical and clinical success rates of 97.9 % and 96.4 % (Manta R et al., End Int open 2018).

In conclusion, the OTSC should be used as primary tool in all high-risk patients (Rockall 7+), high-risk defined e.g. for patients under anticoagulation, for patients with hemoglobin < 10, for patients with cirrhosis and ulcer bleeding, for patients with a spurting ulcer that can easily be faced with an OTSC, and for patients with ulcers with sclerosed ulcer base. Besides, the OTSC should be used in all re-bleeders, if the lesion can be reached with OTSC and clip housing.

Over-the-scope clip for severe upper GI-bleeding – Time for a change in practice?

Hochberger J¹, Wedi E²
¹Berlin, ²Göttingen

The OTSC as first-line single therapy is as safe and effective as combined therapy for the management of high-risk bleeding peptic ulcers

C. Robles-Medrande et al., Guayaquil, Guayas, Ecuador, reported on an analysis of data on consecutive patients who presented with high-risk ulcer GI bleeding between 05/2014 and 09/2018. High-risk upper GI bleeding was considered as those ulcers located in a major arterial territory, if the lesion had an endoscopically visible large-caliber artery (>2 mm), if there was a fibrotic ulcer with high-risk endoscopic stigmata (Forrest classification types I and II).

95 patients were included, 46 received an OTSC as primary therapy for HR bleeding ulcers and 49 matched cases received TTS hemostatic clips in combination with epinephrine injection (combined therapy). The mean age was 60.9 ± 19.1, 32.6 % female. Most lesions were gastric ulcers (71.6 %). The median number of OTSC used was 1 (1-3), whereas for combined therapy was 2 (1-8) TTS clips. Six cases of rebleeding (6.3 %) were observed: two in the OTSC group and four in the combined therapy group (p=0.444). Two cases of the OTSC group (4.3 %) had rebleeding after 48 hours of the procedure; meanwhile, one case of rebleeding was observed in the combined therapy group at the same period and was treated with APC (p=0.520). Three cases in the combined therapy group had rebleeding in less than 48 hours after the procedure (p=0.088), two treated with an OTSC and one with APC. The median procedure time was 11 (10-15) mins for OTSC and 20 (15-40) for combined therapy (p<.001).

The authors concluded that the OTSC as first-line single therapy is as safe and effective as combined therapy for the management of high-risk bleeding peptic ulcers; improving the procedure time.

Over-the-scope clip as first-line therapy in the management of high-risk bleeding peptic ulcers: a case-match control study.

Robles-Medrande C¹, Alcivar-Vasquez JM¹, Baquerizo-Burgos ROJ¹, Ignacio Olmos J¹, Rubio-Cordova M¹, Pitanga Lukashok H¹

¹Guayaquil, Guayas, Ecuador

The OTSC is more effective in obliterating arterial blood flow in severe NVUGIB than standard visually guided endoscopic hemostasis

D.M. Jensen et al., David Geffen School of Medicine at UCLA, Santa Monica, California, USA, held a state-of-the-art lecture on the treatment of patients with severe non-variceal upper gastrointestinal bleeding (NVUGIB). Recurrent NVUGIB after standard visually guided endoscopic hemostasis is common in high-risk patients. A recent randomized controlled trial (Gastro 2017; 152:1310-18) found the 30-day rebleeding rate to be 26.3 % (20/76) with visually applied MPEC and/or standard hemoclips with or without epinephrine. The rebleeding rate was reduced to 11.1 % (8/72) when blood flow was monitored with Doppler endoscopic probe (DEP) and used as a guide for hemostasis. However, when residual arterial blood flow was not obliterated, the rebleeding rate was very high – 88.8 % (8/9 patients). D.M. Jensen et al performed a prospective cohort study with OTSC in 20 patients with severe NVUGIB as primary hemostasis with DEP monitoring before and after hemostasis. 19 patients had bleeding peptic ulcers (12 duodenal, 7 gastric) and 1 Dieulafoys lesion. Results were compared to previously studied patients from the Gastro RCT and to results of a cohort study of DEP in patients with peptic ulcer bleeding before and after visually guided hemostasis (GIE 2016; 83: 129-36).

Residual arterial blood flow detection after OTSC and DEP guided complete hemostasis were similar (5 % - 1/20 vs. 0 % - 0/63 respectively) but were significantly lower than standard visually guided therapy – 24.2 % (23/95) in the cohort study. Low 30 day rebleeding rates were seen after OTSC or successful DEP hemostasis (5 % - 1/20 vs. 0 % - 0/63 respectively) which were significantly lower than standard visually guided hemostasis – 26.3 % (20/76) in the Gastro RCT. Compared to standard through the scope hemoclips, the OTSC was able to imbed in fibrotic based chronic ulcers, grasp a greater volume of tissue with the stigmata of hemorrhage in the center, and more effectively obliterate blood flow underneath NVUGI lesions.

The authors concluded that OTSC was more effective in obliterating arterial blood flow in severe NVUGIB lesions than standard visually guided endoscopic hemostasis. Residual arterial blood flow highly correlated with lesion rebleeding rates. The OTSC as primary treatment of NVUGIB lesions has the potential of significantly reducing rebleeding rates compared to other, standard visually guided hemostasis techniques. A new RCT has been planned by the study group to compare OTSC with standard hemostasis in patients with severe NVUGIB.

Why over-the-scope clip is potentially more effective than standard endoscopic hemostasis as primary treatment of severe non-variceal upper gastrointestinal bleeding

Jensen DM^{2,1,3}, Kovacs TO^{2,1}, Ghassemi KA^{2,1}, Kaneshiro M^{2,1}, Dulai G^{2,1}, Machicado GA^{2,1}

¹Santa Monica, ²Los Angeles, ³Los Angeles, CA, USA

Case control study shows decreased rebleeding rates with OTSC in high-risk peptic ulcer bleeding when compared to conventional endoscopic treatment

G. Ermerak et al., Liverpool Hospital, Sydney, Australia, presented a case control study comparing patient outcomes including risk of re-bleeding and mortality in patients with bleeding peptic ulcer disease (PUD) undergoing conventional endoscopic intervention versus OTSC application at initial endoscopy.

16 cases of bleeding PUD managed with primary OTSC application over a period of 2 years were identified from a prospectively maintained database of GI bleeding at a large tertiary center. Age and sex matched controls undergoing endoscopic intervention with conventional hemostatic treatment were used from the same database. Indications for primary OTSC use included a bleeding vessel >4mm (n=12), concurrent dual antiplatelet or anticoagulant use (n=3), likely concurrent perforation (n=2), difficult access to the bleeding site (n=1) or failure of other interventions at the initial endoscopy (n=1). Cases treated with OTSC were more likely to be hospital inpatients (12 vs 4, P=0.005), hypotensive (Median SBP 100 vs 118, P=0.04), tachycardic (Median HR 101 vs 95,

P=0.02) and have a greater proportion of Forrest I lesions (12 vs 4). There was a trend towards decreased re-bleeding within 30 days in the OTSC group (1 vs 5, P=0.07). The OTSC rebleed case required angioembolisation. All control rebleeds were managed endoscopically. 1 had OTSC salvage therapy. 30-day readmission, angioembolisation or mortality were not significantly different between the two groups.

The authors concluded that despite the presence of more high risk features the patients treated with primary OTSC application in this series had a trend towards reduced rebleeding rates and similar other outcomes when compared to conventional endoscopic therapy.

Over the Scope Clips for primary therapy of bleeding upper gastrointestinal ulcers: a retrospective case control study

Ermerak GG¹, Behary J², Koo JH^{1,3}, Levy MT^{1,3}, Abi-Hanna D¹, Edwards PD¹, Bassan MS^{1,3}
^{1,2,3} Sydney, Australia

After OTSC closure of iatrogenic colonic perforations only 7.5 % of patients required surgical intervention

D. Horenkamp-Sonntag et al., Technicians' Health Insurance, Hamburg, Germany, performed a nationwide evaluation of the use of OTSC in the colon in hospital reality using administrative codes in a large health insurance data base with about 10 Mio insured patients. OTSC in the colon was administered in 500 patients (mean age 66 years, 61 % males) in 212 hospitals. Application in combination with polypectomy was the predominant indication (62.2 %) whereas perforation during colonoscopy (16.0 %) and colonic bleeding (9.8 %) were less common indications. Various clinical settings (e. g. closure of anastomotic leaks, fistulas etc.) were applied in 12.0 % of patients. After closure of iatrogenic perforations by OTSC, only 7.5 % of cases required early surgical intervention. Operative procedures beyond 30 days after OTSC application were predominantly due to underlying diseases (colorectal cancer, diverticulitis etc.).

The authors concluded that OTSC application in the colon is predominantly employed in the context of polypectomy and iatrogenic perforations. OTSC use seems to be safe also in the colon, and in the case of iatrogenic perforation, use of OTSC is an attractive alternative to surgical closure.

Efficacy and safety of over-the-scope clips (OTSC) application in the colon: evidence generated from administrative data for an innovative endoscopic procedure

Horenkamp-Sonntag D¹, Liebertaut J¹, Engel S¹, Koop H²

¹Hamburg, ²Berlin

Successful closure of PEG-related fistula achieved with OTSC

M. A. Al Samman et al., Alpert Medical School of Brown University, Providence, Rhode Island, United States, presented a case of an elderly patient with a non-healing draining fistula from previous feeding tube placement. The 86-year-old female patient was a high-risk surgical candidate, therefore, endoscopic closure of the fistula using an OTSC was planned. OTSC placement was performed uneventful and resulted in immediate closure of the fistula. Follow-up for several months did not reveal any further leak.

A new scope for the management of gastrocutaneous fistulas using over the scope clips

Al Samman MA¹, Gao T¹, Feller E¹, Shah SA¹
¹Providence, Rhode Island, USA

Experimental study suggests that the OTSC is feasible for duodenal C-SEMS fixation

Migration of duodenal covered self-expandable metallic stents (C-SEMS) is the main cause of stent dysfunction in patients with malignant gastric outlet obstruction. The ideal method to prevent migration has not been clarified. Y. Hori et al., Nagoya City University Graduate School of Medical Sciences, Nagoya, Japan, evaluated over-the-scope clip, suture and clip methods for fixation of duodenal C-SEMS based on pathological findings. Gripping force of each device and invasion depth were assessed. The OTSC and suturing systems had a significantly higher mean gripping

force compared with the clipping system (OTSC vs clip: 13.2 vs 1.0 N, p < 0.001; suture vs clip: 8.5 vs 1.0 N, p < 0.001). OTSC compression was stronger when compared with suturing (OTSC vs suture: 13.2 vs 8.5 N, p = 0.006). The submucosal layer, but not the muscle layer, was compressed more widely and deeply by OTSC compared with clips based on pathological findings by hematoxylin and eosin staining. The authors concluded that both OTSC and suturing methods used for duodenal C-SEMS fixation were feasible compared with the clipping method. The pathological evaluation of invasion depth indicated that OTSC may be safe even for preventive use. The study suggests that these methods can be applied clinically for duodenal C-SEMS fixation.

Feasibility and safety of duodenal covered self-expandable metallic stent fixation

Hori Y¹, Hayashi K¹, Naitoh I¹
¹Nagoya, Japan

RESECT+

New demilune shaped device (Coag Dissector) for ESD allows for rapid and safe dissection

H. Neumann et al., University Medical Center Mainz, Mainz, Germany, presented a study evaluating the efficacy and learning curve of a newly developed demilune shaped device for endoscopic submucosal dissection (ESD). The demilune device allows for fast submucosal cutting due to its special design which allows rapid movements above the muscle layer. In addition, the device allows for selective grasping of the vessels thereby enabling ad hoc hemostasis. Ex vivo porcine models were utilized in an advanced endoscopic simulator of interventional endoscopy. Artificial lesions, each 2x2 cm, were created at the fundus, corpus and antrum. ESD was performed after marking of the lesions with the ESD instrument, followed by lifting of the mucosa with submucosal injection of colored saline. Afterwards, circular incision of the lesions was performed with the new ESD instrument. For resection, the submucosa was lifted with a distal clear cap and cut with the new Demilune device. Resection specimens were retrieved to evaluate if all marks were included (R0). Average size of the removed lesions was 30 mm. En-bloc resection rate was 100 % and R0 resection rate was 95 %. Mean total procedure time was 25 minutes and not dependent on the location. No perforations occurred despite the rapid dissection speed through the submucosa. Satisfaction of the endoscopist and the supporting nurse staff was high through all cases.

The authors concluded that the new demilune shaped device for ESD allows for rapid and safe dissection of the submucosa due to its inherent design. Further studies are now focusing on in vivo studies to confirm these initial results.

Evaluation of a new demilune shaped device for endoscopic submucosal dissection

Neumann H¹, Grimminger PP², Zimmermann T¹, Thieringer F¹, Rahman KF¹, Galle PR¹, Kneist W²
^{1,2}Mainz, Germany

EMR+ allows for fast enbloc resection of lesions up to 30 mm

B. Meier and K. Caca, Ludwigsburg Hospital, Ludwigsburg, Germany, reported on the ex-vivo and later in-vivo evaluation of the EMR+ technique in pigs' stomachs. The stomachs were adjusted in a special simulation model (EASIE R) to be accessible to endoscopy. Imaginary lesions were created with a template (circular, 30 mm) by coagulation. An additional working channel (AWC) was mounted on a gastroscope and used for a resection snare. The conventional working channel of the scope was used for an anchor device. For submucosal injection the LiftUp agent was used. After injection, the anchor device was used simultaneously with the snare to facilitate resection. 22 resections were performed and evaluated ex vivo, 13 resections were performed in vivo. Ex-vivo results were the following: Median size of enbloc resection specimen was 30 x 26 x 11 mm. Median procedure time (time from injection to extraction of the resection specimen) was 7 minutes. No perforations occurred. In-vivo results were the following: Median size of enbloc resection specimen was 35 x 35 x 11 mm. Median procedure time was 5 minutes.

92.3 % (12/13) of lesions could be resected macroscopically complete. No major adverse events were observed. In one case (7.3 %) a minor peritumoral bleeding occurred and was managed by coagulation. The authors concluded that EMR+ allows for fast enbloc resection of lesions up to 30 mm. The technique needs to be further evaluated in vivo and in actual mucosal neoplasia.

A novel technique for endoscopic enbloc resection for lesions up to 30 mm (EMR+)

Meier B¹, Caca K¹
¹Ludwigsburg, Germany

For further information:

ovesco endoscopy AG
 Friedrich-Miescher-Strasse 9
 72076 Tuebingen
 Germany
 science@ovesco.com

August 2019 | Conference Report of the German Coloproctology-Congress 2019

The 45th Coloproctology-Congress of the German Society for Coloproctology (Deutsche Gesellschaft für Koloproktologie – DGK) together with the Association of Coloproctology Professionals (BCD) and the Surgical Working Group for Coloproctology of the German Society for General and Visceral Surgery (CACP of DGAV) took place on April 14 – 16, 2019 in Munich, Germany.

Crohn's disease and stoma were central topics of the congress. Almost 1250 participants and 45 exhibiting companies attended. As one highlight of the congress, a prize winners meeting was held on the last day of the congress, during which PD Dr. S. Blaumeiser and colleagues presented a multicentric long-term study with the OTSC Proctology.

The OTSC Proctology is safe, effective and continence preserving

Multicentric data on the OTSC Proctology confirms good healing rates (87 %), low risk of complications and maintenance of continence in all cases

Blaumeiser S. et al. presented a multicenter-analysis evaluating the application of the OTSC Proctology regarding surgical practicability and long-term results. Data of all patients treated with the OTSC Proctology between 11/2012 and 04/2018 in 4 centers were retrospectively analyzed. Permanent fistula closure was defined as primary endpoint. Secondary endpoints were the surgical procedure, intra- and postinterventional complications, relapse rate and incontinence. Overall 62 patients (m:f = 42:20, median age 50 years, range 25 – 84 years) were included in the study. The fistula was cryptoglandular in 87 % of cases. 27 patients were presented with a low transsphincteric fistula, 19 patients with a high transsphincteric fistula, 5 patients with an intersphincteric fistula, 9 patients with a suprasphincteric fistula and 2 patients with an extrasphincteric fistula. 34 % of fistulas were relapse fistulas after other surgical therapy. Pre-interventional drainage was performed in all fistulas, median seton drainage time was 20 weeks. Healing after clip application was observed in 54/62 patients (87 %). In 2 patients, a relapse fistula developed after presumed healing, in 2 patients a fistula persistence was observed after clip removal, in 3 patients the clip detached before fistula healing and in 1 patient the clip had to be removed prematurely due to pain. The average procedure time was 32 minutes; in all hospitals, the same standardized procedure was performed: (1) debridement, (2) rinsing of the fistula tract, (3) cutting out of the anoderm before fistula closure, (4) excision of the outer fistula orifice. Complications occurred in 6 patients: in 1 case pain after clip application, in 5 cases abscess/halt and in 1 case transmural clip migration. The follow-up time was 26 weeks on average. In 39 cases the clip was removed during follow-up, in 8 cases spontaneous clip loss was found during follow-up. Impairment of continence was not observed in any of the patients.

The authors conclude, that the OTSC Proctology is a technically simple, fast and safe option for the treatment of

cryptoglandular anal fistulas. Good healing rates are achieved, the risk of complications is low, and the method does not carry the risk of incontinence.

Clinical experience with the OTSC Proctology: a multicenter-analysis.

(Klinische Erfahrung mit dem OTSC Proctology: eine Multicenter-Analyse)

S. Blaumeiser¹, J. Lutonsky², M. Maurus³, U. Nitsche¹, W. Kauer¹

¹Munich, ²Schneeberg, ³Immenstadt.

For further information:

ovesco endoscopy AG
Friedrich-Miescher-Strasse 9
72076 Tuebingen
Germany
science@ovesco.com

July 2019 | EMR+ achieves significantly higher R0 resection rate and requires significantly lower procedure time when compared to conventional EMR

In lesions of 3 cm, EMR+ reaches its best discriminatory power: R0 resection rate of 86.36 % was achieved compared to R0 resection rate of 18.18 % with conventional EMR

On April 4-6, 2019, ESGE Days (European Society of Gastrointestinal Endoscopy) took place in Prague, Czech Republic. At the conference, Knoop RF et al., University Medicine Goettingen, presented a prospective trial systematically comparing EMR+ to the gold standard of classical EMR. Standardized lesions, measuring 1 cm, 2 cm, 3 cm or 4 cm were set in an ex-vivo animal model with pig stomachs placed into the EASIE-R simulator, a well-established endoscopic model.

Overall 152 procedures were performed. In lesions of 1 cm, both EMR and EMR+ were very reliable with a R0 resection rate of 100 %. In 2 cm lesions, EMR already dropped to 54.55 %. Classical EMR did not provide sufficient resection rates for lesions with 3 cm or even 4 cm (18.18 % and 0 %). EMR+ still presented very satisfying results in 3 cm lesions (86.36 %) but also relevantly decreased at 4 cm (60.00 %). From 3 cm on, EMR+ was significantly faster than conventional EMR. A perforation rate of 15 % was observed in the 4- cm-group treated with EMR+.

The authors summarized that EMR+ enables a grasp-and-snare technique and consequently facilitates en-bloc resection of larger lesions compared to conventional EMR. In lesions > 2 cm, EMR+ outdoes its advantages, especially concerning the rate of R0 resections. At 3 cm, EMR+ reaches its best discriminatory power. At 4 cm, also EMR+ comes to its inherent limits and the risk of perforations rises. Then, ESD or surgery should be considered. The authors concluded that EMR+ could help to close a therapeutic gap in interventional endoscopy with manageable technical complexity, time and costs.

The use of an additional working channel (AWC) in endoscopic mucosal resection (EMR+) compared to conventional EMR

Knoop RF, Wedi E, Ellenrieder V, Neesse A, Kunsch S. Endoscopy 2019; 51(04): S39-S40

July 2019 | Meta-analysis comprising 475 patients demonstrates success of OTSC® as first and second line therapy for non-variceal gastrointestinal bleeding

93 % primary hemostasis was achieved with primary OTSC therapy (288 patients), 91 % with rescue OTSC therapy (187 patients)

Oforu A et al., The Brooklyn Hospital Center, Clinical Affiliate of the Mount Sinai Hospital, Brooklyn, NY, USA, performed a systematic review and meta-analysis to evaluate primary hemostasis rates and re-bleeding rates of the OTSC for primary therapy and rescue endoscopic interventions in patients with non-variceal gastro-intestinal bleeding (NVGIB).

A total of 16 studies which involved 475 patients were included. 288 patients were treated with OTSC as primary

therapy while 187 patients were treated with OTSC as rescue therapy. Primary hemostasis rate achieved with primary endoscopic therapy with OTSC was 0.93 (95 % CI: 0.89 – 0.96). Similarly, primary hemostasis rate achieved with rescue endoscopic therapy with OTSC was 0.91 (95 % CI: 0.84 – 0.95). Re-bleeding rates after primary endoscopic therapy with OTSC was 0.21 (95 % CI: 0.08 – 0.43) and 0.25 (95 % CI: 0.17 – 0.34) with rescue therapy. There was a decreased risk of re-bleeding in patients treated with OTSC as primary therapy versus rescue therapy RR = 0.52 (95 % CI: 0.31-0.89).

The authors concluded that this meta-analysis demonstrates success on the use of OTSC as primary and rescue therapy in the management of NVGIB. Further studies are however needed.

Over-the-scope-clips as primary and rescue therapy for non-variceal gastrointestinal bleeding: a systematic review and meta-analysis

Oforu A, Ramai D, John F, Barakat M, Sunkara T, Sharma S, Gaduputi V, Adler DG, Reddy M.

Minerva Gastroenterol Dietol. 2018 Nov 7. doi: 10.23736/S1121-421X.18.02513-8. [Epub ahead of print]

Ovesco Research Update 31

June 2019 | FTRD®: Multicenter Italian experience confirms effectiveness and safety of the FTRD for several high-risk indications in the colorectum

Technical success (lesion reached and resected) was achieved in 94.4 %. R0 resection was achieved on lateral and deep margins in 90 % and 92 %, respectively. Major complications (requiring surgical treatment) occurred in 1.8 %.

Andrisani G et al., University Campus Bio-Medico, Rome, Italy, reported on 110 patients (61m/49f, mean age 68±11 years, range 20-90) from 12 Italian endoscopy centers who underwent EFTR using the FTRD System. Indications were: residual/recurrent adenoma (39), incomplete resection at histology (26), non-lifting lesion (12), adenoma involving the appendix (2) or diverticulum (2), subepithelial lesions (10), suspected T1 carcinoma (16), and diagnostic resection (3).

Technical success (lesion reached and resected) was achieved in 94.4 % of cases. R0 resection was achieved on lateral and deep margins in 90 % and 92 %, respectively. Mean size of specimens was 20 mm (range 6-42). In residual/recurrent adenomas, final analysis revealed: low-risk T1 (11), adenoma with low-grade dysplasia (LGD) (24) adenoma with high-grade dysplasia (HGD) (3), scar tissue (1). Histology reports of R1 resections were: adenoma with LGD (6), adenoma with HGD (1), low-risk T1 (6) and high-risk T1 (1), scar tissue (12). Non-lifting lesions were diagnosed as: adenoma with HGD (3), low-risk T1 (7) and high risk T1 (2). Adverse events occurred in 12 patients (11 %), 2 patients (1.8 %) required further surgical treatment: In one patient, stenosis of the sigmoid colon occurred after clip deployment and had to be finally treated surgically after post-dilation perforation. One of 2 patients with resection involving the appendix developed acute appendicitis and underwent surgery. Three-months follow-up was available in all patients and revealed residual disease in 7 patients (6.4 %).

The authors concluded that EFTR using the FTRD System is a feasible, effective and safe technique for the treatment of different high-risk colorectal lesions.

Colo-rectal endoscopic full-thickness resection (EFTR) with the over-the-scope device (FTRD): A multicenter Italian experience.

Andrisani G, Soriani P, Manno M, Pizzicannella M, Pugliese F, Mutignani M, Naspetti R, Petruzzello L, Iacopini F, Grossi C, Lagoussis P, Vavassori S, Coppola F, La Terra A, Ghersi S, Cecinato P, De Nucci G, Salerno R, Pandolfi M, Costamagna G, Di Matteo FM. Digestive and Liver Disease 2019; 51:375-381.

June 2019 | Conference of DGE-BV 2019 with FTRD® Lunch Symposium

The 49th Conference of the German Society for

Endoscopy and Imaging Procedures (DGE-BV) 2019 took place together with the learned societies of CAES, CATC, DEGEA, DEGUM, DGBMT, DGD, ÖGGH and bng on March 28 – 30, 2019 in Stuttgart, Germany.

Ovesco products were presented in nine workshops on three different topics: "Hemostasis techniques", "Management of complications: perforations and post-operative leakages" (OTSC®) and "ESD techniques" (RESECT+). Additionally, a workgroup meeting of young endoscopists on OTSC® and RESECT+ took place. Besides the FTRD® Lunch Symposium on the topic "EFTR with the FTRD® – where are we today?", several talks and posters discussed products of Ovesco.

Advancement of ESD and EMR is subject of current endoscopic research – at Ovesco, too!

RESECT+

The novel high viscosity injection solution LiftUp® for ESD leads to a long-lasting lifting effect and could make the ESD procedure more effective

Endoscopic submucosal dissection (ESD) is an established procedure for endoscopic treatment of early-stage neoplasms. A crucial step of the complex procedure is the submucosal injection of an ideally highly viscous injection solution.

E. Wedi et al., University Hospital Goettingen, Germany, presented the results of a prospective, randomized preclinical study comparing the novel copolymer injection fluid LiftUp to normal saline solution (NaCl 0.9 %) and hydroxyethyl starch (HAES 6 %). A total of 60 standardized ESD procedures were performed in artificial lesions, each 3 x 3 cm in size, in an ex vivo porcine model (n=20 per injection solution). All 60 lesions were successfully resected using the standard ESD technique. R0 resection was achieved in 95 % (n=19) with LiftUp, in 100 % (n=20) with HAES and in 80 % (n=16) with NaCl. LiftUp had no procedure related perforations, one perforation occurred in the HAES group, and two perforations in the NaCl group. Adequate lifting was achieved in 16/20 cases (80 %) using LiftUp, in 6/20 cases (30 %) using HAES and in 6/20 cases (30 %) using NaCl. General ESD procedure time was shorter in the LiftUp group than in the other two groups, the difference, however, did not reach statistical significance.

The authors concluded that LiftUp appears to be a safe alternative to established fluids for ESD. It had a significantly improved lifting effect and required significantly less injected volume compared to well-established lifting solutions. With LiftUp, the ESD procedure could become more effective.

A novel high viscosity injection solution (LiftUp) for Endoscopic Submucosal Dissection (ESD). A prospective comparison study with two established lifting solutions

(Eine neue hochvisköse Injektionslösung (LiftUp) für die Endoskopische Submukosa Dissektion (ESD). Eine prospektive Vergleichsstudie mit zwei etablierten Injektionslösungen)

E. Wedi¹, C. Jung¹, J. Hochberger², J. Maiss³, C.-N. Ho⁴, G. Conrad⁴, U. Baulain⁵, V. Ellenrieder¹, P. Koehler⁵

¹Goettingen, ²Berlin, ³Forchheim, ⁴Tuebingen, ⁵Mariensee

The EMR+ technique allows bi-manual operation by triangulation

R.F. Knoop et al., University Hospital Goettingen, Germany, presented the results of a preclinical study comparing efficacy and safety of a modified EMR technique (EMR+) using an Additional Working Channel (AWC®) to the established ESD procedure. In total 40 standardized lesions (3 x 3 cm) were manually placed in ex-vivo porcine models secured in an EASIE-R simulator. 20 lesions were resected using the novel EMR+ technique and 20 lesions by ESD. Median procedure time was significantly shorter in the EMR+ group (median 10.5 min, range 4.4 – 24 min) than in the ESD group (median 32 min, range 14 – 61.6 min; p<0.0001). The median resection area per minute was significantly larger in the EMR+ group (median 53.2 cm²/min, range 15.6 – 370.4 cm²/min) than in the ESD group (median 28.1 cm²/min, range 12.7 – 64.4 cm²/min, p<0.0001). The median size of the resected specimen was smaller in the EMR+ group (median 6 cm², range 2.5 –

30 cm² vs. constantly 9 cm²). The rate of en-bloc resections was significantly lower in the EMR+ group (38 %) vs. the ESD group (95 %, $p < 0.0001$). Two perforations occurred in the EMR+ group, in the ESD group no perforations occurred.

The authors concluded that the EMR+ technique using the AWC enables the transformation of a standard one-channel endoscope into a double-channel tool, which makes bimanual operation by triangulation possible. The technique is a cost-effective option to overcome limitations of the conventional EMR as well as the ESD technique.

Comparison study evaluating Endoscopic Submucosal Dissection and a modified EMR technique using an Additional Working Channel (EMR+)

(Vergleichsstudie zwischen der Endoskopischen Submukosadisektion und einer modifizierten EMR-Technik mit einem zusätzlichen Arbeitskanal (EMR+))

R.F. Knoop¹, C. Jung¹, C.-N. Ho², G. Conrad², J. Maiss³, U. Baulain⁴, V. Ellenrieder¹, P. Koehler⁴, E. Wedi¹

¹Goettingen, ²Tuebingen, ³Forchheim, ⁴Mariensee

The BougieCap leads to high dilation rate in benign esophageal stenosis and to significant improvement of the clinical symptoms of dysphagia

BougieCap

Bougienage with the BougieCap leads to significant improvement of dysphagia-associated symptoms (DHI score regression from 51.6 to 27.9)

B. Walter and colleagues, University Hospital Ulm, Germany, presented the first prospective multicenter-study evaluating application of the BougieCap for endoscopic treatment of benign esophageal stenosis. Use of the BougieCap allows in contrast to conventional Savary-Bougies a direct optical control of the dilation process, x-ray imaging is not necessary with the BougieCap. 50 patients (25 m/25 f, median age 67.1 ± 16.8 years) with benign esophageal strictures and clinically apparent dysphagia symptoms were included in the study. Etiology of the stenoses was peptic (n=23), radiation-induced (n=13), post-surgery (n=6), consequence of alkali burn (n=4), condition after ESD (n=2), eosinophilic esophagitis (n=1) and unknown (n=1). The pre-interventional diameter of the stenoses was median 7.5 mm (± 2.4 mm). Bougienage was successful in 96 % of cases (n=48). The average number of bougienage procedures per patient was 2.3 (± 0.7). A guide wire was used in 10 cases, in 8 of these a pediatric endoscope was in use, in the 2 remaining cases a standard gastroscope. In 2 cases, bougienage had to be discontinued because passage of the endoscope through the pharynx into the esophagus was not possible. In those two cases no guidewire had been used. Dysphagia symptoms were regredient from a median DHI score of 51.6 (± 14.4) before bougienage to a score of 27.9 (± 9.3) 14 days after bougienage (Mann-Whitney: $p < 0.0001$). No major complications occurred. In 2 cases a BougieCap was lost in the stomach, which did not lead to clinical symptoms. **Results of a prospective multicenter study evaluating endoscopic treatment of benign esophageal stenoses (BougieCap study)**

(Ergebnisse einer prospektiven multizentrischen Studie zur endoskopischen Behandlung gutartiger Ösophagusstenosen (BougieCap Studie))

B. Walter¹, S. Schmidbaur¹, I. Rahman², B. Schumacher³, D. Albers³, A. Meining¹

¹Ulm, ²Southampton/UK, ³Essen

The FTRD is used for risk stratification of colorectal carcinomas

FTRD System

Full-thickness resection of colorectal cancer with the FTRD System allows exact histologic risk stratification and spares patients with "low-risk" features (62 % of the cohort) a surgical intervention

A. Kuellmer et al., University Hospital Freiburg, Germany, presented a retrospective multicenter study assessing efficacy, safety and clinical value of the FTRD System for colorectal cancer. Data of 1234 patients treated with the FTRD System for various indications at 96 endoscopic centers were screened for eligibility. 156 patients with

histological evidence of an adenocarcinoma were identified. This cohort comprised 64 patients undergoing FTRD resection after incomplete resection of a malignant polyp (group 1) and 92 non-lifting lesions (group 2). Endpoints of the study were: technical success, R0 resection rate, adverse events, successful discrimination of "low-risk" and "high-risk" tumors, as well as the need for surgical oncological resection.

Technical success was achieved in 144/156 patients (92.3 %). Mean procedural time was 42 min. R0 resection was achieved in 122/156 cases (71.8 %). Subgroup analysis showed a R0 resection rate of 87.5 % in group 1 and of 60.9 % in group 2 ($p < 0.001$). Severe procedure-related adverse events were recorded in 3.9 % of patients. Discrimination between high-risk versus low-risk tumor was successful in 155/156 patients (99.3 %). In group 1, 84.1 % were identified as low-risk lesions, whereas 16.3 % in group 2 had low-risk features. In total 53 patients (34 %) underwent oncologic resection due to high-risk features whereas 98 patients (62 %) were followed endoscopically. The authors concluded, that endoscopic full-thickness resection with the FTRD for colorectal cancer is feasible, effective and safe. It allows exact histological risk stratification and can avoid surgery for patients with "low-risk" tumors.

Endoscopic full-thickness resection of colorectal carcinomas with the FTRD-System - a retrospective multicentric study

(Endoskopische Vollwandresektion von kolorektalen Karzinomen mit dem FTRD-System - eine retrospektive multizentrische Studie)

A. Kuellmer¹, J. Mueller¹, K. Caca², P. Aepli³, D. Albers⁴, B. Schumacher⁴, A. Glitsch⁵, H. Albrecht⁶, I. Wallstabe⁷, C. Hofmann⁸, A. Erhardt⁹, B. Meier², D. Bettinger¹, R. Thimme¹, A. Schmidt¹, and the FTRD study group

¹Freiburg, ²Ludwigsburg, ³Luzern/Switzerland, ⁴Essen, ⁵Greifswald, ⁶Neumarkt i.d.OPf., ⁷Leipzig, ⁸Mainz, ⁹Wuppertal

Increased risk of appendicitis following FTRD resection of adenomas arising from the appendiceal orifice

A. Wannhoff et al., Ludwigsburg Hospital, Germany, presented a study assessing the risk of post-interventional appendicitis after FTRD resection at the appendiceal orifice. Data of all patients at the Ludwigsburg Hospital or the University Hospital Ulm, who had undergone full-thickness resection with the FTRD System at the appendiceal orifice between 2014 and 2018, was retrospectively analysed. Available follow-up data was analysed regarding development of appendicitis. Patients with appendectomy before FTRD application were excluded from the study.

Overall 44 patients (median age 68 years, range 47-85 years, n=25 female) matched the inclusion criteria. EFTR was successfully performed in all patients. During follow-up (median follow-up time 21 weeks, range 0-126 weeks) acute appendicitis occurred in 9 patients (23.7 %). 5 patients developed appendicitis within 10 days after resection, the remaining patients more than a month after the intervention. Six patients received appendectomy and 3 patients conservative treatment. All patients recovered completely.

The authors concluded that acute appendicitis can develop after full-thickness resection with the FTRD at the appendiceal orifice. Development of appendicitis can occur early after the intervention or after a period of longer latency. During informed consent discussion, patients should be informed about the appendicitis risk and a possible need for surgical treatment.

Evaluation of the risk of appendicitis following FTRD resection of adenomas arising from the appendiceal orifice

(Untersuchung des Appendizitis-Risikos nach FTRD-Resektion am Appendixabgang)

A. Wannhoff¹, B. Walter², T. Kreuzer¹, S. Schmidbaur², B. Meier¹, A. Meining¹, K. Caca¹

¹Ludwigsburg, ²Ulm

Application of the FTRD System can avoid surgical revision in many patients with colorectal relapse adenoma

A. Schmidt, University Hospital Freiburg, Germany, held an expert review on the treatment of colorectal relapse adenomas. Non-pedunculated adenomas are mostly resected by EMR. Meta-analyses show that the risk of a relapse adenoma following this resection method is about 15 %. Piece-meal resection is a significant risk factor for recurrence of the adenoma. The risk of recurrence is about 20 % after piece-meal resection, while it is 3 % after en-bloc resection ($p < 0.0001$). The ESD technique carries a risk of recurrence of about 4.8 %. Especially problematic for the follow-up resection is the fact that the previous resection results in scarring of the resection area, the relapse adenoma shows a non-lifting sign in most cases. In non-lifting lesions, EMR resection is usually not possible, ESD resection is time-consuming and technically extremely difficult when the lesion is located beyond the rectum. Besides, the risk of perforation is high when using the ESD technique. Another problem in the therapy of relapse adenomas is the fact that only 62 % of patients undergo endoscopic follow-up examination, this results in non-discovery of residual lesions and relapse adenomas or large size of residual and relapse lesions at a late follow-up date.

The FTRD System is very suitable for therapy of colorectal relapse adenomas. On the one hand it can be used within the entire colon and rectum, on the other hand, also non-lifting lesions can be resected with the FTRD System. Several mono- and multicenter clinical studies, among them the prospective multicenter study WALL RESECT incorporating 181 patients, showed technical success rates of about 90 %, R0 resection rates of about 80 % and rates of major complications of about 2 %. FTRD application is related to short procedural times and is relatively simple. FTRD usage is limited in large lesions, the optimal lesion size for this method is ≤ 2 cm. In larger lesions often hybrid techniques are possible (i.e. EMR + FTRD). In conclusion, application of the FTRD represents for many patients with colorectal relapse adenoma the avoidance of revision surgery.

Colorectal relapse adenomas: always FTRD?

(Rezidiv-Adenome im Kolorektum: Immer FTRD?)

A. Schmidt¹,

¹Freiburg

Lunch Symposium on the FTRD System at the Conference of the DGE-BV 2019

During the 49th conference of the German Society for Endoscopy and Imaging Procedures (DGE-BV), which took place on March 28th – 30th, 2019 in Stuttgart, experts evaluated and discussed in a lunch symposium the current status of the clinical application of the FTRD System.

From the point of its introduction into clinical use until now, the FTRD System has become a well-established instrument for colonic full-thickness resection and can be used safely and effectively in hospitals of all levels of care

Prof. Dr. Karel Caca, Ludwigsburg: EFTR with the FTRD System is considered the standard procedure for resection of colonic lesions, which are not resectable or not safely resectable using conventional techniques

K. Caca, Ludwigsburg hospital, Germany, presented several recent studies on FTRD application. Especially in German-speaking countries the FTRD System is already in wide clinical use and has established as an efficacious and safe procedure. For resection of colorectal lesions, which are not resectable (i.e. non-lifting adenomas) or not safely resectable (i.e. lesions at diverticula) using conventional methods, EFTR with the FTRD System nowadays is considered the method of choice. For large, partly non-lifting lesions, hybrid-techniques in combination with EMR or ESD are used.

The FTRD is also increasingly used in the upper gastrointestinal tract. Here, it is currently preferably used for resection of relapse duodenal adenomas with non-lifting sign and for resection of submucosal tumors. Furthermore, full-thickness resection with the FTRD allows for definite diagnosis of gastric and duodenal lesions.

In conclusion, Caca stated that endoscopic full-thickness resection with the FTRD System has developed into a

standard procedure in interventional endoscopy. Limitations of FTRD application are size and localization of the lesion, as well as availability of endoscopic equipment and expertise. Caca emphasized the importance of expertise and close cooperation with surgical colleagues.

Is the FTRD nowadays a standard in endoscopic resection procedures? A positioning statement (Ist das FTRD heute ein Standard bei endoskopischen Resektionsverfahren? Eine Positionsbestimmung)

K. Caca¹

¹Ludwigsburg

PD Dr. Arthur Schmidt, Freiburg: prospective data of the German colonic FTRD registry confirm efficacy and safety of the colonic FTRD in broad clinical application

A. Schmidt, University Hospital Freiburg, Germany, presented the final evaluation of the German colonic FTRD registry. Between 09/2015 and 12/2018, prospective data from 64 endoscopic centers in Germany have been collected in the registry in order to assess safety and efficacy of the colonic FTRD in routine clinical use in hospitals of all levels of care. Data from overall 1176 patients (median age 69 years, range 11 – 93 years, 39.7 % female) have been evaluated. Primary endpoints were technical success (successful clip application and resection with the FTRD System) and R0 resection rate. Secondary endpoints were rate of full-thickness resection, procedural time, complications, and follow-up findings. Indications for FTRD application were: non-lifting relapse adenomas (n=454), native non-lifting adenomas (n=223), adenoma at diverticula (n=10), adenoma at appendiceal orifice (n=87), follow-up resection of malignant polyp (n=171), early carcinoma (n=64), submucosal tumor (n=74), diagnostic FTRD without lesion (n=20) and others (n=73). Location of lesion and intervention was in 52.1 % rectal, in 14.2 % sigmoidal, in 5.3 % in the descending colon, in 12.8 % in the transverse colon, in 18.3 % in the ascending colon, in 19.2 % coecal, and in 5.2 % in other locations (right flexure, appendix, terminal ileum, recto-sigmoidal transition, left flexure, coecal pole). Average lesion size was 15 mm (range 2 – 56 mm, lesions with very large diameter were resected by hybrid-technique. Median procedure time was 35 min (range 2 – 203 min), technical success rate was 89.5 %, the rate of full-thickness resection was 87.3 % and of R0 resection 77.3 %. Minor or moderate complications (no need for surgery) occurred in 14.5 % of cases, major complications (need for surgery) in 1.8 %. Follow-up data after median 4.8 months was available for 595 patients. Rx/R1 situation was diagnosed in 43 patients (7.2 %), suspected local recurrence of initially diagnosed R0 resected lesion in 38 patients (6.4 %). Most of these patients underwent EFTR once more (33.3 %) or EMR/ESD (30.8 %).

(note: the listed data is preliminary, outstanding information on interventions is still delivered by users. Therefore, some values might still change. Final data can be extracted from the following publication)

For summary, Schmidt presented a pooled-proportion analysis of the registry data in comparison with recent meta-analysis data of EMR/ESD of large adenomas (>2 cm). With regard to technical success rate and R0 resection rate, the EFTR is in the range of the ESD, significantly higher than the success rates of the EMR. In regard to complication rates the following applies: the EFTR corresponds in peri-interventional bleeding rates to the conventional techniques, the perforation rate of EFTR is about as high as with EMR, which is significantly lower than with ESD. The need for surgery because of complications is with EFTR technique approximately as high as with the conventional resection techniques, however, this is frequently owed to a wrong sequence of the interventional steps.

Schmidt concluded, that the German colonic FTRD registry confirms that the colonic FTRD is safe and efficacious; results of the registry equal those of clinical studies, for example of the WALL RESECT study, a prospective multicenter study. Hence, the colonic FTRD has found broad-based appliance with comparable quality to large centers.

German colonic FTRD registry – final evaluation of 1176 cases

(Deutsches colonic FTRD Register – Abschlussauswertung über 1176 Fälle)

A. Schmidt¹

¹Freiburg

Prof. Dr. Peter Bauerfeind, Zuerich: FTRD application in the upper gastrointestinal tract requires more experience, but is feasible, effective and safe

P. Bauerfeind, city hospital Triemli, Zuerich, Switzerland, presented on the basis of exemplary clinical cases the application of the FTRD in the upper gastrointestinal tract. In three cases, the FTRD was used in the duodenum, in one case in the distal esophagus and in another case in the stomach. Histological examination of the full-thickness resectate of patient 1 showed a well differentiated neuroendocrine tumor (NET G1) with infiltration of the muscular layer (ca. 4 mm DM), the resection margins were tumor-free. Patient 2 had a tubulovillous adenoma with focal high-grade dysplasia at the duodeno-jejunal transition, full-thickness resection was successful, no complications occurred, histological examination showed adenoma-free and dysplasia-free lateral and basal margins. In patient 3, a duodenal relapse adenoma with non-lifting sign was found, also in this case full-thickness resection with the FTRD was successful. Due to the lack of time, Bauerfeind then directly mentioned patient 5, this patient presented with a GIST at the gastric basis with Rx/R0 situation. Follow-up resection with FTRD was effective and without complications.

FTRD application in the upper intestinal tract is in some points technically more challenging than application in the lower GI-tract: for one thing, the pylorus represents a difficult-to-pass narrow and has to be dilated, on the other hand the mucosal layer of the stomach is very motile, which can cause difficulty especially in submucosal and small lesions in identifying whether the target tissue is sucked in correctly into the cap of the endoscope. Before application of the FTRD in the duodenum, the papilla should be reliably identified. The risk of post-intervention hemorrhage is increased in the upper gastrointestinal tract when compared to the lower GI-tract, it occurs in about 10 % of cases. This should be addressed in the informed consent discussion. Visible vessel endings should always be obliterated after FTRD application. Close monitoring of the patient after the intervention allows timely endoscopic re-intervention in case of hemorrhage.

Despite these technical specialties, also for the upper gastrointestinal tract the FTRD is a very promising instrument, first clinical applications show very good results. Indications for FTRD application are in a particular way submucosal tumors (NET, GIST...) up to 1.5 cm, adenomas and relapse adenomas with or without non-lifting sign and adenomas in patients at high surgical risk. Besides, the FTRD is very well suitable for follow-up resections of endoscopically pre-treated lesions with R1-finding.

FTRD in the upper gastrointestinal tract: duodenum and stomach

(FTRD im oberen Verdauungstrakt: Duodenum und Magen)

P. Bauerfeind¹

¹Zuerich, Switzerland

For further information:

Ovesco Endoscopy AG
Friedrich-Miescher-Strasse 9
D-72076 Tuebingen
science@ovesco.com

May 2019 | Systematic review and meta-analysis of OTSC® literature data covering 1868 patients confirms high technical and clinical success rates

Pooled durable clinical success proportions of OTSC use were 87.5 % in hemorrhage, 81.4 % in acute leaks/perforations and 63.0 % in chronic leaks/fistulae. Since its market launch in 2007, the OTSC System has been the object of intensive clinical research. Weiland T et al, Tuebingen, performed a systematic collection of all data for post-market clinical follow-up (PMCF) aiming to

systematically assess efficacy and safety of the OTSC.

The PMCF database was systematically searched for clinical data on OTSC therapy of GI hemorrhage (H), acute leaks/perforations (AL) and chronic leaks/fistulae (CL). Major outcomes were successful clip application and durable hemostasis/closure of defects.

457 publications were reviewed. 58 articles comprising 1868 patients fulfilled criteria to be included in the analysis. These consisted of retrospective analyses, prospective observational trials, one randomized-controlled trial (STING) and one quasi-controlled study (FLETROCK). The pooled proportion analysis revealed high overall proportions of technical success: H – mean 93.0 % [95 % CI 90.2 – 95.4], AL – mean 89.7 % [95 % CI 85.9 – 92.9] and CL – mean 83.8 % [95 % CI 76.9 – 89.7]. Pooled durable clinical success proportions were: H – mean 87.5 % [95 % CI 80.5 – 93.2], AL – mean 81.4 % [95 % CI 77.0 – 85.3] and CL – mean 63.0 % [95 % CI 53.0 – 72.3].

By pooling all clinical data gained, the authors concluded that OTSC application in GI hemorrhage and closure of GI lesions is safe and effective in real clinical use.

Efficacy of the OTSC System in the treatment of GI bleeding and wall defects: a PMCF meta-analysis.

Weiland T, Rohrer S, Schmidt A, Wedi E, Bauerfeind P, Caca K, Khoshab MA, Hochberger J, Baur F, Gottwald T, Schurr MO.

Minimally Invasive Therapy & Allied Technologies, DOI: 10.1080/13645706.2019.1590418.

<https://doi.org/10.1080/13645706.2019.1590418>

April 2019 | FTRD® resection in the colorectum ensures local radical excision where other endoscopic techniques do not suffice and reduces the need for surgery

Prospective data from 51 FTRD-procedures shows technical success in 88 % (45/51) and a R0-resection rate of 89 % (40/45). EFTR-specimens, obtained for indeterminate previous T1 colorectal carcinoma resection, were free of residual carcinoma in 25 of 26 cases (96 %)

Van der Spek B et al., Department of Gastroenterology and Hepatology, Northwest Hospital group, Alkmaar, The Netherlands, published a single-center case series evaluating resection of colorectal lesions with the FTRD device.

Between 07/2015 and 10/2017, 51 EFTR procedures were performed in 48 patients (63 % men, median age 69). Indications for EFTR were non-lifting adenoma (n = 19), primary resection of malignant lesion (n = 2), resection of scar tissue after incomplete endoscopic resection of low-risk T1 colorectal carcinoma (n = 26), adenoma involving a diverticulum (n = 2) and neuroendocrine tumor (n = 2). Two lesions were treated by combining endoscopic mucosal resection and EFTR. Technical success was achieved in 45 of 51 procedures (88 %). Five resections were macroscopically incomplete and in one case no specimen could be obtained due to inability to mobilize the lesion into the cap. Mean estimated lesion size was 12.2 mm (range 2 – 30). Mean maximum resection specimen diameter was 23 mm (range 11 – 45). Histopathology confirmed full-thickness resection in 43 of the remaining 45 specimens (96 %) and radical resection (R0) in 40/45 procedures (89 %). EFTR-specimens, obtained for indeterminate previous T1 colorectal carcinoma resection, were free of residual carcinoma in 25 of 26 cases (96 %). In six patients (13 %) a total of eight adverse events occurred within 30 days after EFTR. Four of these patients had minor bleeding not necessitating blood transfusion. One patient suffered major bleeding needing blood transfusion. One perforation occurred, which was corrected endoscopically. No emergency surgery was necessary. Six patients needed additional surgery because of either high risk for lymph node metastases (n = 4; three patients with T2-T3 CRC and one patient with lymphatic invasion in T1 CRC), technical failure of EFTR (n = 1) or endoscopically untreatable adenoma recurrence at surveillance (n = 1).

The authors concluded that this study confirmed safety and efficacy of the FTRD device for the resection of colorectal lesions. They proposed a clinical algorithm for EFTR case selection. EFTR allows en-bloc and transmural resection

where other advanced endoscopic techniques are unsuitable and reduces the need for surgery in selected cases.

Endoscopic full-thickness resection in the colorectum: a single-center case series evaluating indication, efficacy and safety

Van der Spek B, Haasnoot K, Meisli C, Heine D
Endoscopy International Open 2018; 06:E1227-E1234

April 2019 | First case series with the AWC (Additional Working Channel) shows bimanual instrumentation to facilitate endoscopic resection of large GI lesions

The AWC was used in 8 patients with GI lesions 17-45 mm (average 35.9 mm) in size for EMR with modified grasp-and-snare technique (4 cases) or ESD (4 cases). R0-resection was achieved in 6 cases. Complications were acute arterial bleeding post-EMR in two cases treated by endoscopic clipping.

Walter B et al., Department of Gastroenterology, InExEn, University Hospital Ulm, Ulm, Germany reported on the use of the AWC (Additional Working Channel) in eight patients with large, flat lesions or early stages of cancer in the upper or lower gastrointestinal tract.

Endoscopic en-bloc resection of large, flat GI lesions is challenging. No bimanual tasks are possible using standard endoscopes. Dual-channel endoscopes are not available everywhere and have a small distance between the channels. The AWC can be fixed to the tip of a standard gastroscope or pediatric colonoscope, a second endoscopic tool can be inserted through the AWC and used for a distinct tissue traction and counter-traction during endoscopic resection, as the distance between working channels can be adjusted.

Eight patients with large, flat lesions or early stages of cancer in the upper or lower gastrointestinal tract were treated with endoscopic mucosal resection (EMR) with a modified grasp-and-snare technique (4 cases) and endoscopic submucosal dissection (ESD, 4 cases) using the AWC. Mean procedure time (scope-in to scope-out) was 68.5 minutes. R0-resection was achieved in all of the four cases treated by ESD and in two cases treated by EMR. The remaining two cases were treated by piecemeal EMR with positive lateral margins. Complications were acute arterial bleeding post-EMR in two cases treated by endoscopic clipping. No delayed bleeding, no perforation and no further severe adverse events occurred.

The authors concluded that the AWC enables endoscopic resection of large lesions in the upper and lower gastrointestinal tract. Benefits are its suitability for EMR and ESD, no need for a dual-channel endoscope and an adjustable distance of working channels.

Improved endoscopic resection of large flat lesions and early cancers using an external additional working channel (AWC): a case series.

Walter B, Schmidbauer S, Krieger Y, Meining A.
Endoscopy International Open 2019;07: E298-E301.

March 2019 | FTRD® is feasible and safe in T1 early rectal cancer: a case series and video report

6 consecutive T1-ERC patients underwent FTRD application which resulted in complete endoscopic resection in all cases. No complications nor disease recurrence were observed during the 1-year follow-up period.

Soriani P et al., IRCCS Policlinico San Donato, San Donato Milanese, Milan, Italy performed a study aiming to assess the feasibility and long-term clinical impact of endoscopic treatment with the FTRD in patients with high-risk malignant rectal polyps = T1 early rectal cancer (T1-ERC). According to Ueno et al. histopathological criteria for high risk malignant polyps are defined as follows: low tumor differentiation grade (G3), Haggitt's levels (pedunculated polyps) 3-4, Kikuchi's levels (sessile polyps) sm3, width of submucosal invasion $\geq 4000 \mu\text{m}$, depth of submucosal invasion $\geq 2000 \mu\text{m}$, positive tumor budding, distance from the excision margin $< 1 \text{ mm}$, presence of vascular invasion. 6 consecutive patients (5 men, mean age 63 years, range

51-78 years), who had undergone a rectal EMR within the previous 1-3 months and thereby received the histopathological diagnosis of high risk malignant polyp, were treated with the FTRD System. Previous endoscopic mucosal excision had been incomplete in all cases. Three patients were unfit for surgery (ASA IV), the other three patients refused surgery.

The endoscopic full-thickness resection was technically feasible in all cases within 8-15 minutes. No immediate or late complications occurred. All patients could be discharged home within the following 24 hours. The histopathologic analysis performed on the en-bloc-resected specimen demonstrated a complete endoscopic resection in all patients. During follow-up, all patients underwent an oncologic work-up with endoscopy, CT scan, and rectal EUS every 6 months. All patients were in oncologic remission after a median follow-up of 12 months (range 12-18 months) without any radio- or chemotherapy. One patient died from cardiac failure at the 8th month of follow-up after showing no sign of disease recurrence at the 6-month oncological work-up.

The authors concluded that this study provides initial evidence in favor of the FTRD System as therapeutic option for T1 early rectal cancer in patients either unfit for surgery or refusing the standard surgical approach.

A full-length demonstrative video as well as high resolution endoscopic images displaying the original rectal lesions, the EFTR procedures, and the follow-up have been published at *Endoscopy International Open*: <https://doi.org/10.1055/s-0043-118657>.

Endoscopic full-thickness resection for T1 early rectal cancer: a case series and video report.

Soriani P, Tontini GE, Neumann H, de Nucci G, De Toma D, Bruni B, Vavassori S, Pastorelli L, Vecchi M, Lagoussis P.

Endoscopy International Open 2017; 05:E1081-E1086.

Risk factors for an adverse outcome in early invasive colorectal carcinoma.

Ueno H, Mochizuki H, Hashiguchi Y et al.
Gastroenterology 2004; 127: 385 – 394.

Ovesco Research Update 30

March 2019 | Application of the FTRD® in duodenal lesions is feasible, efficacious and safe

EFTR of duodenal lesions in 20 patients showed technical success in 17/20 (85.0 %), a R0 resection rate of 12/17 (70.6 %) and only minor peri-procedural bleeding in 3/20 cases (15.0 %). No major bleedings or perforations occurred.

Bauder M et al., Department of Gastroenterology and Oncology, Ludwigsburg Hospital, Ludwigsburg, Germany, performed a study investigating the safety and efficacy of FTRD application in the duodenum. Between 03/2014 and 06/2017, a total of 20 patients underwent EFTR of a duodenal lesion. Indication for EFTR was: adenomas (n=13, seven treatments naïve, six pretreated), subepithelial tumors (n=5) and T1 adenocarcinoma (n=1). The FTRD could be advanced to the lesion in 19/20 cases (95 %). In one case advancing the FTRD through the pylorus was not possible despite balloon dilatation. Overall technical success was 17/20 (85.0 %). In two cases the FTRD clip was deployed correctly, but the integrated snare could not be closed because of device dysfunction. Both lesions were then resected with a standard snare above the FTRD clip after extraction of the device. In both cases there was no macroscopic evidence of residual adenoma. However, R0 resection could not be confirmed in both cases. R0 resection rate in the technically successful FTRD applications was 12/17 (70.6 %). Minor bleedings occurred at the first post-interventional day in 3/20 patients (15.0 %). No major bleedings or perforations occurred. During follow-up after 3 and 12 months, there were two recurrent adenomas that were successfully re-resected by FTRD. The authors concluded that the FTRD indicates good technical efficacy and safety for resection of duodenal non-ampullary adenomas and subepithelial tumors. It offers the possibility of re-resections at the same site. Especially in pretreated or difficult lesions, such as non-lifting adenomas, EFTR should be considered.

Endoscopic full-thickness resection of duodenal lesions-a retrospective analysis of 20 FTRD cases.

Bauder M, Schmidt A, Caca K.

United European Gastroenterol J 2018 Aug;6(7):1015-1021.

February 2019 | Large multicenter study with 286 patients confirms safety and efficacy of the OTSC® for hemostasis of high-risk lesions

96.4 % hemostasis rate is reported from 286 emergency endoscopies for either upper or lower gastrointestinal bleeding in whom the OTSC was used as first-line therapy.

Manta R et al., Digestive Endoscopy Unit, S. Agostino-Estense Hospital, Modena, Italy, published a study comprising data from all patients with upper or lower gastrointestinal bleeding from eleven Italian tertiary endoscopic referral centers who underwent OTSC placement as first-line treatment between 2014 and 07/2017. The decision to use an OTSC as first-line endoscopic treatment depended on the endoscopist's evaluation. Briefly, all high-risk patients and/or those with high-risk bleeding lesions were considered for the OTSC approach. Patients on antithrombotic therapy and those with relevant comorbidities (heart, kidney, and hepatic impairment) were considered as high-risk. Bleeding lesions were classified as high-risk when at least one of the following conditions was present: lesion with bleeding artery or vessel larger than 2 mm in diameter visible, lesion deeply penetrating, excavated or fibrotic in which the presence of a microperforation could not completely be ruled out or thermal therapy could increase risk of perforation, and lesion not suitable for safe treatment by other endoscopic devices.

Overall 286 patients with either UGIB (n=214) or LGIB (n=72) were included. 112 patients (39.2 %) were receiving antithrombotic therapy at the time of the intervention. Technical success and primary hemostasis rates were 97.9 % and 96.4 %, respectively. Early re-bleeding occurred in 4.4 %, more frequently in those with antithrombotic therapy, and no late re-bleeding was observed. Following a successful primary hemostasis, only 5.2 % of patients needed blood transfusions, the median hospital stay was 4 days (range 3 -11). 18 patients with either technical failure (n=6) or re-bleeding (n=12) underwent radiological or surgical approaches. Overall, bleeding-related deaths occurred in 5 patients (1.7 %), including 3 patients with technical procedural failure and 2 in the re-bleeding group.

The authors conclude that the data from this large multicenter study show that OTSC placement is an effective first-line treatment for hemostasis in high-risk patients with lesions in the upper or lower gastrointestinal tract.

First-line endoscopic treatment with over-the-scope clips in patients with either upper or lower gastrointestinal bleeding: a multicenter study.

Manta R, Mangiafico S, Zullo A, Bertani H, Caruso A, Grande G, Zito FP, Magniavillano B, Pasquale L, Parodi A, Germana B, Bassotti G, Monica F, Zilli M, Pisani A, Mutignani M, Conigliaro R, Galloro G.

Endoscopy International Open 2018;06:E1317-E1321.

February 2019 | OTSC® is safe and effective in pediatric patients for acute GI bleeding throughout the GI tract

Case series of 11 OTSC applications in 10 pediatric patients shows 100 % technical success with immediate hemostasis and no complications.

Tran P et al., UT Southwestern Medical Center, Children's Health – Children's Medical Center Dallas, TX, USA, published the center's experience utilizing OTSCs for nonvariceal gastrointestinal bleeding in pediatric patients. Overall 10 patients (median age 14.7 years, range 3.9 – 16.8 years, median weight 39 kg, range 17.4 – 85.8 kg) underwent 11 endoscopic procedures utilizing the OTSC System for hemostasis. Upper GI bleeding due to stomach or duodenal ulcer was seen in 4 patients, 2 of these had ulcer disease of the stomach and duodenum, respectively,

secondary to nonsteroidal anti-inflammatory drug use. 1 patient had peptic ulcer disease of unknown etiology and 1 had duodenal ulcers secondary to active *Helicobacter pylori* infection. Upper intestinal bleeding was found in 2 other patients, 1 with postpolypectomy bleeding in the stomach and 1 with biliary sphincterotomy. Lower intestinal bleeding was seen in the remaining 4 patients. 1 patient had an ulcer located in the sigmoid colon presumed to be secondary to intestinal ischemia, 1 had postpolypectomy bleeding in the sigmoid colon, and 2 patients had anastomotic ulcers, 1 at an ileoileal and 1 at an ileocolonic anastomosis. 4 patients (40 %) had OTSC placed as first-line intervention.

Placement of the OTSC was technically successful in all patients resulting in immediate hemostasis. No complications occurred. The two patients with anastomotic ulcers have continued to have clinical bleeding resulting in chronic anemia. One of these patients continues to require monthly iron infusions and the other remains transfusion dependent. The remaining 8 patients have had no evidence of recurrent bleeding at follow-up (median follow-up time 32.9 months, range 21.2 – 39.4 months). The authors concluded that the OTSC System is a reliable and effective tool for active GI bleeding or high-risk lesions and should be considered for high risk or urgent/emergent cases of bleeding in children.

Over the Scope Clips for Treatment of Acute Nonvariceal Gastrointestinal Bleeding in Children Are Safe and Effective.

Tran P, Carroll J, Barth BA, Channabasappa N, Troendle DM.
JPGN 2018;67: 458-463.

February 2019 | FTRD® is rated a highly effective, time-sparing and safe alternative to ESD for the treatment of non-lifting colorectal lesions and/or scars from R1 resection

Italian single-center study of 20 consecutive patients with non-lifting lesions, adenoma recurrence/relapse of previous endoscopic resections, and scars from incomplete endoscopic resection reports on 100 % full-thickness resection and no major bleeding or perforation.

Andrisani G et al., University Campus Bio-Medico, Rome, Italy, published a single center study evaluating consecutive patients with superficial colorectal neoplasms, who underwent endoscopic full-thickness resection (EFTR) using the FTRD device. Inclusion criteria were non-lifting lesions, adenoma recurrence/relapse of previous endoscopic resections with a negative lifting sign, and scars from incomplete endoscopic resection with a positive deep margin (R1 resection).

Between 1/2016 and 9/2016, 20 patients (12 m/ 8f, mean age 67, range 51-79 years) met the inclusion criteria. Indications included scars of incomplete resections of polypectomy/mucosectomy (7/20), non-lifting lesions (9/20), and adenoma recurrence/relapse (4/20). The lesions were located in the rectum (n=10), descending colon/sigmoid (n=5), transverse colon (n=1), ascending colon (n=3) and caecum (n=1). Technical success, defined as full-thickness resection, was obtained in 100 % and was subsequently histologically confirmed. Among the 7 resection scars treated with FTRD, all were histologically negative for neoplasia. Among the 9 non-lifting lesions, seven were histologically diagnosed as T1/G1/sm1, one as T1/G1/sm2, and one as T1/G1/sm3; this patient underwent a surgical resection. Mean size of the resected lesions was 26 mm, ranging from 10 to 42 mm. No major immediate or delayed bleeding was observed. There was no perforation or need for emergency surgery. One (5 %) patient developed abdominal pain, fever and leukocytosis after the intervention. Perforation or abscess was excluded by CT scan; the patient was treated conservatively with medical therapy. The endoscopic follow-up after 3 months showed spontaneous OTSC dislocation in 100 % of patients and the scar histology was negative for neoplasia. The authors concluded that EFTR is a feasible and effective technique and a valid alternative to EMR and ESD in the management of recurrent adenomas, non-lifting lesions and scars of R1 resections.

Endoscopic full-thickness resection of superficial colorectal neoplasms using a new over-the-scope clip system: A single-centre study.

Andrisani G, Pizzicannella M, Martino M, Rea R, Pandolfi M, Taffon C, Caricato M, Coppola R, Crescenzi A, Costamagna G, Di Matteo FM.
Digestive and Liver Disease 2017; 49:1009-1013.

January 2019 | OTSC®: A meta-analysis of 1517 cases over 9 years confirms its outstanding clinical effects for rescue therapy

OTSC use in patients with refractory gastrointestinal diseases achieved an overall clinical success rate of 78 %, 85 % for bleeding, 85 % for perforation, 52 % for fistula, 66 % for anastomotic dehiscence, and 95 % for other conditions. Overall OTSC-associated complications were 1.7 %, severe OTSC-associated complications 0.59 %.

Kobara H et al., Departments of Gastroenterology und Neurology, Faculty of Medicine, and Gastroenterological Surgery, Faculty of Medicine, Kagawa University, Takamatsu, Japan, published a meta-analysis clarifying the current status and limitations of OTSC according to different indications of GI refractory disease, including refractory bleeding, perforation, fistula, and anastomotic dehiscence. An extensive literature search identified studies reporting on 10 or more cases, in which the OTSC System was applied. A total of 1517 cases described in 30 articles were retrieved. The clinical success rates and complications were calculated overall and for each indication.

The average clinical success rate was 78.3 % (n = 1517) overall, 84.6 % for hemorrhage (n = 559), 84.6 % (n = 351) for perforation, 51.5 % (n = 388) for fistula, 66 % (n = 97) for anastomotic dehiscence, and 95.1 % (n = 122) for other conditions, respectively. The authors rated these results, despite the lower performance of the OTSC System for fistula, as more than satisfactory when considering that there are no other effective endoscopic methods currently available and these refractory conditions hitherto required surgical interventions. With respect to safety, the overall OTSC-related complication rate was 1.7 % (26/1517 cases), the incidence rate of severe complications that required surgery was 0.59 % (9/1517 cases).

The authors concluded that the OTSC system serves as a safe and effective device for GI refractory diseases, which hitherto required surgical interventions.

Over-the-scope clip system: A review of 1517 cases over 9 years.

Kobara H, Mori H, Nishiyama N, Fujiwara S, Okano K, Suzuki Y, Masaki T.
J Gastroenterol Hepatol 2018 Aug 2 doi:10.1111/jgh.14402

January 2019 | Multicentric analysis: OTSC® is highly effective as first- and second-line treatment for NVUGIB in high-risk patients with cardiovascular disease and complex, large ulcers

In a multicentric cohort of 100 consecutive patients with mean age 72 y, 51 % severe cardiovascular comorbidity, and 73 % on antiplatelet or/and anticoagulation therapy presenting with non-variceal upper GI bleeding from ulcers with median size of 3 cm, OTSC therapy led to 94 % primary hemostasis and 86 % long-term clinical success.

Non-variceal upper gastrointestinal bleeding (NVUGIB) is a common clinical problem with high rates of morbidity and a mortality rate between 5-10 %. An aging patient population with a high prevalence of cardiovascular comorbidity has led to increase in the incidence of NVUGIB in patients older than age 70 years. A high comorbidity rate has also been identified as an independent risk factor for complications and mortality after NVUGIB.

Wedi E et al., Department of Gastroenterology und GI Oncology, University Medical Centre Goettingen, Germany presented prospectively collected multicentric data from 02/2009 to 09/2015 from all patients who underwent emergency endoscopy for high-risk NVUGIB in two

academic centers and were treated with OTSC as first-line (n=81) or second-line therapy (n=19). 100 consecutive patients (mean age 72 years, range 27-97) were included in the study. 51 % had severe cardiovascular comorbidity (ischemic heart disease, congestive heart failure, hypertension, valvular heart disease, peripheral arterial occlusive disease and atrial fibrillation) and 73 % were on antiplatelet or/and anticoagulation therapy, other comorbidities in the patient cohort included kidney disease (n=20), former or recent malignancy (n=23), respiratory disease (n=16), liver disease (n=6) and diabetes mellitus (n=26). The median size of the treated ulcers was 3 cm (range 1-5 cm), the patients in 85 % of cases were classified as ASA 3 to 5 (ASA 3: 40 %, ASA 4: 40 %, and ASA 5: 5 %). In 94 % of patients (n=94) primary hemostasis with OTSC was achieved. The primary endpoint with successful initial hemostasis and no early (≤ 24 h) or delayed (≤ 30 days) re-bleeding was achieved in 86 % of patients. Mean number of OTSCs placed on the initial endoscopic exam was 1 (range 1-3). Long-term hemostasis (6 months) was achieved in 86 %.

The authors concluded that in this cohort the OTSC System was demonstrated to be a safe and effective first-and second-line treatment for NVUGIB in high-risk patients with cardiovascular disease and complex, large ulcers.

Use of the over-the-scope-clip (OTSC) in non-variceal upper gastrointestinal bleeding in patients with severe cardiovascular comorbidities: a retrospective study

Wedi E, von Renteln D, Gonzalez S, Tkachenko O, Jung C, Orkut S, Roth V, Tumay S, Hochberger J
Endoscopy International Open 2017; 05: E875-E882 | <http://dx.doi.org/10.1055/s-0043-105496>

Ovesco Research Update 29

December 2018 | Significant reduction of rebleeding rates in patients with high-risk NVUGIB by OTSC®

OTSC use decreased the rebleeding rate in high-risk (RS ≥ 8) patients with statistical significance compared to the rates reported by the Rockall study (0 % vs. 53 %, $p < 0.01$). Also in intermediate-risk (RS = 4 – 7) patients rebleeding was reduced (0 % vs. 24 %, $p = 0.08$).

Asokkumar et al., Singapore General Hospital, Singapore, studied rebleeding and mortality rates of patients treated with OTSC for high-risk adverse outcome (HR-AO) non-variceal upper gastrointestinal bleeding (NVUGIB).

The Rockall data and a historic cohort of the same institution (52 patients with peptic bleeding) were used for comparison. 18 patients with 19 bleeding lesions were included: 9 (47 %) duodenal ulcers, 4 (21 %) Dieulafoy's lesions, 3 (16 %) gastric ulcers, and 3 (16 %) bleedings after gastric biopsy, gastric polypectomy and endoscopic ultrasound-guided fine-needle aspiration of peri-gastric mass. OTSC was applied as first-line treatment in 10 (53 %) and as second-line treatment in 9 (47 %) lesions. Complete hemostasis was achieved in all patients. There were no complications associated with OTSC placement. OTSC use significantly decreased 0 % vs. 53 %, $p < 0.01$ and reduced (0 % vs. 24 %, $p = 0.08$) the rebleeding rate in high-risk (RS ≥ 8) and intermediate-risk (RS = 4 – 7) Rockall score patients as compared to the rates reported by the Rockall study, respectively. When compared to the institution's prior study, a decrease in the rebleeding rate was found with OTSC (0 % vs. 21 %, $p = 0.06$) in the intermediate-to-high risk Rockall score patients (RS ≥ 4). There was no difference in mortality rates as compared to both control studies.

The authors concluded that use of OTSC is safe, efficacious and appears superior to standard treatment for HR-AO NVUGIB. OTSC should be considered as first-line treatment for HR-AO bleeding.

Use of over-the-scope-clip (OTSC) improves outcomes of high-risk adverse outcome (HR-AO) non-variceal upper gastrointestinal bleeding (NVUGIB).

Asokkumar R, Soethniko R, Sanchez-Yague A, Wie LK, Salazar E, Ngu JH.
Endoscopy International Open 2018; 06: E789-E796 C.

December 2018 | Conference Report United European Gastroenterology Week (UEGW)

2018

- **RESECT+**: additional working channel (AWC®) and temperature-dependent agent for submucosal injection (LiftUp®) enable fast endoscopic en-bloc resection of specimens up to 30 mm
- **BougieCap**: prospective multicenter study shows 96 % successful bougienage and no complications
- **FTRD**®: endoscopic full-thickness resection of rectal neuroendocrine tumors is feasible, safe and effective and allows for definite diagnosis and treatment in the same session
- **OTSC**®: large systematic review (2462 patients) shows 77-96 % clinical success by OTSC in various indications without the need for further intervention

The 26th United European Gastroenterology Week (UEGW) was held on October 20-24, 2018, in Vienna, Austria. Several workshops, talks and posters presented original research with Ovesco technology and procedures. Hands-on training sessions in the ESGE learning area with the OTSC System attracted lively interest.

RESECT+**Additional working channel (AWC) effectively supports endoscopic resection of large lesions in the upper and lower GI tract**

B. Walter et al., Department of Internal Medicine I, University Hospital Ulm, Germany, presented first experiences using the additional working channel (AWC). The device can be fixed at the tip of a standard gastroscope or pediatric colonoscope. The distance of the two working channels can be adjusted by the endoscopist. Via the AWC a second endoscopic tool can be inserted and used for bi-manual handling. ESD and EMR with a modified 'grasp-and-snare' technique was performed, EMR in 4 patients (1 with lesion in the upper GI tract, 3 with lesions in the lower GI tract), and ESD in 4 patients (2 with lesions in the upper GI tract, 2 with lesions in the lower GI tract). Mean procedure time was 68.5 min. Reported complications were acute arterial bleeding post EMR in two cases treated in the same session. No delayed bleeding or perforation were reported. Passage with the AWC-equipped endoscope was possible in all cases. The authors concluded that the AWC effectively supports endoscopic resection of large lesions in the upper and lower GI tract. Potential benefits are its suitability for EMR and ESD, no need for dual-channel endoscope and an adjustable distance or working channels.

EMR+: the new technique allows for fast endoscopic en-bloc resection of lesions up to 30 mm

B. Meier and K. Caca, Department of Internal Medicine, Klinikum Ludwigsburg, Germany presented preclinical data on a new EMR technique (EMR+). This technique allows for en-bloc resection of specimen > 20 mm, which are usually resected in piecemeal EMR or by ESD, which however is time-consuming and associated with a higher risk for complications. EMR+ was developed and evaluated in an ex vivo porcine stomach. The stomach was adjusted in a special simulation model to be accessible to endoscopy. An additional working channel (AWC) was mounted on a standard gastroscope and used for a resection snare. The conventional working channel of the scope was used for an anchor device. For submucosal injection a newly developed agent with a temperature-dependent viscosity (LiftUp) was used. The agent has liquid consistency at room temperature, which allows submucosal injection. At body temperature, the agent gels and forms a stable cushion within seconds, which provides stable resection conditions (no deminishing over time, re-injections are not necessary). The effectiveness and safety of this agent has already been shown in vivo in domestic pigs. Imaginary lesions of 30 mm were marked by coagulation. After injection, the anchor device was used for tissue lifting simultaneously with the snare to facilitate resection. After the resection technique was established, 22 resections were performed and evaluated. The median size of the en-bloc resection specimens was 30 x 26 x 11 mm (max. 40 x 33 x 14 mm). The procedure times were between 6-7 minutes. No perforations occurred. The authors concluded that the EMR+ technique allows for fast

en-bloc resection and obtains resection specimens of 30 mm.

BougieCap**Endoscopic treatment of benign stenosis using the BougieCap enables direct visual control of the bougienage**

B. Walter et al., Department of Internal Medicine I, University Hospital Ulm, Germany, presented a prospective interventional study on patients with a benign oesophageal stenosis and with clinical symptoms of dysphagia treated with the BougieCap at three endoscopy units in Germany and UK. 50 patients (m/f 25/25) underwent the procedure, mean age was 67.1 years (± 16.8). Etiology of strictures was peptic (n=23), radiation (n=13), anastomosis (n=6), caustic ingestion (n=4), post ESD (n=2), EoE (n=1) or unknown (n=1). Successful dilatation with the BougieCap was possible in 96 % (n=48). On average 2.3 (± 0.7) BougieCaps of subsequent sizes were used per patient. A stiff guide-wire was used in 10 cases to aid with bougienage, using a pediatric scope in 8 cases and a standard gastroscope in 2 cases. In two cases with a narrow stricture and no usage of guide wire treatment failed as a result of high resistance at the site of stricture causing buckling of the endoscope in the pharynx. Symptoms of dysphagia (as assessed per Dysphagia Handicap Index score) decreased significantly after bougienage in short-time follow-up (14 days post-interventional). No severe complications were reported. Adverse events were loss of 2 BougieCaps in the stomach causing no symptoms. The authors concluded that endoscopic treatment of benign stenosis using the BougieCap enables direct visual control of the bougienage procedure and therefore of mucosal damage within the area of strictures. This might help to adapt endoscopic treatment even more precisely to the stricture. Symptoms of dysphagia are improved in short-term follow-up. Additional wire guidance is reasonable.

FTRD System**3-year multicenter UK experience: EFTR highly successful in the treatment of colonic lesions not previously amenable to endoscopy**

I. Rahman et al., Department of Gastroenterology University Hospital Southampton, UK, presented data from the UK FTRD registry. Registry data from 04/2015 – 01/2018 comprised 52 cases of FTRD application in 8 centers. Patients had a median age of 72 years (39-93). The target lesion could be reached with the FTRD mounted on top of the endoscope in 51/52 patients (98 %). 1 caecal lesion could not be reached due to sigmoid diverticulosis. Median total procedure time was 45 minutes (10-150). Median FTRD insertion time was 5 minutes (1-100). Median specimen size was 22 mm (10-30). Technical success was achieved in 88 % (45/51). Technical difficulty was experienced in 9 cases: In 6 cases snare closure was not possible, in 3 cases the lesion slipped from the grasper on clip deployment. R0 resection was achieved in 74 % (38/49), for two patients, histological data was incomplete. Residual/recurrent lesions at follow-up were found in 7 % (2/30). Complications occurred in 3 patients; 1 acute appendicitis at day 6 after resection of appendix base adenoma, 1 arterial fibrillation and hypotension, and 1 rectal bleeding. There were no cases of perforation or fistula. The authors concluded that treating colonic lesions with the FTRD shows high success rates and low complication rates, making EFTR a viable alternative to surgery.

Pooled analysis from all studies that report on FTRD use (532 patients): 77.5 % R0-resection rate, 5.4 % complication rate

A. Wannhoff et al., Department of Internal Medicine, Klinikum Ludwigsburg, Germany, reported on a study analyzing all so far published data with the FTRD System (published studies and relevant congress abstracts). A total of 18 studies were included, 9 of them published as a full-text and 9 as congress abstracts, which comprised a total of 532 patients from 7 countries. The target lesion was reached with the FTRD mounted on top of the endoscope in 522 (98.1 %) patients and technical success was achieved in 486 (91.4 %) patients. The full-thickness

resection was histologically confirmed in 326 of 401 (81.3 %) patients, in the remaining 131 no data on this endpoint was reported. The R0 resection rate was 77.5 % and achieved in 383 of 494 patients for which data on resection margins was reported. Technical problems were mostly related to the resection snare, which occurred in 34 cases. In most of this cases a successful resection however was achieved by use of a conventional resection snare following clip application with the FTRD. Complications included minor bleeding and post-polypectomy syndrome in 14 (2.6 %) patients each. Severe bleeding occurred in 2 (0.4 %) patients and perforations were reported in 13 (2.4 %) patients. A surgical intervention due to a FTRD related complication was necessary in 9 (1.7 %) patients. The authors concluded that the FTRD system provides high efficacy in the colorectum. The complication rate is low and most complications can be managed conservatively or endoscopically.

EFTR with the FTRD for rectal NET is feasible, safe and effective and allows for definite diagnosis and therapy at once

B. Meier and K. Caca, Department of Internal Medicine, Klinikum Ludwigsburg, Germany, presented a study evaluating EFTR for rectal neuroendocrine tumors. All cases of rectal NETs in the German FTRD registry, which comprises data of FTRD procedures of 31 German centers, were retrospectively analyzed. 40 patients (19 male, 21 female, median age 58 years, range 28-81) met the inclusion criteria. Lesions were located in the lower (n=13), middle (n=24) and upper rectum (n=3). Median size of the lesions was 8.4 mm (3-25). Biopsies were taken before EFTR in 19 patients and EMR had been performed in 10 patients prior to EFTR, histology had shown well differentiated NET in all cases. However, in all cases resection status was unclear or incomplete. 6 NET (15 %) were recurrent NET and had been treated previously (multiple forceps biopsies or snare resection). Mean procedure time of EFTR was 23 minutes (range 7-60 minutes). A full-thickness resection specimen could be obtained in all cases. R0-resection was achieved in all cases. However, in 7 cases (28 %) a NET could no longer be proven. Adverse events occurred in 5 cases (12.5 %), 4 patients suffered peri-interventional bleeding, which could be managed endoscopically in all cases, in 1 patient a technical problem occurred (rupture of the FTRD snare, resection was performed with a conventional snare). Follow-up data was available for 32/40 patients. Mean follow-up time was 17 weeks (1-45 weeks). Residual or recurrent tumors were not found during follow-up. The authors concluded that EFTR of rectal NET < 20 mm is feasible, safe and effective and allows diagnosis/risk stratification and therapy (R0 resection) at once. The technique should be considered as first-line therapy.

OTSC System**Lively interest in Hands-On Trainings with the OTSC System**

The European Society of Gastrointestinal Endoscopy (ESGE) offered an ESGE Learning Area to all delegates of the UEGW to provide a platform for live encounter and interaction among aspiring endoscopists and renowned experts in the field.

In the ESGE Learning Area, three 90-minute Hands-On Trainings with the OTSC System were offered. All Hands-On Trainings were fully booked.

Besides, a talk on the OTSC System was held in the ESGE Lunch Session (A. Caputo: "Advantages of the OTSC System in the treatment of UGIB") and the exhibition of Ovesco products attracted lively interest.

Large systematic review shows 77-96 % clinical success of OTSC in various indications without the need for further intervention

N. Bartell et al., Department of Gastroenterology and Hepatology, University of Rochester, United States, reported on a systematic review with the OTSC System. The study evaluated a large body of literature to determine the overall efficacy and safety of OTSC. 81 case series/retrospective reviews/prospective studies (Group A with a total of 2285 patients) and 157 case reports (Group B with

a total of 177 patients) were included.

In Group A, technical success of OTSC placement was 95.3 %, with a clinical success of 77.2 %. Indications for OTSC placement were fistula closure (30.6 %), bleeding (28.9 %), perforation closure (16.3 %), leaks (15.1 %), EFTR (8.4 %) and stent fixation (0.7 %). Complete luminal obstruction (n=1) was the only reported adverse event across all studies. 24/81 papers reported the need for surgery despite OTSC placement (90/673 patients, 13.4 %).

Indications for OTSC placement in Group B were fistula closure (37.9 %), perforation closure (33.9 %), bleeding (14.1 %), EFTR (7.9 %) and leaks (6.2 %). Pooled technical success in this group was 99 % and clinical success was 96.0 %. 7/177 (4 %) patients required surgical intervention despite OTSC placement. Complete luminal obstruction in 1/177 patients and small bowel fixation with pneumoperitoneum in 1/177 patients were the only OTSC related adverse events reported.

The authors concluded that the OTSC is a safe and effective, surgery-sparing endoscopic tool in today's GI practice with 77-96 % of patients achieving clinical success without the need for further intervention. Technical success of > 95 % has been reported across all indications.

OTSC for high-risk peptic ulcer bleeding: one and done in 75 %

S. Gölder et al., Department of Internal Medicine III, Klinikum Augsburg, Germany, presented a study evaluating the use of OTSC for the treatment of high-risk peptic ulcer bleeding (HRUB).

Between 4/2014 and 03/2018, 100 patients with peptic ulcer bleeding (Forrest Ia-Ib), in the stomach of the duodenum were treated with OTSC. The OTSC was used as first-line procedure in 66 patients. Successful primary hemostasis could be achieved in 89.4 %. The OTSC was used as secondary treatment after failure of an initial endoscopic treatment in 34 patients. OTSC clipping led to successful primary hemostasis in 94.1 %. Recurrent bleeding occurred in n=9 for primary OTSC (15.3 %) and in n=7 patients with secondary OTSC (21.9 %) (p=0.812).

No treatment beside the single OTSC clip was necessary in 75.8 % (n=50) in the primary-OTSC arm and in 73.5 % (n=25) in the secondary-OTSC arm, respectively.

OTSC failure occurred more often in large ulcers (> 3 cm, p=0.03), in the duodenal bulb (p=0.03) and in ulcers with negative helicobacter test (p=0.045). The patients with OTSC failure received more blood transfusions (p=0.002). No statistical difference was found for the Rockall score (median 7.5), the Glasgow Blatchford score (median 15.5), NSAID use or anticoagulation.

The authors concluded that the OTSC has a high rate of bleeding control in first- and second line treatment of peptic ulcer bleeding. Potential risk factors for treatment failure are location in the duodenal bulb, longer ICU stay, higher amount of transfusions and a higher reimbursement per case.

For questions and further information:

Ovesco Endoscopy AG
Scientific Information Service
Dorfackerstraße 26
D-72074 Tübingen
science@ovesco.com

November 2018 | Korean multicenter study confirms efficacy and safety of OTSC® for GI fistulas, leaks and perforations

All 19 patients were treated successfully with the OTSC System. In 74 % of cases, complete healing of the leakage was achieved using OTSC alone

Lee HK et al., Department of Internal Medicine, Hanyang University College of Medicine, Seoul, Korea, performed a prospective multicenter study at seven centers in Korea aiming to examine the therapeutic success rate of endoscopic treatment of gastrointestinal perforations, fistulas and anastomotic leakages using the OTSC Ssystem in Korean patients.

A total of 19 patients were included, with gastrointestinal leakages from gastrojejunostomy sites (n=3), esophago-

jejunostomy sites (n=3), esophago-gastrostomy sites (n=4), esophago-colonostomy sites (n=1), jejuno-jejunal sites (n=1), endoscopic full thickness resection site closures (n=2), Boerhaave's syndrome (n=1), esophago-bronchial fistulas (n=2), gastro-colonic fistula (n=1), and colono-pseudocyst fistulas (n=1). The size of the leakage ranged from 5 to 30 mm (median diameter 10 mm). The median procedure time was 16 min. All cases were technically successful. Complete healing of the leakage was achieved in 14 of 19 patients (74 %) using OTSC alone. There were no complications associated with the OTSC procedures.

The authors concluded that the OTSC System is a safe and effective method for the management of gastrointestinal defects, especially in cases of anastomotic leakage after surgery.

Efficacy of the Over-the-Scope Clip System for Treatment of Gastrointestinal Fistulas, Leaks, and Perforations: A Korean Multi-Center Study.

Lee HK, Cho JY, Cho JH, Park JJ, Kim CG, Kim SH, Han JH (2017).

Clin Endosc 2017, Aug 29. [Epub ahead of print]
<https://doi.org/10.5946/ce.2017.027> Print ISSN 2234-2400
• On-line ISSN 2234-244.

November 2018 | Visually controlled dilatation with the BougieCap is effective and prevents complications due to overdilatation

Multicenter study shows high technical success rate and significant improvement of dysphagia symptoms.

B Walter et al. presented a multicenter study (Ulm, Southampton and Essen) evaluating dilatation of benign esophageal stenoses with the BougieCap. The BougieCap allows, in contrast to Savary bougies, direct visual control of the process without the need for x-ray.

50 patients (25 f, 25 m, median age 67.1 ± 16.8) with benign stenosis of the esophagus and clinically apparent symptoms of dysphagia were included. Cause of the stenosis was peptic (n=23), radiation (n=13), anastomotic (n=6), caustic ingestion (n=4), Post-ESD (n=2), EoE (n=1) and unknown (n=1). Dilatation was successful in 96 % of all cases (48/50). In eight cases a pediatric gastroscope with guidewire was used. In two cases a standard gastroscope with guidewire was used. In the two cases, passage of the stenosis was not possible, no attempt with guidewire had taken place. BougieCap (median 2.3 ± 0.7) of different sizes were used per session. Dysphagia symptoms were reduced from a median DS value of 3.0 ± 0.6 before dilatation to 1.6 ± 0.7 after dilatation (Mann-Whitney, p < 0.0001). No major complications occurred. In two cases, a BougieCap was lost in the stomach; no clinical discomfort of complications resulted.

The authors stated that endoscopic treatment of benign esophageal stenoses with the BougieCap allows direct visual control of the dilatation process and of beginning mucosal lacerations. Thus, in contrast to the conventional blind method, overdilatation and re-traumatization are reduced and the dilatation process can be performed with better adaptation to the stenosis. Usage of a guidewire is reasonable and necessary in special cases (i.e. very high-grade stenosis, usage of a pediatric gastroscope).

The BougieCap: a new method for endoscopic treatment of esophageal strictures*.

Walter B, Schmidbaur S, Rahman I, Schumacher B, Albers D, Meining A.

*UEG Week, Vienna, Austria, October 20-24, 2018.

November 2018 | Conference Report DGVS / DGAV

The 73th annual conference of the German Society of Gastroenterology, and Digestive and Metabolic Diseases with Endoscopy Section (DGVS) took place together with the 12th autumn conference of the German Society of General Surgery and Visceral Surgery (DGAV) on September 12-15, 2018 in Munich, Germany.

Ovesco products were presented in talks, posters, research, innovation and video forums and hands-on training sessions. Dr med. Edris Wedi (University Hospital Goettingen) received the DGVS endoscopy research award and the award of the Olympus Europe foundation

2018 for his work.

FTRD® is described as effective and safe resection device for lesions otherwise difficult to treat endoscopically.

colonic FTRD

Meta-analysis of all to date published data (777 patients) regarding FTRD application in the colorectum shows 78 % R0 resection rate and < 1 % surgery because of complications

A Wannhoff et al. presented a pooled analysis of all published data (full texts and conference contributions) evaluating FTRD application in the colorectum. 21 studies comprising overall 777 patients were included. The target lesion was reached in 746 (96 %) cases. Resection was technically successful in 684 (88 %) cases. Main reason for technical failure were problems with the snare in 35 cases, in 29 of these cases resection succeeded with subsequently introduced conventional snare. Histological examination confirmed full-thickness resection in 326 of 401 (81.3 %) cases and R0 resection in 383 of 494 (78 %) cases; respective information was not available for the remaining cases. Complications occurred in overall 8 %, thereof 13 post-polypectomy-syndroms (1.7 %), 16 minor haemorrhage (2 %), 2 major haemorrhage (0.3 %), 12 perforations (2 %, partially due to wrong order of operational steps), and 5 appendicitis (0.6 %); surgery because of complications was necessary in < 1 % of patients. The authors concluded that FTRD application in the colorectum is safe and the target lesion can be successfully resected in the majority of cases.

Effektivität und Sicherheit des Full-Thickness Resection Device (FTRD) im Kolorektum: Ergebnisse einer gepoolten Analyse bisher veröffentlichter Daten (Efficacy and safety of the Full-Thickness Resection Device (FTRD) in the colorectum: Results of a pooled analysis of to date published data).

Wannhoff A, Meier B, Caca K, Ludwigsburg.

University Hospital of Erlangen achieves 85 % R0 resection rate of adenomas and early adenocarcinomas with FTRD

T Rath and colleagues presented the experiences gathered at the University Hospital of Erlangen with the FTRD System. Between 06/2015 and 09/2017, the FTRD was applied in 14 patients (7 m, 6 f, median age 64.5 ± 6.1y) with colorectal adenomas and early T1 adenocarcinomas. The lesions had a median size of 16 ± 4.7 mm and were located in the rectum (n=6), caecum (n=2), ascending colon (n=2), left flexure (n=1), and right flexure (n=3). The technical success rate was 100 %. The procedural time was 72 ± 40 min. R0 resection was achieved in 85 % of cases (11/13). Histological examination of the specimen yielded the following findings: adenoma with low-grade intraepithelial neoplasia (n=6), adenoma with high-grade intra-epithelial neoplasia (n=4), fibrotic area without dysplasia (n=2), adenocarcinoma (n=2). No complications occurred. In one patient, a relapse lesion was found at the resection site during follow-up, which could be resected once again with the FTRD. In conclusion, the authors rated the endoscopic full-thickness resection with the FTRD a procedure with few complications and high technical success and R0 resection rate.

Endoskopische Vollwandresektion mittels FTRD für die Resektion von kolorektalen Adenomen und frühen T1 Adenokarzinomen: eine Single Center Erfahrung (Endoscopic full thickness resection using the FTRD for resection of colorectal adenomas and early T1 adenocarcinomas: A single centre experience).

Vitali F, Naegel A, Siebler J, Neurath M, Rath T, Erlangen.

colonic FTRD is suitable for therapy of early colorectal carcinoma

S Herrmann et al. presented the experiences of the Clinical Center in Neuperlach, Munich, with the FTRD System. Between 01/2015 and 04/2018 indication for eFTR was present in 30 patients. In 11 patients malign histology was previously known (8 adenocarcinomas, thereof 2 verified by biopsy, 6 pre-treated with resulting Rx/R1 situation, 3 NETs). The remaining 19 patients showed residual or recurring adenoma, non-lifting sign or difficult localisation.

Technical success was achieved in 25 of 30 patients (83.3%). In 5 patients, the procedure had to be discontinued due to insufficient mobilisation of the lesion into the cap (n=4) or due to failure of the snare (n=1). 4 of the 5 technically unsuccessful procedures took place in the first 12 patients. The R0 resection rate in technically successful procedures was 92 % (23/25). Resected lesions were 27.8 ± 6.4 mm in size. 12 of the 25 resected specimen proved to be malignant, 11 (91.7 %) of those were resected in R0. Oncological surgery was necessary in 3 patients. Thus, colonic FTRD was the curative treatment in 84 % (21/25) of patients. Complications occurred in 3 cases: post-polypectomy-syndrome (n=1), clip failure (n=1), and perforation (n=1). All complications could be managed conservatively. In 1 patient, a relapse polyp was found during follow-up after 3 months (SSA without dysplasia). In conclusion, the authors stated that the future primary application field of the FTRD could be small pretherapeutically verified carcinomas, because eFTR yields a save R0 resection option and enables definitive evaluation of sm-invasion for histologically based therapy stratification.

Endoscopic full-thickness resection (eFTR): Effektivität der eFTR für komplexe kolorektale Läsionen, insbesondere als Therapie des kolorektalen Frühkarzinoms (Efficacy of eFTR for complex colorectal lesions, especially for the treatment of early stage colorectal carcinomas). Herrmann S, Götzberger M, Blöchliger M, Dollhopf M, Ulm.

Retrospective multicenter study shows 71 % R0 resection using the FTRD in different conventionally-not-resectable lesions

I Krutzenbichler and colleagues presented "real life" data gathered in the Clinical Center in Munich evaluating FTRD application in various cases. Overall, data from 61 procedures in 59 consecutive patients undergoing eFTR was retrospectively analysed. Indications for eFTR were: 25 % adenocarcinoma in colon and stomach, 11 % flat adenoma with non-lifting sign, 12 % relapse adenoma and 11 % neuroendocrine tumors. The size of the resected lesions was 20.2 ± 5.5 mm. The primary technical success was 70.5 %. Minor peri-interventional bleeding occurred in 4.9 %. Post-interventionally, further complications occurred in 4 cases (3 bleedings, 1 gangrenous appendicitis treated by emergency ileocecal resection). The R0 resection rate was 71 %, the full thickness resection rate was 80.36 %. The authors concluded, that eFTR shows a high success rate in resecting different lesions across the entire colon and a low rate of procedural complications. The FTRD is regarded as alternative to surgery for lesions that cannot be resected with conventional methods (EMR/ESD).

Anwendung der endoskopischen Vollwandresektion mithilfe des „full thickness resection device“ (FTRD) bezüglich technischer Erfolgsrate sowie prozeduraler Komplikationen – eine Single Centre Studie (Application of the endoscopic full-thickness-resection using the "full thickness resection device" (FTRD) regarding technical success rate and procedural complications – a single-center study).

Krutzenbichler I, Fuchs M, Lewerenz B, Leimbach T, Nehrlich A, Schopp W, Gundling F, München.

FTRD resection at the appendiceal origin can spare over 80 % of patients a surgical procedure

T Kreutzer and colleagues presented a study evaluating the risk of post-interventional appendicitis following FTRD application at the appendiceal origin. All patients of the Clinical Center Ludwigsburg and the University Hospital Ulm undergoing endoscopic full-thickness resection at the appendiceal area using the FTRD between 2014 and 2018 were analysed retrospectively. The available follow-up data was analysed in regard of the development of appendicitis. Patients that had undergone an appendectomy prior to FTRD application were not included in the study. Overall 38 patients (65.8 % female, median age at FTRD application 68 years (47-85)) met the inclusion criteria. FTRD application was successful in all cases. During follow-up (average of 21 weeks, range 0-126 weeks) 9 patients (23.7 %) developed acute appendicitis. In 5 patients the appendicitis occurred within 10 days after FTRD

application, in the 4 remaining cases more than a month after the procedure. In 6 cases, an appendectomy was performed, the remaining 3 patients were treated conservatively. The authors concluded that about a fifth of all patients undergoing FTRD application at the appendiceal origin developed acute appendicitis. The complication may occur early after FTRD application or with greater latency. Patients should be informed about the risk of appendicitis development before FTRD application at the appendiceal origin.

Untersuchung des Appendizitis-Risiko nach endoskopischer Vollwandresektion von Adenomen im Bereich der Appendix mit dem FTRD System (Evaluation of the risk of appendicitis following endoscopic full-thickness resection of adenomas close to the appendix using the FTRD System).

Kreutzer T, Walter B, Schmidt A, Meier B, Wannhoff A, Schmidbaur S, Meining A, Caca K, Ludwigsburg/Ulm.

gastric FTRD

RESET study: reliable dignity determination of gastric SETs using gFTRD

Meier B and colleagues presented a multicenter prospective pilot study evaluating the use of the gFTRD for endoscopic full-thickness resection of sub-epithelial tumors (SETs) of the stomach. Gastric SETs are rare, mostly benign and usually coincidentally found during gastroscopy. Superficial biopsy is often insufficient for reliable histological assessment. Endoscopic resection with standard methods (EMR/ESD) is often not possible and associated with an increased risk for complications. The study assessed feasibility, efficacy and safety of endoscopic full-thickness resection using the gFTRD for resection of gastric SETs in 29 patients. Lesions up to 15 mm in size were included. In 77 % of cases initial histology could not provide a reliable dignity determination of the SET. With full-thickness resection, the dignity of all SETs could be reliably determined. Average lesion size was 11 mm (range 5 – 15 mm). Median procedure time was 36.3 min (24 – 90 min). 76 % (22/29) of the specimen were resected in R0, 65.5 % (19/29) in full-thickness. In 31 % of cases peri-interventional minor bleeding occurred, which could be directly treated endoscopically. In the follow-up examination after 3 months, clips were already dislocated in 81 % of the cases, there was no evidence for relapse or residual lesions in any case. The authors concluded, that endoscopic full-thickness resection with the gFTRD is a safe and effective procedure, which enables in contrast to conventional biopsy a reliable dignity determination of gastric SETs. Sufficient risk stratification (in case of GIST/NET) is possible. Besides, sufficient therapy by R0 resection is achieved in most cases.

Endoskopische Vollwandresektion subepithelialer Tumoren des Magens mit dem gFTRD-System – Eine prospektive Pilotstudie (RESET Studie) (Endoscopic full thickness resection of sub-epithelial tumours of the stomach with the gFTRD-system – A prospective pilot study (RESET study)).

Meier B, Schmidt A, Meining A, Caca K, Ludwigsburg, Freiburg, Ulm.

OTSC® System – presented studies confirm superiority of the OTSC in acute gastrointestinal haemorrhage

Marburg: OTSC highly effective for the treatment of acute ulcer bleeding

A Waldthaler presented retrospective data gathered in the University Hospital of Giessen and Marburg evaluating different endoscopic modes of therapy for non-variceal upper gastro-intestinal bleeding (NV-UGIB). Between 09/2016 and 1/2018, 131 patients (median age 68 years, 77 male) with NV-UGIB were treated. In 68 patients, the bleeding required intervention at the time of examination. Cause of hemorrhage was a peptic ulcer in 47 cases (69.1 %; 31 duodenum, 13 stomach, 1 cardia, 2 anastomosis), a Mallory-Weiss syndrome in 7 cases (10.3 %), tumor bleeding in 6 cases (8.8 %), angiodysplasia in 5 cases (7.4 %), and other causes in 3 cases (4.4 %). Primary endoscopic therapy consisted of a combination approach using injections and hemoclipping (n=15), injections (n=10), hemoclipping (n=9), OTSC

(n=12, thereof 8 for duodenal ulcer) and thermal coagulation (n=1). 9 of the 68 treated patients suffered from recurrent ulcer bleeding (6 from a duodenal ulcer, 2 from anastomosis, 1 patient with Mallory-Weiss syndrome), none of these had received OTSC as primary therapy (rebleeding rate primary OTSC vs primary other treatment 0 % vs 8 %; p=0.001). 4 of the 6 patients suffering rebleeding from duodenal ulcer were treated with OTSC. The two remaining patients received a combination therapy consisting of injection and hemoclipping, both patients developed a second rebleeding which in turn was treated using an OTSC Clip. The authors concluded that therapy of acute ulcer bleeding with the OTSC proves to be highly efficient as primary and secondary therapy. They enhanced the fact that none of the patients included in the present study, which received an OTSC, developed recurrent bleeding. Advantages of OTSC treatment especially arose in the therapy of duodenal ulcer not only in cases of recurrent bleeding but also as primary therapy.

OTSC Therapie der nicht varikösen oberen gastrointestinalen Blutung im Klinikalltag – eine retrospektive Analyse (OTSC treatment of non-variceal upper gastrointestinal bleeding in hospital routine – a retrospective analysis).

Lerner C, Waldthaler A, Wißniewski TT, Bauer C, Grote T, Gallmeier E, Gress TM, Denzer U, Marburg.

Augsburg: closure of ulcer bleedings with high risk of recurrence: one and done in 75 %

S Gölder et al. presented a retrospective study comprising all patients with high-bleeding-risk ulcers (Forrest Ia-Ib), treated with OTSC at the Augsburg Hospital. A total of 100 patients with peptic ulcer, primarily or secondarily treated with OTSC, were included (n=25 with gastric ulcer, n=75 with duodenal ulcer, primary OTSC treatment n=66, secondary OTSC treatment n=34). Primary hemostasis by OTSC without further endoscopic treatment was achieved in 92 patients (92 %, n=60 primary therapy, n=32 secondary therapy). In 8 cases hemostasis could not be achieved with one single OTSC clip. In 17 cases recurrent bleeding occurred 1-12 days after initially successful hemostasis (n=10 primary therapy, n=7 secondary therapy). The group of patients with unsuccessful OTSC treatment showed significantly larger ulcers (median size 3 cm, IQR 2 – 3, 13; p=0.03), more frequent bleeding in the duodenal bulb (22 vs. 2, p=0.033), more frequent negative H.p. status (p=0.045) and significantly higher number of transfused ECs (p=0.002). No significance was reached regarding the Rockall score (median 7.5, p=0.69) nor regarding the Glasgow-Blatchford score (median 15.5, p=0.15). Also, NSAID or anticoagulant treatment was not significantly different between the groups (p=0.53 and p=0.44, respectively). The authors concluded, that OTSC Clip application for peptic ulcer bleeding shows high clinical success rates as primary and secondary therapy. Possible risk factors for therapy failure are ulcer size, localization of the bleeding source in the duodenal bulb, negative H.p. status and increased demand for transfusion.

Over the Scope Clip (OTSC) bei Magen- und Duodenalulcera mit hohem Blutungsrisiko – One and done? (Over-the-Scope Clip (OTSC) for gastric and duodenal ulcers with high bleeding risk – one and done?)

Gölder S, Neuhaus L, Stücker J, Ebigo A, Braun G, Probst A, Weber T, Freuer D, Messmann H, Augsburg, Deutschland.

Analysis of the STING treatment cases: haemorrhage treatment with OTSC in comparison to standard therapy not only cost-effective, but cost-cutting

A Küllmer et al. presented results of a study based on data gathered during a prospective randomized study (STING), exploring whether OTSC treatment is more cost-effective than conventional clips due to the higher success rate, despite of the higher price per clip. Two parameters for cost effectivity were calculated: (1) ICER (Incremental Cost Effectiveness Ratio): defines additional expenses for additional clinical results, meaning Δ costs of both alternatives divided by Δ clinical effect. (2) ACER (Average Cost Effectiveness Ratio): costs arising from a specific clinical result. The clinical status that had to be achieved

was similar to the primary outcome of the STING study: successful hemostasis without any recurrent bleeding. The parameters for the total procedure, including costs for accommodation etc. were calculated as well as the costs for the endoscopic treatment only. The overall costs of standard treatment approaches were 13,025.95 €, versus 12,776.19 € for OTSC treatment; costs for the endoscopic procedure alone were 2,100.03 € (standard therapy) versus 1,960.17 € (OTSC-therapy). The ICER regarding the overall treatment was -589.01 € and -329.86 € for the endoscopic treatment. The ACER for the overall costs was 30,721.58 € for standard therapy and 15,066.26 € for OTSC therapy. ACER for the endoscopic procedure showed 4,952.90 € and 2,311.52 € for standard and OTSC treatment respectively. As a conclusion, OTSC therapy of recurrent ulcer bleeding was rated cost-effective and cost-cutting when compared to standard approaches.

OTSC- versus Standard-Therapie der Rezidiv-Ulkusblutung: eine Kosteneffektivitätsanalyse (OTSC versus standard treatment of recurrent ulcer bleeding: an analysis of cost effectiveness).

Küllmer A, Behn J, Glaser N, Thimme R, Caca K, Schmidt A, Freiburg Ludwigsburg, Deutschland.

Cross-sector routine data from social health insurance confirms safety and efficacy of colonic OTSC

D Horenkamp-Sonntag et al., German Technicians' Health Insurance, Hamburg, presented a study based on cross-sector routine data gathered by social health insurance (>10 million insured parties), examining OTSC application in the colon. Indication, patient characteristics, outcome and complications were assessed in the actual care setting. 348 patients (median age 67 years, 60 % male) were subject to colonic OTSC (OPS-Code 5460s3). Using further codes from different performance sectors, suspected indications were identified: (iatrogenic) perforation (n=58), polypectomy (n=210), bleeding (n=34) and others (n=46). A total of 16 patients (4.6 %) underwent an additional endoscopic intervention within 10 days of the initial procedure, 43 patients (12.4 %) within 100 days of the initial procedure. 12 patients (3.4 %) received abdominal surgery within 10 days after OTSC procedure, 41 patients (11.8 %) within 100 days of the procedure. Surgery after more than 30 days after OTSC application was mostly due to treatment of the underlying disease (carcinoma, diverticulitis etc.). Overall 9 patients (2.6 %) deceased within 100 days after the intervention. The authors concluded that, in the actual care setting, OTSC is mostly applied for polypectomies and iatrogenic perforations. The presented data supports first findings indicating that OTSC application in the colon is safe and helps to prevent surgery due to iatrogenic complications.

Sind OTS-Clips am Kolon effektiv und sicher? Evidenz-Generierung von endoskopischen Innovationen durch GKV-Routinedaten (Are OTS-Clips in the colon effective and safe? Evidence generation of endoscopic innovations with health insurance routine data).

Horenkamp-Sonntag D, Liebentraut J, Engel S, Knoop H, Hamburg bzw. Berlin, Deutschland.

OTSC as part of combination therapy of esophageal perforations and anastomotic insufficiencies following oncological resections

C Jung et al. presented a retrospective evaluation of all patients, that had been treated since 2014 at the University Hospital Goettingen for iatrogenic esophageal perforation (IEP) or post-surgical anastomotic insufficiency (PAI) with the EndoVac system, with esophageal stents and OTSCs. A total of 21 patients were recorded, 4 out of these with iatrogenic esophageal perforation and 17 with PAI. 12/17 PAI patients had received a preoperative radio/chemotherapy (5 CROSS, 1ICF, 1 FLOT+RT, 2 FLOT, 1 RTC, 1 GASTRIPEC, 1 unknown). Overall 8 patients received a fully-covered esophagus stent as primary therapy whereas 13 patients received an EndoVac as primary therapy. Complementary therapy was necessary in 6 patients (28.6 %) (2 stent + EndoVac, 1 EndoVac + Stent, 1 EndoVac + stent + fibrin, 1 stent + EndoVac + OTSC, 1 stent + OTSC). In overall 16/21 patients (76.2 %) complete restoration of the anastomosis was achieved. In 5 cases, continuity could not be restored,

2 of the patients died, 3 patients received a cervical drainage. The authors concluded that the group of patients examined was heterogenic and showed complex disease courses. The concept of combination therapy using EndoVac, esophageal stent, OTSC and endoscopic debridement seems to be promising. Further large scale studies are necessary to reliably describe the efficacy of this approach.

Multimodale endoskopische Behandlung von Ösophagusperforationen und postoperativen Anastomoseninsuffizienzen nach onkologischen Resektionen. Was ist die richtige Strategie? (Multimodal endoscopic treatment of esophageal perforation and post-surgical anastomosis insufficiency following oncological resection. Which is the correct strategy?)

Jung C, Kunsch S, Müller-Domieden A, Gaedcke J, Schüller P, Seif Amir Hosseini A, Ghadimi M, Ellenrieder V, Wedi E, Göttingen, Deutschland.

remOVE – registry study on clip removal and first multicenter case series evaluating the BougieCap

remOVE System – endoscopic removal of OTSC and FTRD clips is effective and safe

M Bauder and colleagues presented multicentric prospective registry data regarding application of the remOVE System. Data on 119 patients from 63 centers were submitted. Main indications for clip removal were: necessity of local re-therapy (62/119), local clip-associated complications (27/119), and ineffective clip placement (16/119). Cutting of the clip through both bows was successful in 89.1 % of cases, endoscopic retrieval of both clip fragments was possible in 82.4 %. Uncovering the clip from granulation tissue before application of the remOVE System was necessary in 23 cases. Average procedure time was 62 min, whereby a correlation to the thickness of the nitinol scaffold of the clip was seen (statistically significant between OTSC 11 and FTRD). Complications occurred in 3.4 % (4/119). These were in all cases minor bleedings, which could be managed endoscopically. The authors concluded that removal of OTSC and FTRD clips using the remOVE System is effective and safe.

Prospektive multizentrische Registerstudie zur Entfernung von OTSC und FTRD-Clips mit einem Gleichstrom-Schneideinstrument (Prospective multicenter registry study evaluating removal of OTSC and FTRD clips with a DC cutting instrument)

Bauder M, Wannhoff A, Meier B, Caca K, Ludwigsburg.

BougieCap – visually controlled dilatation with the BougieCap is effective and prevents complications due to overdistention

B Walter et al. presented a multicenter study (Ulm, Southampton and Essen) evaluating dilatation of benign esophageal stenoses with the BougieCap. The BougieCap allows, in contrast to Savary bougies, direct visual control of the process without the need for x-ray. 50 patients (25 f, 25 m, median age 67.1 ± 16.8) with benign stenosis of the esophagus and clinically apparent symptoms of dysphagia were included. Genesis of the stenosis was peptic (n=23), radiation (n=13), anastomotic (n=6), caustic ingestion (n=4), Post-ESD (n=2), EoE (n=1) and unknown (n=1). Dilatation was successful in 96 % of all cases (48/50). In eight cases a pediatric gastroscop with guidewire was used. In two cases a standard gastroscop with guidewire was used. In the two cases, passage of the stenosis was not possible, no attempt with guidewire had taken place. The number of subsequent endoscopic bougienages was median 2.3 ± 0.7. Dysphagia symptoms were regressive from a median DS value of 3.0 ± 0.6 before dilatation to 1.6 ± 0.7 after dilatation (Mann-Whitney, p < 0.0001). Severe complications did not occur. In two cases, a BougieCap was lost in the stomach; no clinical discomfort or complications resulted. The authors stated that endoscopic treatment of benign esophageal stenoses with the BougieCap allows direct visual control of the dilatation process and of beginning mucosal lacerations. Thus, in contrast to the conventional blind method, overdistention and re-traumatization are reduced and the dilatation process can be performed with better adaptation to the

stenosis. Usage of a guidewire is reasonable and necessary in special cases (i.e. very high-grade stenosis, usage of a pediatric gastroscop).

The BougieCap: a new method for endoscopic treatment of esophageal stenosis

Walter B, Schmidbaur S, Rahman I, Schumacher B, Albers D, Meining A.

RESECT+ – presented studies on first (pre-)clinical data confirm beneficial effort of the additional working channel (AWC®) and a new injection solution (LiftUp®) for optimized endoscopic resection

EMR+ – new endoscopic resection technique for en-bloc resection of lesions up to 30 mm

B Meier und K Caca presented a preclinical study evaluating a new endoscopic resection technique, which allows en-bloc resection of large lesions. The endoscopic mucosal resection (EMR) is regarded as standard procedure for endoscopic resection of mucosal intestinal neoplasms. However, when the lesion size surpasses 20 mm, the en-bloc resection rate is below 40 %. For lesions > 20 mm with urgent need of en-bloc/R0 resection, the endoscopic submucosal dissection is on hand, which albeit proves to be very sophisticated and time-consuming and associated to a higher complication rate. Submucosal injection is a crucial part of both techniques, with longer duration of the intervention, however, the so created cushion is more and more absorbed. The authors reported on a new procedure technique (EMR+). For the EMR+ technique, a standard endoscope with additional working channel (AWC) is used. Through the working channel of the endoscope, a grasping anchor is conducted, a resection snare through the additional working channel. The grasping anchor is led through the resection snare. For submucosal injection of the lesion, a new polymer injection solution (LiftUp) is used. After injection, the target lesion is lifted with the anchor, the snare is conducted over the lesion and the tissue below the lesion is cut. For the preclinical study, EMR+ was performed in two sessions with 11 resections each in ex-vivo porcine models. The en-bloc specimen had an average size of approx. 30 x 26 x 11 mm (maximal 40 x 33 x 14 mm). The overall procedure time was in average 6-7 minutes. Perforations did not occur. The authors rated the EMR+ as technically easy and fast technique for the en-bloc resection of lesions up to 30 mm in size.

Neue endoskopische Resektionstechnik zur en-bloc Resektion für Läsionen bis 30 mm (EMR+) (New endoscopic resection technique for lesions up to 30 mm (EMR+))

Meier B, Caca K, Ludwigsburg.

First clinical data shows that the AWC (additional working channel) makes en-bloc resection of large polyps in the upper and lower GI tract possible

B Walter and colleagues presented first clinical results with the AWC (additional working channel). The AWC, through which an additional instrument can be introduced to allow bi-manual working, can be fixed at the tip of a standard endoscope. In contrast to the double-channel endoscope, the distance between the working channels can be adjusted individually. In 4 patients, an ESD was performed and in 4 patients an EMR in a modified "grasp-and-snare-technique". For this technique, the target lesion is lifted by submucosal injection and subsequently with an OTSC anchor, which was previously conducted through a snare (EMR+) or the target lesion is held in tension with a grasper during the cutting process (ESD+). Lesions for ESD+ were located in the stomach (n=2, size 17 and 37 mm, respectively) and rectum (n=2, size 33 and 37 mm, respectively). Lesions for EMR+ were located in the stomach (n=1, size 31 mm) and colon (n=3, size 42 mm and two times 45 mm). En-bloc resection was successful in 6 cases, 1 EMR in the colon was achieved in two fragments and 1 EMR in the colon in 3 fragments. Median procedure time was 68.5 min. Complications occurred in the form of acute arterial bleeding directly after EMR in 2 cases. No perforations or latent haemorrhage occurred. The authors concluded, that the AWC makes en-bloc resection of large lesions possible. Advantages are the usability for ESD as well as EMR, the possibility of bi-manual working without 2-

channel endoscope and the individual adjustment of the channel distance.

Verwendung eines zusätzlichen externen Arbeitskanals (AWC) zur verbesserten endoskopischen Großflächenresektion (Usage of an additional external working channel (AWC) for improved endoscopic resection of large areas).

Walter B, Schmidbaur S, Hann A, Meining A, Ulm.

LiftUp – first preclinical data shows that this high-viscosity injection solution creates a stable and long-lasting cushion thereby increasing safety and efficacy of EMR and ESD

LiftUp is a polymer injection solution consisting of surface-active mass polymers, which is used for submucosal injection of early neoplasms in the course of endoscopic submucosal dissection (ESD). E Wedi et al. presented a prospective randomised study comparing LiftUp, NaCl 0.9 % and hydroxyl-ethyl-starch (HES 6 %) in an EASIE-R model. Overall 60 standardised ESD procedures were performed (n=20 per injection solution) in artificial lesions of 3 x 3 cm size in the corpus of a pig's stomach. ESD technique led to successful resection of all 60 lesions. R0 resection was achieved with LiftUp in 95 % (n=19), with HES in 100 % (n=20), and with NaCl in 80 % (n=16). Adequate mucosal lifting was reached in 80 % (n=16) with LiftUp, in 30 % (n=6) with HES and in 30 % (n=6) with NaCl (p<0.0002). Three perforations occurred, one in the HES-group and 2 in the NaCl-group. The authors rated the LiftUp injection solution a safe alternative for HES and NaCl. A particularly advantageous characteristic of LiftUp is the creation of a stable submucosal cushion, which retains for hours. Thereby, the ESD procedure as well as presumably the EMR procedure could become more safe and effective.

Evaluierung einer neuen submukosalen, hochviskösen Injektionslösung (LiftUp) für die Endoskopische Submukosa Dissektion (ESD) am EASIE-R Model: Eine prospektiv randomisierte Vergleichsstudie (Evaluation of a new high-viscosity solution (LiftUp) for submucosal injection in the course of endoscopic submucosal dissection (ESD) in the EASIE-R model: a prospective randomised comparison study).

Wedi E, Köhler P, Hochberger J, Dammer SS, Maiss J, Kunsch S, Ho N, Conrad G, Baulain U, Ellenrieder V, Jung C, Göttingen.

New demilune ESD-device (Coag Dissector) allows for rapid, effective and safe endoscopic submucosal dissection

Endoscopic submucosal dissection (ESD) has been established as an effective treatment option for early gastrointestinal cancer. To date, various devices for ESD are available. H Neumann and colleagues presented a prospective preclinical study evaluating the efficacy and learning curve of a new demilune device for ESD, which potentially allows for fast submucosal cutting above the muscular layer due to its special design. In addition, the device can be opened like scissors therefore also acting for hemostasis. The study was performed in two steps. First, ex vivo porcine models were utilized in an advanced endoscopic simulator or interventional endoscopy. After the initial learning curve, the study was repeated in living pigs under general anesthesia. For both study arms, artificial lesions, each 25 x 25 mm in size, were created in the fundus, corpus and antrum of the stomach. ESD was performed after marking of the lesions with the ESD instrument, followed by lifting of the mucosa with submucosal injection of colored saline. Afterwards, circular incision of the lesions was performed with the new ESD instrument. For resection, the submucosa was lifted with a distal clear cap and cut with the new demilune device. Resection specimen were retrieved to evaluate if all marks were included (R0). Average size of removed lesions was 30 mm. En-bloc resection rate was 100 % and R0 resection rate was 95 %. Mean total procedure time was 25 minutes and not dependent on the location or if the resection was performed in ex vivo models or in vivo. No perforations occurred during the study despite the rapid dissection speed through the submucosa. Satisfaction of the endoscopist and the supporting nurse staff was high throughout all cases. The authors concluded that the new

demilune device for ESD is safe and efficient and allows for rapid dissection of the submucosa due to its inherent design.

Preliminary report of a new demilune device for rapid endoscopic submucosal dissection (ESD).

Neumann H, Zimmermann T, Grimminger P, Rahman F, Thieringer F, Galle PR, Kneist W, Mainz.

For questions and further information:

Ovesco Endoscopy AG
Scientific Information Service
Dorfackerstraße 26
D-72074 Tübingen
science@ovesco.com

October 2018 | Successful application of OTSC® in GI bleeding under antithrombotic / anticoagulant therapy

100 % primary hemostasis rate and improved management of rebleeding with OTSC used as first line therapy

Lamberts R and colleagues, HELIOS Park-Hospital, Department II for Internal Medicine, Leipzig, Germany, conducted a retrospective study examining success rates in hemostasis of acute upper and lower gastrointestinal bleeding with the OTSC System as first or second line therapy in patients taking different regimens of antithrombotic and/or anticoagulant therapy. Overall 75 consecutive patients (mean age 71.7, 55 men, 20 women) with active gastrointestinal bleeding were analysed. 34 patients (45.3 %) were under antiplatelet monotherapy, 10 patients (13.3 %) under dual antiplatelet therapy, 13 patients (17.3 %) under inhibitors of plasminic coagulation, and 18 patients (24.0 %) no antithrombotic/anticoagulant therapy. OTSC was the first-line treatment in 45 (60 %) patients, in 30 patients (40 %) it was used in second line after preceding hemostasis attempts with conventional clips, adrenalin injection, fibrin glue and/or APC therapy. Key outcomes measured were: success rate with the OTSC therapy, rebleeding episodes, their management and the influence of antithrombotic or anticoagulant therapy.

Application of the OTSC resulted in immediate hemostasis in all 75 patients (100 % primary success rate). In 26 patients (34.7 %) a rebleeding episode was noted. In the group of first-line OTSC treatment the rebleeding rate was 28.9 % (13/45) compared to 43.3 % (13/30) in the group of second line OTSC treatment. In 23 patients rebleeding could be treated by further endoscopic interventions. Only 3 patients had to undergo radiological or surgical treatment because of final failure of endoscopic therapy attempts. In the rebleeding group the use of antiplatelet therapies was higher (73.1 % vs 48.9 %).

The authors concluded that primary OTSC application should be the treatment of choice in this high-risk patient population. Repeated endoscopic treatments to achieve definitive hemostasis may be justified and show promising results.

Use of over-the-scope clips (OTSC) for hemostasis in gastro-intestinal bleeding in patients under antithrombotic therapy.

Lamberts R, Koch A, Binner C, Zachaeus M, Knigge I, Bernhardt M, Halm U (2017).

Endoscopy International Open 2017; 05: E324-E330.

October 2018 | FTRD® offers endoscopic approach in the management of non-lifting and submucosal colorectal lesions and avoids surgical interventions

87.9 % R0-resection of non-lifting and submucosal colonic lesions with FTRD

Aeppli P et al., Gastroenterology and Hepatology Unit, Luzerner Kantonsspital, Lucerne, Switzerland, and Division of Gastroenterology/Hepatology, Kantonsspital St. Gallen, St. Gallen, Switzerland, reported on the clinical experience of the two tertiary referral centers with FTRD procedures. 33 consecutive patients with colonic neoplasms (21 colon, 12 rectum) were scheduled for EFTR using the FTRD device between 05/2015 and 11/2016. Indications were

residual adenoma with non-lifting sign after previous polypectomy (n=18), non-lifting adenoma without previous polypectomy (n=4), staging following resection of a malignant polyp (n=4), adenoma at appendiceal orifice (n=2), primary EFTR of polyps suspected to be malignant (n=2), adenoma involving a diverticulum (n=1), non-lifting adenoma recurrence after EFTR (n=1), and incomplete resection of neuroendocrine tumor G1 (n=1).

31 resections were successfully performed. In one case the target lesion could not be reached because of sigmoid stenosis due to diverticulosis, the other failure was due to snare malfunction.

Resection was en bloc and histologically complete (R0) in 87.9 % (29/33) of patients. The mean diameter of resected specimen as assessed by the pathologist was 2.7 cm (range 18-43 mm). Three post-procedure bleedings and one perforation were seen.

The authors concluded that the FTRD System offers an additional endoscopic approach in the management of non-lifting colorectal lesions and helps to avoid surgical interventions.

Endoscopic full thickness resection (EFTR) of colorectal neoplasms with the Full Thickness Resection Device (FTRD): Clinical experience from two tertiary referral centers in Switzerland.

Aeppli P, Cribblez D, Baumeler S, Borovicka J, Frei R (2017.) United European Gastroenterology Journal 0(0) 1-8.

October 2018 | OTSC® safe and effective for treatment of leak at the tip of the “J” ileal pouch

66.6 % of patients with leaking from the “J” of the tip of an ileo-anal pouch anastomosis were spared surgery by use of the OTSC

The tip of the “J” ileal pouch is the vulnerable location for leak after restorative proctocolectomy, which has normally been treated with surgery.

Lian L and Shen B, Interventional IBD Center, Digestive Disease and Surgery Institute, Cleveland Clinic, Cleveland, OH, USA, described the first case of endoscopic treatment of the leak at the tip of the “J” with OTSC in 2014. Since then, OTSC therapy has become the first-line approach for this lesion in the Cleveland Clinic.

Recently, Kochhar GS and Shen B, same affiliation, published a cohort study comprising 12 consecutive patients with a leak at the tip of the “J” from the Center's prospectively maintained Pouch Registry. In all patients, OTSC was used for leak closure.

All 12 patients had successful deployment of OTSC during endoscopy. No excessive bleeding or perforation was observed. Eight patients (66.6 %) achieved complete closure of the leak documented by endoscopy confirmed with guidewire and/or contrasted pouchogram, with 6 requiring a single endoscopic session and 2 undergoing a repeat session. Four patients (33.3 %) had a persistent leak and required surgical intervention, of whom 1 developed abscess in the pre-spine region 14 days after the endoscopic procedure and underwent pouch revision surgery.

The authors concluded that leaking from the tip of the “J” in patients with ileo-anal pouch anastomosis can be effectively and safely treated with the over-the-scope clipping system.

Endoscopic treatment of leak at the tip of the “J” ileal pouch.

Kochhar GS, Shen B (2017).

Endoscopy International Open 2017; 05: E64-E66.

September 2018 | Clinical experience with the remOVE System, a bipolar cutting device for OTSC® removal

Bauder M et al., Department of Gastroenterology and Oncology, Ludwigsburg Hospital, Ludwigsburg, Germany, presented a clinical study evaluating OTSC removal with the remOVE System.

OTSC (or FTRD) removal may be indicated when OTSC/FTRD-associated complications occur, when the clip was misplaced, re-therapy after full-thickness resection is needed or an OTSC-fixed stent has to be removed.

The remOVE System is a bipolar grasping device with

which short direct current (DC) impulses can be applied to cut the OTSC at two opposing sites. The DC impulses are delivered by a special electric generator connected to the grasping device. The bipolar grasper can be advanced through a 2.8 mm endoscope working channel. Its tip consists of three electrodes that are brought in contact with the thinnest parts of the nitinol clip. Application of DC impulses then selectively heat up and melt the nitinol. As soon as the clip is cut or the contact to the nitinol is lost during the cutting process, an integrated safety feature automatically stops the current flow.

In the next step, OTSC fragments are extracted with a standard forceps. For extraction, a plastic cap at the tip of the endoscope is used to avoid tissue damage.

Data of all consecutive patients with indication for OTSC removal were collected and analysed retrospectively. Between 12/2012 and 02/2016, a total of 42 OTSC removals in the upper (n=25) and lower (n=17) gastrointestinal tract have been performed. Overall technical success, defined as cutting the OTSC at two opposing sites and extraction of both fragments, was achieved in 92.9 % (39/42) of all cases. Successful fragmentation of the OTSC was achieved in 97.6 % (41/42). Minor bleedings were rare and could be managed endoscopically in all cases. There were no perforations and no major or delayed bleedings.

The authors concluded that endoscopic OTSC removal with the remOVE System is feasible, safe and effective. The technique can be applied in the upper and lower gastrointestinal tract.

Endoscopic removal of over-the-scope clips: Clinical experience with a bipolar cutting device.

Bauder M, Meier B, Caca K, Schmid A (2017).

United European Gastroenterol J 2017 Jun; 5(4):479-484.

August 2018 | Systematic review: leaks and fistulae after laparoscopic sleeve gastrectomy successfully closed by OTSC® in 86.3 %

A systematic literature review evaluating efficacy and safety of the OTSC System in the management of leak and fistula after laparoscopic sleeve gastrectomy (LSG) was published by Shoar S et al., Bariatric and Metabolic Institute, department of surgery, The Brooklyn Hospital Center, Icahn School of Medicine at Mount Sinai, Brooklyn, NY, USA.

A total of 10 eligible studies including 195 patients with post-LSG leaks/fistulae were identified. The time interval between LSG and leak/fistula ranged from 1 day to 803 days. Most of the leaks/fistulae were located at the proximal staple line, and had a size from 3 to 20 mm. Time between leak diagnosis and OTSC clipping ranged from 0 to 271 days.

Details for endoscopic management of post-LSG leak/fistula by OTSC were available for nine studies (73 patients). Of the 73 patients with post-LSG leak treated with OTSC, 63 patients had an overall successful closure (86.3 %). Number of the deployed OTSC was reported by six studies (53 patients). Of these, 33/53 patients (63.5 %) required one clip for closure of the lesion, 14 patients (36.9 %) required one or more clips and 5 patients (9.6 %) required two clips.

Regarding OTSC-related complications, OTSC migration was reported in one patient (1.4 %), stenosis in one patient (1.4 %), and tear in one patient (1.4 %).

In conclusion, the authors stated that the OTSC System is a promising endoscopic approach for management of post-LSG leaks in appropriately selected patients. Unfortunately, most studies are series with a small sample size, short-term follow-up, and mixed data of concomitant procedures with OTSC. Further studies should distinguish the net efficacy of the OTSC system from other concomitant procedures in treatment of post-LSG leak.

Efficacy and Safety of the Over-the-Scope Clip (OTSC) System in the Management of Leak and Fistula After Laparoscopic Sleeve Gastrectomy: a Systematic Review.

Shoar S, Poliakin L, Khorgami Z, Rubenstein R, El-Matbouly M, Levin JL, Saber AA (2017).

Obes Surg 2017 Sep;27(9):2410-2418.

August 2018 | Multicenter experience comparing simple suction and OTSC® Twin Grasper®

Kobara H and colleagues, Department of Gastroenterology and Neurology, Faculty of Medicine, Kagawa University, Kagawa, Japan performed a retrospective study examining 58 consecutive patients undergoing OTSC placement for gastrointestinal defect closure in 5 medical centers.

The overall rates of technical success, clinical success, complications and procedure time were analysed as major outcomes. Subsequently, 56 patients, excluding 2 cases that used the Anchor device, were divided into two groups: 14 cases of simple suction (SS-group) and 42 cases using the OTSC Twin Grasper (TG-group). Secondary evaluation was performed to clarify the predictors of OTSC success.

Overall clinical outcomes demonstrated efficacy and safety of the OTSC System and were as follows: technical success rate (TSR) 89.7 %, clinical success rate (CSR) 84.5 %, complications 1.8 %, and median procedure time 8 minutes (range 1-36 min).

Significant differences were observed between the two groups in terms of the mean procedure time (5.9 min vs 14.1 min). The clinical success rate of the SS- and TG-groups among cases with a maximum defect size ≤ 10 mm and immediate or acute refractory bleeding was 100 %, which suggests that SS is a better method than TG in terms of time efficacy. The clinical success rate in the SS-group (78.6 %), despite the technical success of the SS method (100 %), tended to decrease due to delayed leakage compared to that in the TG group (TSR 88.1 %, CSR 88.1 %), indicating that the OTSC Twin Grasper may be desirable for leaks and fistulae with defects of the entire layer.

The authors concluded, that the OTSC System is a safe and effective therapeutic option for gastrointestinal defects. Individualized selection of the suction method based on particular clinical conditions may contribute to the improvement of OTSC success.

Outcomes of gastrointestinal defect closure with an over-the-scope clip system in a multicenter experience: An analysis of a successful suction method.

Kobara H, Mori H, Fujihara S, Nishiyama N, Chiyo T, Yamada T, Fujiwara M, Okano K, Suzuki Y, Murota M, Ikeda Y, Oryu M, AboEllail A, Masaki T (2017).

World J Gastroenterol 2017 March 7; 23(9): 1645-1656.

OTSC® Update 28

July 2018 | Closure of acute GI defects with OTSC® successful in more than 75 % of patients of an unselected cohort

Raithel M and colleagues, Department of Medicine II, University Erlangen, Germany presented data from three tertiary referral centers with 24-h emergency endoscopy (Erlangen, Wuerzburg, Fuerth) on patients receiving OTSC therapy for acute gastrointestinal wall defects.

Unselected consecutive patients presenting with acute non-surgical perforations, postoperative anastomotic leaks or inadvertent postoperative perforations underwent attempted OTSC placement as primary closure method after interdisciplinary consensus. Their clinical data and intervention characteristics were evaluated in an intention to treat analysis during a 24-month period to assess closure rates, 30-day mortality, hospitalization and comorbidity.

In total, 34 patients were included with 22 non-surgical perforations and 12 postoperative leaks or perforations. 5 GI defects were located in the oesophagus, 14 in the stomach, 4 in the duodenum, 2 in the jejunum, 1 in the ileum, 2 in the colon, and 6 in the rectum. Definitive closure of the perforations and leaks was achieved in 26 patients (76.5 %). The closure rate among non-surgically caused perforations was 72.7 % and among acute postoperative GI wall defects 83.3 % in total. Successful closure of the GI wall defect resulted in a significantly shorter hospital stay (8 vs 18 days, p = 0.03). In the group with OTSC failure, 6 of 8 patients (75 %) required immediate surgery. In the group with successful OTSC closure, comorbidity rate was significantly higher (19/26 patients vs 4/8 patients in the group with OTSC failure; p = 0.005). Three deaths occurred in the group with successful OTSC closure due to

comorbidity, while one death in the OTSC failure group was related to a refractory perforation. Favourable indications and locations for a successful OTSC procedure were identified as PEG complications, and endoscopic or postoperative leaks of stomach, colon or rectum.

The authors concluded, that OTSC was effective for closure of acute GI wall defects in more than 75 % of patients in an unselected cohort.

Outcome, comorbidity, hospitalization and 30-day mortality after closure of acute perforations and postoperative anastomotic leaks by the over-the-scope clip (OTSC) in an unselected cohort of patients.

Raithel M, Albrecht H, Scheppach W, Farnbacher M, Haupt W, Hagel AF, Schellner V, Vitali F, Neurath MF, Schneider. HT (2016) Surg Endosc 2017 Jun; 31 (6): 2411-2425.

July 2018 | Single center experience: OTSC® especially valuable in treating defects of the upper gastrointestinal tract

Mizrahi I and colleagues, Digestive Disease Institute, Cleveland Clinic Florida, Weston, FL, USA, reported their institutional experience with OTSC for patients with GI defects.

Gastrointestinal tract defects are associated with high patient morbidity and pose a clinical and technical challenge for surgeons and gastroenterologists. Defects such as anastomotic leak and fistula are reported in up to 13 % of patients following upper GI tract surgery for weight loss. Even higher rates are reported after colorectal surgery for malignant and inflammatory indications. The conventional treatment of the defects mentioned above is surgery, often including stoma creation or percutaneous drainage with the obvious related morbidity.

For the study, prospective data from all patients treated with OTSC in the institution were analyzed. Primary outcome was the clinical success of the OTSC for the individual indication.

During the study period, 51 patients (28 females, mean age 54.9 ± 19.4 years) were treated with OTSC, 21 patients had defects in the lower GI tract (LGI) and 30 patients in the upper GI tract (UGI). The most common indication for OTSC was anastomotic leak (n=24; UGI = 12, LGI = 12), followed by chronic fistula (n=17; UGI = 8, LGI = 9), acute perforation (n=4), acute bleeding (n=4) and stent anchoring (n=2). OTSC was the primary therapy in all patients with bleeding and stent anchoring, in all other patients alternative endoscopic or surgical attempts had been made prior to OTSC.

Technical success was achieved in 98 % of all patients. Clinical success rates for UGI perforation, bleeding, and stent anchoring indications were 75, 75, and 50 %, respectively. Clinical success rate for the treatment of anastomotic leaks was 59 % (UGI 66 % vs LGI 33 %). A lower success rate was noted for chronic fistulae (UGI 62 % vs LGI 0 %).

The authors conclude that the OTSC is a reasonable minimally invasive option for managing patients with various GI defects before a more invasive operative approach is attempted. It appears that the OTSC is especially valuable in the treatment of acute defects of the upper gastrointestinal tract.

The Clinical Utility of Over-the-Scope Clip for the Treatment of Gastrointestinal Defects.

Mizrahi I, Ellawil R, Haim N, Chadi SA, Shen B, Erim T, DaSilva G, Wexner SD (2016).

J Gastrointest Surg 2016; 20: 1942-49.

July 2018 | Large single center experience presented: Establishment of the OTSC® clip in daily endoscopic routine

Honegger C and colleagues, Division of Gastroenterology and Hepatology, University Hospital Zurich, Zurich, Switzerland, presented data on 262 OTSC placements in a total of 233 interventions. Since 2009, the placement of OTSC has been established at the University Hospital Zurich for the entire spectrum of indications. OTSC has become a device of daily practice. A retrospective study now presents data of all patients treated with the OTSC device at the institution, focussing on indications, anatomic site of OTSC deployment, complications, and immediate

and 30-day success rates.

Patient age ranged from 14 to 93 years with a median of 61 years. 51.5 % were male. Immediate success of OTSC treatment was observed in 87.1 % of all sessions (203/233). The success rates per indication were as follows: spontaneous bleeding 84.8 % (28/33); iatrogenic bleeding 100 % (20/20); acute perforation 90.3 % (65/72); prophylaxis for perforation 100 % (24/24); anastomotic leakage 61.1 % (11/18); fistulae 80.7 % (46/57); diameter reduction of the gastro-jejunal anastomosis 100 % (6/6); and stent fixation 100 % (3/3).

At 30-day follow-up, the overall success rate was 67.4 % (157/233). The success rates per indication were as follows: spontaneous bleeding 69.7 % (23/33); iatrogenic bleeding 90 % (18/20); acute perforation 86.1 % (62/72); prophylaxis for perforation 100 % (24/24); anastomotic leakage 33.3 % (6/18); fistulae 29.8 % (17/57); diameter reduction of the gastro-jejunal anastomosis 83.3 % (5/6); and stent fixation 66 % (2/3).

The authors concluded that the treatment with an OTSC is safe and feasible in clinical routine, with high immediate success rates with sustained clinical success at 30-day follow-up.

Establishment of Over-The-Scope-Clips (OTSC) in final endoscopic routine.

Honegger C, Valli PV, Wiegand N, Bauerfeind P, Gubler C (2016).

United European Gastroenterol J. 2017 Mar;5(2):247-254.

June 2018 | 100 % long-term success with OTSC® in acute GI defects, 73 % in chronic defects

Gianfranco G. et al., Interventional Endoscopy Unit, Hospital Privé des Peupliers, Paris, France, presented a retrospective study from a prospectively-maintained database evaluating immediate and long-term success rates of OTSC deployment in acute and chronic gastrointestinal pathologies.

Between 01/2012 and 12/2015 a total of 51 OTSCs were delivered in 45 patients (35 female, average age 56 years, range 24-90) due to GI defects resulting from a diagnostic or interventional endoscopic procedure (acute setting group; n=15) or due to fistula following abdominal surgery (chronic setting group; n=30). All procedures were carried out in a private endoscopic service.

Technical success was always achieved in the acute setting group with an excellent clip adherence and a clinical long-term success rate of 100 % (15/15, median follow-up 9 months, range 1-24 months).

Considering the chronic setting group (OTSC treatment after an average period of 146.6 days (range 5-880 days) after primary surgery), technical success was achieved in 50 % of patients (15/30). Long-term clinical success in patients with succeeded primary fistula closure by OTSC was 73.3 % (11/15; mean follow-up 23 months, range 1-34 months). Two minor complications occurred. A total of three patients died due to causes not directly related to clip deployment.

The authors concluded that OTSC deployment is an effective and minimally-invasive procedure for GI defects in acute settings. It avoids emergency surgical repair and it allows, in most cases, completion of the primary endoscopic procedure. OTSC should be incorporated as an essential technique of today's modern endoscopic armamentarium in the management of GI defects in acute settings. OTSCs were less effective in cases of chronic defects.

Closure of gastrointestinal defects with Ovesco clip: long-term results and clinical implications.

Donatelli G, Cereatti F, Dhumane P, Vergeau BM, Tuszyński T, Marie C, Dumont JL, Meduri B (2016). Therap Adv Gastroenterol. 2016 Sep;9(5):713-21.

June 2018 | Breaking news: Ovesco OTSC® Clip superior to standard hemostatic therapy in randomized-controlled trial

OTSC has long been described in the scientific literature as a highly effective device for the treatment of upper GI hemorrhage. Now a randomized-controlled trial at 9 academic referral centers (in Germany, Switzerland, and

Hong Kong) has proven OTSC to be superior to standard methods. The trial, published by Dr. Arthur Schmidt, Ludwigsburg, Germany, enrolled 66 patients with recurrent bleeding and randomized them to receive either OTSC therapy or standard techniques (a combination of 2 methods from through the scope clipping, injection or electrical coagulation).

Persistent bleeding after per-protocol hemostasis was observed in 42.4 % of patients in the standard therapy group and 6.0 % in the OTSC group (P=.001). Further bleeding occurred in 57.6 % in the standard therapy group and 15.2 % in the OTSC group (absolute difference, 42.4 %; 95 % CI 21.6-63.2; P=.001).

Over the Scope Clips are More Effective Than Standard Endoscopic Therapy for Patients With Recurrent Bleeding of Peptic Ulcers.

Schmidt A¹, Gölder S², Goetz M³, Meining A⁴, Lau J⁵, von Delius S⁶, Escher M⁷, Hoffmann A⁸, Wiest R⁹, Messmann H², Kratt T², Walter B⁴, Bettinger D¹⁰, Caca K¹¹.

Author information

¹Department of Gastroenterology, Klinikum Ludwigsburg, Ludwigsburg, Germany; ²Department of Medicine II, Medical Center, Faculty of Medicine, University of Freiburg, Germany.

³Department of Gastroenterology, Klinikum Augsburg, Augsburg, Germany.

⁴Interdisciplinary Endoscopy, University of Tübingen, Tübingen, Germany.

⁵Department of Gastroenterology, University of Ulm, Ulm, Germany.

⁶Department of Surgery, University of Hong Kong.

⁷Department of Gastroenterology, Klinikum Rechts der Isar, TU München, München, Germany.

⁸Department of Gastroenterology, Robert Bosch Krankenhaus Stuttgart, Stuttgart, Germany.

⁹Department of Gastroenterology, Horst Schmidt Kliniken Wiesbaden, Wiesbaden, Germany.

¹⁰Department of Gastroenterology, Inselspital Bern, Bern, Switzerland.

¹¹Department of Medicine II, Medical Center, Faculty of Medicine, University of Freiburg, Germany; Berta-Ottenstein-Programme, Faculty of Medicine, University of Freiburg.

¹²Department of Gastroenterology, Klinikum Ludwigsburg, Ludwigsburg, Germany.

Gastroenterology. 2018 May 24. pii: S0016-5085(18)34570-0. doi: 10.1053/j.gastro.2018.05.037. [Epub ahead of print].

May 2018 | Single-center study confirms safety and efficacy of the FTRD® in the colorectum

Valli and colleagues, Division of Gastroenterology and Hepatology, Zurich University Hospital, Switzerland, reported about their data on the colonic FTRD® between June 2012 and October 2016. Full-thickness resection in the colorectum (52) and off-label use in the upper GI (8) were performed in 60 patients with the following indications: recurrent and primary non-lifting adenomas, combined procedure of EMR and EFTR in large polyps (EFTR for non-lifting area), primary and re-resection of T1 tumors not suitable for surgical resection, submucosal lesions and adenomas at difficult locations (appendix base, diverticulum).

The overall technical success rate was 97 % (58/60), the overall RO and full-thickness resection rates were 79 % and 88 % respectively. Clinical success was shown to be even better (88 %) based on follow-up histology. In median a procedure took 60 minutes with a mean resection size of 24 mm. The adverse event rate was 7 %. One patient (2 %) developed appendicitis after resection at the appendix base and needed surgical appendectomy. All other complications (minor bleeding, perforation due to accidental lack of clip deployment) were treated endoscopically. The data corresponds with the results of the Wall Resect trial.

The authors conclude that EFTR with the colonic FTRD is safe and feasible with respective prior training. It enables endoscopic resection of all gut layers with low risk of severe adverse events. EFTR with the FTRD shows also an

alternative procedure to surgery for lesions that were previously not endoscopically resectable.

Safe and successful resection of difficult GI lesions using a novel single-step full-thickness resection device (FTRD®).

Valli P.V., Mertens J., Bauerfeind P.

Surg Endosc. 2017 Jun 29. doi: 10.1007/s00464-017-5676-9.

May 2018 | Conference Report

The 48th Conference of the German Society for Endoscopy and Imaging Procedures (DGE-BV) took place together with the learned societies of CAES, CATC, DEGEA, DEGUM, DGBMT, DGD, ÖGGH and bng on March 15 – 17, 2018 in Munich, Germany.

Ovesco products were presented in six workshops on two different topics (hemostasis techniques held by M. Mühleck and S. Loeffler, respectively, and management of complications held by C. Hamperl and S. Loeffler, respectively). Additionally, several talks and posters discussed products of Ovesco.

FTRD® System

Diagnosis of amyloidosis with FTRD full-thickness rectal tissue sampling

A. Braun and H. Dawson, Gastroenterology and Endoscopy, SRO Langenthal, Switzerland and Institute of Pathology, University of Bern, Switzerland, presented their study on amyloidosis diagnosis with the FTRD System. This is the first description of FTRD use in this indication.

Amyloidosis is a heterogeneous group of diseases with accumulation of abnormal protein, known as amyloid fibrils, which build up in interstitial tissue, leading to manifold clinical problems. The GI tract is also affected. Diagnosis must be confirmed by biopsy and histological examination and samples must contain submucosal vessels and muscularis propria. There is currently no reliable minimally invasive sampling technique. The study investigated feasibility, performance and safety of endoscopic full-thickness resection with the FTRD System in the rectum for gastrointestinal amyloidosis diagnosis.

Between 2015 and 2017, full-thickness excision of rectal wall with the FTRD System was performed in 12 patients (5 female, median age 73 years (29-81)) with suspected amyloidosis. Sigmoidoscopy was performed under light sedation in all patients. Biopsies were taken from the upper third of the rectum 14 – 18 cm ab ano. FTRD application and full-thickness resection were successful in all cases. Maximal procedure time was 20 minutes. No adverse events occurred and the clinical course was uneventful in all cases. Clear diagnosis was possible in all histological examinations. In 7 of the 12 patients (m=5, f=2) amyloid fibril accumulation was found in small submucosal vessels and in the muscularis propria layer, confirming the diagnosis of amyloidosis.

The authors conclude that FTRD application for diagnosis of gastrointestinal amyloidosis is a safe and very effective method and mitigates the diagnostic challenges that amyloidosis can pose.

Amyloidose-Diagnostik durch Full-Thickness Resection Device (FTRD).

Braun A, Langenthal and Dawson H, Bern.

Neoplasms, which were not accessible for EMR, in over 90 % resected with the FTRD System

A. von Helden and colleagues, Community Hospital Bonn, presented a case series on endoscopic full-thickness resection (EFTR) with the FTRD System in colorectal neoplasms with difficult localisation or extensive fibrosis. In 27 patients minor adenomas were found, which were not accessible for endoscopic mucosal resection (EMR) due to localisation (5 x appendix, 1 x diverticulum) or extensive fibrosis in consequence of prior attempts of endoscopic resection or chronic inflammatory bowel disease. 74 % of the reported lesions were located in the right colon, 15 % in the left colon, and 11 % in the rectum. All interventions were performed with the FTRD System on standard colonoscope with peri-interventional antibiotics (2g Ceftriaxon) and under analgo-sedation. Technical success was achieved in 25/27 patients (92.59 %), the target lesion could not be

reached with FTRD cap in two cases due to distal stenosis. The size of the resected specimen was median 27 mm (12 – 33 mm). Histopathological examination showed 2 carcinomas (8%; 1 x curative, 1 x elective oncological resection due to high-risk histological findings), 6 high-grade intra-epithelial neoplasms (24%), 11 low-grade intra-epithelial neoplasms (44%) and 6 serrated adenomas (24%). R0-resection-rate was 76% (19/25). Two perforations occurred (8%), one was treated by emergency surgery. In one patient with native appendix, appendicitis developed after two days and required ileocecal resection. Relevant haemorrhage was not observed. 30-days mortality rate was 0%.

In summary this study confirms that application of the FTRD System makes endoscopic resection of neoplasms with difficult localisation and extensive fibrosis possible. Technical problems arise from stenosis distal to the target lesion. The authors warn against using the device in cases with native appendix.

Endoskopische Vollwandresektion mit dem FTRD System: Effektivität der Methode bei 27 Patienten mit kolorektalen Neoplasien, die einer Mukosaresektion nicht zugänglich waren. Von Helden A, Sido B, Hildenbrand R, Dumoulin FL, Bonn.

FTRD application in adenomas with non-lifting sign, submucosal localisation or progressed histology: analysis from 3 hospitals

H. Albrecht and colleagues presented an analysis of data from patients from 3 hospitals, who had been admitted for EFTR due to adenomas in the lower gastrointestinal tract. The respective adenomas showed non-lifting sign, submucosal localisation or suspected (pre-) malign histology. Aim of the study was to describe the histological findings of the resected lesions, the resulting proceeding and the effects on the clinical course of the patient.

Between 11/2014 and 02/2017, a total of 55 patients were admitted to the three centers for FTRD application due to above-mentioned indication. Two of the adenomas initially deemed suitable for FTRD therapy presented endoscopically too big and were subsequently treated by primary surgery. Four lesions could not be resected because of inaccessible localisation or lacking retrieval. The other 49 lesions were diagnosed as relapse adenoma or adenoma with non-lifting sign (n=21), high-grade intra-epithelial neoplasm and/or intra-mucosal adenocarcinoma (n=21), submucosal lesion / NET (n=6) and metastasis of malign melanoma (n=1). The resected specimen had a median size of 2.5 ± 2 cm. R0-resection was achieved in 38/49 (77.6%), 11 lesions were incompletely resected (R1 or R2). In 8 cases (8/49 = 16.3%), surgical revision was necessary. In three cases (3/49 = 6.1%) the full-thickness specimen showed early infiltration of lymphatic vessels, these patients also underwent surgery. The following complications occurred (n=4): one Hb-relevant haemorrhage, which could be managed endoscopically, one perforation, which could be managed without surgery by OTSC application, in one case not enough tissue could be drawn into the cap, and in one case resection was incomplete because of snare dislocation during resection. In summary these data show that surgery could be avoided by FTRD application in more than two thirds of the patients. The authors recommend, however, that the indication for EFTR for big lesions (>4 cm) and in case of evidence of mucosal carcinoma in biopsies should be rather strict.

Endoskopische transmurale Vollwandresektion (EFTR) am unteren Gastrointestinaltrakt (GIT): Welche Patienten profitieren?

Albrecht H, Neumarkt i.d.OPf., Raithel M, Nagel A und Braun A, Erlangen, Stegmaier A, Schwabach, und Schaefer C, Neumarkt i.d.OPf.

Successful eFTR at the appendiceal basis with the FTRD System

C. Schaefer and colleagues reported on a 66-year old female patient, who was admitted for EMR attempt of a coecal polyp at the appendiceal basis. Colonoscopy showed the aforesaid adenomatous polyp located at the appendiceal cavity of the caecum. Histological examination yielded the diagnosis of a tubular adenoma without

evidence of dysplasia. Submucosal injection and following resection by EMR was not successful due to the location of the polyp. A second attempt of resection by eFTR with the FTRD System was performed and succeeded. The resected specimen was sized 3x1x1 cm. Histological examination showed a partial appendix, at the appendiceal basis a 9 x 8 mm sized polyp with R0 resection status. The patient received peri-interventional antibiotic prophylaxis with Cefuroxim and Metronidazol. The post-interventional clinical course was uneventful, step-wise return to normal diet was tolerated without problems and the patient could be discharged after two days.

Erfolgreiche endoskopische Appendektomie mittels EFTR.

Schaefer C, Michalek M, Albrecht H, Hemmel M, Reitingner S, Neumarkt i.d.OPf., und Evert M, Universität Regensburg.

OTSC® System

100 % hemostasis with OTSC reported for first-line emergency treatment of acute hemorrhage

A. Braun and S. Peter, SRO Langenthal, Switzerland, presented personal data on OTSC application for emergency treatment of acute hemorrhage.

Between 2011 and 2017, 48 patients (29 female, median age 75.5 years (61-92)) each received one OTSC clip for first-line treatment of acute gastrointestinal hemorrhage. Bleeding was located in the upper GI tract in 34 cases (14 Forrest Ia, 15 Forrest Ib, and 5 Forrest IIa) and in the lower GI tract in 14 cases (4 Forrest Ia, 7 Forrest Ib and 3 Forrest IIa). Patients with upper GI bleeding received peri-interventional PPI medication (80mg i.v. bolus, 320 mg i.v. / 24 h).

OTSC application and primary hemostasis were successful in all cases. Maximal procedure time was 20 minutes. No relapse hemorrhage occurred. 26 patients (15 F Ia, 9 F Ib) received follow-up endoscopy on day 1 to 4, which showed the clip in situ and no bleeding stigmata. The other 22 patients received no follow-up examination. The clinical course was uneventful in all cases.

The authors concluded that OTSC application for emergency endoscopic treatment of acute hemorrhage is safe and very effective, and related to short procedural time.

Endoskopische Behandlung von akuten Blutungen mit einem over-the-scope clip (OTSC).

Braun A, Peter S, Langenthal, Schweiz.

April 2018 | Diagnosis of amyloidosis with FTRD® full-thickness rectal tissue sampling

A. Braun and H. Dawson, Gastroenterology and Endoscopy, SRO Langenthal, Switzerland and Institute of Pathology, University of Bern, Switzerland, presented their study on amyloidosis diagnosis with the FTRD System at the 48th DGE-BV Conference (DGE-BV: German Society for Endoscopy and Imaging Procedures) in March 2018. This is the first description of FTRD use in this indication. Amyloidosis is a heterogeneous group of diseases with accumulation of abnormal protein, known as amyloid fibrils, which build up in interstitial tissue, leading to manifold clinical problems. The GI tract is also affected. Diagnosis must be confirmed by biopsy and histological examination and samples must contain submucosal vessels and muscularis propria. There is currently no reliable minimally invasive sampling technique. The study investigated feasibility, performance and safety of endoscopic full-thickness resection with the FTRD System in the rectum for gastrointestinal amyloidosis diagnosis.

Between 2015 and 2017, full-thickness excision of rectal wall with the FTRD System was performed in 12 patients (5 female, median age 73 years (29-81)) with suspected amyloidosis. Sigmoidoscopy was performed under light sedation in all patients. Biopsies were taken from the upper third of the rectum 14 – 18 cm ab ano.

FTRD application and full-thickness resection were successful in all cases. Maximal procedure time was 20 minutes. No adverse events occurred and the clinical course was uneventful in all cases. Clear diagnosis was possible in all histological examinations. In 7 of the 12 patients (m=5, f=2) amyloid fibril accumulation was found in small submucosal vessels and in the muscularis propria layer, confirming the diagnosis of amyloidosis.

The authors conclude that FTRD application for diagnosis of gastrointestinal amyloidosis is a safe and very effective method and mitigates the diagnostic challenges that amyloidosis can pose.

Amyloidose-Diagnostik durch FTRD (Full-Thickness Resection Device).

Braun A and Dawson H (2018).

April 2018 | 100 % hemostasis with OTSC® reported for first-line emergency treatment of acute hemorrhage

The 48th Conference of the German Society for Endoscopy and Imaging Procedures (Deutsche Gesellschaft für Endoskopie und Bildgebende Verfahren, DGE-BV) took place on March 15-17, 2018 in Munich, Germany. Dr. A. Braun, SRO Langenthal, Switzerland, presented personal data on OTSC application for emergency treatment of acute hemorrhage.

Between 2011 and 2017, 48 patients (29 female, median age 75.5 years (61-92)) each received one OTSC clip for first-line treatment of acute gastrointestinal hemorrhage. All patients had shown acute hemoglobin decrease and secure bleeding signs such as hematemesis, melena or hematochezia. Bleeding was located in the upper GI tract in 34 cases (14 Forrest Ia, 15 Forrest Ib, and 5 Forrest IIa) and in the lower GI tract in 14 cases (4 Forrest Ia, 7 Forrest Ib and 3 Forrest IIa). Patients with upper GI bleeding received peri-interventional PPI medication (80mg i.v. bolus, 320 mg i.v. / 24 h). For placement of the OTSC in some cases an endoscopic forceps was used to grasp tissue. No further local therapies were applied. All OTSC applications were performed by one single endoscopist.

OTSC application and primary hemostasis were successful in all cases. Maximal procedure time was 20 minutes. No relapse hemorrhage occurred. 26 patients (15 F Ia, 9 F Ib) received follow-up endoscopy on day 1 to 4, which showed the clip in situ and no bleeding stigmata. The other 22 patients received no follow-up examination. The clinical course was uneventful in all cases.

The author concludes that OTSC application for emergency endoscopic treatment of acute hemorrhage is safe and very effective. Primary hemostasis is achieved in a large fraction of patients, which makes reduction of the mortality rate of acute gastrointestinal hemorrhage possible. OTSC application is related to short procedural time.

Endoskopische Behandlung von akuten Blutungen mit einem over-the-scope clip (OTSC).

Braun A (2018).

March 2018 | One step application of OTSC® for salvage hemostasis and simultaneous perforation closure

EI Douaihy Y et al., Department of Internal Medicine, Staten Island University Hospital - Northwell Health System, Staten Island, New York, USA, reported on a case of active bleeding from a gastroduodenal artery pseudoaneurysm and simultaneous perforation, which was treated by deployment of a single OTSC-Clip.

The 61-year-old male patient with history of duodenal ulcer and angiographic embolization of a gastroduodenal artery pseudoaneurysm 6 months before presented to the emergency room for bright red blood per rectum with signs of upper GI bleed. In esophagogastroduodenoscopy two vessels were identified at the base of an ischemic ulcer correlating with the previous location. The endoscopist at the time elected to inject epinephrine and apply electrocautery which resulted in an arterial pulsatile bleed and a perforation. The field of vision was extremely compromised in addition to the difficult location. Mounting of a cap on the gastroscope to improve stability of the scope and applying point pressure to decrease bleeding, as well as vigorous water irrigation permitted the identification of the exact bleeding site. Then, an OTSC was deployed in a single attempt which resulted in immediate adequate hemostasis and closure of the perforation.

The authors emphasize that deployment of the OTSC requires pin-point precision to achieve satisfactory hemostasis. They rate the use of the over-the-scope clip simple yet very effective. The device was not only a rescue tool for hemostasis from a recurrent actively bleeding GDA

pseudoaneurysm, but also for simultaneous perforation closure.

The video can be viewed directly from the GIE website or by using the QR code below.



Over-the-scope clip to the rescue of a bleeding gastroduodenal artery pseudoaneurysm

El Douaihy Y, Kesavan M, Deeb L, Abergel J, Andrawes S (2016).

Gastrointest Endosc. 2016 Jun 12. pii: S0016-5107(16)30236-X. doi: 10.1016/j.gie.2016.05.043.

OTSC® Update 27

February 2018 | OTSC® prevents rebleeding in over 70 % of high-risk GI bleeding cases

J Brandler and colleagues, Department of Internal Medicine, Mayo Clinic, Rochester, Minnesota, USA, performed a study on 67 patients with gastrointestinal bleeding from high-risk lesions who were treated with the OTSC System.

The definition of high-risk lesions was lesions situated in the area of a major artery and larger than 2 mm in diameter and /or a deep penetrating, excavated fibrotic ulcer with high-risk stigmata, in which perforation could not be ruled out or thermal therapy would cause perforation, or lesions that could not be treated by standard endoscopy (epinephrine injections, hemoclips, coagulation). Between 12/2011 and 02/2015, data from 67 patients with high risk non-variceal gastrointestinal bleeding, of which 49 received OTSCs as primary and 18 as rescue therapy, was prospectively collected and retrospectively analysed. Clinical severity was determined based on the Rockall score and a modified Blatchford score.

Out of 67 patients, 47 (70.1 %) remained free of rebleeding at 30 days after OTSC placement. No difference was found in the proportion of patients with rebleeding who received primary or rescue therapy (hazard ratio .639; 95 %CI .084 – 4.860; P=.6653). Only 9 rebleeding events were linked clearly to OTSCs and required intervention, indicating an OTSC success rate of 81.3 %.

The authors concluded that OTSCs have a valuable role in managing and averting high-risk radiologic or surgical interventions for non-variceal gastrointestinal bleeding, despite the presence of high risk of adverse outcome and severe prognostic scores.

Efficacy of Over-the-Scope Clips in Management of High-Risk Gastrointestinal Bleeding

Brandler J, Baruah A, Zeb M, Mehfooz A, Pophali P, Wong Kee Song L, AbuDayyeh B, Gostout C, Mara K, Dierkhising R, Buttar N (2017)

Clin Gastroenterol Hepatol. 2017 Jul 26. pii: S1542-3565(17)30857-1. doi: 10.1016/j.cgh.2017.07.020.

November 2017 | 96 % hemostasis with OTSC® as first-line treatment in patients with gastrointestinal bleeding: an Italian multicentric study comprising 201 consecutive patients

Mangiafico S et al., Azienda Ospedaliero, University of Modena, Italy, presented at the 25th UEG week (October 28 – November 1, 2017, Barcelona) data from 9 Italian tertiary referral centers comprising a large series of patients with non-variceal upper and lower gastrointestinal bleeding lesions in whom OTSC was used as first-line endoscopic treatment.

Over a period of three years (01/2014 - 01/2017), data on 201 consecutive patients (mean age 68 years, range 28-89 years), who underwent emergency endoscopy for severe acute nonvariceal gastrointestinal bleeding and were treated with OTSC as primary first-line therapy, was prospectively collected and analyzed.

106/201 patients were treated with the a version of the OTSC system while in 95/201 patients the t clip was preferred. Indications for OTSC treatment included duodeno-jejunal ulcer Forrest 1a (n=29) and Forrest 1 b (n=35), gastric ulcer Forrest 1a (n=19) and Forrest 1b (n=28), Mallory Weiss (n=19), Dieulafoy's lesion (n=9), post

gastric- ESD bleeding (n=14), post EMR bleeding (n=15), post ESD bleeding (n=12), traumatic rectal ulcer (n=2), colonic diverticulum (n=4), and surgical anastomosis bleeding (n=15).

Technical success was achieved in all cases (100 %). Primary hemostasis was achieved in 193/201 patients (96 %). In the remaining 8 patients hemostasis was obtained with radiological vascular embolization (n=5) or surgery (n=3).

Early rebleeding (within the first 24 hours) occurred in 9/201 patients (4 %) and it was treated with epinephrine injection with or without use of through the scope clips or radiological vascular embolization. No late rebleeding was observed in the series.

The authors concluded that the use of OTSC as first-line therapy in acute high-risk gastrointestinal bleeding is safe and highly effective.

High efficacy of OTSC as first-line endoscopic treatment in patients with gastrointestinal bleeding: an Italian multicentric experience in a large cohort of patients.

Mangiafico S, Russo S, Lupo M, Caruso A, Grande G, Zito F, Bertani H, Conigliaro R, Pisani A, Germaná B, Galloro G, Pasquale L, Mangiavillano B, Bassotti G, Mutignani M, Manta R (2017)

November 2017 | Recommendation for OTSC® as first-line therapy in non-variceal upper gastrointestinal bleeding

Chan SM and Lau JYW, Prince of Wales Hospital, The Chinese University of Hong Kong, Hong Kong, China, published an editorial in Endoscopy International Open on the question: "Can we recommend OTSC as first-line therapy in case of non-variceal upper gastrointestinal bleeding?"

The authors explicate that 8 to 15 % of patients with non-variceal upper GI bleeding (NVUGIB) continue to bleed after endoscopic hemostasis and acid suppression therapy. Further bleeding remains one of the most important predictors of mortality. These facts make research on methods to improve endoscopic hemostasis so important.

The authors list several limitations to conventional hemostatic methods such as the impossibility to consistently seal larger vessels with thermocoagulation, the difficulty of tangential application of hemostatic clips, the frequent dislodgement of the clips and the difficulty of clip application in chronic ulcers with a fibrotic base. The authors argue that the Over-the-Scope-Clip, with a wider jaw and greater strength, has the advantages of a firm grip over a larger amount of tissue. Clip retention is almost universal. The editorial names the study from Wedi et al with 100 patients with NVUGIB and first-line OTSC management and a reported 94 % success rate for primary hemostasis. Besides, the study of Richter-Schrag et al is cited, including 100 patients with both NVUGIB and lower GI bleeding and showing similar results.

However, the paper also names problems that can lower the success of OTSC hemostasis, namely tangential application or OTSC deployment with scope in retroflexion (when ulcers are located in the lesser curve or the posterior wall of the duodenal bulb). The text offers a solution to this problem: usage of a smaller OTSC and an anchoring device to puncture near the bleeding site to guide the OTSC. Second, pretreatment with adrenaline injection is recommended to improve visualization in case of actively bleeding ulcers.

The authors narrate to eagerly await the publication of the STING trial, which randomized patients with refractory bleeding to OTSC or conventional treatment. They propose an RCT comparing OTSC as primary treatment to current standards.

In summary, the editorial recommends the application of OTSC in patients with hemodynamic instability, comorbid illness, with active bleeding ulcers, large ulcers and ulcers at posterior duodenum and lesser curve. The authors speculate that the added cost in managing further bleeding after standard treatment likely outweighs the cost of OTSC.

Can we now recommend OTSC as first-line therapy in case of non-variceal upper gastrointestinal bleeding?

Chan SM, Lau JYW

Endoscopy International Open 2017; 05: E883–E885

September 2017 | Prospective multicenter study confirms efficacy and safety of the FTRD® device for difficult-to-resect colorectal lesions

Schmidt A, Beyna T, Schumacher B et al., affiliated to different German hospitals (Ludwigsburg, Freiburg, Düsseldorf, Essen, Ulm, Augsburg, Dortmund, Krefeld and Tuebingen) published a prospective multicenter study which evaluates efficacy and safety of the FTRD device for resection of colorectal difficult adenomas (non-lifting and/or at difficult locations), early cancers and subepithelial tumours (SET). Hitherto, small retrospective studies had shown feasibility of the EFTR technique in the before-mentioned indications.

In total, 181 patients were recruited. The overall technical success rate was 89.5 %, the overall R0 resection rate 76.9 %. In 127 patients with difficult adenomas and benign histology, R0 resection rate was 77.7 %. In 14 cases, lesions harboured unsuspected cancer, another 15 lesions were primarily known as cancers. Of these 29 cases, R0 resection was achieved in 72.4 %; 8 further cases had deep submucosal infiltration >1000 µm. Therefore, curative resection could only be achieved in 13/29 (44.8 %). In the subgroup with SET (n=23), R0 resection rate was 87.0 %. In general, R0 resection rate was higher with lesions equal to or smaller than 2 cm vs >2 cm (81.2 % vs 58.1 %, p=0.0038). Adverse event rate was 9.9 % with a 2.2 % rate of emergency surgery. Three month follow-up was available from 154 cases and recurrent/residual tumour was evident in 15.3 %.

In summary, this first prospective multicenter study on the FTRD System demonstrates its efficacy for colorectal lesions that otherwise would have required more invasive techniques (mainly surgery), especially in lesions equal to or smaller than 2 cm with acceptable complication rates. The authors state that the curative resection rate for early cancers was too low though to recommend its primary use in this indication for now. They state that further comparative studies will show the clinical value and the longterm outcome in such lesions.

Colonoscopic full-thickness resection using an over-the-scope device: a prospective multicentre study in various indications

Schmidt A, Beyna T, Schumacher B, Meining A, Richter-Schrag HJ, Messmann H, Neuhaus H, Albers D, Birk M, Thimme R, Probst A, Faehndrich M, Frieling T, Goetz M, Riecken B, Caca K (2017)

Gut 2017;0:1–10. doi:10.1136/gutjnl-2016-313677

OTSC® Update 26

July 2017 | Video case report: OTSC® hemostasis in patients with refractory bleeding due to chronic peptic ulcer

Xiao X and Lau JY, Department of Surgery, Chinese University of Hong Kong, Hong Kong, published an article on VideoGIE, the official video journal of the American Society of Gastrointestinal Endoscopy, showing OTSC treatment in two patients with refractory peptic ulcer bleeding.

The first patient was an 89-year-old woman admitted with fresh hematemesis and a haemoglobin of 4.8g/dl. Endoscopy revealed bleeding from a 2- cm chronic bulbar ulcer. She was treated by angiographic embolization to her right gastroduodenal artery (GDA). Three days later, she again experienced massive bleeding. A pulsatile vessel at the ulcer base was discovered and treatment with an OTSC clip induced. The cap was adjusted to encompass the vessel, and a trip string was pulled to deploy the OTSC. The patient was discharged 4 days later without further bleeding.

Patient two was a 76 year old man presenting with fresh melena and a haemoglobin of 7.5 g/dl. He reported on a history of recurrent bleeding from a chronic gastric ulcer. Additionally, he had previously been on warfarin therapy for the treatment of deep vein thrombosis complicated by pulmonary embolism. Endoscopy revealed bleeding from a chronic ulcer at the ankle incisura of the stomach. The first

attempt to stop the bleeding with heaterprobe and hemoclips failed. Then an OTSC anchor device was used to target the ulcer base and deploy an OTSC clip without suction. Complete hemostasis was achieved and the patient had an uneventful recovery.

The authors concluded that OTSC is useful in the treatment of chronic peptic ulcerations with refractory bleeding. The anchor device was rated a helpful tool, which allows accurate targeting of the bleeding artery.

Over-the-scope clip treatment of refractory peptic ulcer bleeding

Xiao X, Lau JY (2016)

Gastrointest Endosc. 2016 Feb;83(2):458-9. doi: 10.1016/j.gie.2015.05.040.

https://www.youtube.com/watch?v=G6u_szn_Yqc&feature=youtu.be 309

June 2017 | 80 % success in endoscopic closure of post-surgical gastrointestinal leaks

R Manta et al., Niguarda-Ca Granda Hospital, Milan, Nuovo S. Agostino Hospital, Modena, Nuovo Regina Margherita Hospital, Rome, Baggiovara Hospital Modena and Federico II University of Naples, Naples, all Italy, published a large case series on patients with post-surgical gastrointestinal leaks managed with endoscopy as initial approach.

A total of 76 patients underwent endoscopic treatment for a leak either in the upper (47 cases) or lower (29 cases) gastrointestinal tract. The first attempt for leak closure was the application of one or more OTSC clips. Fibrin glue was used as an adjuvant treatment to close the gap between two OTSCs where needed. A covered self-expanding metal-stent (SEMS) was applied when the closure was considered incomplete at endoscopy. When dehiscence characteristics were not fitting for OTSC positioning, a SEMS was directly used. Endosponge was the first line therapy, when an abscess cavity was present beyond the anastomotic leak.

Leak closure was successful in 39 patients with upper GI-leaks (83 %) and 22 patients with lower GI leaks (75.9 %), accounting for an overall 80.3 % success rate. Leak closure failed in 15 (19.7 %) patients, and the surgical approach was successful in all 14 patients who underwent re-intervention, whilst one patient died due to sepsis at day 7 post-op.

The authors conclude that an endoscopic approach is successful and safe in the majority of patients with anastomotic gastrointestinal leaks. Therefore, endoscopic treatment should be attempted before resorting to more invasive, costly and risky re-interventions.

Endoscopic management of patients with post-surgical leaks involving the gastrointestinal tract: A large case series.

Manta R, Caruso A, Cellini C, Sica M, Zullo A, Mirante VG, Bertani H, Frazzoni M, Mutignani M, Galloro G, Conigliaro R (2015)

United European Gastroenterology Journal 0(0) 1-8 DOI: 10.1177/2050640615626051

OTSC® Update 25

March 2017 | High-risk GI bleeding: primary hemostasis in first-line OTSC® treatment in 95 %

HJ Richter-Schrag and colleagues, Center of Interdisciplinary Gastrointestinal Endoscopy and Department of General and Visceral Surgery, University of Freiburg, Germany, performed a retrospective study evaluating rebleeding, primary failure and mortality of patients, in whom OTSCs were used as first-line and second-line endoscopic treatment (FLET, SLET) of upper and lower gastrointestinal bleeding (GIB). All patients with upper and lower GIB who underwent FLET and SLET with OTSCs between 04/2012 and 05/2016 were included. In addition, patients with upper GIB were categorized by complete Rockall risk score, and the data were used to calculate predictors of OTSC success and mortality.

A total of 93 patients (58 males, median age 72, range 19-98) with altogether 100 severe acute GIB lesions fulfilled the inclusion criteria. One patient had 3 OTSC applications,

and five other patients had 2 OTSCs on different lesions. First-line OTSC treatment was performed in 61 cases and second line OTSC treatment in 42 cases. The mean hospital stay was 19.8 d (range 1-79). Primary hemostasis was achieved in 88 % of cases (88/100). Clinical success (no in-hospital rebleeding) was achieved in 78 % of cases (78/100). Primary failure was significantly lower when OTSCs were applied as FLET compared to SLET (4.9 % vs 23 %, $P=0.0008$). Patients with Rockall scores ≥ 7 had a significantly higher in-hospital mortality compared to those with scores < 7 (35 % vs 10 %, $P=0.034$). No significant differences were observed in patients with scores < 7 in rebleeding and rebleeding-associated mortality. The authors concluded that the reduction of primary failure in endoscopic treatment of severe acute gastrointestinal bleeding was best achieved when OTSC was used for first line treatment. In this series, first line OTSC treatment seemed to be a predictor of successful reduction of rebleeding rates.

First-line endoscopic treatment with over-the-scope clips significantly improves the primary failure and rebleeding rates in high risk gastrointestinal bleeding: A single-center experience with 100 cases.

Richter-Schrag HJ, Glatz T, Walker C, Fischer A, Thimme R (2016) World J Gastroenterol 2016 Nov 7; 22(41): 0000-0000. ISSN 1007-9327 (print) ISSN 2219-2840 (online)

December 2016 | Experience with the FTRD® System in Halle, Germany: FTRD® broadens endoscopic therapeutic spectrum and reduces surgery rate

The 25th annual meeting of the Society for Internal Medicine in Sachsen-Anhalt (25. Jahrestagung der Gesellschaft für Innere Medizin Sachsen-Anhalt) took place on November 18-19, 2016 in Halle (Saale), Germany. Ohse C et al., hospital Martha-Maria Halle-Doelau, Germany, presented in a poster their experience with the FTRD device. This poster received the Poster Champ Award at the meeting. The FTRD System is used at the hospital since January 2015. Until October 2016, 31 endoscopic full-thickness resections with the FTRD in 30 patients were performed. Indications were: non-lifting adenoma ($n=9$), subepithelial tumor ($n=7$), adenoma at appendiceal orifice ($n=1$), adenoma at appendiceal orifice after Stapler-appendectomy ($n=1$), hybrid EMR-EFTR in large adenoma with non-lifting parts ($n=2$), follow-up resection of a carcinoma in situ after piece-meal-EMR ($n=1$), Tis-carcinoma ($n=2$), T1sm1 carcinoma ($n=2$), T1sm2 carcinoma ($n=1$), T1sm3 carcinoma ($n=3$), T2 carcinoma ($n=2$). The lesions were located in the rectum ($n=10$), sigmoidal colon ($n=3$), left flexure ($n=2$), transverse colon ($n=3$), right flexure ($n=3$), ascendant colon ($n=4$), and cecum ($n=6$). The median size of the resected specimen was 19 mm (range 13-35 mm) Resection of the lesions with the FTRD succeeded in 100 % (31/31). En-bloc-R0-resection was achieved in 94 % (29/31). One patient (3.2 %) suffered a post-polypectomy syndrome after resection of an adenoma at the appendiceal orifice, this complication could be managed conservatively. No other complications occurred. In 5 patients, correct oncological follow-up resection was undertaken due to the T-stage (3xT1sm3, 2xT2). During 3-months follow-up, which could be carried out in 26 patients, a relapse adenoma was found in one patient and treated with a second EFTR.

The authors concluded that EFTR with the FTRD System broadens the therapeutic spectrum of lesions in the lower gastrointestinal tract and helps to avoid surgery in selected patients.

Endoskopische Vollwandresektion (EFTR - 'endoscopic full-thickness resection') mit dem FTRD-System (Full-Thickness-Resection-Device): Dölauer Daten

Ohse, Geelhaar H, Hopf K, Buchmann J, Krummenerl P (2016)

Poster auf der Tagung der Gesellschaft für Innere Medizin Sachsen-Anhalt, 18/19.11.16 in Halle (Saale)

November 2016 | Large single-center experience with 101 OTSC® applications in

patients with severe hemorrhage, perforations and fistulae: 89 % overall primary clinical success

Wedi E and colleagues, Strasbourg University Hospitals, Strasbourg, France and St. Bernward Academic Teaching Hospital, Hildesheim, Germany, and Icahn School of Medicine at Mount Sinai, New York, United States, and Boston Children's Hospital, Boston, United States, conducted a retrospective study to investigate efficacy and clinical outcome of patients treated with an OTSC clip for gastrointestinal (GI) emergencies and complications. From 02/2009 to 10/2012, 84 patients were treated with 101 OTSC clips. 41 patients (48.8 %) presented with severe upper-GI bleeding, 3 (3.6 %) patients with lower-GI bleeding, 7 patients (8.3 %) underwent perforation closure, 18 patients (21.4 %) had prevention of secondary perforation, 12 patients (14.3 %) had control of secondary bleeding after endoscopic mucosal resection or endoscopic submucosal dissection (ESD) and 3 patients (3.6 %) had an intervention on a chronic fistula. In 78/84 patients (92.8 %), primary treatment with the OTSC was technically successful. Clinical primary success was achieved in 75/84 patients (89.28 %). In detail OTSC application lead to a clinical success in 35/41 (85.36 %) patients with upper GI bleeding and in 3/3 patients with lower GI bleeding. Technical success of perforation closure was 100 % while clinical success was seen in 4/7 cases (57.14 %) due to attendant circumstances unrelated to the OTSC. Technical and clinic success was achieved in 18/18 (100 %) patients for the prevention of bleeding or perforation after endoscopic mucosal resection and ESD and in 3/3 cases of fistula closure. Two application-related complications were seen (2 %). In conclusion, this study confirms the high value of the OTSC for the treatment of severe gastrointestinal bleeding, fistula closure and the non-surgical management of perforations.

One hundred and one over-the-scope-clip applications for severe gastrointestinal bleeding, leaks and fistulas

Wedi E, Gonzalez S, Menke D, Kruse E, Matthes K, Hochberger J (2016) World J Gastroenterol. 2016 Feb 7; 22(5): 1844-1853.

October 2016 | New patient series with the FTRD® System in Switzerland: EFTR efficacious and safe

The annual congress of the Swiss Society for Gastroenterology (SSG) this year took place in Interlaken, Switzerland on September 22-23, 2016. Clinical experience with the FTRD from two tertiary referral centers in Switzerland was presented by P. Aeppli and colleagues, canton hospitals of Lucerne and St. Gallen, Switzerland. Nineteen consecutive patients with colonic polyps were treated with the FTRD during a one-year period (05/15 - 05/16). Thirteen procedures were performed in the colon and six in the rectum. Indications were adenoma recurrence or residual adenoma with non-lifting sign after previous polypectomy ($n=7$), staging resection following presumed incomplete polypectomy of early carcinoma ($n=6$), treatment-naïve adenoma with non-lifting sign ($n=5$), and one adenoma located at the appendiceal orifice ($n=1$). In only one case (polyp at appendiceal orifice) the lesion could not be reached due to diverticular narrowing of the sigmoid. En bloc and histologically complete (R0) resection was achieved in 17 of 18 cases (94.4 %), complete full-thickness resection succeeded in 15 of 18 cases (83.3 %), i.e. 11/12 (91.7 %) in the colon and 4/6 (66.7 %) in the rectum. The mean diameter of resection specimen was 2.6 cm (range 1.8-3.2 cm). Two post-procedure minor bleedings were seen (one requiring re-colonoscopy and adrenalin injection). During a one-month follow-up no further complications were seen. The authors concluded, that treatment of non-lifting polyps ≤ 30 mm with the FTRD device was efficacious and sufficiently safe. Full thickness resection was evaluated an adjunctive colonoscopic technique, which offers minimally invasive treatment to a group of patients that would otherwise undergo surgery.

Full thickness resection device (FTRD): A novel tool for colonoscopic adenoma resection. First clinical experience from two tertiary referral centers in Switzerland.

Aeppli P, Frei R, Borovicka J, Criblez D (2016) Swiss Med Wkly. 2016 Sep; 146 (Suppl 218): 2S

September 2016 | Complete endoscopic resection of early colorectal cancer with FTRD®: case study of a high-risk anticoagulated patient

A case report on the use of the FTRD System for R0-resection of early colorectal cancer was published by P. Lagoussis et al., Polyclinic San Donato, San Donato Milanese, Italy.

Colonoscopy was performed due to hematochezia in a 78-year-old man with a history of coronary artery disease and recent pulmonary embolism, therefore under anticoagulant therapy. A 3-cm non-pedunculated colorectal polyp was observed and resected by en-bloc endoscopic mucosal resection. Histology revealed an adenocarcinoma (pT1 G2 Sm3) with a positive resection margin (0.7 mm) and deep submucosal invasion (1.4 mm). Total body computed tomography and rectal endoscopic ultrasound showed no lymphatic or metastatic disease. Because of the patient's comorbidities and anticoagulant treatment, endoscopic resection with the FTRD System was considered the appropriate therapeutic option to achieve R0-resection. After antibiotic prophylaxis with an intravenous cephalosporin and the last dose of low-molecular-weight heparin being administered 12 hours before, the procedure was carried out. The lateral margins of the scarred resection site were marked with argon plasma coagulation (APC). The FTRD was mounted on the tip of an operative gastroscope. Through a tissue anchor, the whole scarred lesion was pulled into the cap and the clip was deployed. The pseudopolyp thereby created was resected using the preloaded snare and a standard electrosurgical setting. The procedure took about 8 minutes and there were no complications. Low-molecular-weight heparin was re-introduced 24 hours thereafter and the patient was discharged.

Histological examination of the full-thickness specimen 15 mm in size revealed no remnant dysplasia. This outcome was confirmed in the biopsy samples taken from the rectal scar 3 months later. Endoscopic ultrasound and CT further confirmed the absence of lymphatic or metastatic disease, therefore a chemotherapy was relinquished. In summary, full-thickness resection with the FTRD System was a feasible and safe treatment for early colorectal cancer in this high risk patient, where standard surgery would carry considerable risks and require aggressive strategies.

Over-the-scope clip-assisted endoscopic full-thickness resection after incomplete resection of rectal adenocarcinoma

Lagoussis P, Soriani P, Tontini GE, Neumann H, Pastorelli L, de Nucci G, Vecchi M (2016) Endoscopy. 2016; 48: E59-E60. doi: 10.1055/s-0042-100197.

September 2016 | Endoscopic full-thickness resection of gastric submucosal tumors: Use of OTSC® significantly associated with shorter hospital stays

Yang F et al., Endoscopic Center, Shengjing Hospital of China Medical University, Shenyang, China, conducted a retrospective cohort study to identify factors that impact the procedure and treatment outcomes for endoscopic full-thickness resection (EFTR) of gastric submucosal tumors (SMTs).

Gastric SMTs include gastrointestinal stromal tumors leiomyomas, schwannomas, malignant lymphomas, lipomas, carcinoids, lymphangiomas, and hemangiomas. They are usually detected incidentally during upper gastrointestinal endoscopy, and have an estimated prevalence of 0.4 %.

For the study, the medical records of all patients with gastric SMTs who underwent EFTR procedures in Shengjing Hospital between June 2012 and April 2014 were reviewed. EFTR procedures had been performed as follows: the tumor was completely removed endoscopically, including surrounding mucosa, muscularis propria, and serosa, without injury to the tumor capsule. The post-resection gastric defect was closed immediately using metallic clips or an OTSC system.

In total, 41 patients (13 males and 28 females, mean age 53.9 ± 14.1 years) were treated. All patients underwent endoscopic ultrasound (EUS) before EFTR, the mean tumor size based on EUS was 16.3 ± 5.9 mm. Of the 41 SMTs, 1 was located in the anterior wall of the antrum, 2 in the greater curvature of the antrum, 6 in the anterior wall of the corpus, 6 in the greater curvature of the corpus, 3 in the lesser curvature of the corpus, 9 in the posterior wall of the corpus and 13 in the fundus of the stomach. EFTR of the gastric SMTs was successfully performed in all 41 patients. The final pathologic analyses revealed R0 Resection in all cases. An OTSC system was used in 6 patients, metal clips were used in 35 patients. Maximum tumor size on EUS and tumor location in the greater curvature were significantly associated with the length of the procedure ($P=0.000$ and $P=0.026$, respectively). A pneumoperitoneum occurred in 26 patients during EFTR. There were no cases of bleeding, peritonitis, or abdominal abscesses. A pneumoperitoneum was more likely to occur during EFTR in tumors with a larger EUS size ($P=0.017$). The use of the OTSC system was significantly associated with shorter hospital stays ($P=0.047$) and a higher cost of defect closure ($P=0.001$). The authors conclude that endoscopic full-thickness resection is an effective and safe treatment for patients with gastric submucosal tumors.

Factors associated with endoscopic full-thickness resection of gastric submucosal tumors

Yang F, Wang S, Sun S, Liu X, Ge N, Wang G, Guo J, Liu W, Feng L, Ma W (2015)

Surg Endosc. 2015; 29(12): 3588-3593.

August 2016 | Post-ESD duodenal ulcer closure with OTSC® vs conventional clips: Significantly less adverse events and shorter closure times with OTSC®

Mori H and colleagues, Departments of Gastroenterology and Neurology, Faculty of Medicine, Kagawa University, Japan, compared in a retrospective study safety of closure methods, closure times and medical costs between two groups of patients who had post-endoscopic resection (ER) artificial ulcers in the duodenum.

Nineteen patients with duodenal adenoma, early duodenal cancer, and subepithelial tumors underwent resection of the lesion by endoscopic submucosal dissection (ESD) between 09/2009 and 09/2014. Ulcer floor closure was achieved either by conventional clips (in 12 patients) or over-the-scope clips (in 7 patients). The closure method was selected at the discretion of each operator (i.e. three endoscopists, each having performed > 100 cases of gastric ESD).

Delayed bleeding was observed in three patients from the conventional clip group, in the OTSC group no delayed bleeding occurred ($p=0.049$). No perforations were observed in either group. The mean procedure time for ulcer closure was 33.26 ± 12.57 min with conventional clips and 9.71 ± 2.92 min with OTSC clips ($p=0.0001$). The resection diameters were 18.8 ± 1.30 mm and 22.9 ± 1.21 mm for the conventional clip group and the OTSC group, respectively ($p=0.039$).

In conclusion, the authors state that if the post-endoscopic resection ulcer is over 20 mm, the OTSC closure should be selected with regard to safety and reliable closure.

Suitable closure for post-duodenal endoscopic resection taking medical costs into consideration

Mori H, Ayaki M, Kobara H, Fujihara S, Nishiyama N, Matsunaga T, Yachida T, Masaki T (2015) World J Gastroenterol 2015; 21 (17): 5281-5286 ISSN 1007-9327 (print) ISSN 2219-2840 (online)

August 2016 | OTSC® safe and effective for closure of gastric access in NOTES appendectomy and other procedures

In NOTES (natural orifice transluminal endoscopic surgery), surgical procedures in the abdominal cavity are performed with an endoscope passed through a natural orifice, which allows to decrease the number of trocars placed through the abdominal wall or eliminating them completely. However, obtaining reliable closure of access points at the completion of a NOTES procedure is currently the most significant

limitation and one of the main factors that will determine the future of NOTES in clinical use. Magdeburg R and Kaehler G, Department of Surgery, Medical Faculty Mannheim, University of Heidelberg, Germany, conducted a retrospective study on prospectively collected data of patients who underwent transgastric flexible endoscopic NOTES with final gastric access closure by OTSC clip application.

Between 04/2010 and 03/2014, a total of 43 patients (mean age 34 years) had received an OTSC clip for gastric access closure after a transgastric NOTES procedure. Indications have been acute appendicitis ($n=36$), prophylactic bilateral salpingo-oophorectomy ($n=6$), and uterus myomatosis ($n=1$). In all 43 cases, the endoscopic access to the abdominal cavity succeeded without any difficulty. After performing the operation, closure of the access by OTSC application was also possible in all cases. Overall, three adverse events occurred: in two patients (4.7 %) clinical signs of acute gastrointestinal bleeding appeared, gastroscopically there was no demand for action in one case and in the other case the bleeding could be stopped by hemoclip application. The third patient (2.3 %) suffered insufficiency of the gastric closure with local peritonitis; this patient received laparoscopic suturing and the abdominal cavity was cleaned with liquid solution. With antibiotic therapy, no further problems occurred.

In summary, this retrospective study showed in more than 40 consecutive patients that OTSC application is safe and effective for closure of gastrotomy after NOTES procedures, however, further investigation in clinical settings is absolutely necessary to establish clear indications and guidelines for the use of transgastric NOTES. **Natural orifice transluminal endoscopic surgery in humans: feasibility and safety of transgastric closure using the OTSC system**

Magdeburg R, Kaehler G (2015)

Surg Endosc DOI 10.1007/s00464-015-4163-4, Epub March 24, 2015

July 2016 | >90 % success rate in closure of post-ESD gastric ulcers by combined use of a single OTSC® and supplemental through the scope clips

A five to seven day hospital stay is usually needed after endoscopic submucosal dissection (ESD) of gastric tumor because of the risk of delayed perforation or bleeding. S. Maekawa et al., Department of Gastroenterological Medicine, Niigata Rosai Hospital Japan, evaluated in a prospective study the efficacy of a new technique for post-ESD ulcer closure by combined use of a single over-the-scope clip (OTSC) and additional through-the-scope clips (TTSCs).

Twelve patients (average age 71 years) with early gastric cancer or gastric adenoma underwent the surgical procedure. En-bloc and margin-free resection of the lesion by ESD was achieved in 100 %. Mean size of post-ESD ulcer was 54.6 mm. The mean time for the closure procedure was 15.2 min, and the success rate was 91.7 % (11/12). Mean numbers of OTSC and TTSCs were 1 and 5.8 per patient, respectively. No complications (ex. delayed perforation, postoperative bleeding and gastric stenosis) occurred. The patients with successful ulcer closure could be discharged the day after ESD.

In conclusion, the authors considered the method of post-ESD ulcer closure using a combination of a single OTSC and TTSCs a safe and effective technique for shortening the period of hospitalization and reducing treatment cost. Further investigations are needed to establish the appropriate role of the closure method.

Complete closure of artificial gastric ulcer after endoscopic submucosal dissection by combined use of a single over-the-scope clip and through-the-scope clips (with videos)

Maekawa S, Nomura R, Murase T, Ann Y, Harada M (2015) Surg Endosc 29:500-504 DOI 10.1007/s00464-014-3725-1

July 2016 | OTSC® Proctology: 70 % closure rate in recurrent complex anal fistulae in retrospective analysis

R. Mennigen, et al., Department of General and Visceral Surgery, University Hospital Muenster, Germany, conducted a retrospective study to evaluate the efficacy of

the OTSC Proctology in patients with multiple previous fistula operations, Crohn's disease, or anovaginal fistulae. Only patients with refractory anal fistulae were included, in which alternative surgical approaches had failed or were not feasible any more, e.g. due to scarring. By this policy, only the most problematic and refractory anal fistulae were selected.

Between 10/2012 and 06/2014, five male and five female patients, with a median age of 41 years (range 26-69 years) met the inclusion criteria. The etiology of the fistula was cryptoglandular in four patients, and perianal Crohn's disease in six patients (including one patient with an anovaginal fistula).

The surgical procedure was technically successful in all patients. There were no intraoperative or postoperative complications. Permanent fistula closure was achieved in seven out of ten patients (70 %) within a median time of 72 days (range 31-109 days). Median total follow-up time was 230.5 days (range 156-523 days). There were three failures (30 %), including two cryptoglandular and one Crohn's disease-associated fistula. In all three cases, the OTSC was lost spontaneously on days 22, 23, and 40, respectively, and persistence of the fistula was diagnosed thereafter. In three of the seven patients with successful closure, the OTSC was removed after complete healing of the fistula, because of slight anal discomfort or soiling. These symptoms disappeared completely after clip removal. There was no postoperative incontinence.

The authors conclude that the OTSC Proctology system is a safe and effective method for the closure of even complex and recurrent fistulae. It will compete with established surgical procedures in the management of transsphincteric or suprasphincteric fistulae.

The OTSC proctology clip system for the closure of refractory anal fistulas

Mennigen R, Laukötter M, Senninger N, and Rijcken E (2015)

Endoscopy. 2015 0;47(S 01):E115-E116. Epub 2015 Mar 11

June 2016 | OTSC® vs. cSEMS for intestinal leakage: clipping associated with significantly higher clinical success rates

A retrospective study comparing over-the-scope clips (OTSC) and covered self-expanding metal stents (cSEMS) for upper gastrointestinal perforation or leakage was conducted by Prof. Dr. med. H. Farnik, university hospital Frankfurt am Main, and colleagues in four German tertiary endoscopic centers (Frankfurt, Tübingen, Jena, Dortmund). Technical success, outcome (e.g. duration of hospitalization, in-hospital mortality), and complications were assessed and analyzed with respect to etiology, size and location of leakage.

Between 2006 and 2013, overall 106 patients underwent endoscopic treatment for postoperative leakage, endoscopic perforation or spontaneous rupture of the upper gastrointestinal tract. Of these, 72 (69 %) were treated by cSEMS and 34 (31 %) by OTSC.

OTSC was preferred in small-sized lesions and in perforation caused by endoscopic interventions, cSEMS in patients with concomitant local infection or abscess. For cSEMS vs. OTSC, mean treatment duration was 41.1 vs. 25 days ($p < 0.001$), median leakage size was 10 mm (range 1-50 mm) vs. 5 mm (range 1-30 mm), and complications were observed in 68 % vs. 8.8 % ($p < 0.001$), respectively.

Clinical success for primary interventional treatment was observed in 29/72 (40 %) vs. 24/34 (70 %, $p = 0.006$), and clinical success at the end of follow-up was 46/72 (64 %) vs. 29/34 (85 %) for patients treated by cSEMS vs. OTSC; $p = 0.04$.

In conclusion, cSEMS and OTSC are rather complementary means than to be mutually exchangeable. The authors suggest, that, due to its low complication profile and high effectivity rates, the OTSC should be the first choice in all cases when it is technically feasible and the diameter of the lesion is not too large. In patients with larger defects and already infection accompanying the leak, cSEMS placement might be preferred.

Indication for 'Over the Scope' (OTS)-Clip vs. Covered

Self-Expanding Metal Stent (cSEMS) Is Unequal in Upper Gastrointestinal Leakage: Results from a Retrospective Head-to-Head Comparison

Farnik H, Driller M, Kratt T, Schmidt C, Fährndrich M, Filmann N, Königsrainer A, Stallmach A, Heike M, Bechstein WO, Zeuzem S, Albert JG (2015)

PLoS One. 2015 Jan 28;10(1):e0117483. doi: 10.1371/journal.pone.0117483. eCollection 2015.

May 2016 | 100 % defect closure rates with OTSC® after EFTR of gastric tumors and no complications

Guo J and colleagues, Shengjing Hospital of China Medical University, Shenyang, China, reported on their experience of defect closure with OTSC after endoscopic full-thickness resection (EFTR) of gastric subepithelial tumors. Between October 2013 and March 2014, 23 patients underwent EFTR of a gastric subepithelial tumor ≤ 2 cm originating from the muscularis propria by intentional transection of the gastric wall. Defect closure was achieved by tissue approximation with an OTSC clip. Endoscopic follow-up was performed at 1 week, 1 month and 6 months after operation to check OTSC closure.

The full-thickness resection rate was 100 % (23/23), the success rate of defect closure was also 100 % (23/23), and the average time of defect closure was 4.9 min (range 2-12 min). No post-operative complications such as bleeding and perforation were seen, and no premature OTSC detachment was found.

The authors conclude that OTSC placement is a simple, convenient, safe and effective way of defect closure after EFTR of gastric subepithelial tumors.

Endoscopic full-thickness resection with defect closure using an over-the-scope clip for gastric subepithelial tumors originating from the muscularis propria

Guo J, Liu Z, Sun S, Liu X, Wang S, Ge N, Wang G, Qi Y (2015)

Surg Endosc. doi: 10.1007/s00464-015-4076-2

April 2016 | Technical success rates and long-term clinical outcomes of fistula closure with OTSC®

Dr. R. Law and Dr. L.M. Wong Kee Song, Division of Gastroenterology and Hepatology, Mayo Clinic, Rochester, MN, USA together with Dr. S. Irani, Digestive Disease Institute, Virginia Mason Medical Center, Seattle, WA, USA and Dr. T.H. Baron, Division of Gastroenterology and Hepatology, University of North Carolina at Chapel Hill, NC, USA investigated long-term clinical outcomes of fistula closure with OTSC.

The authors note that most of the literature focuses on short-term success, so they conducted a two-center, review of all patients at Mayo Clinic, Rochester and Virginia Mason Medical Center, Seattle, who were treated with OTSC for closure of chronic fistulas from October 2011 to September 2012. The retrospective study includes 47 unique patients (24 men, 23 women) at a mean age of 57 years. Fistula locations were the small bowel (18), stomach (16), colon-rectum (10) and esophagus (3). Previous percutaneous endoscopy gastrostomy/jejunostomy and prior bariatric procedure were the cause of these fistulas in 10 cases respectively. Patients were treated with the gc (gastric closure) OTSC clip for gastric indications and t (traumatic) clips for all other indications. A guidewire, tissue graspers and tract ablative therapies (i.e., argon plasma coagulation or cytology brush) were used where deemed beneficial by the endoscopist.

Initial technical success rate, measured by lack of contrast extravasation after OTSC placement, was 89 % (42/47). At a median of 39 days (range: 26-86 days), however, fistulas recurred in 19 patients (46 %). The OTSC was still present adjacent to these fistulas in 16 cases (84 %) at repeat intervention. The authors cite current literature reporting similar technical success rates to their own findings, but note that reported long-term success varies widely, from 74 % in a systematic meta-analysis by Weiland et. al. to three studies reporting success rates of 38, 50 and 67 % respectively.

Since nearly 2/3 of fistula recurrence in their study was discovered over 4 weeks after initial technical success, the authors conclude that longer term follow-up is necessary to

properly assess success rates. They call for future studies to determine which specific types of fistula are most likely to respond well to OTSC closure and which assistant therapies might facilitate long-term treatment success. **Immediate technical and delayed clinical outcome of fistula closure using an over-the-scope clip device.** Law R, Wong Kee Song LM, Irani S, Baron TH. Surg Endosc. 2015 Jul;29(7):1781-6. doi: 10.1007/s00464-014-3860-8.

March 2016 | OTSC® helps to avoid emergency laparotomy, according to review on endoscopic treatment of iatrogenic gastrointestinal perforations

The "Deutsches Aerzteblatt" (German journal of physicians) reports on a review conducted by Schmidt and colleagues assessing endoscopic treatment of iatrogenic gastrointestinal perforations. Emergency laparotomy after iatrogenic gastrointestinal perforation can often be avoided by application of modern endoscopic occlusion techniques. The review presents different strategies for management of iatrogenic perforations. These include supportive therapies such as broad-spectrum antibiotic treatment or decompression of high intraabdominal pressure. For endoscopic closure of the perforations different techniques are displayed, namely standard and OTSC-clips and, for oesophageal perforations, covered stents. Standard through-the-scope (TTS-) clips are used for small colonic perforations and slit-shaped gastric perforations < 10 mm. However, these clips normally only grasp mucosa and submucosa, closure of larger lesions especially in oesophagus and stomach cannot be accomplished by standard clips.

Over-the-scope clips (OTSC) on the contrary also grasp deeper layers of the gastrointestinal wall. This enables closure of 20 -30 mm lesions. Besides, deployment of the OTSC is a time saving one-step action in contrast to the sequential occlusion with standard clips. Mainly retrospective clinical studies can be found assessing the use of OTSC-clips for iatrogenic gastrointestinal perforations. These studies examined between 3 and 48 patients with a follow-up of one week to 92 weeks. Technical success rates reported were between 50 and 100 %. Clinical success rates vary because of heterogeneous size and localisation of the perforations and time of diagnosis. The prospective, multicentric CLIPPER - study comprises 36 patients with iatrogenic gastrointestinal perforations < 30 mm. Technical and clinical success rates reached 92 and 89 %, respectively. A systematic review conducted by Weiland et al. in 2013 reports on a technical success rate of 80-100 % and a clinical success rate of 60-100 %. Authors of recently published retrospective multicentre studies with 106 and 48 patients recount similar results. Based on the named results the occlusion with OTSC-Clips in stomach and colon found entry into the recommendations of the current position paper of the European Society of Gastrointestinal Endoscopy (ESGE).

The endoscopic treatment of iatrogenic gastrointestinal perforation

Schmidt A, Fuchs KH, Caca K, Küllmer A, Meining (2016) Dtsch Arztebl Int 2016; 113: 121-8. DOI: 10.3238/arztebl.2016.0121

March 2016 | Use of OTSC® system is safe and efficient in endoscopic full-thickness-resection in the upper and lower digestive tract

Use of OTSC system is safe and efficient in endoscopic full-thickness-resection in the upper and lower digestive tract Fährndrich M and Sandmann M, Department of Gastroenterology, Klinikum Dortmund, report about a study to evaluate the efficacy and safety of the over-the-scope clip (OTSC) system for endoscopic full-thickness-resection (EFTR) of endoluminal gastrointestinal tumors. The retrospective, observational, open-label, single-arm, consecutive case study from May 2010 to May 2014 included 17 patients (8 men and 9 women with a median age of 57.65 years). Six patients suffered from carcinoids, located in stomach, duodenum or rectum, seven patients from R1 situations after conventional endoscopic

polypectomy with low risk for colorectal cancer of the colon, three from adenoma relapses in the colon and one patient from a submucosal lesion of the stomach. The mean diameter of the lesions was 22.7 mm with a range of 10 - 25 mm.

In 16 cases, a combination of the OTSC System and the Inoue Cap for EFTR was used and in one case a new, dedicated full-thickness resection device (FTRD). All cases were performed using OTSCs (12 /6gc or 14/6t) or FTRD clips.

The application of the clips was successful in 16/17 cases, only in one case the clip did not deploy correctly. An endoscopic successful resection was accomplished in all remaining patients (16/16, 100 %).

On final pathology, a full-thickness resection was achieved in 69 % (11/16) and a deep muscle margin of resected specimen (DMR) was accomplished in the remaining 31 % (5/16). Also a Complete resection (R0) was achieved in all patients.

The authors discuss that the OTSC System is safe and feasible for EFTR and eDMR of gastrointestinal lesions up to 25 mm in diameter. In their opinion the technique does not replace existing methods but may become a useful addition of the interventional endoscopy and an effective and valid alternative to surgical resection. As another advantage of OTSC, the article describes that the dynamic clip does not induce strangulation necrosis (like endoloops may do).

The report concludes that the use of the over-the-scope clip system in endoscopic full-thickness-resection is a safe and efficient addition for patients and endoscopists. **Endoscopic full-thickness resection for gastrointestinal lesions using the over-the-scope clip system: a case series**

Fährndrich M, Sandmann M (2014)

Endoscopy. 2015 Jan;47(1):76-9. doi: 10.1055/s-0034-1377975. Epub 2014 Sep 15.

February 2016 | Complete closure of artificial gastric ulcer after endoscopic submucosal dissection by combined use of a single over-the-scope clip and through the scope clips

Maekawa S, Nomura R, Murase T, Ann Y and Harada M from the Department of Gastroenterological Medicine, Japan Labour Health and Welfare Organization Niigata Rosai Hospital, Japan, reported about a prospective study to evaluate the combined use of a single OTSC and multiple TTSCs in closure of artificial gastric ulcer after ESD. From June 2013 to March 2014 nine patients were included with early gastric cancer and three patients with gastric adenoma with a mean age of 71 ± 8 years. The tumors were 3 cm or less in diameter and did not involve the cardiac or pyloric region. Mean size of the artificial ulcer after ESD was 54.6 mm. The average operating time of post-ESD artificial gastric ulcer was 15.1 minutes. 1 OTSC and a mean of 5.8 TTSCs per patient were used.

The success rate of complete defect closure was 91.7 % (11/12). No complications like delayed perforation, postoperative bleeding or gastric stenosis occurred. The OTSC dropout rate was 0 % on the day following ESD and 8.3 % two months later.

The authors argue that ESD-associated complications are common in the treatment of gastrointestinal tract tumors. Studies have shown that postoperative bleeding occurs in 6-16 % of patients after gastric ESD. Maekawa et al. see an advantage of the OTSC System at this point. Beside the closure of fistulas, perforation sites, leaks and severe bleeding of the GI tract the OTSC could improve therapeutic options in endoscopy. Often it has been difficult to close large artificial ulcers only with TTSCs. But the combined use of a single OTSC with its closure power and easy handling showed excellent results in this study. The new closure technique allowed the authors to discharge their patients from hospital after only 2 days instead of the traditional 5-7 days.

In conclusion, Maekawa and colleagues consider this new closure method of artificial gastric ulcers after ESD safe and feasible. It is useful for shortening the period of hospitalization and reducing treatment cost compared to traditional therapy with through-the-scope clips only.

Complete closure of artificial gastric ulcer after endoscopic submucosal dissection by combined use of a single over-the-scope clip and through the scope clips

Maekawa S, Nomura R, Murase T, Ann Y, Harada M. (2014) Surg Endosc. 2015 Feb; 29(2):500-4. doi: 10.1007/s00464-014-3725-1. Epub 2014 Jul 23.

February 2016 | OTSC® as successful last resort before surgery for challenging bleeding lesions

Endoscopedia, the official blog of "GIE: Gastrointestinal Endoscopy" recently published a video on OTSC use by Dr. James Y. W. Lau. The video, called "Over-the-scope clip treatment of refractory peptic ulcer bleeding", presented two cases in which OTSC was used to treat chronic peptic ulcerations with refractory bleeding after failed angiographic embolization and endoscopic treatment respectively. Dr. Lau concluded that, "Endoscopists should consider the use of OTSC when tackling challenging bleeding lesions especially when other standard treatments have failed and certainly before referring your patients to surgery."

February 2016 | Use of Over-The-Scope Clip for treatment of refractory upper gastrointestinal bleeding: high technical and clinical success rates reported

Chan SM, Chiu PWY, Teoh YB and Lau JYW from the Department of Surgery, Institute of Digestive Disease, Prince of Wales Hospital, Chinese University of Hong Kong, China, reported about a prospective case series to evaluate the safety and efficacy of the Over-The-Scope Clip in patients with refractory GI bleeding.

The case series from included nine patients (4 men, 5 women) with a median age of 72.5 years (range 39 - 91 years), suffering from bleeding gastric ulcers (n=2), bleeding duodenal ulcers (n=5), bleeding gastrointestinal stromal tumor in the stomach (n=1), and bleeding from ulcerative carcinoma of the pancreas (n=1). Median diameter of the ulcers was 2.5 cm (1-4 cm). Six of the nine patients underwent previous endoscopic hemostasis attempts. A total of 10 OTSCs were applied in the nine patients. The technical success rate of OTSC was 100 % (10/10). Endoscopic hemostasis was achieved in all patients. No local complications occurred. The clinical effectiveness was 77.8 % (7/9), while two patients with specific conditions developed rebleeding after OTSC application due to chronically fibrotic ulcers because of residual tumor infiltration and previous radiotherapy.

Chan and colleagues discuss that in 8/10 patients, the bleeding was located in difficult positions, where application of conventional clips would have been complicated as the endoscope approach to the ulcer would have been at a steep angle. The OTSCs allowed a larger amount of tissue to be captured for compression compared to common clips while avoiding the possibility of thermal injury with its high risk of perforation, as can happen with thermal hemostasis methods.

As numerous methods of endoscopic hemostasis have been developed, the authors recommend considering the OTSC System in refractory gastrointestinal bleeding before conventional clips, surgery or angiographic embolization.

Use of Over-The-Scope Clip for treatment of refractory upper gastrointestinal bleeding: a case series
Chan SM, Chiu PWY, Teoh YB, Lau JYW
Endoscopy. 2014 May;46(5):428-31. doi: 10.1055/s-0034-1364932. Epub 2014 Feb 6.

January 2016 | OTSC® as successful treatment of massively bleeding jejunal varix, which had resisted previous interventions

S Kothari, T Kothari and V Kaul of the Center for Advanced Therapeutic Endoscopy, Division of Gastroenterology and Hepatology at the University of Rochester/Strong Memorial Hospital in Rochester, NY, USA presented a case of successful treatment of massive gastrointestinal bleeding from a jejunal varix with OTSC after several other treatment options had failed. The 67-year old male patient had a medical history of coronary artery disease, chronic renal

insufficiency and Laennec's cirrhosis before he was admitted for a laparoscopic left radical nephrectomy for renal cell carcinoma. Afterwards, the patient suffered from several complications, including superior mesenteric vein thrombosis, melena with a significant drop in hematocrit and clinical signs of bleeding, which led to identification and unsuccessful treatment of several possible bleeding sites. The patient underwent anticoagulation, a tagged red blood cell scan, angiography, coil-embolization, repeat mesenteric angiography and repeat (push) enteroscopy. The patient also received a total of 38 units of packed red cells, 13 units of thawed plasma, 9 units of fresh frozen plasma, 3 units of platelets and 2 units of cryoprecipitate. Due to multiple comorbidities, he was deemed as a high-risk patient unfit for surgery.

Finally, a tortuous, varix-like, prominent blood vessel with a central small ulceration, bleeding actively, was identified in the proximal jejunum. Ethanolamine injection into the varix did not achieve hemostasis. Finally, a size 12/6t OTSC clip was applied over the actively bleeding jejunal varix using a pediatric colonoscope. Instant and complete hemostasis was achieved with this single clip. No additional transfusions were required and his hematocrit stabilized over the next few days. Due to his overall poor prognosis and multiple comorbidities, the patient's family opted for "comfort measures only" and he passed away several days later. The authors emphasize the fact that they were able to quickly and effectively treat a massively bleeding jejunal varix, which had resisted multiple evaluations and courses of treatment. They deem the OTSC device a major advance in endoscopic management of high-risk patients in a variety of challenging clinical settings, especially in case of poor candidates for surgical intervention. They also note that endoscopic perforation management with the OTSC clip may avoid the cost and morbidity of surgery and other interventions. Statement by Ovesco Endoscopy: the treatment of jejunal varix hemorrhage is not a common indication for the use of OTSC and there is limited experience with such application.

The Over-The-Scope-Clip Device: An Indispensable Tool in Interventional Endoscopy: A Case Series.
Kothari S, Granato cm; Sharma S, Kothari T, Fagan N, Adamciewicz M, Wang G, Ullah A, Kaul V. Program No. P234. ACG 2015 Annual Scientific Meeting Abstracts. Honolulu, HI: American College of Gastroenterology.

January 2016 | Novel remove DC ClipCutter for endoscopic removal of OTSC®s: a retrospective case series

Dr. Schmidt and colleagues from the Department of Gastroenterology and Oncology, Klinikum Ludwigsburg, Germany, report about a new method of removing OTSCs using a bipolar DC impulse cutting device. The retrospective study includes 11 consecutive patients who underwent endoscopic OTSC removal between December 2012 and November 2013. The clips were located all along the digestive tract, in the esophagus, gastric antrum, pylorus, pyloric-jejunal anastomosis, descending duodenum, sigmoid colon and rectum. By that time, the DC ClipCutter was applied under compassionate use conditions, given the lack of suitable other methods. The clip time in situ ranged between 31 and 469 days (mean clip time in situ: 138 days) before extraction. The DC ClipCutter is an endoscopic, bipolar instrument device connected to an electrical generator (DC Impulse) producing direct current impulses which are applied at two opposite sides of the clip. The maximum strength of 140 A selectively heats up the nitinol and separates it thermally. The mean procedure time for clip removal in this case series was 47 minutes (range 35 - 75 minutes), given the specific circumstances of compassionate use with a pre-series system. Cutting of the clip was successful in 100 % of cases (11/11). All clip fragments were successfully removed in 10/11 patients (91 %). In one patient a single clip fragment could not be removed due to deep ingrowth into the duodenal wall. The fragment was left in place and had not caused any complications in follow-up visits three months later. During clip cutting, no local or systemic complications were observed. Minor bleeding occurred in two patients after clip

removal and was treated with epinephrine injection. One patient had a superficial mucosal tear with minor bleeding which presented during extraction of a clip fragment. The bleeding could be managed by application of two standard clips. OTSCs have been demonstrated to be effective tools for endoscopic treatment of gastrointestinal perforation, leakages and fistulas and are also used for full-thickness resection. They usually dislodge from the tissue and are expelled after weeks or months. Although OTSCs are made of biocompatible material and suitable as implants, active removal can be indicated in special situations. The remOVE technique is described as safe and effective with a high overall success of 91 %. In conclusion, Dr. Schmidt et al. expect the remOVE System to be a valuable tool for OTSC removal for emergency situations, e.g. inadvertent misplacement of the clip as well as for elective OTSC removal, where clip extraction is indicated.

Endoscopic removal of over-the-scope clips using a novel cutting device: a retrospective case series

Schmidt A, Riecken B, Damm M, Cahyadi O, Bauder M and Caca K Endoscopy 2014, DOI: 10.1055/s-0034-136549

December 2015 | OTSC® is a safe and efficient technique in treatment of colorectal postsurgical leaks and fistulas

Bonino MA, Verra M, Salvai A, Bullano A, Rapetti L, Arezzo A and Morino M, Department of Surgical Sciences, University of Turin, Italy, reported at the 20th National Congress of Digestive Diseases in Napoli, Italy, March 19th-22nd 2014, about a prospective study including 26 patients treated with OTSC for postsurgical leaks and fistulas of the colorectum. Anastomotic leakage is a feared and serious complication in colorectal surgery associated with increased morbidity and mortality. The prevalence in literature ranges from 1 to 39 %, whereas clinically relevant leaks commonly occur in 3-6 %.

Out of 26 consecutive patients treated with OTSC in this study, 10 patients suffered from acute and 16 patients from chronic leaks (fistulas). The mean defect diameter was 8.7 mm. 14 cases were complicated by recto-vaginal, rectovesical or colo-cutaneous fistula. In 3 cases OTSC was used to complete endoscopic vacuum-assisted closure of large defects.

The overall success rate was 77 % (20/26): 90 % (9/10) in acute and 69 % (11/16) in chronic cases. No OTSC-related complications occurred in the study. Further surgery was required in two cases.

The authors conclude that the endoscopic OTSC closure of colorectal postsurgical leaks is a safe technique with a high success rate in acute and chronic cases, including the treatment of fistulas.

Efficacy of the Over-the-Scope Clip (OTSC) for treatment of colorectal postchirurgical leaks and fistulas.

Bonino MA, Verra M, Salvai A, Bullano A, Rapetti L, Arezzo A and Morino M (2014)

20th National Congress of Digestive Diseases, Napoli, Italy, March 19-22 2014

November 2015 | Promising case series on novel OTSC® removal device

Dr. A. Schmidt, Dr. B. Riecken, Dr. M. Damm, Dr. O. Cahyadi, Dr. M. Bauder and Dr. K. Caca, Department of Gastroenterology and Oncology, Ludwigsburg Hospital, Germany, reported the results of a case series of OTSC clip removal with a specifically developed cutting device (future trade name: remOVE, courtesy of Ovesco).

The authors note that over-the-scope clips (OTSC) have proven their efficacy in an ever-growing variety of indications. Still, there is a lack of follow-up studies on the course of the clip once it is deployed in the gastrointestinal tract. Clinical experience shows that clips usually fall off after several weeks or months, depending on the amount of tissue grasped. Since OTSC clips are fully biocompatible, they may stay in place indefinitely. However, there are a few situations which call for active removal. In these situations, clinicians have to rely on techniques that lack proof of safety and efficacy. In the case series, the prototype of a specifically developed bipolar cutting device for OTSC removal was used.

A total of 11 patients (male: 7, female: 4) between 43 and 73 (median age: 62) were treated with the device under the "compassionate use" statute. Average procedure time was 47 minutes (range: 35-75 minutes) and cutting of the clip was successful in all cases (100 % success rate). Fragments were removed successfully in all but one case (91 %), where a fragment was deeply grown into the duodenal wall. In a follow-up after 3 months, this fragment had not caused any complications. Indications for clip removal included intermittent epigastric/abdominal pain, the need for a repeat biopsy at the resection site and patients' wishes.

The cutting device consists of a grasper connected to a generator. Two electrodes in the grasper tips conduct a direct current impulse, which heats up the clip and thermally cuts the material. The current has a maximum strength of 140 ampere and maximum duration of 100 milliseconds. After successful cutting, the device stops automatically, and there is no relevant current pathway through the patient's tissue.

During the procedures, the clip was grasped at its thinnest visible part and cut by applying 1 to 4 direct current impulses. Then the clip was cut again at the opposite site, creating two clip fragments, which were then removed from the tissue with standard forceps. To avoid tissue damage during extraction of sharp fragments, a soft clear plastic cap was used (prototype accessory, courtesy of Ovesco). After removal of all fragments, the OTSC site was inspected for bleeding or perforation. Minor bleeding was observed in two patients and treated immediately. One patient had a superficial mucosal tear with minor bleeding in the lesser curvature, caused by fragment removal without the use of a cap. In all subsequent patients, the cap was used for fragment extraction.

The authors stress that the removal of an OTSC clip is only necessary in a minority of cases. Even if it does not fall off on its own, the clip is designed as a durable implant. The few indications for endoscopic clip removal include wrongly placed clips, accidental clipping of the retraction device to the tissue, local complications due to the clip (e.g. luminal obstruction, ulcers), scheduled removal of esophageal stents fixed with an OTSC, and the need for a repeat biopsy of the OTSC site. Additionally, patients may insist that symptoms are caused by the clip's presence.

According to the group, recommended removal times have yet to be established. However, they assume that for repeat biopsy, a good time span is 3 months, while for esophageal stent removal they favor 6-8 weeks. Longer periods of time are correlated with deep ingrowth into the intestinal wall, which impedes removal. They also caution that clips placed for perforation closure may not be removed too early to prevent re-perforation. In the case series, clips had been placed at least 4 weeks prior, while the clip in-situ for the longest time (15 months) was difficult to remove and had the longest procedure time (75 minutes).

Overall, the authors see the new device for OTSC clip removal as an indispensable addition to the clinician's repertoire since it proved safe and effective in all treated patients.

Endoscopic removal of over-the-scope clips using a novel cutting device: a retrospective case series.

Schmidt A, Riecken B, Damm M, Cahyadi O, Bauder M, Caca K.

Endoscopy. 2014 Sep;46:9:762-6.

doi: 10.1055/s-0034-1365493

182

September 2015 | OTSC® as effective treatment of GI fistulae: abscess drainage increases healing rates to 88 %

Dr. P. Mercky, Dr. J.-M. Gonzalez, Dr. E. Aimore Bonin, Dr. O. Emungania, Dr. J. Brunet, Dr. J.-C. Grimaud and Dr. M. Barthet of the Departments of Gastroenterology and Digestive Surgery, North Hospital, Méditerranée University, Marseille, France, presented the results of a retrospective study in two teaching hospital centers.

The study encompassed 30 patients (12 male, 18 female) of 23 to 75 years of age (mean age: 48) suffering from GI fistulae (upper GI: 24; lower GI: 6). 60 % of fistulae (18 patients) occurred after laparoscopic sleeve gastrectomy (LSG), the other 12 patients suffered from fistula at mixed

locations (rectovaginal, urethrorrectal, rectovesical, gastrogastic, gastrocutaneous, esophagojejunal fistulae and one colorectal anastomotic leak). Fistula orifice sizes ranged from 3 mm to 20 mm (mean: 7.2 mm) and mean time between fistula diagnosis and OTSC placement was 12.4 months (8 days to 10 years). 18 patients (60 %) had previously undergone endoscopic or surgical treatment attempts of their fistula in another center.

The OTSC clip was placed successfully in all cases, and 16 patients (53 %) recovered without further intervention. Others required secondary treatment. Overall final success rate in the whole group was 70 %.

Regarding efficacy, the paper notes several points: In patients with a previously drained abscess, the success rate was significantly higher (88.2 % healing vs. 53.8 %). The highest primary efficacy, however, could be reported for fistulae related solely to LSG (88.9 %).

Efficacy was related to operator experience since the overall success rate of the high-volume center was higher than that of the low-volume center (74 % vs. 60 %).

The authors note that OTSC placement allows closure of much larger fistulae (<30 mm) in a single procedure than standard clips, and that the procedure is similar to the well-established technique of band ligation, which helps the operators' learning curve. They also state that it enables treatment of difficult-to-treat fistulae, e.g. those at the lower extremity of the staple line after LSG. Judging from the eight patients who did not benefit from OTSC placement, the authors speculate that previous radiotherapy (and resulting microvascular lesions) as well as cardiovascular risk factors might impede fistula closure with OTSC. Overall, they deem the OTSC system a safe and effective treatment for GI fistulae, both as a primary and secondary option. They also emphasize that the treatment of serious and difficult-to-treat leakages after LSG especially benefitted from the use of OTSC.

Usefulness of over-the-scope clipping system for closing digestive fistulas.

Mercky, P, Gonzalez, JM, Aimore Bonin E, Emungania O, Brunet J, Grimaud JC, Barthet M.

Dig Endosc. 2015 Jan;27(1):18-24.

doi: 10.1111/den.12295.

185 / 238

OTSC® Update 20

August 2015 | Clinical experience with OTSC® shows high success rate for recurrent bleeding and complex resections

E. Wedi and J. Hochberger of the Department of Hepato-Gastroenterology at the University Hospital of Strasbourg, France, reported on clinical experiences with the over-the-scope clip system and its application aids such as the OTSC Twin Grasper and OTSC anchor for coarse tissue.

They present a review of 14 clinical studies on OTSC use and add their own experiences with 84 patients (101 OTSC applications). All Strasbourg patients suffered from recurrent bleeding/lesions of perforations and fistulae or post-operative leakage. 78 out of 84 cases (92.85 %) could be treated successfully.

The report notes that misapplication and complications are rare (<3 % according to the literature) when using OTSC. If they occur, they include narrowing of the organ lumen in case of small passageways and comprehensive aspiration of tissue. The authors remark that this should be kept in mind, especially in narrow sections of the esophagus or bowel.

They also maintain, however, that OTSC is often a good option for achieving quick closure in case of acute perforation or severe bleeding. Possible complications might then be compensated in a subsequent procedure, e.g. through dilation of a stenosis. They also report isolated cases of a grasper getting caught in the OTSC clip. However, if the clip is deployed before the grasper is pulled back fully into the cylinder, this complication is extremely unlikely in clinical practice.

The authors report that a device for cutting the clip for later removal as well as other instruments based on the OTSC concept show promising results in experiments.

The paper concludes that OTSC is an asset in interventional endoscopy, especially in case of complex endoluminal resections. According to the literature, OTSC is especially

useful for closure of perforations of up to 1.3 cm (and much larger in individual cases) and bleeding lesions with a high risk of recessive bleeding, e.g. in anti-coagulated patients or treatment of acute Forrest Ia/b hemorrhage.

Chronic fistulae, which have limited chances of successful treatment, regardless of treatment method, due to insufficient circulation in scarred and calloused tissue, remain a challenge, even with OTSC. Caution is also in order when closing no longer fresh postoperative leakage or perforations since these may require sufficient external drainage.

Klinische Erfahrungen mit dem Over-the-Scope Clip (OTSC)

Wedi E, Hochberger J (2014)

Endo-Praxis, 30.1, 14–17

dx.doi.org/10.1055/s-0034-1370894

179

June 2015 | Prophylactic OTSC® application for prevention of complications after colorectal ESD

A team from Kagawa University, Japan, around Dr. S. Fujihara, Department of Gastroenterology and Neurology, together with Dr. Kazi Rafiq, Department of Pharmacology, conducted a study exploring prophylactic closure after colorectal endoscopic submucosal dissection (ESD) using the OTSC clip as a preventative measure.

The authors note that in spite of ESD's increasing prevalence, the technique is not widely used in treatment of large superficial colorectal neoplasms. They assume that this is due to the technical skill required and higher incidence of complications, like inflammation and especially perforation (reported at a 1.4–10.4 % incidence in previous studies) as well as postpolypectomy syndrome and transmural burn syndrome. They also note the risk of delayed perforation after ESD, possibly linked to excessive coagulation in the muscular layer, which sometimes even requires emergency surgery (incidence at 0.3–0.7 % in previous studies).

The study included 68 patients (39 male, 29 female) whose colorectal tumors (mean tumor size 35.4 mm) were treated with ESD. Prophylactic closure was performed on patients with excessive coagulation in the muscularis propria or larger resection size (n=27). Closure was performed either with conventional clips (n=18) or OTSC (n=9). OTSC was used for large mucosal defects (>30 mm), in case of flexure of the colon, excessive coagulation in the muscularis propria and when closure could not be achieved with conventional clips. OTSC closure required more time than conventional closure since in 5 cases tumors were located at a sharp bend in the sigmoid colon and required the use of the Twin Grasper. The median of clips needed was 8 for conventional clips (range 4–12) and 1 (range 1–3) for OTSC.

The study showed the efficacy of endoscopic closure after ESD in preventing local peritoneal inflammation and abdominal symptoms without any adverse effects. It did, however, not demonstrate a reduction of perforation and post-operative bleeding. The authors point out the limited sample size of their study and lack of randomization. They call for more randomized prospective studies with more patients to confirm their encouraging results.

Efficacy and safety of over-the-scope clip: including complications after endoscopic submucosal dissection

Nishiyama N, Mori H, Kobara H, Rafiq K, Fujihara S, Kobayashi M, Oryu M, Masaki T

World J Gastroenterol. 2013 May 14;19(18):2752-60

doi: 10.3748/wjg.v19.i18.2752

242

May 2015 | Conference Report: German Society for Endoscopy and Imaging Procedures (DGE-BV)

The conference, held from March 26–28, 2015 in Munich, featured papers on fistula closure and perforation management (even after EVT or stent therapy) with OTSC®, full-thickness resection with FTRD® at success rates of well over 80 %, treatment of therapy-resistant complicated fistulae with OTSC® Proctology, and one presentation on the new remOVE System for OTSC clip removal (currently in development).

Conference Report | talks, posters and videos

45th Conference of the German Society for Endoscopy and

Imaging Procedures (DGE-BV)

Munich, Germany, March 26 – 28, 2015

45. DGE-BV-Kongress, München, 26. – 28. März 2015 Deutschen Gesellschaft für Endoskopie und Bildgebende Verfahren e.V.

Chairman: Prof. Dr. Hans-Dieter Allescher, Garmisch-Partenkirchen

Ovesco products were presented in four workshops on two different topics (hemostasis techniques, held by E. Wedi and A. Nägel respectively, and management of complications led by T. Lankisch and J. Bernhardt respectively). Additionally, several talks, posters and videos discussed products by Ovesco.

OTSC System

Improved fistula closure with OTSC: treatment of fistulae with highly fibrotic openings

A Meining, University Hospital Ulm, together with M. Bajbouj, Technical University Munich, well as H. Feußner und D. Wilhelm, Technical University Munich, presented a video of a 41-year-old female patient who had developed a large fistula into the bronchial system after resection of the esophagus with gastric interposition due to lye ingestion. Because of coarse, scarred mucosa at the fistula opening in the gastro-esophageal junction it was impossible to permanently and securely anchor an endoscopic clip. Overstenting treatment was also unsuccessful, and the patient had to be tube fed for several weeks.

In order to ensure secure anchoring of the OTSC clip, the mucosa around the fistula opening was incised in a diameter of about 15 mm (corresponding to the opening of the OTSC applicator) with an HF knife. In the uncovered submucosa, the clip could be securely fastened and applied around the fistula with the mobilized mucosal edges covering the opening of the fistula. Administration of a contrast agent confirmed fistula closure, the patient suffered no more aspirations, and endoscopic follow-up after three months confirmed treatment success.

The group recommends incision of the mucosa before application of the OTSC clip in case of coarse and chronically indurated mucosa tissue in order to improve clip anchoring.

Erfolgreicher Verschluss einer großen ösophago-bronchialen Fistel durch mukosale Inzision vor OTSC-Klipp-Platzierung

A. Meining, Ulm; M. Bajbouj, H. Feußner, D. Wilhelm, München

OTSC as treatment option for fistulae, even after EVT or stent therapy

M. Laukötter, T. Vowinkel, D. Palmes, N. Senninger und R. Mennigen, University Hospital Münster, presented two cases in which patients suffered from leaking anastomoses after sleeve gastrectomy.

In the first case, three separate insufficiencies with abscess formation manifested along the staple line in a 50-year-old male. He was treated with endoscopic vacuum therapy (EVT) and a total of 25 polyurethane sponges. The group contrasted this with the case of a 45-year-old woman, who had been unsuccessfully treated with stents after an anastomosis fistula. After the stent had been extracted, one OTSC clip was applied. Treatment duration was eight days. In both cases, no complications occurred during or after and treatment was successful.

The group recommends the OTSC system for fistulae in fresh lesions free from inflammation as well as for rescue treatment after EVT or stent therapy failed.

Endoskopisches Komplikationsmanagement nach Sleeve-Magen-Resektion.

M. Laukötter, T. Vowinkel, D. Palmes, N. Senninger und R. Mennigen, Münster.

OTSC as important treatment option for management of perforations

A. Meining, University Hospital Ulm, emphasized the important role of timing (acute or chronic), pathogenesis and etiology, size and location of a perforation for treatment. He favors the OTSC clip as the best treatment option for perforations during ESD since it is easy to place and offers a higher success rate than other, smaller clips. All in all, he sees OTSC as a useful therapy tool, even for late complications (e.g. after unsuccessful stent placement) and for fistula closure.

Rationales Handeln nach Perforation: Clippen, Stenten

oder Sponge?

A. Meining, Ulm.

“Live Demonstration on Stage” on correcting OTSC clip placement

During a live demonstration using a porcine stomach model, K. Caca, Ludwigsburg hospital, showed possible mistakes in placing OTSC clips and offered suggestions for preventing or correcting these errors. He also demonstrated OTSC removal with the remOVE System.

OTSC-Clip/Clip Fehlplatzierung.

K. Caca, Ludwigsburg; R. Landschoof, Düsseldorf.

FTRD (Full-Thickness Resection Device)

FTRD for endoscopic full-thickness resection in case of NET

P. Klare, B. Neu, M. Bajbouj, R.M. Schmid and S. von Delius, TU Munich, together with R. Burlefinger, Maria Theresia Hospital, Munich, und K. Specht, TU Munich, presented a video the treatment of a 50-year-old male with a neuroendocrine tumor (NET). During a screening procedure, the asymptomatic patient had a rectal polyp removed with a snare. Histological examination identified a NET with a proliferation rate < 2 percent, resection was deemed incomplete, a rectal endosonography and PET-CT revealed no pathologic findings.

Six weeks later, an endoscopy revealed a scar at the resection site, and it was decided to perform an endoscopic full-thickness resection using the FTRD System. The edges of the target area were marked and the tissue was pulled completely into the 23 mm long FTRD cap with the FTRD grasper. Then the FTRD Clip was deployed and full-thickness resection was performed using the integrated HF snare. No complications occurred during the procedure and the patient was fit to leave the hospital the next day.

The group stated that NET < 1 cm in size with a low or moderate proliferation rate (G1 and G2) may be resected endoscopically, but that complete resection is a prerequisite for successful treatment. In case of incomplete resection, they deem endoscopic full-thickness resection a minimally invasive treatment option, which enables definitive histological examination.

Endoskopische Vollandresektion nach R1-Resektion eines neuroendokrinen Tumors im Rektum unter Verwendung eines neuen Over-the-scope Device.

P. Klare, R. Burlefinger, B. Neu, M. Bajbouj, K. Specht, RM Schmid und S. von Delius, München.

Experiences and Recommendations regarding full-thickness resection using the FTRD System

In his presentation, K. Caca, Ludwigsburg hospital, presented clinical data from Ludwigsburg and Zurich of 38 patients treated using FTRD.

In 36 patients, the target lesion could be reached with the FTRD, in 31 of these cases, treatment was technically successful. Average size of the resected tissue was 23.6 mm (12–40 mm) and histological examination confirmed full resection (R0) in 80.5 % of cases and full-thickness resection 86.1 % of cases. Average duration of hospital stay was 4 days (1–12).

Caca recommends identifying and marking the target lesion without mounted FTRD, marking the lesion with a clip next to the lesion site if necessary, and then rinsing before using the FTRD system. If finding the lesion proves difficult, he recommends fluoroscopy using a TTS balloon (15–20 mm) if appropriate, verifying air insufflation, repeatedly pulling back the endoscope and repositioning the patient. After the resection is performed, the endoscope should be reintroduced without mounted attachment to inspect the resection site. Routine second-look endoscopy the next day is unnecessary if not contraindicated. After 8 to 12 weeks, however, another endoscopy, possibly including clip removal, is advised, although the clip is fully biocompatible and, based on experience, 2/3 of all clips will already have fallen off by then.

Limitations to the procedure are a maximum lesion size of 2.5–3 cm (although the largest lesion treated in this trial was 4 cm) and impaired vision because of the size of the cap. Since the rectum is closely attached to the para rectal tissue in the small pelvis, the procedure often results in “deep resection” rather than full-thickness resection. Caca assesses the FTRD system as valuable and currently the only instrument for full-thickness resection in the lower

gastrointestinal tract. He sees FTRD as an effective therapy option for the main indication of non-lifting neoplasia regarding both diag-nosis and treatment.

Endoskopische Vollwandresektion im Kolon mit FTRD. *K. Caca, Ludwigsburg.*

„Live Demonstration on Stage“ on FTRD usage

Conference chairman H.-D. Allescher, Garmisch-Partenkirchen hospital, presented the use of FTRD for full-thickness resection on a porcine gastric model in a special 3D demonstration.

3D Bildgebung im klinischen Einsatz. Live-Demonstration on Stage. Sonographie, Laparoskopie und Endoskopie.

Vorsitz: H.-D. Allescher, Garmisch-Partenkirchen. F. Hagemüller, Hamburg, U. Beilenhoff, Ulm, H. Feußner, München und A. Melzer, Dundee/Schottland.

OTSC Proctology

OTSC Proctology deemed effective for closure of therapy-resistant complicated anal fistulae

R. Mennigen, University Hospital Münster, presented his center's experience in using OTSC Proctology. Technical success rate for fistula closure with OTSC Proctology was 100 percent in 10 patients (5 male, 5 female) with a median age of 41 years (26 to 69). All patients had been treated unsuccessfully before (e.g. mucosa flap, Anal Fistula Plug), and some had undergone several unsuccessful procedures. 4 patients suffered from crypto-glandular fistulae, and in 6 patients fistulae were associated with Crohn's disease.

In 70 percent of patients permanently successful closure and healing was achieved. In 3 patients fistulae recurred (2 cryptoglandular, one patient with Crohn's disease), and the clip dislodged spontaneously after over three weeks. In 3 patients the clip was surgically removed after successful healing of the fistula. Dr. Mennigen deemed the OTSC Proctology system a safe and effective treatment option, even for recurrent complicated anal fistulae.

Verschluss von komplizierten therapierefraktären Analfisteln mit dem OTSC Proctology System.

R. Mennigen, M. Laukötter, T. Vowinkel, N. Senninger und E. Rijcken, Münster.

remOVE (System for OTSC clip removal, in development)

remOVE prototype for treatment of severe dysphagia due to magnetic implant

A new, promising treatment approach for reflux disease is to augment the lower esophageal sphincter with a magnetic system (LINX® Reflux Management System) consisting of a band of interlinked magnetic beads implanted laparoscopically. A video was presented by M. Bauer, M. Kranzfelder, D. Wilhelm, R. Schirren, A. Jell, H. Friess and H. Feußner, TU Munich, together with A. Meining, University Hospital Ulm.

The video showed a patient suffering from esophageal erosion and resulting dysphagia caused by the LINX® system.

Since there were no other therapy options available, the team was able to treat the patient with the Ovesco clip cutter (trade name remOVE system, market launch planned for 2015), which is not yet available for routine treatment, within "compassionate use". They were able to perform a full endoscopic removal, leading to complete remission of symptoms.

Perforation eines magnetischen Antirefluxsystems nach intraösophageal-Ösophaguserosion durch ein LINX® system.

M. Bauer, M. Kranzfelder, D. Wilhelm, R. Schirren, A. Jell, H. Friess, H. Feußner, München und A. Meining, Ulm.

May 2015 | OTSC® Proctology evaluated favorably at two medical conferences in March 2015

Clinicians presented new data on fistula closure with the OTSC Proctology at the 41st Conference of Coloproctologists as well as at the 45th Conference of the German Society for Endoscopy and Imaging Procedures (DGE-BV), both held in Munich, Germany:

OTSC Proctology proves to be effective, especially in initial intervention

Dr. Lothar Duschka, department of surgery/coloproctology at DKD Helios hospital Wiesbaden, Germany, presented his

experience in treating 44 patients using OTSC Proctology. Patient median age was 47 years (24 to 72) with 14 female and 30 male patients. Operating times were between 16 and 95 minutes. Apart from 25 trans-sphincter fistulae, inter-sphincter (3), supra-sphincter (2), recto-vaginal (4) and anal fistulae (8) were included as well as one recto-vesical and one pouch fistula.

During initial intervention, successful closure was achieved in 78 percent of cases. With recurrent fistulae, the success rate was 42 percent. In patients with chronic inflammatory bowel diseases, 64 percent of fistulae were closed successfully. Dr. Duschka praised minimal preparation time at the site of the fistulae when using OTSC Proctology as well as the fact that special variants could often be treated successfully (e.g. 50 percent of recto-vaginal fistulae and the one pouch fistula). He also emphasized a positive learning curve for the procedure, increasing success rates over time.

Das OTSC-Verfahren und seine (Miss)-Erfolge im klinischen Alltag

Al-Haidary J, Zieker D, Borschitz T, Rimpel J, Duschka L Wiesbaden (41. Koloproktologen-Kongress)

OTSC Proctology deemed effective for closure of therapy-resistant complicated anal fistulae

Dr. Mennigen, department of general and visceral surgery, university hospital Münster, Germany, presented his center's experience in using OTSC Proctology. Technical success rate for fistula closure with OTSC Proctology was 100 percent in 10 patients (5 male, 5 female) with a median age of 41 years (26 to 69). All patients had been treated unsuccessfully before (e.g. mucosa flap, Anal Fistula Plug), and some had undergone several unsuccessful procedures. 4 patients suffered from cryptoglandular fistulae, and in 6 patients fistulae were associated with Crohn's disease.

In 70 percent of patients permanently successful closure and healing was achieved. In 3 patients fistulae recurred (2 cryptoglandular, one patient with Crohn's disease), and the clip dislodged spontaneously after over three weeks. In 3 patients the clip was surgically removed after successful healing of the fistula. Dr. Mennigen deemed the OTSC Proctology system a safe and effective treatment option, even for recurrent complicated anal fistulae.

Verschluss von komplizierten therapierefraktären Analfisteln mit dem OTSC Proctology System

Mennigen R, Laukötter M, Vowinkel T, Senninger N, Rijcken E Münster. (45. DGE-BV-Kongress)

April 2015 | Study identifies OTSC® as effective and safe endoscopic therapy for acute gastrointestinal bleeding

In an observational retrospective case series, Dr. Matthew Skinner, Dr. Juan P. Gutierrez, Dr. Helmut Neumann, Dr. C. Mel Wilcox, Dr. Chad Burski and Dr. Klaus Mönkemüller of the Basil I. Hirschowitz Endoscopic Center of Excellence, Department of Gastroenterology, University of Alabama, Birmingham, USA, evaluated the efficacy and safety of OTSC clip placement in patients with upper gastrointestinal bleeding after traditional endoscopic methods had failed.

The study was conducted at a large tertiary care hospital, comprising 12 patients (8 male, 4 female), mean age of 59 (range: 29–86) with ongoing upper gastrointestinal bleeding despite two or three previous sessions of endoscopic therapy. Patients had a mean ASA score of 3 (range: 2–4), mean hemoglobin of 7.2 g/dL (range: 5.2–9.1), and shock was present in 75 % of patients. They had all received packed red blood cells (mean 5.1 units, range 2–12). Bleeding was due to duodenal ulcer (6), gastric ulcer (2), Dieulafoy lesion (2), anastomotic ulceration (1) and Mallory-Weiss tear (1). Hemostasis was achieved in all patients, but rebleeding occurred in two patients 1 day and 7 days* after OTSC placement. There were no complications associated with OTSC application.

The OTSC System was loaded onto a standard gastro-scope and introduced into the upper gastrointestinal tract under standard direct visualization. The bleeding lesion was located and the gastroscope was maneuvered towards it. Once the OTSC cap was positioned on top of the lesion, full suction was applied to tissue, so that the lesion was fully engulfed inside the transparent cap before the clip was

released. In two patients with post-bulbar ulcers, a wire or a wire placed inside a feeding tube was advanced into the distal duodenum and the scope advanced alongside the wire. This measure helped to prevent small-bowel luminal occlusion, which has previously been reported as a major adverse event.

The authors conclude that the OTSC System provides tissue apposition far superior to traditional clipping and can function as a "rescue therapy" in patients with severe non-variceal upper gastrointestinal bleeding in whom prior endoscopic therapies have failed, avoiding more invasive procedures such as embolization or surgery. They point out that the OTSC System appears promising for the treatment of bleeding lesions with large-diameter visible vessels or those located in awkward positions, such as the greater curvature of the stomach or the posterior duodenal wall, which may not always be amenable to treatment with standard endoscopes and endoscopic devices.

*Any bleeding occurring more than 72 hours after therapy is usually regarded as a new incident.

Over-the-scope clip placement is effective rescue therapy for severe acute upper gastrointestinal bleeding

Skinner M, Gutierrez JP, Neumann H, Wilcox C, Burski C, Mönkemüller K

Endoscopy International Open 2014; 02: E37–E40 171 March 2015 | JFHOD Congress, France: OTSC® in hemostasis – high success rate in anti-coagulated patients

Paris, March 19–22: The JFHOD congress – Journées Francophones d'Hépatogastroentérologie et d'Oncologie Digestive 2015 – the major French clinical congress on gastroenterology, hepatology and GI oncology was held under the presidency of Prof. Jean-Christophe SAURIN.

The group around Prof. J. Hochberger presented their data on OTSC in the treatment of severe gastrointestinal hemorrhage.

Dr. E. Wedi, Dept. Of Gastroenterology and Digestive Endoscopy, University Hospital, Strasbourg, France summarized cases admitted to emergency care due to severe GI bleeding. All patients had Hb <7g/dl upon admission. Median age 73 yrs (29–97). 80 % were under concomitant anti-coagulation or antiplatelet therapy. In 87.2 % (48/55) bleeding was from a gastroduodenal peptic ulcer, and the remaining from various other causes.

OTSC clips were placed. The overall clinical success rate without recurrence was 83.6 % (46/55). In 32.7 % OTSC was placed as a salvage treatment due to prior ineffectiveness of other techniques. 12.7 % (n=7) had to undergo surgical treatment. Out of these 7 patients 4 died, giving a total mortality rate of 7.2 % in this highly challenging case series. The report by Wedi et al. underlines that patients with a high Rockall score can profit from OTSC treatment, especially patients with anticoagulative or antiplatelet therapy.

March 2015 | German Congress of Coloproctology: further clinical data on OTSC® Proctology

Munich, March 12–14, 2015:

The annual congress of the German Society for Coloproctology took place in Munich under the presidency of Prof. Dr. Dr. h. c. W. Hohenberger, Erlangen.

Dr. L. Duschka and colleagues from the department of colorectal surgery and proctology, DKD Helios hospital, Wiesbaden reported in a plenary lecture about their clinical experience in using OTSC Proctology.

In their abstract (Coloproctology (2015); 1:76) they summarize the data of 22 patients, treated between March and August 2014. The majority had trans-sphincteric fistulas (n=18), followed by inter-sphincteric, rectal, recto-vaginal and pouch fistula (one case each). 13 patients had prior fistula surgery and 8 patients suffered from IBD. Post-surgical follow-up was 3–9 months.

68 % of the patients had healing of their fistula, 32 % had recurrence. The authors found that patients without prior history of fistula surgery had a higher probability for healing. They conclude that selection of patients is important to optimize the clinical result.

Das OTSC-Verfahren und seine (Miss)-Erfolge im klinischen Alltag

Al-Haidary J, Zieker D, Borschitz T, Rimpel J, Duschka L
Deutsche Klinik für Diagnostik, Wiesbaden Coloproctology,
February 2015; 37(1):76
41. Deutscher Koloproktologen-Kongress

OTSC® Update 19

February 2014 | Two Studies show efficacy of anchoring esophageal SEMS (self-expanding metal stents) with OTSC®

In a retrospective study of 13 patients, Dr. Irani, Dr. Gluck, Dr. Gan, Dr. Ross and Dr. Kozarek of the Department of Gastroenterology, Virginia Mason Medical Center, Seattle, Washington, together with Dr. Baron, Division of Gastroenterology & Hepatology, Mayo Clinic Rochester, Rochester, Minnesota, explored the efficacy of OTSC for securing self-expanding metal stents (SEMS) in patients who had experienced stent migration before. Indications included post-operative leaks in 4 patients (2 esophagogastric anastomotic leaks, 2 fistulae after bariatric surgery), 1 perforation, 3 benign esophageal strictures (2 peptic, 1 anastomotic), and 5 malignant esophageal strictures. The patients (8 male, 4 female) had a median age of 67 years with a range from 40 to 89 years. Before anchoring the stent with OTSC, three patients had a new SEMS deployed, in the other ten cases the migrated stent was merely repositioned. A standard upper endoscope was advanced to the proximal end of the stent and suction was applied, attempting to position tissue and stent in such a way that upon deployment half of the OTSC would grasp the stent and the other half the esophageal wall. Clip placement was successful in all cases with a median procedure time of 3.5 minutes (range of 2 to 5.5 minutes). Cutting and OTSC/SEMS removal was 6 minutes on average. Migration occurred in two patients, both with benign esophageal strictures at a mean of 32 days compared with a mean of 3.5 days for prior migrations without OTSC use. In 4 of 5 patients with malignant disease, the SEMSs remained in place indefinitely. Successful healing occurred in 11 of 13 patients (85 %) with a median stent dwell time of 57 days (2–226 days range). Of three patients with refractory benign esophageal strictures, one recurred after stent removal.

The authors discuss several uses of esophageal stents, noting that fully covered SEMSs have high migration rates, which call for external or internal fixation. They note that in spite of improved stent fixation times from 26.4 to 12.5 minutes when using a suturing device, the use of OTSC is much more time efficient at a median of 3.5 minutes. Additionally, OTSC placement does not require the use of an overtube or double-channel upper endoscope and prevents additional costs of 700 to 800 USD for suturing. Taking into account the fact that all patients in this study had previously experienced stent migration, the reduction of stent migration from 100 % to 15 % through OTSC use seems promising. The authors note that prospective data is needed to define an optimal approach for OTSC use in SEMS placement.

Dr. Mudumbi, Dr. Velazquez-Aviña of the Basil I. Hirschowitz Endoscopic Center of Excellence, Division of Gastroenterology and Hepatology, University of Alabama at Birmingham, United States cooperated with Dr. Baig and Dr. Mönkemüller, Department of Medicine, University of Erlangen, Nürnberg, Germany as well as Dr. Neumann, affiliated with both institutions, in a single-center, retrospective cohort study of SEMS anchoring with OTSC and subsequent removal of clip and stent with an inject-and-resect technique. The study covered a total of 12 patients (8 male, 4 female) at a median age of 57 years (range: 45–72 years). Indications included different nonstricturing benign or malignant esophageal diseases (tracheoesophageal fistula, postoperative leaks, and esophageal perforation). Application and initial anchoring of the OTSC clip was successful in all cases, in two cases clip dislodgment and subsequent stent migration was documented during follow-up. After complete healing, stent and clip were removed in six patients, while the stent was left indefinitely in four patients to treat their underlying condition.

The authors emphasize the potential to significantly reduce stent migration rates currently at up to 40 % in esophageal stent placement by using OTSC as an anchoring device. The inject-and-resect technique used for removal was

successful in all cases in which the underlying condition had been healed, although the authors caution that a possible risk of perforation is connected with anchoring the clip to deeper tissue. While the group believes benefit to be greater for nonstenosing diseases, they point out that patients with malignant stricture receiving chemo-radiation may also benefit from clipping, as a reduction in the size of the tumor may also lead to stent migration. They also mention ex-vivo trials showing that sutures are most resistant to tensile forces (average of 20.4 newton needed for stent removal) compared to clip-anchored stents (16 newton on average) and unanchored stents at 4.0 to 5.6 newton. The paper identifies clips as the least expensive device, but advises that costs associated with the possible use of multiple clips or stent migration should also be taken into account. In conclusion, the authors evaluate the use of OTSC for anchoring of fully covered SEMSs as an easy and safe avenue of treatment with the potential to significantly reduce stent migration rates and call for further studies to assess and refine the technique.

* The fixation of stents is not a common indication for the OTSC System and there is only very limited experience.

Preventing migration of fully covered esophageal stents with an over-the-scope clip device (with videos).

Irani S, Baron TH, Gluck M, Gan I, Ross AS, Kozarek RA
Gastrointest Endosc. 2014 May;79(5):844-51

Epub 2014 Jan 25

Anchoring of self-expanding metal stents using the over-the-scope clip, and a technique for subsequent removal

Mudumbi S, Velazquez-Aviña J, Neumann H, Kyanam Kabir Baig KR, Mönkemüller K

Endoscopy. 2014 Dec;46(12):1106-9

Epub 2014 Sep 30

Dr. Irani et al. videos: To watch videos of the procedures performed by the group around Dr. Irani, please visit the website of Gastrointestinal Endoscopy at <http://dx.doi.org/10.1016/j.gie.2013.12.012> (behind a paywall)

January 2015 | OTSC® clip for closure of pancreatocolonic fistulas – new case studies

K. Ito, Y. Igarashi, T. Mimura et al., Division of Gastroenterology and Hepatology, Toho University Omori Medical Center, Tokyo, Japan, published a case study on the successful OTSC closure of a colonic fistula complicating severe pancreatitis. Dr. E.C. Gorospe, Dr. S. Desai, Dr. B. Al-Bawardy et al., Division of Gastroenterology and Hepatology, Mayo Clinic, Rochester, Minnesota, USA, describe the clip closure of a pancreatocolonic fistula caused by severe necrotizing pancreatitis. Y. Koike, T. Kudo, T. Shigesawa et al., Department of Gastroenterology, Sapporo City General Hospital, Sapporo, Japan, presented the closure of a colonic fistula complicating a pancreatic pseudocyst.

The Tokyo case was a 44-year-old male with hyperlipidemic acute pancreatitis including an abscess in the left abdominal cavity and immense peripancreatic fluid collection. Conservative treatment resulted in middling success. A disruption on the tail of the main pancreatic duct was suspected and confirmed after 90 days and treated with an ENPD tube and a pancreatic stent, which proved ineffective. After conservative management options had been exhausted, surgical therapy was considered, but postponed due to presence of E. coli and MRSA in abscess culture.

Finally, an OTSC clip was used on day 148 to endoscopically seal the fistula. A follow-up indicated complete sealing of the leak and improved healing. After the patient had been upgraded to a full diet, an endoscopic pancreatic stent was placed to deal with increased percutaneous drainage. Several follow-ups showed the success of the procedure with improvement of the abscess, clip in situ and healed perforation site with no signs of inflammation, ulceration or pancreatic duct disruption. Patient is well and now receives outpatient care for hyperlipidemia. The Mayo Clinic reported the case of a 69-year-old female with necrotizing pancreatitis, who had had 3 previous transgastric necrosectomies. When she was hospitalized 3 weeks after her last necrosectomy, there was reflux of fecal-like material into the debrided cavity as well as a sigmoid structure, likely caused by pancreatic necrosis

and pancreatic secretions directly into the colon through fistulae. Two fistulae were located and closed from a colonic approach using OTSC clips. Closure was confirmed fluoroscopically and endoscopically. After 7 months, the patient remained asymptomatic.

The team from Sapporo City General Hospital reported the case of a 53-year-old man with a history of alcohol-induced chronic pancreatitis. A fistula between a pancreatic pseudocyst and descending colon did not respond well to traditional endoscopic drainage and was reinfected. It was also believed to be the cause of repeated cyst infections. The fistula was finally closed using the OTSC clip and closure was confirmed via radiographic imaging. The patient's pancreatic pseudocyst has decreased in size with no signs of reinfection. In cases where endoscopic drainage alone cannot ensure colonic fistula closure, the OTSC System is an interesting treatment option since it is less invasive than surgery.

Severe Acute Pancreatitis with Complicating Colonic Fistula Successfully Closed Using the Over-the-Scope Clip System

Ito K, Igarashi Y, Mimura T, Kishimoto Y, Kamata I, Kobayashi S, Yoshimoto K, Okano N

Case Rep Gastroenterol. 2013 Jul 23;7(2):314-21

Print 2013 May

Over-the-scope clip closure of pancreatocolonic fistula caused by severe necrotizing pancreatitis

Gorospe EC, Desai S, Al-Bawardy B, Baron TH, Buttar NS, Wong Kee Song LM

Gastrointest Endosc. 2014 May;79(5):71

Epub 2013 Dec 12

Pancreatic pseudocyst with complicating colonic fistula successfully closed using the over-the-scope clip system

Koike Y, Kudo T, Shigesawa T, Fujita Y, Endo A, Ono Y, Nakamura M, Nagasaka A, Nishikawa S

Endoscopy. 2014;46 Suppl 1 UCTN:E178-9

Epub 2014 Apr 22

January 2015 | ASGE: Over-The-Scope Clipping device is safe and effective for management of GI defects

ASGE – The American Society for Gastrointestinal Endoscopy issued a press release concerning a publication in its GIE-Gastrointestinal Endoscopy journal:

„An international multicenter study reports that over-the-scope clip (OTSC) placement is a safe and effective therapy for the closure of gastrointestinal (GI) defects, which includes anastomotic leaks, fistulae and perforations. Clinical success was best achieved in patients undergoing closure of perforations or leaks when OTSC placement was used for primary or rescue therapy. The overall clinical success for the closure of perforations and leaks ranged between 90 percent and 73 percent; however, successful closure of fistulae was achieved in less than half of the patients. The type of defect (i.e. perforation or leak) is the best predictor of successful long-term closure. The study appears in the October issue of GIE: Gastrointestinal Endoscopy, the monthly peer-reviewed scientific journal of the American Society for Gastrointestinal Endoscopy (ASGE).“

For Immediate Release:

Media Contact: Anne Brownsey abrownsey@asge.org

American Society for Gastrointestinal Endoscopy

www.asge.org | www.screen4coloncancer.org

AN OVER-THE-SCOPE CLIPPING DEVICE FOR ENDOSCOPIC MANAGEMENT OF GASTROINTESTINAL DEFECTS IS SAFE AND EFFECTIVE

DOWNERS GROVE, Ill – October 23, 2014 – An international multicenter study reports that over-the-scope clip (OTSC) placement is a safe and effective therapy for the closure of gastrointestinal (GI) defects, which includes anastomotic leaks, fistulae and perforations. Clinical success was best achieved in patients undergoing closure of perforations or leaks when OTSC placement was used for primary or rescue therapy. The overall clinical success for the closure of perforations and leaks ranged between 90 and 73 percent; however, successful closure of fistulae was achieved in less than half of the patients. The type of defect (i.e. perforation or leak) is the best predictor of successful long-term closure. The study appears in the October issue of GIE: Gastrointestinal Endoscopy, the monthly peer-

reviewed scientific journal of the American Society for Gastrointestinal Endoscopy (ASGE). Conventional treatment of GI defects is with surgical management, which is associated with significant morbidity and mortality.

Technological advances in endoscopic devices have allowed for endoscopic closure of GI defects. Endoscopic therapies include placement of self-expandable metal stents (SEMSs) and application of clips and sealants, all of which have proven their utility in different clinical scenarios with varying degrees of success.

The most common endoscopic approach for treatment of GI defects before the over-the-scope clip was the use of temporary SEMSs.

A large case series reported a success rate as high as 75 percent with SEMSs; however, this practice was associated with a high rate of adverse events (46 percent). The OTSC provides more durable closure than standard clips because of its wider mouth and ability to grasp larger amounts of tissue.

In addition, full-thickness closure is achievable because of greater compressive force. The current study is the largest to-date to assess outcomes of OTSCs in the management of GI fistulae, perforations and leaks. "The primary goal of this study was to describe a large, international, multi-center experience with OTSCs for the management of GI perforations, fistulae and anastomotic leaks and to determine the overall success of GI defect closure. Secondary goals were to determine success rate by type of defect and type of therapy, primary vs. rescue, and to determine predictors of OTSC success," said study lead author Yamile Haito-Chavez, MD, Johns Hopkins University, Baltimore, Maryland. "Our study found that long-term success was achieved in 60.2 percent of patients. The rate of successful closure of perforations was 90 percent, closure of leaks was 73.3 percent and closure of fistulae was 42.9 percent. Long-term success was significantly higher when OTSCs were applied as a primary therapy."

Methods A retrospective review of consecutive patients who underwent attempted OTSC placement (either as primary or rescue therapy) for the indication of GI leak, fistula or perforation at 16 academic centers in the United States, The Netherlands, Germany, Italy, and Chile was conducted between May 2006 and November 2012. Patients were identified by using endoscopic databases at each institution. Anastomotic leak was defined as disruption at a surgical anastomosis resulting in a fluid collection with or without evidence of extravasation of contrast medium on radiologic evaluation. Fistula was defined as abnormal communication between two epithelialized surfaces. Perforation was defined as an unintentional, acute iatrogenic, full-thickness defect in the GI tract. The main outcome measurement was the long-term success of the procedure.

Results A total of 188 patients (108 fistulae, 48 perforations, 32 leaks) were included. Long-term success was achieved in 60.2 percent of patients during a median follow-up of 146 days. The rate of successful closure of perforations (90 percent) and leaks (73.3 percent) was significantly higher than that of fistulae (42.9 percent). Long-term success was significantly higher when OTSCs were applied as primary therapy (primary 69.1 percent vs. rescue 46.9 percent). On multivariate analysis, patients who had OTSC placement for perforations and leaks had significantly higher long-term success compared with those who had fistulae. In an accompanying editorial, Danny Cheriyan, MB, BCh, MRCP and Robert Enns, MD, FRCP, Division of Gastroenterology, St. Paul's Hospital, University of British Columbia, Vancouver, Canada, state "The future of the OTSC in mainstream endoscopy is promising. It is arguable that these clips should be available in every well-stocked unit, and because application and deployment are similar to those of a standard ligation bander, it would be appropriate for implementation into training and clinical practice in gastroenterology."

December 2014 | Full-thickness resection of adenoma in colonic diverticulum using Ovesco FTRD® System

PV Valli, M Kaufmann and P Bauerfeind, Dept of Gastroenterology and Hepatology, University Hospital

Zurich and B Vrugt, Institute of Pathology, University Hospital Zurich in Switzerland published the first case where colonic adenoma located from a diverticulum, a rare finding, was treated using the FTRD full-thickness resection device. The patient was a 66-year old woman with extensive diverticulosis in the entire colon. Colonoscopy revealed a 10 mm lesion (pathological size: 13 mm) inside a diverticulum in the ascending colon. Using a standard colonoscope with Indian ink injection and a hemoclip, the diverticulum was marked before a therapeutic colonoscope, fitted with the FTRD System, was introduced and advanced to the adenoma, located 10 cm proximal to the hepatic flexure. Adenoma and inverted diverticulum were mobilized into the cap of the FTRD with grasping forceps and additional suction before the FTRD clip was placed. Then resection was performed above the clip with the electrical snare integrated in the FTRD device. Histopathology showed successful full-thickness resection of a tubular adenoma with low-grade dysplasia and the resected diverticulum. The patient received single-shot peri-interventional antibiotic prophylaxis and was kept overnight for observation. Free of pain, the patient was discharged the next day, and no signs of complication arose over a 3 month follow-up period. The authors consider the new FTRD System as a secure treatment option for the resection of high-risk polyps without the risk of leakage of bowel content into the peritoneal cavity and see potential for use in an outpatient setting.

To watch a video of the procedure, please visit the website of Gastroenterology journal at:

x.doi.org/10.1053/j.gastro.2014.07.053

Endoscopic Resection of a Diverticulum-Arisen Colonic Adenoma Using a Full-Thickness Resection Device

Valli PV, Kaufmann M, Vrugt B, Bauerfeind P

Gastroenterology 2014, 147.5: 969-71

221

November 2014 | Clinical presentations confirm efficacy of OTSC® clipping at German Congress of Visceral Medicine in Leipzig

Clinical presentations at German Congress on Visceral Medicine confirm efficacy of OTSC clipping and show clinical data on novel Ovesco products FTRD and remOVE Leipzig, September 17–20, 2014.

The 69th annual congress of the German society for gastroenterology, digestive and metabolic diseases, DGVS, was held under the presidency of Prof. Dr. med. Peter R. Galle. A significant number of presentations provided clinical data on OTSC clipping and confirmed clinical efficacy and safety in the primary indications of the product, hemostasis, closure of acute lesions/perforations and closure of chronic lesions/fistulae (source: www.viszeralmedizin.com).

Conference report

Clinical presentations at German Congress on Visceral Medicine confirm efficacy of OTSC clipping and show clinical data on novel Ovesco products FTRD and remOVE.

Leipzig, September 17–20, 2014: The 69th annual congress of the German society for gastroenterology, digestive and metabolic diseases, DGVS, was held under the presidency of Prof. Dr. med. Peter R. Galle.

A significant number of presentations provided clinical data on OTSC clipping and confirmed clinical efficacy and safety in the primary indications of the product, hemostasis, closure of acute lesions/perforations and closure of chronic lesions/fistulae (source: www.viszeralmedizin.com).

Reports on Ovesco OTSC

OTSC proves to be preferable treatment option for several indications.

Glitsch A, Schreiber A, Boldt J, Keßler W, and Mayerle J, Greifswald, reported about a cohort of 46 patients treated with OTSCs clips. The indications include postoperative anastomotic insufficiency (n=13), bleeding (n=7), perforations (n=15), pancreatic fistulae in the colon (n=4), fistulae in patients with inflammatory bowel disease (n=5) and OTSC use after ESD (n=2). Successful closure and complete healing was achieved for all but one indication. In the case of postoperative anastomotic insufficiencies 3 out

of 13 insufficiencies could not be closed successfully (76.93 % success rate). These patients had to undergo further laparoscopic treatment. In all other cases (n=43), no complications were observed and thus no further treatment was necessary. The authors declare OTSC to be a procedure with significantly lower morbidity and mortality in comparison to conventional treatment options and emphasize that it spares patients elaborate and more complication-prone methods of treatment.

Retrospektive Auswertung der OTSC Anwendung an der Universitätsmedizin Greifswald

Glitsch A, Schreiber A, Boldt J, Keßler W, Mayerle J

OTSC for stopping acute bleeding in the gastrointestinal tract.

Braun A, Freiburg, and Kirschniak A, Tübingen, presented data about a total of 16 patients (median age=75.5 years, R=61–92 years; m=9, f=7) over three years with acute bleeding, who were treated with OTSC application during emergency endoscopy. 8 procedures were performed in the upper and 8 in the lower gastrointestinal tract. Patients with upper GI bleeding were given a highly dosed proton pump inhibitor (80mg i.v. Bolus, 320 mg i.v./ 24h). Hemorrhages were classified as follows: Forrest Ia (n=7), Forrest Ib (n=7), and Forrest IIa (n=2). All patients suffered an acute drop in hemoglobin and showed definite signs of bleeding. No further local therapies were administered. All OTSC applications were performed by the same clinician and took 20 minutes or less.

OTSC application and thus primary hemostasis was successful in all cases. None of the patient suffered recurrent gastrointestinal bleeding. 6 patients (4 Fl a, 1 Fl b, 1 Fl a) had a follow-up endoscopy between day 1 and 7; all clips were in-situ with no signs of bleeding. The other 7 patients were not reexamined due to good response to treatment. The authors see OTSC as a safe and very effective treatment option in emergency endoscopy. Primary hemostasis is possible for a large percentage of patients, which improves lethality, and examination time is low.

Endoskopische Behandlung von akuten Blutungen mit einem Over-The-Scope Clip (OTSC)

Braun A, Kirschniak A

Interdisciplinary treatment regime for thoracic anastomosis insufficiencies.

Pauthner M, May A, Lorenz D, and Ell C, Offenbach, introduced the complication management regime for thoracic anastomosis insufficiencies (AI) at HSK Wiesbaden hospital. From 07/2000 to 12/2013, they counted 632 cases of resections in the esophagus, 557 of which included transthoracic esophageal resections with intrathoracic anastomosis of a gastric sleeve. 49 of these 557 patients (8.8 %) suffered from confirmed AI.

Of these cases, 13 (26.5 %) were treated conservatively with a triple-lumen jejunal feeding and gastric decompression tube (TLT), 14 (28.6 %) were treated with a primary stent, in 7 patients an OTSC clip was placed endoscopically (14.3 %) and 2 (4.1 %) received transluminal vacuum therapy. 12 patients (26.5 %) had to undergo repeat thoracotomy, 7 of which had a stent placed during the procedure (14.3 %). Hospital lethality after all esophageal resections was 3.6 %, with only 2.3 % in the last 4 years. If gastric sleeve shows good blood circulation and the AI is small, either a clip or a TLT is placed. Routine check after 36–48h. Larger AIs are treated with stents; routine check after 24 h and placement of TLT to protect stent from bile. If pleural empyema occurs, a repeat thoracotomy is performed, including decortication and sewing-over of the AI as well as stent placement; stent is fixated with an absorbable suture. The authors report that this indicationspecific, standardized complication management regime reduced lethality after occurrence of thoracic AI from 14.3 % to 3.1 %.

Therapie der Anastomoseninsuffizienz nach Ösophagus-Resektion – die viszeralmedizinische Herausforderung

Pauthner M, May A, Lorenz D, Ell C

New treatment option for chronic, therapy-resistant esophageal-bronchial fistulae.

Wedi E, Sportes A, and Hochberger J, Strasbourg, France, presented the case of a 68-year-old patient with a chronic esophageal-bronchial fistula. In early 2010 he presented

with haemoptysis and an unidentified pulmonary lesion. Examination of the mediastinum and subsequent exploratory thoracotomy showed giant-cell granuloma with no indication of malignancy. In December 2011 patient reported repeated incidents of coughing during food intake with recurrent bronchio-pulmonary infections, and an esophageal-bronchial fistula was diagnosed. Initial treatment included a fully covered Nitinol stent (23/18 mm wide, 12 cm long), which dislocated two days after food was reintroduced. In the following 1.5 years, a variety of endoscopic treatment options were explored (partially covered stent, fibrin glue, standard hemoclips, etc.), but all proved unsuccessful while a 5–7 mm wide fistula tract had formed. In collaboration with pulmonologists and thoracic surgeons, an experimental course of treatment was employed. Deep tissue in the fistula tract was excised using endoscopic submucosal dissection (ESD). Then a flexible bronchoscopy with APC and chafing of the fistula tract with a brush was performed. Finally, the fistula tract was closed with a 17.5 mm OTSC macro clip. Preliminary endoscopic-radiological follow-up over a 4 month period has shown no recurrence. The authors conclude that excision of the fistula and subsequent closure with an OTSC macro clip is a promising new treatment option, which should be further evaluated.

Therapie-refraktäre ösophago-bronchiale Fistel – Was tun, wenn alle Therapieoptionen versagen?

Wedi E, Sportes A, Hochberger J

Report on Ovesco FTRD (Full-Thickness Resection Device), newly launched by Ovesco Preliminary clinical experience with the FTRD system in the lower gastrointestinal tract.

Schmidt A, Damm M and Caca K, Ludwigsburg, together with Gubler C and Bauerfeind P, Zurich, Switzerland, reported their experience with endoscopic full-thickness resection in the lower GI tract of 21 patients from July 2012 to March 2014. Resection was always performed using the FTRD system mounted onto a standard endoscope. Indications included recurrent or incompletely resected adenoma with non-lifting sign (n=9), adenoma with high-grade prostatic intraepithelial neoplasia (HGPIN) (n=1), adenoma on base of appendix (n=3), broad-based adenoma in patient with coagulation disorder (n=1), diagnostic (re-)resection in patients with T1 carcinoma (n=3), adenoma on diverticulum (n=1), a submucosal tumor (n=2), a diagnostic FTR for a patient with suspected Hirschsprung's disease (n=1). Lesions were located in the cecum (3), ascending colon (4), transverse colon (2), descending colon (4), sigma (2), recto-sigmoid (3), and rectum (3). Navigation to target lesion with FTRD mounted onto endoscope was possible in all but one case (95.2 % success rate).

Once the lesion was reached, resection was technically successful in all cases and macroscopically complete in 19 out of 20. Histological findings confirmed complete full-thickness resection in 17/20 cases (85 %). No perforations or relevant bleeding occurred. Two patients developed postpolypectomy syndrome (PPCS) after resection in the cecum, which was treated conservatively and with success. The authors conclude that endoscopic full-thickness resection in the lower GI tract with the FTRD System is technically feasible, effective and safe. Larger studies are necessary for further evaluation of this technique.

Endoskopische Vollwandresektion im unteren GI-Trakt mit dem FTRD® System: eine retrospektive Studie Schmidt A, Damm M, Gubler C, Caca K, Bauerfeind P

Report on OTSC Proctology

Prospective study deems OTSC Proctology a very promising method for anorectal fistula closure.

Probst R presented the experiences of a prospective pilot study at St. Anna hospital, Stuttgart, and edz center of excellence in proctology, Mannheim, regarding the use of OTSC Proctology. The study included 20 patients (14 male, 6 female), aged 56.1 years on average (R 25–73 years). There were 14 transsphincteric and 6 suprasphincteric anorectal fistulae. Average procedure time for clipping of fistulae was 32 minutes (R 17 to 66 minutes). There were no intraoperative complications. Follow-up endoscopy after

six months or more showed proper healing in 18 of 20 cases (90 % success rate). The clip had remained in-situ in 13 patients. The clips fell off spontaneously (10 to 4 days post-op) in three patients. The clip was surgically removed in two patients (clip dislocation, severely impaired wound healing). The fistula persisted/recurred in two of 20 cases (10 %). Reasons were spontaneous clip displacement on day 3 post-op and failure to heal. In the second case, the clip was removed and another loop was placed.

The author pointed out that, since this is a new method, clinical experience is limited and efficiency with regards to ano-/recto-vaginal fistulae is unclear. The report claimed OTSC Proctology to be a very promising new treatment option with lower morbidity and fewer complications than established methods. It prevents traumatizing the sphincter, eliminates the risk of postoperative incontinence, and does not impede further treatment options. Additionally, patient satisfaction is high and the instrument is easy to use.

OTSC ohne Endoskop: Anorektaler Fistelverschluss mittels OTSC Proctology

Probst R

(cf. Probst RL, Joos AK, Ehni W, Bussen D, Herold A Prospective pilot study of anorectal fistula closure with the OTSC | Proctology. Colorectal Dis. 2014 Sep 1. Epub)

Report on remOVE System

(product approval in Europe pending)

Early clinical data on new bipolar DC cutting instrument for OTSC clips show efficacy and safety.

Rische S, Schmidt A, Damm M, Cahyadi O, Bauder M, and Caca K, Ludwigsburg, summarized a retrospective study of compassionate use cases, which used a prototype of the remOVE System to remove OTSC clips in 13 patients. Clips remained in the upper and lower gastrointestinal tract for an average of 70 days (R 7–469 days). Fragmentation of clip was successful in all cases. Mean procedure time was 45 minutes (R 35–75 minutes). Endoscopic removal of clip fragments was possible in all but one case (92.31 % success rate). No serious complications were observed. The authors conclude that OTSC clip removal with the remOVE System prototype is easy, fast and safe, and is thus suited for elective procedures as well as endoscopic emergencies.

Endoskopische Entfernung von Over-The-Scope Clips (OTSC) mit Hilfe eines bipolaren Schneideinstruments (DC ClipCutter): retrospektive Fallserie von 13 Patienten Rische S, Schmidt A, Damm M, Cahyadi O, Bauder M, Caca K

October 2014 | ESGE Position Statement on the treatment of endoscopic perforations: OTSC® endorsed as preferential device.

Under the lead of GA Paspatis, Gastroenterology Department, Benizelion General Hospital, Heraklion, Greece, a group of experts analysed best clinical practice in the treatment of iatrogenic endoscopic perforations.

Iatrogenic perforation of the gastrointestinal tract related to diagnostic or therapeutic endoscopy is a severe adverse event, associated with significant morbidity and mortality. In the Position Statement, ESGE aimed to define the main risk factors for iatrogenic perforations as well as clear diagnostic and therapeutic algorithms for their management.

The authors allude to the OTSC Clip by Ovesco Endoscopy as a first line therapy for the use in different GI organs. For esophageal perforations that are associated with particularly high risk of secondary complications OTSC is recommended for lesions between 1 and 2 cm. In gastric perforations the authors cite the high therapeutic efficacy of OTSC published in the literature. OTSC is the device of choice for closing gastric defects larger than 10 mm. The use of conventional clips is not recommended in such lesions. Also combinations of such clips and endoloops is only appropriate in case no OTSC clip should be available. For colonic perforations the guideline states that OTSC is the preferential device for larger lesions. Only small holes can be treated with conventional clips.

For the complete list of recommendations by ESGE please check the resources linked at the bottom of this news message.

The position paper is an official statement of the European

Society of Gastrointestinal Endoscopy. ESGE guidelines are intended to be an educational device to provide information that may assist endoscopists in providing best care to patients.

Diagnosis and management of iatrogenic endoscopic perforations: European Society of Gastrointestinal Endoscopy (ESGE) Position Statement

Paspatis GA, Dumonceau JM, Barthet M, Meisner S, Repici A, Saunders BP, Vezakis A, Gonzalez JM, Turino SY, Tsiamoulos ZP, Fockens P, Hassan C

Endoscopy. 2014 Aug;46(8):693-711 doi: 10.1055/s-0034-1377531

Epub 2014 Jul 21. Free paper (download German version) ESGE: www.esge.com/esge-guidelines.html

September 2014 | OTSC® clip for closure – new case studies on alimentary tract fistulas

Dr. T. H. Kothari and Dr. G. Haber, Division of Gastroenterology, Lenox Hill Hospital, New York, USA published a three-case series on closure of fistulas in the alimentary tract, while Dr. T. Meister, Dr. J. Kuhlitz and Dr. M. Floer, Helios-Albert-Schweitzer-Klinik, Northeim, Germany published a case report on closure of a postoperative enterocutaneous fistula.

The first New York case was an 80-year-old female with dysphagia, history of aortic valve regurgitation, esophageal perforation due to TEE, repair with feeding jejunostomy and venting gastrostomy. Healing resulted in a gastrocutaneous fistula, which persisted for 9 months despite PPI therapy. Tissue was grasped with OTSC Anchor and pulled into the cap, then an OTSC clip was successfully deployed. Follow-up after 3 months showed that the fistula had fully healed.

The second case was a 36-year-old female with a history of diverticulosis, presenting with passage of feculent material from her vaginal canal and recurrent urinary tract infections. CT scan revealed a fistula between sigmoid colon and vagina, but locating it gastroscopically was difficult. The OTSC clip was deployed, resulting in good tissue entrapment. The patient was symptom-free for several weeks. When symptoms recurred, surgery revealed an abscess communicating with the colovaginal fistula, which prevented healing.

The third case was a 41-year-old female with colonic interposition after lye ingestion and PEG tube placement. After the tube was removed, gastrostomy tract did not close for several months and conventional methods of closure failed. The fistula tract presented with some exudate at the gastric orifice. After several attempts to draw sufficient tissue to the cap with the OTSC Anchor, the OTSC clip was successfully deployed. After a few weeks, patient started having secretions. It was hypothesized that the diameter of the fistula (> 1 cm) was to blame for inefficient healing.

The German case report was about a 48-year-old female suffering from an enterocutaneous fistula for four months, leading to malnourishment. Prior attempts to close the fistula with fibrin glue had failed. The OTSC clip was placed onto the fistula opening under continuous suction and success fully deployed. Follow-up after 12 months showed continued success.

The over-the-scope clip system – a novel technique for gastrocutaneous fistula closure: the first North American experience

Kothari TH, Haber G, Sonpal N, Karanth N

Can J Gastroenterol. 2012 Apr;26(4):193-5

75

Over-the-Scope Clip (OTSC) application as rescue therapy for postoperative enterocutaneous fistula closure

Meister T, Kuhlitz J, Floer M

Acta Chir Belg. 2014 Jan-Feb;114(1):87-9

186

August 2014 | Management of postoperative anastomotic defects: OTSC® System as preferred treatment option

Two case reports recently published by Dr. Tontini and colleagues, Dept. of Medicine I, University of Erlangen-Nuernberg, Germany and by Dr. Chen and colleagues, Dept. of Gastroenterology, Riverside, University of California, USA illustrate the complete closure of larger anastomotic leaks with the OTSC System when other techniques have failed or deemed unsuitable. A 69-year old man had an Ivor-Lewis esophagectomy due to esophageal

carcinoma. After the development of chest pain 6 days later a CT scan revealed an anastomotic dehiscence. During an EGD a metal stent was placed but a persistent defect was found 10 days later by a CT. An additional stent was placed overlapping the first stent. Because of a continued leakage another EGD was performed demonstrating a gastric conduit fistula. After application of conventional endoscopic clips the patient was discharged but presented with worsening symptoms. Another EGD showed a persistent fistula that was finally closed with an OTSC clip resulting in a complete healing.

The other patient (71-year old woman) presented with hypotension, melena and low hemoglobin level 3 weeks after a Billroth I gastroenteral anastomosis. EGD showed an oozing bleeding and as well a defect at the anastomosis. In this case the dehiscence extended over half the circumference of the anastomosis. Since other techniques seemed inappropriate due to large leak, massive bleeding and difficult target position it was decided to use the Twin Grasper and an OTSC clip to close the defect. Complete closure was confirmed by a subsequent endoscopic examination. According to the authors, the OTSC clip should be considered as the first choice for sealing of intermediate leaks.

Successful over-the-scope clip (OTSC) treatment for severe bleeding due to anastomotic dehiscence

Tortini GE, Naegel A, Albrecht H, Vieth M, Vecchi M, Neurath MF, Neumann H

Endoscopy. 2013;45 Suppl 2 UCTN:E343-4 146

Over-the-scope clip for closure of persistent post-esophagectomy gastric conduit fistula

Chen AI, Lim BS, Ma JS, Chaya CT

Gastrointest Endosc. 2014 Apr;79(4):546 147

July 2014 | Retrospective study confirms safety and effectiveness of OTSC® in the endoscopic treatment of GI bleeding, perforation and fistula

Dr. Vijay Jayaraman and colleagues, Cedars Sinai Medical Center, Los Angeles, recently presented a retrospective study on their experience with the OTSC System in the treatment of GI bleeding, fistula and perforation. Their case series consisted of 24 consecutive patients treated between January 2011 and April 2012 (mean age 70 years) included the following indications for OTSC placement (28 clips): postsurgical enterocutaneous fistula (n=10), spontaneous perforation (n=1), anastomotic leak (n=4), perforation after mucosal resection (n=3), prophylactic closure of mucosal defect after EMR (n=1), postpolypectomy bleeding (n=2), postendoscopic perforation (n=2), tracheoesophageal fistula (n=1) and leakage from a percutaneous jejunostomy site (n=1). Instruments or modalities used to grasp the tissue were dedicated devices (OTSC Twin Grasper and OTSC Anchor) in 16 and nondedicated devices (rat tooth/alligator forceps or suction alone) in 15.

Median follow-up time was 2.9 months; mean defect size 10 mm (range 5–2 mm). The overall success rate was 61 %. In their experience the success rate of closure of an acute defect is higher compared to chronic fistula. 9 out of 24 lesions were chronic (>1 month) in this series which might explain the lower overall success rate in comparison to the literature (72–100 %).

Furthermore, a trend towards higher success rate was noted in defects <10 mm compared to defects >10 mm. No patient reported any complications associated with OTSC placement.

Endoscopic therapy is still the initial choice before any surgical intervention to manage GI bleedings, fistulae, perforations and leaks. As through the scope clips are limited by their smaller wing span and low force of closure leading to suboptimal results, the OTSC clip provides a safe and effective endoscopic alternative.

Clinical Application and Outcomes of Over the Scope Clip Device: Initial US Experience in Humans

Jayaraman V, Hammerle C, Lo SK, Jamil L, Gupta K

Diagn Ther Endosc. 2013;2013:381873 137

OTSC® Update 17

June 2014 | Spanish researchers receive award for successful OTSC® case presented at the National Digestive Congress Spain, June 14–16, 2014 in Valencia

D. López Peñas and colleagues, Servicios de A. Digestivo, Otorrinolaringología y Oncología médica, Hospital de Llerena, Spain received an award for their presentation at the Semana de las Enfermedades Digestivas (SED 2014) in Valencia.

The researchers report on a successful closure of a pharyngo-cutaneous fistula after total laryngectomy in a 58-year old patient. The 3-stage treatment of dilatation, subsequent percutaneous gastrostomy and closure with an OTSC clip implicated a substantial improvement of food supply and quality of life.

Here you can see the presentation (in Spanish):

http://www.ovesco.com/fileadmin/Downloads/Flyers_for_events/2014_National_Digestive_Congress_Spain.pdf

June 2014 | Three case reports on surgery-sparing uses of the OTSC® clip in multiple indications

Three different case reports lately published by Dr. V. Gómez et al., Dept. of Gastroenterology and Hepatology, Mayo Clinic, Jacksonville, USA, Dr. S. Singhal et al., Div. of Gastroenterology, The Brooklyn Hospital Center, New York, USA and Dr. J. Albert, Center of Internal Medicine, Johann Wolfgang Goethe University Hospital, Frankfurt/Main, Germany illustrate the broad spectrum of indications for which placement of OTSC clips can be useful.

The first case report describes the use of the OTSC System in the management of a Dieulafoy lesion. A 74-year old man suffered from a recurrent, obscure, life-threatening gastrointestinal bleeding. EGD revealed a non-bleeding Dieulafoy lesion at the lesser gastric curvature. Due to the large size and difficult position of the lesion, conventional through-the-scope clips were not used, but an OTSC clip was successfully deployed.

In another case report a 61-year old woman presented for EGD for evaluation of dysphagia. Four arteriovenous malformations were found in the duodenum, which were cauterized. On withdrawing the endoscope, a 2-cm gastric perforation was identified on the lesser curvature. Using the suction technique an OTSC clip was applied to close the defect.

The third paper presents the case of a patient with severe bleeding from a duodenal ulcer that could not be controlled by conventional clips and injection of fibrin glue. Angiographic placement of coils into the afferent vessel then stopped the bleeding. After 3 days a fistula penetrated into the dorsal duodenum leading to a peritoneal leakage. Successful closure of the fistula was achieved with an OTSC clip. All the authors agree that the OTSC System is an effective tool for endoscopic control of bleedings, perforations and fistulas.

Novel treatment of a gastric Dieulafoy lesion with an over-the-scope clip

Gómez V, Kyanam Kabir Baig KR, Lukens FJ, Woodward T

Endoscopy. 2013;45 Suppl 2 UCTN:E71. 129

Endoscopic closure of gastric perforation using over-the-scope clip: a surgery-sparing approach

Singhal S, Atluri S, Changela K, Gupta SS, Krishnaiah M, Anand S

Gastrointest Endosc. 2014 Jan;79(1):23 150

Closure of an Ischemic Duodenal Fistula with an Over-The-Scope Clip

Albert JG

Video Journal and Encyclopedia of GI Endoscopy. 2013 June, pages 219–20 155

June 2014 | Prophylactic closure of large mucosal defects after colorectal ESD significantly reduces the inflammatory reaction and abdominal symptoms of patients with neoplasms

Fujiyama et al., Department of Gastroenterology and Neurology, Faculty of Medicine, Kagawa University,

Kagawa, Japan, assessed the efficacy and safety of a prophylactic closure for large mucosal defects after colorectal ESD.

From April 2010 to December 2012, 68 patients with colorectal tumors were treated with ESD. The prophylactic closure using a conventional clip and the over-the-scope clip (OTSC) system was indicated for patients with excessive coagulation in the muscularis propria or larger resection size. The closure group reduced the peritoneal inflammatory reaction and abdominal symptoms without increasing complications. The closure group also had a significantly lower WBC count (post operative day 1), CRP (post operative day 4) and abdominal pain after colorectal ESD compared to the non-closure group. Perforation occurred in 1 case, and postoperative bleeding in 2 cases, with only 1 bleeding case needing an emergency endoscopy in the non-closure group. One perforation case needed emergency surgery because the endoscopic treatment was ineffective. Without increasing adverse effects, the prophylactic closure efficiently reduced the inflammatory reaction and abdominal symptoms of colorectal ESD in patients with large superficial colorectal neoplasms.

The efficacy and safety of prophylactic closure for a large mucosal defect after colorectal endoscopic submucosal dissection

Fujiyama S, Mori H, Kobara H, Nishiyama N, Kobayashi M, Rafiq K, Masaki T

Oncol Rep. 2013 Jul;30(1):85-90 | [Epub 2013 May 14] 129

May 2014 | Sleeve gastrectomy leaks: Closure with the OTSC® System

Sleeve gastrectomy is increasing in popularity for the treatment of morbid obesity. The most serious and dreaded complication of this procedure is an anastomotic leak typically at the gastroesophageal junction.

Dr. Ahmad Aly and colleague, Upper GI & Bariatric Unit, Austin Hospital, Heidelberg, Australia present two case reports on managing a sleeve leak with the OTSC System. A 58-year old woman with a BMI of 45 underwent sleeve gastrectomy without intraoperative incident. In the case of a 44-year old woman a conversion from laparoscopic adjustable gastric band to a sleeve gastrectomy was performed.

To prevent leakage from the resection line, Seamguard®, a staple line reinforcement product was used in both cases. After initial recovery both patients presented with abdominal pain and fever (8th/30th postoperative day) and a CT scan confirmed leaks at the gastroesophageal junction. Intravenous antibiotics and nutritional support were instituted and fluid collections drained percutaneously and laparoscopically. In the case of the 58-year old woman conservative management was continued for 6 weeks, but the leak persisted. Therefore, it was decided to use the OTSC System. By applying an OTSC clip complete closure was achieved in both patients. After 6 and 8 months respectively, there was no evidence of further leaks and inflammatory markers remained normal. As spontaneous closure of a gastric staple line fistula is rare, many patients require further complex surgery for definitive closure. The OTSC System has the potential to significantly simplify the management of leaks after sleeve gastrectomy by offering a simple endoscopic solution.

The use of over the scope clip (OTSC) device for sleeve gastrectomy leak

Aly A, Lim HK

J Gastrointest Surg. 2013 Mar;17(3):606-8 101

April 2014 | Multipurpose use of the OTSC® System to treat endoluminal gastrointestinal disorders

Recently Mönkemüller et al. from Birmingham, AL, USA report the analysis of an observational retrospective case series of 16 patients (median age 65.8 years) with mixed indications for the treatment with the OTSC System. The overall success rate of 75 % is well in line with other reports and with the meta-analyses of Weiland et al. with a 71 % success rate in fistulas and anastomotic leaks, 79 % in acute perforations, and 88 % in acute GI hemorrhages.

The range of indications included gastrointestinal bleeding (n=6), gastrocutaneous fistulas (n=3), esophagotracheal

and/or esophagopleural fistulae (n=3), resection of submucosal tumor (n=2), stent fixation (n=1), and anastomotic leak after esophagectomy (n=1). The overall per case success rate was 70 % (14 of 20 applications). Mean follow-up was 10 months (range 1–10). There were no complications (0 %) related to endoscopy, sedation or application of the clipping device.

The authors pointed out in the discussion that OTSC allows for the entrapment of a larger amount of tissue, allowing closure of fistula holes and, as shown in these cases, hemostasis superior to other devices. In their critical remarks they also discuss situations where they experienced certain limitations to the system such as the tubular structure of the esophagus which at times might impede an adequate apposition of the device.

Comment Ovesco: especially in cases where the apposition of the OTSC System might seem difficult, the OTSC Anchor is usually a very useful device to facilitate the successful application of a clip with the Anchor functioning as guide wire for both scope and System (e.g. esophagus, cardia, postpyloric duodenum).

In essence the authors draw a very positive conclusion stating "that the OTSC device is ideally suited to treat soft tissue leaks or fistulizing lesions and high-risk bleeding lesions such as ulcers in the posterior duodenum or Dieulafoy's lesions" with the main underlying mechanism being compressing the surrounding tissue around the vessel. They continue "...The OTSC device may become a better device to treat bleeding ulcers located in difficult positions because of its barrel-shaped transparent cap design which allows it to suction the bleeding lesion. It is well known that these bleeding ulcers and lesions are of a higher risk and also more difficult to treat because of their awkward location and/or position...". This statement is followed by an elaborate discussion of the shortcomings of alternative devices. It is important to underline also that the authors support "... multiple OTSC applications in a single session..." as sometimes being useful and allowing approximation of tissue to facilitate subsequent closure. "Interestingly, the device does not tear tissue, as it snaps it together. So far, there have been no reports of GI wall tearing..."

Finally the authors discuss the issue that once OTSC is deployed it cannot be removed easily, and report of two methods they have been using in this case: the "wire technique" as described by Mönkemüller et al., and the use of an Nd-YAG laser, as described by Fährndrich et al. earlier.

Comment Ovesco: we are aware of this issue and are currently finalizing the development of an own, easy to use clip cutter.

Multipurpose use of the 'bear claw' (over-the-scope clip system) to treat endoluminal gastrointestinal disorders
Mönkemüller K, Peter S, Toshniwal J, Popa D, Zabielski M, Stahl RD, Ramesh J, Wilcox AM
Dig Endosc. 2014 May;26(3):350-7

133

April 2014 | Conference report | 44th DGE-BV Congress, Hamburg

The 44th DGE-BV Congress of the German Society for Endoscopy and Imaging Procedures/Diagnostics was held in Hamburg, April 3–5, 2014 under the presidency of Prof. Dr. Thomas Rösch.

Again a significant number of both oral presentations and posters have been featured at this year's event. In summary they all reported their mostly positive experiences with the OTSC® System in all main indications. In addition, our products were featured in several hands-on courses alongside the conference (Chairs: Hochberger J., Maiss J., Kraus F.). Ovesco presented their new products, the DC Clip Cutter and the FTRD® device which are both due to be launched later this year. The reaction of the medical world was more than promising.

• Neue Clips für Blutung und Verschlussstechniken
Caca K, Ludwigsburg, Germany
K. Caca gave a talk on "New tools for the treatment of GI-hemorrhage and perforation". Even though also mentioning other devices he mainly elaborated on the OTSC System. In his summary of clinical cases his take home message was: "the OTSC device achieves hemostasis more quickly

than all other devices and is more effective particularly regarding acute, difficult and heavy bleedings." For the treatment of perforation OTSC was the standard choice. Also, he showed first experiences with the all new DC Clip Cutter device as an important tool for removing the OTSC which will be launched later this year.

• Update Endoskopie – meine Toppapers

Häfner M, Vienna, Austria

M. Häfner updated the plenary session on important recent papers on GI hemorrhage. There he cited two papers by Manta et al. (2013) and Chan et al. (2014) where OTSC had proven to be safe, effective and efficient also in severe bleeding when other procedures had already failed.

Over-the-scope clip (OTSC) represents an effective endoscopic treatment for acute GI bleeding after failure of conventional techniques

Manta R, Galloro G, Mangiavillano B, Conigliaro R, Pasquale L, Arrezzo A, Masci E, Bassotti G, Frazzoni M

Surg Endosc. 2013 Sep;27(9):3162-4
doi: 10.1007/s00464-013-2871-1 [Epub 2013 Feb 23]

114

Use of Over-The-Scope Clip for treatment of refractory upper gastrointestinal bleeding: a case series

Chan SM, Chiu PW, Teoh AY, Lau JY

Endoscopy. 2014 May;46(5):428-31
doi: 10.1055/s-0034-1364932 [Epub 2014 Feb 6]

162

• Clip-Karussell

Groth S, Hamburg, Germany

S. Groth elaborated on the endoscopist's option once it comes to use clips. Interesting enough he exempted the OTSC from the rest of all products stating that OTSC is playing in a different league.

(Comment by Ovesco: the comparator of OTSC is surgery!)

• Techniken zum Perforationsverschluss

Fritscher-Ravens A, Kiel

A. Fritscher-Ravens underlined in her talk on techniques of perforation closure the importance and advantages of the use of the Twin Grasper. Other than that she referred to OTSC as standard treatment.

Altogether five posters were dealing with OTSC:

• Over-the-Scope Clip System (OTSC) – One Therapy for Safety Closure

Leonhardt K, Ohse A, Bauer B, Repp M, Altenburg, Germany

Germany report their 3.5-year experience with our system regarding the three major indications: hemorrhage, acute perforation, and chronic fistula/anastomotic leakage where they achieved a 85.7 %, 84.6 %, and 60 % success rate.

33 patients were included in this retrospective analysis. Average age was 69 years (41–92 ys). Three patients received two clips at once. Across the GI tract the number of patients was equally distributed, except for Jejunum and ileum with only one patient each. The authors conclude that OTSC is a useful and effective tool for the endoscopist sparing the surgeon in many cases.

• The OTSC System in the treatment of a perforation of the colon complicated by adhesion of small bowel

H. Albrecht et al., Erlangen, Germany presented a case report. During diagnostic colonoscopy they experienced an acute perforation which was as usual closed with an OTSC clip. After a few days the patient developed acute symptoms with free sub-diaphragmatic air which led to a sigmoidal resection and overstitch of the small bowel. The authors point out the possible risks of the use of OTSC.

Comment by Ovesco: the authors used suction only for getting colonic tissue into the cap. Already during this process small bowel was trapped in the cap and could clearly be seen between the teeth of the closed clip. The IFU of the product recommend the use of the Twin Grasper for fresh perforations. Suction may be used after the edges of a fresh perforation have clearly been identified and pulled into the cap.

Einsatz des Over-the-Scope Clips (OTSC) zur Behandlung einer Colon-Perforation verursachte eine Dünn darmfixation mit nachfolgender chirurgischer Resektion

Albrecht H, Nügel A, Hagel A, Rösler W, Förtsch T, Neurath MF, Raitheil M

• Comparison of the OTSC and cSEMS in the treatment of gastrointestinal leakages: results of retrospective

multicenter analysis

H. Farnik et al. reported in a very important paper of the interventional therapy of postoperative (73, 69 %), postinterventional leakages (24, 23 %), and 9 (8 %) spontaneous perforations. Primary closure was done with an cSEMS in 72 patients (69 %), and with OTSC in 31 pts. (29 %). Average duration of the treatment needed 45.6 days with cSEMS versus only 19.8 days in patients of the OTSC group. Treatment was complicated in 66.7 % of cSEMS patients, and only in 5.9 % in the OTSC group. 1.12 interventions (1.00-1.23) were needed with OTSC, 2.44 (2.12 – 2.76) in the cSEMS group. The diameter of the defect was larger in the cSEMS group (12.6 mm, 10.9-15.2) than in the OTSC group (7.1, 4.4-9.7).

	OTSC	cSEMS	P
Age (mean, years)	60.82	63.3	n.s.
Gender (m/f)	22/12	53/19	
Diameter of leakage (mean mm, 95 % ci)	7.1 (4.4-9.7)	12.6 (10.9-15.2)	
Duration of defect (days)	10.97	7.96	n.s.
Local infection (n, %)	15(44 %)	51(71 %)	
Number of interventions (mean, 95 % ci)	1.12 (1.00-1.23)	2.44 (2.12-2.76)	
Technical success (n, %)	33 (97.1 %)	71 (98.6 %)	n.s.
Duration of treatment	24.8	38.9	
Dislocation rate (n, %)	1 (3 %)	17 (24 %)	

The authors conclude in their discussion that OTSC is the preferred primary therapy of smaller post-interventional leakages. It might also be used in combination where cSEMS treatment was incomplete. They conclude that the longer treatment period with cSEMS and the higher complication rate might be due to sicker patients, but also due to the relevant dislocation rate of cSEMS.

Vergleich zwischen OTS-Klipp und cSEMS zur Indikationsstellung bei der Behandlung gastrointestinaler Leckagen: Ergebnisse einer retrospektiven, multizentrischen Analyse

Farnik H¹, Driller M¹, Kratt T², Schmidt C³, Fährndrich M⁴, Friedrich-Rust M¹, Filmann N¹, Königgrainer A², Stallmach A³, Heike M⁴, Zeuzem S¹, Albert J¹

¹Frankfurt a. M., ²Tübingen, ³Jena, ⁴Dortmund, Germany

• Endoscopic treatment of acute bleedings with an Over-The-Scope Clip (OTSC)

A. Braun et al. investigated the role of OTSC in the treatment of acute GI hemorrhage in an emergency.

Between 2011 and 2013 they treated 16 patients (median 75.5 y/o (61-92), m=9, f=7) with OTSC for upper and lower-GI bleeding (8 each). Patients with upper-GI bleeding received high PPI-medication simultaneously. 7 patients were classified F Ia, 2 F Ib, and 2 F IIa. All patients presented with an acute decrease of hemoglobin, with hematemesis, melena, and hematochezia. The clip was applied by using a standard forceps. Technical success was achieved in all 16 patients (100 %) with immediate primary hemostasis. None of the interventions took longer than 20 minutes. Only 6 patients underwent follow-up endoscopy between day 1 and 7 after clip application. All control endoscopies were uneventful and showed clinical success. 9 patients did not need any further endoscopy. None of the patients needed any further therapy for bleeding. All patients started normal oral intake from day 2. The authors conclude that OTSC is safe and effective for the treatment of hemorrhage which reduces mortality, with short intervention times.

Endoskopische Behandlung von akuten Blutungen mit einem Over-The-Scope Clip (OTSC)

Braun A, Richter-Schrag HJ, Fischer A, Freiburg, Germany

• Clinical experience in the treatment of perforations, leakages, and fistulas in the GI tract with the Over-the-scope clip (OTSC)

J. Stücker et al. report their retrospective results in the standard indications of OTSC.

21 patients (Median 69 years (30–87), m=11, f=10) were treated for leakages and fistulas (n=11, 52 %) due to anastomotic leaks. 5 patients had complications due to

diagnostic or therapeutic endoscopy. 2 patients had fistulas due to necrotizing pancreatitis. 2 patients suffered from a persistent PEG fistula. Technical success was reported in 20/21 cases (95 %). All 5 endoscopic complications could successfully be treated with OTSC. 7/11 anastomotic leaks could successfully be treated as well. One patient with duodenal leak due to acute necrotizing pancreatitis and a patient with perforated antrum died due to sepsis. The treatment of persistent PEG fistula was clinically not successful in this series.

The authors conclude that altogether the treatment of perforations, leakages, and fistulas with OTSC is very promising, and point out that this is especially true for the management of complications during endoscopy and surgical complications like anastomotic leakage.

Klinische Erfahrungen bei der Behandlung von Perforationen, Leckagen und Fistelungen im Gastrointestinaltrakt mit dem Over the scope Clip (OTSC®)

Stückle J, Probst A, Bittinger M, Scheubel R, Ebigo A, Messmann H, Gölder S, Augsburg, Germany

April 2014 | Efficacy and safety of OTSC® Proctology confirmed by clinical data at two major German conferences

40th Congress of the German Society for Coloproctology, April 3–5, 2014 in Munich under the presidency of Prof. Dr. Dr. h.c. W. Hohenberger

Prospective multicentric trial shows 90 % success rate for OTSC Proctology in anal fistula closure

R. Probst and co-authors, Stuttgart and Mannheim, Germany, presented data from a prospective multicentric trial including 20 patients treated with OTSC Proctology for anal fistula. 18 of 20 patients (90 %) reached the treatment success defined as clinical healing of the anal fistula and absence of recurrence at 6-month postoperative period. The authors conclude that OTSC Proctology is a new minimally invasive device for the treatment of anorectal fistula which is procedurally simple and time efficient. The risk profile is favourable, without relevant risk of fecal incontinence.

Anorektaler Fistelverschluss mittels OTSC Proctology: Ergebnisse einer prospektiven Beobachtungsstudie

R. Probst, A. Joos, A. Herold, D. Bussen, W. Ehni
Proktologisches Institut Stuttgart & Enddarmzentrum Mannheim

The 131st Congress of the German Society for Surgery, March 25–28, 2014 in Berlin under the presidency of Prof. Dr. J. Jaehne

OTSC® Proctology in retrospective analysis of mixed case series: efficacy and safety confirmed

S. Dango and colleagues, Kassel and Goettingen, Germany presented their experience using OTSC Proctology in the treatment of transsphincteric anal fistula. They conclude that OTSC placement is a promising sphincter-preserving minimally invasive method with considerably less complications than in more invasive types of surgical fistula treatment.

Efficacy and safety of the over-the-scope clip in the treatment of anal trans-sphincteric fistula

S. Dango, D. Schrader, M. Ghadimi, F. Antonakis, R. Hesterberg

Departments of General and Visceral Surgery, Rotes Kreuz Krankenhaus, Kassel and University Hospital, Goettingen
R. Menningen et al., Muenster, Germany report about their first experience with OTSC Proctology fistula closure in patients who had recurrence after fistula surgery. 9 consecutive patients were included into the trial. The authors conclude that OTSC is a safe and effective procedure for closing recurrent anal fistula even in more complex cases with Crohn's disease or multiple surgical pretreatments.

Verschluss analer Rezidivfisteln mit dem OTSC Proctology System

R. Menningen, M. Laukoetter, N. Senninger, E. Rijken
Klinik für Allgemein- und Viszeralchirurgie, University Hospital, Muenster

For more detailed information on the studies see reports in a pdf file on:
www.ovesco.com.

April 2014 | The OTSC® System: a surgery-sparing device for the management of iatrogenic duodenal perforation during endoscopic ultrasound

Duodenal perforations are a rare but serious complication during endoscopic ultrasound examinations. The closure of these perforations with the OTSC System can be a surgery-sparing approach. Three case studies published by Dr. Silvia Salord et al., Dept. of Digestive Disease, University Hospital, Barcelona, Spain and by Dr. Gianfranco Donatelli and colleagues, Endoscopy Unit, Hôpital Privé des Peupliers, Paris, France demonstrate the successful use of the OTSC System in case of iatrogenic duodenal perforation.

Two patients (aged 88 and 67) presented with cholangitis, one 74-year old woman with obstructive jaundice. In all three cases perforations occurred during endoscopic ultrasound procedures. Two perforations were located in the duodenal bulb, one at the superior duodenal flexure.

By deploying an OTSC clip successful closure was achieved in all cases, no further surgical interventions were required. The two patients with cholangitis underwent therapeutic endoscopic retrograde cholangiography (ERC) afterwards without any complications. Oral food intake was restarted after 2 or 5 days, respectively.

Endoscopic closure of duodenal perforation with an over-the-scope clip during endoscopic ultrasound guided cholangiopancreatography

Salord S, Gornals JB, Maisterra S, Pons C, Busquets J, Fabregat J

Rev Esp Enferm Dig. 2012 Sep;104(9):489-90

Closure with an over-the-scope clip allows therapeutic ERCP to be safely performed after acute duodenal perforation during diagnostic endoscopic ultrasound

Donatelli G, Vergeau BM, Dritsas S, Dumont JL, Tuszyński T, Meduri B

Endoscopy. 2013 Nov;45 Suppl 2 UCTN:E392-3

March 2014 | OTSC®@FISMAD, Naples, Italy:

77 % success in anastomotic leak treatment

At the 20th National Congress of Digestive Diseases, Napoli, Italy, March 19-22, MA Bonino and colleagues, Department of Surgery, Turin University reported about a consecutive series of 26 patients treated with OTSC for postsurgical colorectal leaks.

The mean defect size was 8.7 mm, in 10 cases there were acute and in 16 cases chronic leaks (fistula). 4 cases were complicated by recto-vaginal, 3 by recto-vesical and 7 by colo-cutaneous fistula. In 3 cases OTSC was used to complete earlier vacuum sponge therapy. The overall success rate was 77 % (20/26), 90 % in acute (9/19) and 69 % (11/16) in chronic cases. There were no OTSC-related complications, additional surgery was needed in 2 cases.

Anastomotic leakage is a serious and non infrequent complication in colorectal surgery. Incidence rates in the literature range from 1 to 39 %. Clinically relevant leaks are commonly seen in 3-6 % of the cases. OTSC closure of colorectal post-surgical leaks and fistula is a safe technique with a high success rate.

Efficacia della clip OTSC per il trattamento di deiscenze e fistole chirurgiche del colon-retto

Efficacy of the Over-The-Scope Clip (OTSC) for treatment of colorectal postsurgical leaks and fistulas

Bonino MA, Verra M, Salva A, Bullano A, Rapetti L, Arezzo A, Morino M

March 2014 | Management of esophageal perforation with the OTSC® System – four new case studies by different authors report favourable results

Spontaneous or iatrogenic esophageal perforation is a life-threatening condition that can lead to severe mediastinitis, sepsis and multiple organ failure. Endoscopic management has contributed to the decrease of morbidity and mortality associated with surgical repair. Four different case reports lately published by Dr. Alexander Braun et al., Div. of General Surgery, University of Freiburg, Germany, Dr.

Daive Bona et al., Div. of General Surgery, University of Milan, Italy and Dr. Alexandre Ferreira, Dept. of Gastroenterology and Hepatology, Hospital de Santa Maria, Lisbon, Portugal illustrate the successful closure of esophageal perforations with the OTSC System.

Two patients presented with Boerhaave's syndrome, one patient had an iatrogenic perforation and one patient suffered from a perforation caused by a fishbone. In all four cases a minimally invasive approach with the OTSC System was chosen. Two patients were treated with the OTSC clip within 12 hours. Although the two other patients were admitted to hospital not until after 48 h after an episode of vomiting, late management of the esophageal perforation with the OTSC System was successful. After complete closure of the defect, all patients were kept on antibiotic therapy and were discharged in stable condition after 10 days (patient with iatrogenic perforation), 21 days (patient with perforation caused by a fishbone) and 20 or 28 days respectively (patients with Boerhaave's syndrome).

3-month follow-up revealed a free esophageal passage and correct placement of the OTSC clip. The OTSC clip is a new, safe and effective treatment alternative for the management of esophageal perforation. Due to the endoscopic approach and shorter hospital stay, the procedure is more cost effective than conventional surgical procedures.

Endoskopischer Verschluss von distalen Ösophagus-Perforationen mit einem Over-The-Scope Clip (OTSC)

Braun A, Hopt UT, Richter-Schrag HJ

Endo heute 2013; 131

Management of Boerhaave's syndrome with an over-the-scope clip

Bona D, Aiolfi A, Rausa E, Bonavina L

Eur J Cardiothorac Surg. 2013 Jul 18 [Epub ahead of print]

Snapper fishbone esophageal perforation closed with an over-the-scope clip

Ferreira AO, Lopes J, Velosa J

BMJ Case Rep. 2013 Oct 25;2013 144

March 2014 | OTSC: easy to use with good results, decreasing morbidity and mortality in diagnostic and therapeutic endoscopy

In the quest to describe the use and the clinical applications of OTSC System in an environment where endoscopic and surgical techniques are increasingly more complex and frequent Singhal et al. have searched and analysed the literature using the key words „endoscopy“ and „over-the-scope clip“ in order to identify human studies evaluating the application of OTSC from January 2001 to August 2012.

The indication, efficacy, complications, and limitations were recorded. The overall success rates of OTSC based on the current literature are in the range of 75 % to 100 % for closure of iatrogenic gastrointestinal perforations, 38 % to 100 % for closure of gastrointestinal fistulas, 50 % to 100 % for anastomotic leaks, and 71 % to 100 % for bleeding lesions. OTSC has shown 100 % success rates in managing postbariatric surgery weight gain secondary to dilation of the gastro-jejunoanastomosis.

The authors conclude that OTSC is easy to use with good results, thus decreasing the morbidity and mortality associated with complications secondary to both diagnostic and therapeutic endoscopy and avoiding surgery in many situations.

Over-the-Scope Clip: Technique and Expanding Clinical Applications

Singhal S, Changela K, Papafragkakis H, Anand S, Krishniah M, Duddempudi S

J Clin Gastroenterol. 2013 Oct;47(9):749-56

Further reviews by Weiland et al.:

Performance of the OTSC System in the endoscopic closure of iatrogenic gastrointestinal perforations: a systematic review

Weiland T, Fehker M, Gottwald T, Schurr MO

Surg Endoscop. 2013; 27:2258-74

Performance of the OTSC System in the endoscopic closure of gastrointestinal fistulae – a meta-analysis

Weiland T, Fehker M, Gottwald T, Schurr MO

Minim Invasive Ther Allied Technol. 2012 Jul;21(4):249-58

Epub 2012 Jun 14 113

February 2014 | New case series on use of OTSC® for treatment of refractory upper GI bleeding

Apart from using the OTSC System in acute and chronic perforations (i.e. perforations, anastomotic leakage, fistulae) the authors of the renowned Institute of Digestive Disease, Department of Surgery, Chinese University of Hong Kong are reporting of patients in whom OTSC was used for endoscopic control of refractory or major upper gastrointestinal bleeding from lesions in the gastroduodenal tract between 1 July and 31 December 2012. Nine patients were included (median age 72.5 years, range 39–91) with bleeding gastric ulcers (n=2), bleeding duodenal ulcers (n=5), bleeding gastrointestinal stromal tumor in the stomach (n=1), and bleeding from ulcerative carcinoma of the pancreas (n=1). The median size of the ulcers was 2.5 cm (range 1–4). Six of the nine patients had undergone previous endoscopic hemostasis. Technical success (defined as hemostasis achieved at index endoscopy) was achieved in all patients and the clinical effectiveness was 77.8 % (defined as technical success with no rebleeding). All procedures were carried out by two experienced endoscopists. Those two patients that experienced rebleeding suffered from complex duodenal ulcer. One of them had been treated with radiotherapy for residual disease after resection of common bile duct cholangiocarcinoma. After several additional EGDs, transarterial embolization, and one surgical intervention which all failed to stop the bleeding, the patient died eventually. The second patient bled from the inferior pancreaticoduodenal artery and needed arterial embolization as well.

The authors discuss a meta-analysis of 1156 patients in 15 randomized trials where endoclips were shown to be superior to injection alone, and as effective as heater probe treatment. The overall rate of rebleeding in those conventionally treated patients ranged between 7.1 % and 9.5 % though. Since rebleeding correlates with the adverse outcome of this indication they speculate that control of bleeding would have a positive impact on patient outcome. Even though the study was carried out in patients with complex duodenal ulcer and underlying malignancies the technical success rate of OTSC was 100 %. They also point out that usually in cases like these the application of conventional clips is difficult; the repeated application of heater probe being associated with a higher risk of perforation. Whereas the application of OTSC allows for larger amounts of tissue and constitutes a quite durable treatment (OTSC in situ after a median of 28 days in this study). The authors conclude that the use of OTSC is a safe and effective method of endoscopic hemostasis for major bleeding from miscellaneous upper gastrointestinal causes and should be considered in refractory bleeding after conventional endoscopic hemostasis, before surgery or angiographic embolization.

Comment by Ovesco: a prospective controlled randomized multicenter trial with 64 patients with recurrent upper GI bleeding is recruiting in Germany (Endoscopic Treatment of Recurrent Upper GI Bleeding: OTSC [Over the Scope Clip] Versus Standard Therapy (STING). ClinicalTrials.gov Identifier: NCT01836900)

Use of Over-The-Scope Clip for treatment of refractory upper gastrointestinal bleeding: a case series

Chan SM, Chiu PW, Teoh AY, Lau JY

Endoscopy. 2014 Feb 6.

[Epub ahead of print]

162

February 2014 | Retrospective study on efficacy and safety of the OTSC® System in the treatment of GI bleeding, fistula and perforation: primary technical success rate 91.3 %, durable clinical success rate 82.6 %

Dr. Noriko Nishiyama and colleagues, Dept. of Gastroenterology and Neurology, Kagawa University, Japan, recently presented their retrospective study on efficacy and safety of the OTSC System in endoscopic closure of gastrointestinal bleeding, fistulas and perforations, concluding that the OTSC System is a highly useful device that can safely be utilized for these indications.

Their case series consisted of 23 consecutive patients treated between November 2011 and September 2012 (mean age 77 years) included the following indications for OTSC placement: stopping GI bleeding (n=9), closing perforation (n=10), closing chronic fistula (n=4) and prevention of post endoscopic submucosal dissection (ESD) duodenal artificial ulcer perforation (n=1). One patient had a perforation that formed a fistula. Lesions were located in the esophagus (n=1), the stomach (n=10), the duodenum (n=5), the sigmoid colon (n=3) and in the rectum (n=4). In 8 patients other therapies preceded OTSC application (e.g. conventional hemostatic clips, local injections, hemostatic coagulation forceps). Median follow-up time was 67 days. The primary technical success rate was 91.3 % (21/23). In two cases application of the OTSC clip was not possible due to stiff, fibrotic lesion edges. The overall clinical success rate (complete closure by using only OTSC clips) was 82.6 %. Major contributing factors for OTSC failure were a large lesion size (greater than 20 mm) and a delayed diagnosis (more than 1 week). No patient reported any complications associated with OTSC placement. In conclusion, the OTSC is an interesting and novel device that enhances the armamentarium of therapeutic gastroenterologists.

Efficacy and safety of over-the-scope clip: including complications after endoscopic submucosal dissection
Nishiyama N, Mori H, Kobara H, Rafiq K, Fujihara S, Kobayashi M, Oryu M, Masaki T

World J Gastroenterol. 2013 May 14;19(18):2752-60

126

OTSC® Update 16

January 2014 | OTSC® in mucosal flap closure after peroral endoscopic myotomy (POEM)

Maintaining the integrity of the mucosal flap and the reliable closure of mucosal entry during peroral endoscopic myotomy (POEM) is paramount in preventing leakage of esophageal contents into the mediastinal space. In a recently published case series (n=2) Payal Saxena, MD and colleagues, Dept. of Medicine and Div. of Gastroenterology and Hepatology, Johns Hopkins Medical Institutions, Baltimore, Maryland, USA describe their positive experience with the application of the OTSC System for reliable and easy flap closure after POEM.

Both patients presented with dysphagia and regurgitation and were diagnosed with achalasia. It was decided to proceed with POEM. After myotomy of the inner circular muscle bundles it was noted that the mucosal incision had elongated from 2 cm to 4 cm in one case. Whereas the distal part of the mucosal entry was successfully closed with conventional hemostatic clips (Resolution Clip, Boston Scientific) in both cases, closure of the proximal half was not possible even with different clips. As the clips were noted to slip to one side of the mucosal incision, there was a risk of displacing clips into the submucosal tunnel. Hence, all partially attached clips were removed with biopsy forceps. Finally, complete closure of the mucosal incision was performed with the OTSC clip and the OTSC Twin Grasper in both cases. Contrast swallow of the esophagus the following day revealed no leaks in either patient.

The authors state that the OTSC clip provides more durable closure than standard hemostatic clips and full-thickness closure is achievable due to greater compressive force. Considering that failure of closure risks serious adverse events, like mediastinitis and sepsis, these features of the OTSC clip appear even more attractive.

An alternative method for mucosal flap closure during peroral endoscopic myotomy using an over-the-scope clipping device

Saxena P, Chavez YH, Kord Valeshabad A, Kalloo AN, Khachab MA

Endoscopy. 2013 Jul;45(7):579-81

122

January 2014 | Avoiding Surgery: Minimally invasive endoscopic management of an iatrogenic colon perforation

Iatrogenic lesions of GI organs are a significant complication of diagnostic or interventional endoscopic procedures. Dr. Pilar Díez-Redondo and colleagues, Dept. of Gastroenterology, Hospital Universitario, Río Hortega, Valladolid, Spain present a case report on OTSC clipping for colon perforation closure: For assessment of iron

deficiency an 82-year old woman was referred to the endoscopic unit. A gastroscopy confirmed a hiatal hernia. Colonoscopy revealed no abnormalities. 18 cm proximal to the anus an iatrogenic perforation with a size of 12 mm occurred. To close the perforation endoscopically an 11/a OTSC clip was chosen. The target tissue and a piece of omentum were pulled into the applicator cap by suction and the clip released successfully, approximating the edges of the lesion. A small residual recess was closed with two conventional, endoscopic clips. The patient was discharged 10 days after the intervention. A 7-month follow-up confirmed the correct placement of the OTSC.

Iatrogenic colon perforations can cause severe complications and often require surgery, as the major drawback of an endoscopic approach with conventional clips is the limited ability of these clips to achieve sufficient apposition of the mucosa and submucosa to ensure tight sealing of the perforation. With the advent of the larger and more powerful OTSC clips, surgery can be avoided and perforations managed in a minimally invasive, endoscopic way. For that reason the authors suggest that the OTSC System should be available to all endoscopy units as a bail-out device.

A novel system for endoscopic closure of iatrogenic colon perforations using the Ovesco® clip and omental patch

Díez-Redondo P, Blanco JI, Lorenzo-Pelayo S, De-la-Serna-Higuera C, Gil-Simón P, Alcaide-Suárez N, Pérez-Miranda M

Rev Esp Enferm Dig. 2012 Oct-Nov;104(10):550-2

103

January 2014 | Closure of gastric fistulas after bariatric surgery with the OTSC® System – two case studies

Iatrogenic gastric fistulas after bariatric surgery are a potentially dangerous situation as they can lead to severe complications, such as peritonitis and abscess formation. Two case reports recently published by Dr. Victoria Gómez and colleagues, Dept. of Gastroenterology and Hepatology, Mayo Clinic, Jacksonville, USA, and Dr. Hany Shehab et al., Dept. of Gastroenterology, Dar Al Fouad Hospital Giza, Egypt, respectively, describe the closure of gastric fistulas with the OTSC System after laparoscopic bariatric surgery. Dr. Gómez reports on a 45-year old woman who was hospitalized for management of complications from a prior sleeve gastrectomy. Postoperatively the patient developed fever and abdominal pain. A CT scan showed a fluid collection in the region of the right liver lobe, free intraperitoneal air and an abscess in the postsurgical bed of the stomach. The results of a barium contrast study were consistent with a significant leak in the proximal third of the gastric sleeve. EGD revealed a gastric fistula 4 centimeters below the esophago-gastric junction. As an initial treatment with an esophageal stent and abdominal drains had failed to seal the leak, the stent was removed and a fully covered esophageal stent was applied. A second stent had to be deployed to bridge the prior stent. Since there was no improvement of the fistula, the stents were again removed. As next treatment approach the fistula was grasped with the OTSC Twin Grasper and closed by application of an OTSC clip. A follow-up radiograph showed no extravasation of contrast.

Dr. Shehab presents the case of a 36-year old man who had undergone a Roux-en-Y gastric bypass for morbid obesity. Postoperatively an anastomotic leak was found. Two attempts of surgical repair failed as well as a conservative approach with drainage and insertion of a feeding jejunostomy. 5 months after the first surgery an EGD revealed a well-epithelialized fistula with a wide lumen. It was decided to close the fistula by OTSC clipping. To remove the epithelium at the fistula orifice, argon plasma coagulation was applied to the proximal lumen of the fistula. Then the OTSC Twin Grasper was used to approximate the edges of the fistula orifice followed by the application of the OTSC clip. After 10 months there was no evidence of a fistula recurrence.

Since a surgical intervention for postoperative fistulas in an obese patient with recent bariatric surgery is most often not desirable, a minimally invasive, endoscopic approach with the OTSC System is an attractive treatment option. In comparison to conventional clips that are only suitable for

small fistulas and only attach to the superficial mucosal layer, the OTSC clip offers a deeper grasp of the tissue and a sturdier closure.

Closure of an iatrogenic bariatric gastric fistula with an over-the-scope clip

Gómez V, Lukens FJ, Woodward TA

Surg Obes Relat Dis. 2013 Mar-Apr;9(2):e31-3

doi: 10.1016/j.soard.2012.09.004

99

Combined endoscopic techniques for closure of a chronic post-surgical gastrocutaneous fistula: case report and review of the literature (with video)

Shehab HM, Elasmr HM

Surg Endosc. 2013 Aug;27(8):2967-70

doi: 10.1007/s00464-013-2839-1 Epub 2013 Feb 23

119

December 2013 | First report on successful management of delayed presentation of Boerhaave's syndrome

Current guidance has advocated surgery for delayed presentations of Boerhaave's syndrome with evidence of mediastinal contamination. However, Dr. Eamon Ramhamadany and colleagues, Dept. of General Surgery, University Hospital Coventry and Warwickshire, UK, present the successful management of Boerhaave's syndrome in a 69-year-old man by means of the OTSC System, sparing the patient surgery and possible associated complications. The man presented to hospital with an episode of forceful vomiting. A chest radiograph was performed revealing a pleural effusion. After several days without improvement a CT chest showed an oesophageal perforation with mediastinitis. Because of the size of the defect and the delay in presentation, it was decided not to perform surgery, but to apply the OTSC clip for endoscopic repair. A contrast swallow confirmed the correct placement of the clip and the successful closure of the leak. After a total parenteral nutrition for 3 days, the patient was fed via a naso-jejunal tube. Intravenous antibiotics and bilateral chest drains led to a resolving mediastinitis. The whole procedure resulted in a favourable outcome without the need for surgery.

The authors conclude that the OTSC can be used to manage patients with delayed presentation of Boerhaave's and that further evaluation is needed to define the indications for minimally invasive techniques like the OTSC System.

A delayed presentation of Boerhaave's syndrome with mediastinitis managed using the over-the-scope clip.

Ramhamadany E, Mohamed S, Jaunoo S, Baker T, Mannath J, Harding J, Menon V

J. surg. case rep. (2013) 2013 (5): rjt020.

120

December 2013 | Management of postoperative gastrointestinal leakages and fistulas with the OTSC® System: long-term success rate of 79 %

Dr. Rudolf Mennigen and colleagues, Dept. of General and Visceral Surgery, University Hospital of Muenster, Germany, recently presented a study on efficacy and safety of the OTSC System in endoscopic closure of postoperative gastrointestinal leakages and fistulas, concluding that the OTSC System dramatically increases the possibilities of defect closure by endoscopic clipping as opposed to conventional endoclips.

Their case series of 14 consecutive patients (May 2011–November 2012) included patients with anastomotic leakage (n=6) e.g. after gastrectomy, perforation after fundoplication (n=1) and post-operative fistulas (n=7, colocutaneous, enterocutaneous, gastrobronchial, rectourethral, rectocutaneous, gastrophleural). 11 of the 14 lesions were chronic (treated by OTSC later than postoperative day 14) and in 9 patients other therapies preceded OTSC application (e.g. covered stent application, fibrin glue injection). Median follow-up time was 5.5 months. The primary procedural success rate was 100 %. 3/14 patients (21 %) required further treatment during follow-up. Reasons for OTSC failure were massive fibrosis of the fistula and application in an actively inflamed bowel segment in Crohn's disease. However, unsuccessful OTSC treatment did not impair subsequent surgical therapies. Complete and clinically durable closure of the defects was achieved in 79 %, indicating from the authors' point of view that the OTSC will play an important role in the therapy of

postoperative leakages.

Endoscopic closure of postoperative gastrointestinal leakages and fistulas with the Over-the-Scope Clip (OTSC)

Mennigen R, Colombo-Benkmann M, Senninger N, Laukoetter M

Gastrointest Surg. 2013 Jun;17(6):1058-65

115

November 2013 | OTSC® in endoscopic treatment of acute GI bleeding after failure of conventional techniques: primary hemostasis of 97 %

The OTSC System can overcome the limitations of conventional clips in the treatment of patients with acute GI bleeding by providing compression of large amounts of tissue, leading to a more efficient hemostasis. Dr. R. Manta and colleagues, Gastroenterology and Endoscopy Unit, New S. Agostino Hospital, Modena, Italy draw this conclusion on the basis of a retrospective analysis of a consecutive case series of 30 patients with severe acute GI bleeding treated with the OTSC System after failure of conventional techniques.

Data were collected from six high-volume endoscopy units in a period between December 2011 and September 2012. All 30 patients suffered from bleeding lesions unresponsive to saline/adrenaline injection and through-the-scope clipping located in the upper and lower GI tract in 23 and 7 cases, respectively. Bleeding lesions included duodenal ulcer (n=12), gastric ulcer (n=6), Mallory-Weiss (n=2), Dieulafoy (n=2) and surgical anastomosis (n=1) in the upper GI tract and endoscopic mucosal resection (n=5), endoscopic submucosal dissection (n=1) and colonic diverticulum (n=1) in the lower GI tract.

Primary hemostasis with OTSC was achieved in 29 of 30 cases (97 %). Rebleeding in two cases was successfully treated with injection of saline and adrenaline. Endoscopic follow-up after 2–4 days and after 1 month revealed correct placement of the OTSC clip and no procedure-related complications. Thus, the OTSC is an effective and safe device for treatment of acute GI bleeding and represents a useful adjunct to the therapeutic armamentarium in endoscopic emergencies.

Over-the-scope clip (OTSC) represents an effective endoscopic treatment for acute GI bleeding after failure of conventional techniques

Manta R, Galloro G, Mangiavillano B, Conigliaro R, Pasquale L, Arezzo A, Masci E, Bassotti G, Frazzoni M

Surg Endosc. 2013 Sep;27(9):3162-4

114

October 2013 | OTSC® System: Effective closure of esophageal fistula following total gastrectomy

Postoperative leaks after total gastrectomy are among the most common early complications. Dr. C. N. Ferreira and colleagues, Serviço de Gastreenterologia e Hepatologia, Hospital de Santa Maria, Lisbon, Portugal report on a 78-year old woman presenting with melena. She was diagnosed with gastric adenocarcinoma and treated with total gastrectomy and esophagojejunal Roux-en-Y anastomosis. On the fifth postoperative day she developed a septic condition caused by a fistulous orifice just above the intact anastomosis. Due to her poor general condition a surgical intervention was unfeasible. Thus, it was decided to treat the fistula endoscopically by means of the OTSC clip. By using the OTSC Twin Grasper to approximate the edges of the fistula and application of an OTSC clip the orifice was effectively closed. The patient was discharged in stable condition two weeks later. In a commentary to this publication Dr. David Robbins, Assistant Editor of the Journal Gastrointestinal Endoscopy emphasizes the significantly higher strength of the OTSC clip for hemostasis and closure of GI tract wall in comparison to conventional endoscopic clips.

Total gastrectomy in an elderly patient complicated by esophageal fistula: rescue by the over-the-scope clip

Noronha Ferreira C, Ribeiro LC, Velosa J, Ferreira J, Ferreira C, Freire JP, Marques J, Ruivo A, Bicha Castelo H

Gastrointest Endosc. 2013 Mar;77(3):497-8 ;

[Epub 2013 Jan 4]

112

October 2013 | Efficacious OTSC® hemostasis in Dieulafoy's gastric lesion resistant to conventional endoscopic treatment

Dr. B. Mangiavillano and colleagues, Gastrointestinal Endoscopy, Azienda Ospedaliera San Paolo University, Hospital-University of Milan, Italy, present a case study of a 69-year old woman with an episode of melena. EGD showed a Dieulafoy's bleeding lesion in the proximal third of the posterior wall of the stomach. The lesion was treated with an epinephrine injection and application of two conventional working-channel delivered metallic clips and the patient was discharged two days later. After three days the patient again presented with melena. Blood transfusions were necessary. An EGD was performed, showing no sign of an actively bleeding ulcer. The patient was admitted to hospital and suffered from another episode of melena with hemorrhagic shock. The now actively bleeding Dieulafoy's lesion was then treated with an OTSC clip, stopping the hemorrhage completely and persistently. Endoscopic follow-up after 30 days displayed correct placement to the OTSC and no signs of further bleeding.

Successful treatment with an over-the-scope clip of Dieulafoy's gastric lesion resistant to conventional endoscopic treatment

Mangiavillano B, Arena M, Morandi E, Viaggi P, Masci E

Endoscopy. 2012;44 Suppl 2 UCTN:E387

105

OTSC® Update 15 EXTRA

October 2013 | OTSC® successful in providing hemostasis in posterior duodenal ulcer bleeding after failure of conventional clips

Ulcer bleeding is one of the key indications for the OTSC System. In a recently published case series (n=4), Prof. Klaus Mönkemüller and colleagues, Dept. of Internal Medicine, Gastroenterology and Infectious Diseases, Marienhospital Bottrop, Germany add to the growing clinical experience in using the OTSC System to control massive gastrointestinal bleedings and achieve life-saving hemostasis. All four patients (mean age 84.5) presented with hypotension and mean hemoglobin of 9 g/dL. After initial fluid resuscitation an emergent EGD displayed actively oozing ulcers in the posterior duodenum. As an initial therapy with injection of epinephrine-saline solution and standard clip placement failed and all patients suffered from rebleeding, the decision to apply the OTSC System was made. Hemostasis was attained successfully and all patients discharged in stable conditions. Even in difficult located ulcers in the posterior duodenum the placement of the OTSC is easy and effective to obliterate bleeding vessels resulting in life-saving hemostasis.

Utility of the „bear claw“, or over-the-scope clip (OTSC) system, to provide endoscopic hemostasis for bleeding posterior duodenal ulcers

Mönkemüller K, Toshniwal J, Zabielski M, Vormbrock K, Neumann H

Endoscopy. 2012;44 Suppl 2 UCTN:E412-3.

[Epub 2012 Nov 20]

108

September 2013 | OTSC® in post-surgical complications: retrospective case review confirms high clinical efficacy

Dr. Alisa Coker and colleagues, Dept. of Surgery, University of California San Diego, USA, report on their experience with the OTSC System in a retrospective review of all cases treated between August 2011 and March 2012.

All 10 patients had clinically significant gastrointestinal post-surgical complications. Indications included: gastric leaks after sleeve gastrectomy (n=4), post-operative colonic leak following extended hemicolectomy and palliative debulking (n=1), gastro-gastric fistulas following Roux-en-Y gastric bypass (n=2), esophageal perforation (n=3).

Three of the four patients with gastric leaks had undergone previous unsuccessful attempts at endoscopic repair (stenting, fibrin glue application, traditional clipping, endoscopic suturing). The overall clinical success rate was 70 %. Re-surgery was needed in the two cases of gastro-gastric fistulas. In the colonic leak patient the clip placement procedure was aborted due to a fixed tortuous sigmoid colon as a result of the metastatic disease and adhesions, limiting

endoscope passage.

For the subgroup of seven patients treated for leaks and perforations a success rate of 87.5 % with complete resolution was achieved. The mean follow-up period was 83 days. No complications occurred.

The authors conclude that the OTSC System is simple to use, safe and effective with a great potential for success in a broad number of applications. For the treatment of gastric leaks following sleeve gastrectomy the OTSC System is their first-line treatment.

Initial Experience with an Innovative Endoscopic Clipping System

Coker AM, Jacobsen GR, Acosta G, Talamini MA, Savides TJ, Horgan S

Surg Technol Int. 2012 Dec;22:39-43

[Epub ahead of print]

111

September 2013 | Preventive closure of duodenal lesion after endoscopic submucosal dissection with the OTSC® System to obviate delayed perforation

The two case reports published in the journal Digestive Endoscopy by Dr. Hirohito Mori and his colleagues, Dept. of Gastroenterology and Neurology, Faculty of Medicine, Kagawa University, Japan illustrate the complete closure of secondary duodenal ulcers after endoscopic submucosal dissection (ESD) with the OTSC System without any complications.

Two elderly patients were diagnosed with early duodenal cancer. ESD was carried out successfully removing the lesions en bloc. In one case the muscle layer was slightly injured but not perforated. Because of the exposure to bile and pancreatic juices the risk of post-ESD delayed perforation is much higher in the duodenum than in other parts of the gastrointestinal tract. As conventional clips are less suitable due to small size and insufficient grasping power, Dr. Mori and his team used the OTSC System to close the lesion completely without any complications. The ulcer closure procedure time was 7 resp. 10 min. In both cases control endoscopy revealed a complete healing of the ulcer after 30 days.

Dr. Mori and his colleagues consider the OTSC System to be one of the most effective devices to prevent delayed perforations in post-ESD ulcer.

Successful closing of duodenal ulcer after endoscopic submucosal dissection with over-the-scope clip to prevent delayed perforation

Mori H, Shintaro F, Kobara H, Nishiyama N, Rafiq K, Kobayashi M, Nakatsu T, Miichi N, Suzuki Y, Masaki T

Dig Endosc., 2013 Jul;25(4):459-61 | Epub 2012 Aug 7

89

September 2013 | First two publications of endoscopic closure of gastrocolic fistula using the OTSC® System

Dr. Alberto Murino, Wolfson Unit for Endoscopy, St Mark's Hospital, London, UK, and his colleagues report on a successfully treated gastrocolic fistula using the OTSC System. A migrated PEG tube caused a gastrocolic fistula in the transverse colon in a 41 y/o male with cerebral palsy. The fistula led to extensive diarrhea and mouth odor. The CT showed an involvement of the greater curvature of the stomach. By using the OTSC Anchor to approximate the tissue the OTSC clip was released precisely closing the fistula orifice completely. Diarrhea and mouth odor were stopped. The 3 months' follow-up revealed a complete healing of the fistula.

The second case report published in the World Journal of Gastrointestinal Endoscopy by Prof. Klaus Mönkemüller and colleagues, Division of Gastroenterology and Hepatology, Basil Hirschowitz Endoscopic Center of Excellence, University of Alabama, Birmingham, USA, describes the effective endoscopic closure of a large gastrocolic fistula with the OTSC System in an extremely malnourished patient with complex post-surgical upper GI anatomy. The 47 y/o man presented with chronic diarrhea and severe weight loss of 32 kg in a 1-year period. He had a history of chronic pancreatitis, alcoholism and Billroth II gastrectomy due to a perforated peptic ulcer.

Endoscopy showed a clean based ulceration at the anastomosis and a second orifice that represented the fistula. Connecting stomach and colon, the fistula measured about 10–12 mm. Because of the patient's poor clinical status he could not benefit from a surgical intervention so an endoscopic procedure using the OTSC System was chosen. To ensure a definitive closure of the fistula the OTSC Twin Grasper was used to approximate the edges of the fistula. The application of the OTSC led to a complete closure of the gastrocolic fistula which was confirmed by an endoscopy.

For Prof. Mönkemüller this case "adds to the growing evidence that the OTSC System is a useful device to treat clinically significant endoluminal GI defects." He believes that "this device is a major breakthrough for the management of various types of discontinuity defects or fistulas of the GI tract (...)" and that "the OTSC System should be incorporated into the therapeutic armamentarium of the advanced endoscopist."

First report of endoscopic closure of a gastrocolic fistula using an over-the-scope clip system (with video)

Murino A, Despot EJ, Vaizy C, Bashir G, Hansmann A, Gupta A, Konieczko K, Fraser C

Gastrointest Endosc. 2012 Apr;75(4):893; discussion 894

72

Endoscopic closure of gastrocolic fistula using the over-the-scope clip system

Mönkemüller K, Peter S, Alkurd B, Ramesh J, Popa D, Wilox C

World J Gastrointest Endosc. 2013 Aug 16;5(8):402-6

136

August 2013 | The interesting case: OTSC® closure of esophagobronchial fistula

Dr. E. Zolotarevsky and colleagues from the Department of Gastroenterology and Nutrition Service at Memorial Sloan-Kettering Cancer Center, New York City report about an interesting case in which OTSC clipping was used for closing an esophagobronchial fistula.

An 83 y/o woman presented with a symptomatic fistula arising from an esophageal diverticulum with recurrent pulmonary infections. Placing a covered self-expanding metal stent was not believed to result in adequate seal of the chronic lesion. The placement of a percutaneous gastrostomy tube was refused by the patient. Also bronchial stenting and surgery were not considered as good options in this case.

In this situation closure of the fistula with the OTSC clip was decided. A 12/6/t clip was placed under endoscopic control and with the aid of the OTSC Anchor for better manipulation and targeting of the fistula orifice. Immediate technical success was achieved and verified by barium esophagogram 2 days later. The patient was discharged from the hospital after 1 week in stable condition. The clip was still found in place at 1 month follow-up by chest X-ray but passed spontaneously and uneventfully as seen in CAT scan 45 days after the procedure. Final follow-up at 3 months revealed no recurrence or postprandial cough.

Esophagobronchial fistula closure using a novel endoscopic over-the-scope clip

Zolotarevsky E, Kwon Y, Bains M, Schattner M

Ann Thorac Surg. 2012 Sep;94(3):e69-70.

96

July 2013 | OTSC® effective in emergency closure of iatrogenic GI perforations instead of abdominal surgery

Dr. Hagel and colleagues, Dept. of Gastroenterology, University of Erlangen-Nuremberg, Germany reported about a consecutive series of 17 cases with perforations of the digestive tract, treated with OTSC clipping. All cases were considered as being candidates for abdominal surgery for closing the perforation. In 11 cases perforation closure with OTSC was immediately successful, thus avoiding surgery in 64.7 %. In 6 cases surgical closure was done. The area size of perforation in the successful cases was 21.1 ± 9.1 sqmm; in the unsuccessful group the area size was 97.6 ± 149 sqmm. Unsuccessful cases had on average a larger size, necrotic margins and required more OTSC clips during closure attempts (2.3 ± 0.5 , $p=0.018$).

The authors conclude: "OTSC application yields a high rate of endoscopic perforation closure in patients with macroscopic gastrointestinal perforation, even in an emergency setting, representing an alternative to surgery, especially when the size of the lesion is not too large and when vital or solid perforation margins are expected."

OTSC application yields a high rate of endoscopic perforation closure in patients with macroscopic gastrointestinal perforation, even in an emergency setting, representing an alternative to surgery, especially when the size of the lesion is not too large and when vital or solid perforation margins are expected."

Over-the-Scope Clip Application Yields a High Rate of Closure in Gastrointestinal Perforations and May Reduce Emergency Surgery

Hagel AF, Naegel A, Lindner AS, Kessler H, Matzel K, Dauth W, Neurath MF, Raitchel M

J Gastrointest Surg. 2012 Nov;16(11):2132-8

[Epub 2012 Aug 18]

88

OTSC® Update 14

July 2013 | OTSC® System in transgastric appendicectomy

Kaehler et al. report the results of their first 15 patients in a prospective trial on "Transgastric appendectomy" which now already recruited 30 patients who are currently under follow-up.

From April 2010 the Mannheim group offered to their patients a transgastric appendectomy. Patients with generalised peritonitis and/or local contraindications were not recruited. Out of 111 eligible candidates 15 agreed to undergo the proposed NOTES procedure.

14 out of 15 were actually operated through NOTES, whereas 1 patient was switched to laparoscopic procedure due to severe inflammation and adhesions. In each case the gastrotomy was closed by a single OTSC System using Twin Grasper and 12/6/t clip. All closures were tight primarily and uneventful throughout the follow-up.

This is the first series of transgastric appendectomy using the OTSC System (and the second series overall). All 30 patients who have been recruited altogether will be reported in a separate publication.

Transgastric appendicectomy

Kaehler G, Schoenberg MB, Kienle P, Post S, Magdeburg R Brit J Surg. 2013; 100: 911-915

130

July 2013 | Recommendation of OTSC® System in complex GI bleeding

In an overview article the authors are referring to the current guideline therapies available and new developments. They report that other new three-dimensional clips seem to be even less efficacious than normal hemoclips. Thus, the authors conclude that obviously one of the key elements to successful hemostasis is the strength of the jaws of a clip and the amount of tissue captured. They state that this is obviously fulfilled by the design of the OTSC System which allows for the capture of a large amount of tissue and is more secure than other clips in the experimental setting. Thus the OTSC System is being recommended and used in complex GI bleeding. According to Leung & Lau a single clip suffices for most circumstances and therefore the procedure is shorter when compared to multiple applications of hemoclips.

Comment by Ovesco: In a recently published series of 83 patients with severe and complicated GI bleedings (e.g. relapses after conventional endoscopic hemostasis or indication for a surgical intervention) the success rate was close to 93 % with OTSC (Kratt T et al., Poster DGE-BV meeting, Munich 3/2013)

New endoscopic hemostasis methods

Leung Ki EL, Lau JY

Clin Endosc. 2012 Sep;45(3):224-9 | Epub 2012 Aug 22

87

June 2013 | Report on successful removal of an OTSC® Clip

Prof. Mönkemüller and colleagues presented a clip removal case in a letter to the editor of Gastrointestinal Endoscopy. Ten days after treating an anastomotic leak with the OTSC System, there was still a leak due to misplacement of the clip. The clip had to be removed to place another OTSC onto the leak. At first clip rising was accomplished by injecting saline solution below the OTSC. A snare was positioned around the clip, slowly closed and retracted. The clip dislodged and was retrieved carefully without injury by catching it with the snare and keeping it close to the distal end of the endoscope. The anastomotic leak was thereafter closed successfully with a new OTSC.

Endoscopic removal of an over-the-scope clip ("bear claw")

Mönkemüller K, Toshniwal J, Zabielski M
Gastrointest Endosc. 2012 Nov;76(5):1077-8 107

June 2013 | German surgical periodical alludes to OTSC® Proctology as a novel therapy for anal fistula

In the German surgical periodical "Chirurgische Allgemeine" Prof. Dr. A. Herold, German Center for the Anorectum (EDZ), Mannheim, Germany, gave an overview on new treatments and devices for anorectal fistula. Prof. Herold is the General Secretary of the German Society for Coloproctology (DGK). In his paper he refers to OTSC Proctology as a new therapeutic alternative.

Neue Techniken bei der Therapie der Analfistel

Herold A.
Chirurgische Allgemeine (2013); 14: 99–102 117

May 2013 | Ovesco's Full Thickness Resection Device (FTRD®) presented in live endoscopy at Endo-Update meeting

During clinical live demonstrations at endo-update which took place under the presidency of Prof. Dr. H. Messmann and Prof. Dr. H.-D. Allescher in Augsburg, Germany, a neuro-endocrine tumor (NET) in the rectum was resected with the new Full-Thickness Resection Device of Ovesco Endoscopy: the FTRD. A 62 year old patient showed a submucosal tumor of about 9 mm diameter. Biopsy revealed a neuroendocrine tumor. Prof. Dr Thomas Rösch (University Hospital Hamburg-Eppendorf) used the FTRD to resect the lesion. The FTRD consists of an elongated OTSC cap premounted with a specially designed, derivative OTSC clip and the cap incorporates a resection snare.

Prof. Rösch grasped the lesion with a grasping forceps and pulled the target tissue into the cap in a full thickness fashion. After mobilizing the tissue into the cap, the clip was released to seal the invaginated tissue before resection. Right afterwards the snare was closed and the tissue resected with HF current.

The resection specimen included the full thickness of the wall carrying the NET, with a safety margin. The serosa was seen in histology, confirming that the specimen was a full-thickness resection.

The FTRD device is not yet commercially available.

endo-update 2012 | 30.11.-01.12.2012

Venue: Klinikum Augsburg, Augsburg, Germany

May 2013 | Iatrogenic digestive tract perforations: OTSC® closure as preferred method

Dr. C. Gubler and Prof. P. Bauerfeind, Dept of Gastroenterology, Zurich University Hospital, Switzerland, report about the use of the OTSC clip for endoscopic closure of iatrogenic organ perforations. In a consecutive patient series (n=14) they investigated technically successful closure of perforations that occurred as a result of an endoscopic intervention. All patients were followed clinically for 24 hrs. Endoscopic closure was achieved in 13 of the 14 cases (92.8%). In 3 patients abdominal pain led to evaluation of the closure site by laparoscopy as a precaution. All 3 OTSC closure sites were found intact and no segmental resection of the bowel was needed. One OTSC gastric closure patient had gastric resection after histology revealed gastric adenocarcinoma after endoscopic mucosal resection. The authors conclude that GI perforations up to 30 mm diameter, observed during endoscopic should be treated with endo- scopical OTSC clip closure.

Endoscopic closure of iatrogenic gastrointestinal tract perforations with the over-the-scope clip

Gubler C, Bauerfeind P

Digestion. 2012;85(4):302-7

Epub 2012 May 17

81

OTSC® Update 13

April 2013 | OTSC® System found safe and appropriate for closure of acute perforations in the stomach

In this first trial from China (after compassionate use cases

in patients earlier on) the authors investigated the feasibility of the OTSC System for the closure of gastric perforations in the fundus. This location is of special interest since the handling of a flexible scope in the retroflex position is sometimes quite challenging. The investigation was done in a dog model. The perforation was performed with electrocautery and a needle knife in seven dogs. Closure was performed with one OTSC clip each. The closure was performed in 18.5 +/- 6.4 minutes (team without prior experience). The following leak pressure test with maximum air insufflation and 500 ml methylen blue solution resulted in one minor leak (laparoscopic control) without clinical consequences though. The authors conclude that the OTSC System is safe and appropriate for the closure of acute perforations in the stomach despite the well known difficulties with the J-maneuver.

Feasibility study of secure closure of gastric fundus perforation using over-the-scope clips (OTSC) in a dog model

Zhang XL, Qu JH, Sun G, Tang P, Yang YS

J Gastroenterol Hepatol. 2012 Jul;27(7):1200-4 86

April 2013 | Conference report | OTSC® at German Endoscopy Conference (DGE-BV 2013 in Munich)

OTSC was well-covered in the scientific programme of this year's German Endoscopy Conference in Munich.

Clinical presentations confirm efficacy of OTSC clipping in a range of indications

Munich, March 14–16, 2013. The 43rd German Endoscopy Congress, DGE-BV, was held under the presidency of Prof. Dr. Christoph F. Dietrich.

A significant number of presentations had clinical data of OTSC clipping as their topic and confirmed clinical efficacy and safety in the primary indications of the product, hemostasis, closure of acute lesions/perforations and closure of chronic lesions/fistula (source: www.dge-bv.de).

Large single center OTSC cohort with hemostatic and organ wall closure indications

Wedi E, Menke D, and Hochberger J, Strasbourg (France) reported about a cohort of 84 patients with OTSC clipping for GI bleeding, fistula and GI wall insufficiency. 101 OTSC clips have been used in this cohort, or 1.2 clips per patient. Indications included mainly severe upper GI peptic ulcer hemorrhage (n=38) and preventive clipping to avoid rebleeding (n=12) or secondary perforation (n=18) after large area ESD. The clinical success rate in peptic ulcer bleeding was 79 %, most patients had already been treated unsuccessfully with other hemostatic techniques before OTSC clipping or had been candidates for surgical hemostasis.

2 complications were encountered: 1 inadvertent clipping of an instrument with OTSC and fixation of the instrument to the tissue and one perforation of the sigmoid with the OTSC cap. The authors state that OTSC application is an effective procedure to deal with endoscopic situations that otherwise would require a surgical approach.

Der Einsatz des OTSC-Makroclips bei 84 Patienten mit schwerer GI-Blutung, Fisteln und Insuffizienzen – ein Résumé

E. Wedi, D. Menke, and J. Hochberger, Strasbourg

Large single center cohort on OTSC hemostasis in severe GI bleeding

Kratt T, Stüker D, Gräpler F, Küper M, Wichmann D, Königsrainer A, Tübingen, showed data from their cohort on OTSC in endoscopic hemostasis (n=85). The bleeding location was in the upper GI tract in 63 % (21 % peptic gastric ulcers and 40 % peptic duodenal ulcers) and in the lower GI tract in 37 % (mostly bleeding after polypectomy in the rectum).

The characteristics of the cohort underline the severity of bleeding: life-threatening bleeding (28.4 %), patient in hemorrhagic shock (31.1 %), immediate blood transfusion (33.8 %), patient under anti-coagulation (21.6 %), Forrest I bleeding (72.3 %). OTSC placement was achieved with cap suction in 72 cases and with an OTSC Anchor in 2 cases. Technically successful hemostasis for 72 hrs was achieved in 92.8 % of cases, a persistent bleeding and an early relapse bleeding (<72 hrs) were seen in 3.6 %, respectively. Late relapse bleeding (>72 hrs) was observed in 3.6 %. No severe complications were observed; in 3 cases mucosal

esophageal lesions from device introduction were seen. In 14.5 % OTSC clipping was done for recurrence of an initially successful other endoscopic therapy and in 13.3 % for failure of other methods in the same treatment session. In 35.1 % OTSC clipping was seen as an ultima ratio and as an alternative to surgical therapy otherwise becoming necessary. The summary of the authors is that the simple and easy to handle OTSC System is an effective treatment in severe GI bleeding and can avoid surgery in several cases.

Das OTSC-Clip-System: Klinische Erfahrungen zur Therapie der schweren GI-Blutung bei 85 Patienten

T. Kratt et al., Tübingen

OTSC to prevent migration of covered self-expanding stents

Fähndrich M, Pohl T, Rolfs S, Sandmann M, and Heike M, Dortmund, presented their technique of using OTSC to avoid migration of covered, self-expandable stents.

Stent migration has an incidence of up to 30 % and represents a significant clinical challenge. To prevent stent migration, the authors used OTSC to fix the stent permanently to the neighboring GI wall. In 24 cases with benign indication for stent placement OTSC fixation was carried out in the following locations: esophagus, small bowel and colon. After 5–8 weeks the OTSC clips were removed by Nd:YAG laser cutting to intentionally remove the stent. In all 24 patients the procedure was technically successful. In 1 patient an undesired stent migration before intentional removal was observed. In another case the stent had to be removed after a few days due to intolerance by the patient in a location close to the upper esophageal sphincter. The authors conclude that OTSC clipping was found to be a safe and practical technique and has prevented stent migration in 96 % of the cases studied.

Verwendung des Ovesco-Clips zur Verhinderung der Migration bei vollgecoverten selbstexpandierenden Stents

M. Fähndrich, T. Pohl, S. Rolfs, M. Sandmann, and M. Heike

Hospitalisation time and 30-days mortality in GI perforations after technically successful and unsuccessful OTSC closure

Hagel A, Nägel A, Raithel S, Diebel H, Neurath M, and Raithel M, Erlangen, showed data on the management of GI perforations with OTSC clips. They studied 19 patients with apparent perforation of a digestive organ wall in various anatomical locations. In 13 patients the perforation could be closed with OTSC ("O+") to avoid emergency surgery. In 6 patients OTSC closure was technically unsuccessful and emergency surgery was needed ("O-").

In the O+ group the duration of hospitalisation was 10.7 +/- 10 days, no mortality. 2 patients in this group had comorbidities unrelated to clip closure, leading to a prolonged hospital stay; excluding these 2 patients, hospitalisation was 5.8 +/- 2 days. In the O-group hospital stay was 12.1 +/- 7 days, one patient with esophageal perforation died after emergency surgery was not able to prevent fatal mediastinitis. The authors draw the conclusion that OTSC treatment can significantly reduce morbidity and mortality in GI perforations.

OTSC-Anwendung bei manifester GI-Perforation: 30-Tages-Mortalität, Hospitalisationsdauer und Out-come nach endoskopisch erfolgreichem und nicht-erfolgreichem Perforationsverschluss

A. Hagel, A. Nägel, S. Raithel, H. Diebel, M. Neurath, and M. Raithel, Erlangen

Monocentric case experience with OTSC in a broad range of wall closure indication: safe transmural closure

Nietsch H, Hammelmann F, and Asperger W, Halle, summarized their initial experience with OTSC in endoscopic closure of the GI organ wall in 10 consecutive applications. Indications included: postsurgical rectal anastomotic leak (n=2), rectal ESD perforation (n=1), gastric ESD perforation (n=2), esophageal perforation after balloon dilation (n=1), Mallory-Weiss tear (n=1), perforated gastric ulcer (n=1), post-surgical duodenal leak (n=1) and post-surgical bariatric suture line leak. All cases were successful. The authors conclude: OTSC enables a safe transmural closure of spontaneous and iatrogenic

perforations. In a majority of cases target tissue handling is possible with suction only and does not require additional instruments. In well-trained endoscopy centers the learning curve for OTSC is short.

Erfahrungsbericht der ersten 10 Anwendungen des endoskopischen OTSC-Clipsystems

H. Nietsch, F. Hammelmann, and W. Asperger, Halle

OTSC for closure of distal esophageal perforation

Braun A, Richter-Schrag H, Hopt U, Fischer A, Freiburg, showed data on OTSC in the treatment of distal esophageal perforation after vomiting (Boerhaave, n=1) and iatrogenic injury (n=1). Esophageal perforation is a life-threatening situation with a high complication and mortality rate. In both cases endoscopic closure of the esophagus was achieved within 12 hrs after the lesion. Both patients received bilateral thorax drainage and antibiotic therapy. No patient developed sepsis. Starting oral intake was without problems. Control endoscopy after 3 months revealed no stenosis and both clips were found in place. The authors summarize that the closure of esophageal perforations with OTSC is a safe and effective method and is significantly more economic than common surgical therapy requiring longer hospital stays.

Endoskopischer Verschluss von distalen Ösophagus-Perforationen mit einem Over-The-Scope Clip (OTSC)

A. Braun, H. Richter-Schrag, U. Hopt, A. Fischer, Freiburg

Consecutive case series of OTSC application in the endoscopic management of complications and emergencies

Thomsen T, Berthold B, Khiabanchain M, and Trabandt I, Neubrandenburg, presented data of a case series (n=11). Indications included upper and lower GI bleeding, PEG-fistula closure, rectal-pelvic fistula closure, sigmoid anastomosis leak, bleeding from diverticulum (Hartmann situation), arterial bleeding from colon anastomosis. The overall clinical success rate in the mixed case series was 82 %. No procedure took more than 30 min. As complications 1 fistula recurrence (required second OTSC procedure), 1 rebleeding and 1 remaining perforation were seen. The authors summarize that OTSC clipping is a fast procedure with a high primary success rate and is quick to learn.

Endoskopische Interventionen mit dem OTSC-System am Klinikum Neubrandenburg

T. Thomsen, B. Berthold, M. Khiabanchain, and I. Trabandt, Neubrandenburg

OTSC for stopping gastroduodenal artery bleeding in duodenal ulcer

Kratt T, Stüker D, Kirschniak A, Heininger A, Wietek B, Königsrainer A, Tübingen, showed a case series (n=7) in which OTSC was applied in upper GI emergency hemostasis to stop bleeding from the gastroduodenal artery. Gastroduodenal artery bleeding is besides aortoduodenal fistula considered the most severe bleeding complication in the digestive tract, associated with high morbidity and mortality. In many cases surgical emergency hemostasis is inevitable.

In all cases reported here the gastroduodenal artery was verified as the bleeding source by angiography after successful endoscopic treatment. In all 7 patients the acute bleeding from an ulcer at the posterior duodenal wall was successfully controlled with OTSC, in 4 cases fibrin glue was additionally applied. After the initial 72 hrs, 3 patients suffered from rebleeding, which was then controlled surgically. No mortality was encountered in this case series. The authors draw the conclusion that OTSC is effective in emergency management of gastroduodenal artery bleeding. In more than half of the cases endoscopic management was the only therapy. In the other patients OTSC was a successful "bridge to surgery" and allowed stabilizing the patient before the operation.

OTSC-basierte Notfall-Hämostase der lebensbedrohlichen A. gastroduodenalis Ulkus-Arrosionsblutung: alleinige endoskopische Therapie oder „bridge-to-surgery“

T. Kratt, D. Stüker, A. Kirschniak, A. Heininger, B. Wietek, A. Königsrainer, Tübingen

Report on Ovesco FTRD (pre-commercial device)

Kratt T, Stüker D, Gräpler F, Schnek M, Adam P, and Königsrainer A, Tübingen, presented data of their first 8

cases with FTRD, a device of Ovesco Endoscopy, not yet commercially available. It combines modified OTSC clipping with tissue resection. In 7 of the 8 cases the procedure was technically feasible; in 1 case the target lesion could not be reached. The cases treated included various indications in which FTRD was used as a device for full-thickness tissue retrieval with the primary purpose of enhanced histological examination of an in-toto full thickness specimen. The target lesions were in the upper GI tract, melanoma metastasis (n=3), GIST (n=1) or in the lower GI tract (adenoma or early colorectal cancer, low risk histology; n=4). The presentation gave a detailed case history of an elderly patient with recurrent adenoma (high grade dysplasia, partially adenocarcinoma) of the rectum. The patient had full thickness resection with FTRD under single-shot antibiotics and was discharged the following day. As histology demonstrated complete removal of the lesion, no further therapy was done. Follow-up was uncomplicated. After 14 weeks control endoscopy revealed that the clip had detached from the tissue, normal scar formation was seen at the resection site and no signs of residual lesion or new recurrence were found.

Klinische Evaluation eines neuen endoskopischen GI-Trakt-Vollwandresektionssystems: das OTSC-basierte „full thickness resection device“ (FTRD)

T. Kratt, D. Stüker, F. Gräpler, M. Schnek, P. Adam, and A. Königsrainer, Tübingen

FTRD is not yet commercially available.

March 2013 | Dr. Thomas Kratt, University of Tuebingen, Germany, wins award for clinical research with Ovesco's FTRD®

Dr. Thomas Kratt, Interdisciplinary Endoscopy, University Hospital, Tuebingen, Germany, received an award for this presentation of clinical research in the field of full-thickness resection at the 43rd Congress of the German Society for Endoscopy and Imaging (DGE-BV), held in Munich, March 14–16, 2013.

Dr. Kratt presented data of his first 8 cases with FTRD, a device of Ovesco Endoscopy, not yet commercially available. It combines modified OTSC clipping with tissue resection. In 7 of the 8 cases the procedure was technically feasible; in 1 case the target lesion could not be reached.

The cases treated included various indications in which FTRD was used as a device for full-thickness tissue retrieval with the primary purpose of enhanced histological examination of an in toto full-thickness specimen. The target lesions were in the upper GI tract, melanoma metastasis (n=3), GIST (n=1) or in the lower GI tract (adenoma or early colorectal cancer, low risk histology; n=4).

The presentation of Dr. Kratt gave a detailed case history of an elderly patient with recurrent adenoma (high grade dysplasia, partially adenocarcinoma) of the rectum. The patient had full-thickness resection with FTRD under single-shot antibiotics and was discharged the following day. As histology demonstrated complete removal of the lesion, no further therapy was done. Follow-up was uncomplicated. After 14 weeks control endoscopy revealed that the clip had detached from the tissue, normal scar formation was seen at the resection site and no signs of residual lesion or new recurrence were found.

www.dge-bv.de/german/home.php

March 2013 | Prospective trial on OTSC® Proctology in anal fistula treatment presents first data

Munich, March 8, 2013. The annual conference of the German Society for Coloproctology (DGK) was held in Munich, March 8 and 9, 2013. At the conference first data were presented from an investigator initiated multicentric prospective observational clinical trial on the use of OTSC Proctology in the treatment of anal fistula. The two participating trial sites are the Stuttgart Institute of Proctology (PD. Dr. R. Probst, Dr. W. Ehni), Stuttgart and the German Anorectal Center (EDZ) (Dr. A. Joos, Prof. Dr. A. Herold, PD Dr. D. Bussen), Mannheim.

The trial presented an interim analysis on the first 15 patients. Inclusion criteria are supra-, extra- or high trans-sphincteric anal fistula, including first recurrence but excluding patients with IBD.

Mean follow-up was 6.9 months (1–15 months) after OTSC placement. 8 patients had already completed follow-up (6 months), 7 patients were still followed. In patients who had already completed the trial, mean follow-up was 10.8 months (6–15 months).

In these patients the healing rate, defined as post-surgical closure of the fistula, absence of drainage from the fistula and absence of recurrence after 6 months was 88 %.

In his presentation PD Dr. R. Probst, Stuttgart, coordinator of the trial, summarized that data were encouraging but completion of the trial had to be awaited. The trial is expected to close in 2013.

www.mcn-nuernberg.de/DGK2013/programm-08032013.php

March 2013 | EndoResect study – Endoscopic full-thickness resection of gastric subepithelial tumors

Meining et al. report of 20 patients with gastric subepithelial tumors (SET) up to 3 cm in diameter. Patients were prospectively enrolled and 14 of them treated by endoscopic resection using the OTSC Anchor and a monofilament snare. In cases where perforation occurred the defect was closed with Twin Grasper and OTSC System. The authors conclude that this method seems to be faster and easier than other endoscopic techniques such as ESD or submucosal tunneling.

Perforations could be adequately managed by the OTSC System (100 % closure). Thus, endoscopic resection without laparoscopic control seems possible in selected patients with purely intraluminal tumors. The authors discuss the malignant potential of SETs, especially GISTs which cannot be reliably determined by either endoscopic or endo-sonographic surveillance. According to guidelines GISTs larger than 2 cm should be resected. However, since also smaller tumors have malignant potential complete resection of all suspected lesions seems advisable according to the authors. They argue that GISTs rarely develop lymph node metastases, and thus local resection with large negative margins and without lymph node resection are considered curative approaches.

Comment by Ovesco: since only tumors without connection to the muscularis propria layer have a 80–100 % resection rate in literature, it might be feasible to perform full-wall resections in SETs and similar tumors. Ovesco is currently completing the development of a new Full Thickness Resection Device (FTRD) for the lower GI tract to start with.

EndoResect study – Endoscopic full-thickness resection of gastric subepithelial tumors

Schlag C, Wilhelm D, von Delius S, Feussner H, Meining A Endoscopy. 2013 Jan;45(1):4-11. Epub 2012 Dec 19

February 2013 | OTSC® Proctology – description of operative technique in MITAT

PD Dr. R. Probst and Dr. W. Ehni, Stuttgart Proctology Institute, Stuttgart, Germany, pioneers in the application of Ovesco's OTSC Proctology system, recently described their preferred technique for anal fistula closure with the device. The procedure consists of 3 steps: local removal of the anoderm around the inner orifice of the fistula, debridement of the fistula tract and clip closure of the fistula.

They also present an indicative case study of a 54-year old female patient suffering from a high transsphincteric anal fistula and recurrence after unsuccessful prior surgery. After transanal clip release from the OTSC Proctology applicator, the internal fistula opening was adequately closed by the clip. Eight months after clip closure the fistula had permanently healed. The authors conclude fistula closure using the OTSC Proctology system represents a promising sphincter-preserving minimally invasive procedure.

The OTSC Proctology clip system for anorectal fistula closure: the „Anal Fistula Claw“: Case report

Probst RL, Ehni W

Minim Invasive Ther Allied Technol. 2012 Jul;21(4):307-12 [Epub 2012 Jun 4]

February 2013 | Retrospective multicentric

review of early OTSC® patients in the US: overall clinical success rate of 71 %

Dr. Todd H. Baron and colleagues, Division of Gastroenterology & Hepatology, Mayo Clinic, Rochester MN, USA report about their experience with 45 patients and 48 OTSC clip placements from March 2011 to January 2012. Median follow-up time in this mixed cohort was 77 days (30–330 days). Indication break-down included hemostasis (n=7), closure of chronic fistula (n=28), closure of iatrogenic perforations (n=5), closure of post-esophagectomy anastomotic leakage (n=3) and miscellaneous (n=2).

Before OTSC placement 49 % of the patients had undergone other therapies for their condition that had failed. The overall clinical success rate was 71 %. Hemostasis was achieved in 100 % of cases. Anastomotic leakage and fistula was closed in 65 %. Also one case of OTSC clip removal by means of APC-cutting of a clip hinge is described.

The authors conclude that the OTSC clip appears clinically effective and is a welcome addition to the therapeutic armamentarium in the closure of leaks, fistula, perforations and non-variceal bleeding.

Use of an over-the-scope clipping device: multicenter retrospective results of the first U.S. experience

Baron TH, Song LM, Ross A, Tokar JL, Irani S, Kozarek RA *Gastrointest Endosc.* 2012 Jul;76(1):202-8

January 2013 | Combined use of OTSC® System and stent to close large EMR-related perforations

Treatment of large EMR-caused perforations with a combined use of OTSC and stenting is reported by Hadj Amor et al.

One patient with a 20-mm esophageal perforation was treated with an OTSC, several other clips and an endoloop. A fully covered stent was placed on top to bypass the perforation. The large duodenal perforation in the other patient was initially unsuccessfully treated with a fully covered stent and several clips to avoid migration. After removal of the stent an OTSC and two other clips were used to close the perforation completely. The perforation was bridged by another fully covered stent that was placed over the closed perforation without fixation. In both patients the stents were removed after several weeks and both sites showed healing of the perforation.

Successful endoscopic management of large upper gastrointestinal perforations following EMR using over-the-scope clipping combined with stenting

Hadj Amor WB, Bonin EA, Vitton V, Desjeux A, Grimaud JC, Barthel M

Endoscopy. 2012;44 Suppl 2 UCTN:E277-8 [Epub 2012 Aug 29]

OTSC® Update 12

January 2013 | OTSC® used to prevent stent migration in the treatment of anastomotic leak

Toshniwal J et al. report about the use of the OTSC System to anchor a fully covered self-expandable metal stent to prevent stent migration. The patient underwent distal esophagectomy with gastric pull-up. The stent was placed to a post-operative anastomotic leak in the esophagus. However, the stent partially migrated into the stomach. The stent was then repositioned onto the leak. The OTSC System was placed using the OTSC Twin Grasper to grasp the stent edge and suction. After application the OTSC clip fixed the stent to the esophageal wall. Follow-up showed successful closure of the anastomotic leak.

The authors conclude that the placement of the clip was easy, fast and prevented stent migration effectively. The fixation of stents is not a common indication for the OTSC System and there is only very limited experience.

Combination of the “bear claw” (over-the-scope clip system) and fully covered stent for the treatment of post-operative anastomotic leak

Toshniwal J, Zabielski M, Fry LC, Mönkemüller K *Endoscopy.* 2012;44 Suppl 2 UCTN:E288-9

December 2012 | Closure of anastomotic leaks and chronic fistulas in the digestive tract: best

results in earlier treatment cases

Dr. Selcuk Dişibeyaz and co-authors, Department of Gastroenterology of Türkiye İhtisas Hospital, Ankara, report about their case series of 9 patients (age 22–65 years). Anastomotic leakage from GI surgical anastomosis was present in 5, fistula in 3 and acute perforation in 1 patient. Type “a” clips were placed in all cases. In 4 cases clip deployment was not undertaken, due to strong tissue fibrosis. In the other 5 patients the clip was successfully deployed and closed the defect without the need of further treatment. The median time between diagnosis of the defect and OTSC clip placement was 35 (20–80 days) days in the cases with successful placement and 70 days (38–94 days) in the unsuccessful cases. The median defect size was 15 mm (5–20 mm). In 4 cases clip deployment was not undertaken, due to strong tissue fibrosis. No clip-related complications were encountered.

Endoscopic closure of gastrointestinal defects with an over-the-scope clip device. A case series and review of the literature

Dişibeyaz S, Köksal AŞ, Parlak E, Torun S, Sağmaz N *Clin Res Hepatol Gastroenterol.* 2012 Dec;36(6):614-21 [Epub 2012 Jun 14]

December 2012 | OTSC® effective in closure of chronic esophago-jejunal anastomotic leaks after total gastrectomy

Prof. Dr. Gennaro Galizia and co-workers from the Second University of Naples, Italy, recently described the application of OTSC clips in the treatment of postsurgical anastomotic failure after total gastrectomy.

In a case series of 3, patients that developed anastomotic leaks after gastrectomy and Roux-en-Y jejunal transposition and esophago-jejunostomy were endoscopically treated with OTSC clipping. The case series was published in the *Journal of Gastrointestinal Surgery*.

In all patients clip closure of the leak was technically simple, clinically effective and did not result in complications. The authors conclude that the OTSC System may represent a new option in the management of postoperative esophago-jejunal leaks. The incidence of anastomotic leaks ranges from 4 to 27 % after total gastrectomy and is a not infrequent challenge in such patients.

The Over-The-Scope Clip (OTSC) System is effective in the treatment of chronic esophago-jejunal anastomotic leakage

Galizia G, Napolitano V, Castellano P, Pinto M, Zamboli A, Schettino P, Orditura M, De Vita F, Auricchio A, Mabilia A, Pezzullo A, Lieto E *J Gastrointest Surg.* 2012 Aug;16(8):1585-9 [Epub 2012 Mar 7]

December 2012 | The interesting case: ERCP-related jejunal perforation managed by OTSC® clipping

In a recent issue of *Gastrointestinal Endoscopy* Dr. F. Buffoli and colleagues, Digestive Endoscopy and Gastroenterology Unit of the Hospital Institutes, Cremona, Italy, presented an interesting case report on OTSC clipping for jejunal perforation closure:

An 85 y/o woman with bile obstruction due to pancreatic cancer presented with jaundice. The patient had Billroth II anatomy from gastric resection due to a peptic ulcer 35 years in the past.

Cholangiography showed a bile duct stricture. An endoscope-related perforation of the jejunum with a size of about 20 mm was visualized distally of the papilla. After placing a covered self-expanding stent through the biliary duct stricture it was decided to close the perforation of the bowel with an OTSC clip. The patient was considered inoperable due to age and comorbidities.

Closure of the jejunum was successful. Retroperitoneal fatty tissue was additionally pulled into the cap by suction and created a “retroperitoneal fat patch”. Abdominal CT revealed retroperitoneal air but no free liquids. The patient received parenteral nutrition and antibiotic treatment. The patient remained symptom-free and the jaundice disappeared. Control CT after 20 days demonstrated complete absorption of the air and the patient was

discharged.

ERCP has a perforation rate of approx. 0.3 to 1.3 %, as described in the clinical literature. The authors conclude that for the endoscopic closure of large ERCP-related perforations OTSC may be considered as a possible treatment.

Endoscopic “retroperitoneal fatpexy” of a large ERCP-related jejunal perforation by using a new over-the-scope clip device in Billroth II anatomy (with video)

Buffoli F, Grassia R, Iiritano E, Bianchi G, Dizioli P, Staiano T *Gastrointest Endosc.* 2012 May;75(5):1115-7 [Epub 2011 Aug 5]

November 2012 | First publication of Japanese experience with OTSC®

In the recent issue of the *World Journal of Gastroenterology* Dr. Hirohito Mori published first Japanese experiences with the OTSC System.

Two elderly patients who had suffered iatrogenic lesions in the rectum (one large rectal perforation with abscess formation and one recto-vesical fistula). Both patients were not subject to a surgical intervention for poor general condition, and thus were successfully treated with one OTSC clip each. Both interventions resulted in a dramatic improvement of the patients' status. It should be noted that both patients underwent direct endoscopic lavage before closure. This is noteworthy especially in the case with the abscess where no pararectal drainage was inserted.

The authors state: “The endoscopic closure of perforations and fistulae with OTSC is a simple and minimally invasive technique. Given the complete closure and healing of large fistulae with OTSC in our two cases, this approach may be less expensive and more advantageous than surgical closure.”

Rectal perforations and fistulae secondary to a glycerin enema: Closure by over-the-scope clip

Mori H, Kobara H, Fujiwara S, Nishiyama N, Kobayashi M, Masaki T, Izuishi K, Suzuki Y *World J Gastroenterol.* 2012 Jun 28;18(24):3177-80

November 2012 | Efficacy of OTSC® for the treatment of colorectal postsurgical leaks and fistulas: 86 % overall success rate

Anastomotic leaks and fistulas are a severe complication in colorectal surgery. The incidence of clinically relevant leaks is in the range of 3–6 % of cases.

Prof. Dr. Alberto Arezzo and colleagues, Dept of Digestive, Colorectal and Minimal Invasive Surgery, University of Turin, Italy report about a prospective case series covering 14 consecutive patients, treated between April 2008 and September 2011. Criteria for treatment with OTSC were a wall opening of <15 mm with no extraluminal abscess and absence of stenosis. The mean defect size treated was 9.1 mm in diameter. One OTSC clip of either size 11 or 12 was sufficient in all defects. In one case two separate defects were treated in the same patient. In 8 cases the leak was a fresh, acute lesion, in 6 cases a chronic fistula.

The overall success rate of durable defect closure in this prospective case series was 86 %; for acute cases it was 87 % and for chronic cases 83 %. No OTSC-related complications were reported. Re-surgery was needed in 1 case, in a second failure case the patient refused re-surgery and was left untreated.

The authors conclude that endoscopic closure of colorectal postsurgical leaks are a safe technique with a high success rate, including rectovaginal and colocolic fistula.

Efficacy of the over-the-scope clip (OTSC) for treatment of colorectal postsurgical leaks and fistulas

Arezzo A, Verra M, Reddavid R, Cravero F, Bonino MA, Morino M *Surg Endosc.* 2012 Nov;26(11):3330-3

November 2012 | Performance of the OTSC® System in the endoscopic closure of gastrointestinal fistulae – a meta-analysis

The recent issue of “Minimally Invasive Therapy & Allied Technologies” publishes a systematic review and meta-analysis on the challenging field of closing gastrointestinal fistulae by means of the OTSC System. The paper provides an extensive overview of relevant primary clinical research,

case reports and conference abstracts published on this topic. The statistical evaluation of, in total, 19 examined articles revealed a high rate of procedural success (mean 84.6 %; 95 % confidence interval 66.6 % to 93.8 %) and durable clinical success (mean 69.0 %; 95 % confidence interval 51.8 % to 82.2 %) in OTSC-mediated closing of GI fistulae.

In summary, the authors rate endoscopic closure of gastrointestinal fistulae by means of the OTSC System as a safe and effective method.

Performance of the OTSC System in the endoscopic closure of gastrointestinal fistulae – a meta-analysis

Weiland T, Fehker M, Gottwald T, Schurr MO

Minim Invasive Ther Allied Technol. 2012; 21(4):249-58

October 2012 | The success rates for hemostasis in severe GI bleeding, perforation closure and chronic fistula closure are 88 %, 79 % and 73 %, respectively

The OTSC System has been described in more than 40 clinical papers in the scientific literature covering a range of indications. In order to summarize the clinical data published so far and to evaluate the overall clinical efficacy, Ovesco Endoscopy has commissioned systematic literature research on the OTSC System.

The study was limited to clinical publications and covered the key applications of the OTSC System, hemostasis, closure of acute GI lesions (perforations) and chronic GI lesions (fistula). Only clinical reports with >4 patients were included into the survey, that was carried out by Dr. Timo Weiland, novineon CRO, a specialized contract research organization for the medical device industry (www.novineon.com).

The success rates defined as permanent achievement of the therapeutic goal for hemostasis in severe GI bleeding, perforation closure (including acute anastomotic suture line failure) and chronic fistula closure are 88 %, 79 % and 73 %, respectively. The OTSC System compares to the effectiveness of a surgical intervention in the respective indications or offers a new therapeutic option in situations where surgery is not feasible.

http://www.ovesco.com/fileadmin/Downloads/OTSC_System_clinical_data_eng_Rev01_2012-10-22.pdf
(English)

http://www.ovesco.com/fileadmin/Downloads/OTSC_System_KlinischeDaten_deu_Rev01_2012-10-22.pdf
(German)

October 2012 | Hemostasis in large gastric ulcer with the OTSC® System

Vormbrock et al. report a successful treatment of gastric ulcer bleeding with the OTSC System.

In an emergency EGD removal of clots and fresh blood revealed an ulcer with a 2-mm thick pulsating vessel. Injection therapy was difficult due to the fibrotic tissue. Thus OTSC placement was decided. To mobilize the target tissue into the cap, two edges of the ulcer were grasped by each of the two jaws of the OTSC Twin Grasper. After retraction of the grasper and additional suction the OTSC was applied and immediate hemostasis achieved.

The authors conclude that the OTSC was effective for hemostasis in this fibrotic ulcer which was very hard to treat with other endoscopic methods. They state that the placement of OTSC was quick and easy resulting in potentially life-saving hemostasis.

Use of the “bear claw” (over-the-scope clip) to achieve hemostasis of a large gastric ulcer with bleeding visible vessel

Vormbrock K, Zabielski M, Mönkemüller K

Gastrointest Endosc. 2012 Oct;76(4):917-8

October 2012 | Postsurgical colorectal anastomotic leaks: OTSC® clip recommended as treatment of choice at SMIT conference

Barcelona, September 21st 2012: The 24th conference of the Society for Minimally Invasive Therapy (SMIT) was held in Barcelona, Spain, under the presidency of Dr. Enric Laporte.

Prof. Dr. Alberto Arezzo and colleagues, 2nd Dept of General

Surgery, University of Turin, Italy, presented latest data of 25 clinical cases with postsurgical anastomotic leaks or fistula after colorectal surgery.

In the general literature anastomotic leaks have an incidence of about 7–9 % after laparoscopic or open colorectal surgery. In the 25 cases prospectively collected in Turin, 21 were successfully treated with endoscopic OTSC clipping alone. This is a success rate of 84 %. In 3 patients the fistula did not heal, and in 1 patient additional surgery was needed to close the defect.

In conclusion the authors recommend the use of endoscopic OTSC clipping for lesions up to 12 mm in size as the primary treatment for patients with postsurgical leaks and fistula after colorectal surgery.

Efficacy of the over-the-scope clip (OTSC®) for treatment of colorectal postsurgical leaks and fistula

ArezzoA, Verra M, Reddavid R, Cravero F, Bonino MA, Morino M