



## FBR Set.....

for safe foreign body retrieval

- Alligator and rat tooth in one grasping device
- Long and elastic retrieval cap
- Safe removal of sharp and foreign objects





## ..... FBR Set Ref. no. 400.05


The FBR Set is a foreign body retrieval set for safe removal of various foreign body objects.

The package contains two sets (ref. no. 400.05); each set includes a remOVE SecureCap 12 and a remOVE Grasper.


The remOVE SecureCap is an elastic and transparent retrieval cap. It allows safe extraction of ingested foreign objects in the digestive tract.

The remOVE Grasper is a grasping forceps used to grasp and retrieve foreign bodies in the digestive tract. It includes alligator and rat teeth to ensure a strong hold and easy retrieval of foreign bodies into the remOVE SecureCap.

### Dimensions




remOVE SecureCap 12



remOVE Grasper

## ..... Application

Extraction of a clip fragment by using the FBR Set.



Source: Prof. K. Caca, Klinikum Ludwigsburg, Germany

### Ovesco Endoscopy AG

Friedrich-Miescher-Strasse 9  
72076 Tuebingen  
Germany

Phone +49 (0) 7071 96528-160

Fax +49 (0) 7071 96528-260

E-Mail [service@ovesco.com](mailto:service@ovesco.com)

### Ovesco Endoscopy USA Inc.

15300 Weston Parkway  
Suite 101  
Cary, NC 27513,  
USA

Phone +1 919 651 9449

Fax +1 408 608 2077

[customerservice@ovesco-usa.com](mailto:customerservice@ovesco-usa.com)

Scan here for more  
product information



**Ovesco Endoscopy AG** is a medical device company specialising in the fields of flexible endoscopy, endoluminal surgery and telemetric diagnostics. The company develops, manufactures and markets innovative products for the detection and treatment of gastrointestinal disease. Ovesco products are synonymous with efficacy and efficient application.

[www.ovesco.com](http://www.ovesco.com)

CE<sup>0297</sup>

English version 8 | 2024-10-15